DIUIA

MANUEL D'UTILISATION USER MANUAL

 **BMB**
MEDICAL

FRANÇAIS

BMB MEDICAL

100 Avenue du Formans
CS 11001
01604 TRÉVOUX CEDEX
Tél : +33(0)4 74 08 71 71
Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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www.bmb-medical.com

This manual is for the purpose of familiarization with the use, maintenance and safe handling of your new medical device. We strongly advise you to carefully read this document and to follow the recommendations for safe, optimal use.

BMB MEDICAL offers a wide variety of options for both models featured in this user manual. Consequently, some features described in this manual may not apply to your couch.

Some illustrations are included for reference only and do not accurately represent the device you own.

Stick the
product reference
label here

The information in this document may be changed without notice, for which BMB MEDICAL can in no way be held responsible. Any representation or reproduction, whether in whole or in part, without permission from BMB MEDICAL, is illegal. Such illegal representation or reproduction, by any process whatsoever, shall constitute an infringement. The illustrations contained in this document are not contractually binding.

CONTENTS

ENGLISH	22
1 INTRODUCTION SPECIFICATIONS	25
1.1 PRODUCT PRESENTATION	25
1.1.1 DIVIA HEIGHT ADJUSTABLE EXAMINATION COUCH.....	25
1.1.2 DIVIA MIXTE HEIGHT ADJUSTABLE EXAMINATION COUCH.....	26
1.2 AREA OF APPLICATION	27
1.3 PURPOSE OF THE USER MANUAL.....	27
1.4 SAFEKEEPING.....	28
1.5 SYMBOLS AND STANDARDS USED	28
2 OPERATING INSTRUCTIONS.....	29
2.1 BATTERY	29
2.2 ADJUSTABLE HEIGHT AND LOWERING	29
2.3 BACKREST	30
2.4 INTERMEDIATE BACKREST (DIVIA MIXTE).....	30
2.5 TRENDelenburg POSITION.....	31
2.6 INDEPENDANT WHEELS BRAKES.....	31
2.7 OPTIONS	32
2.7.1 HEAD REST.....	32
2.7.2 CENTRALLY CONTROLLED BRAKE AND STEERING PEDAL	32
2.7.3 SIDERAiLS.....	32
2.7.4 CLAMPS FOR FASTENING STIRRUPS, CRADLES, LEG SUPPORTS OR EXAMINATION LAMP	33
2.7.5 LATERAL TECHNICAL RAILS	34
2.7.6 RECOVERY TANKS.....	35
3 CLEANING AND DISINFECTION	36
3.1 INTRODUCTION	36
3.2 FREQUENCY	36
3.3 RECOMMENDED MATERIAL AND PRODUCTS	36
3.4 DIRECTIONS FOR CLEANING AND DISINFECTION	37
3.4.1 CLEANING OPERATION BASED ON THE COMPONENTS	37
3.4.2 DIRECTIONS TO FOLLOW	37
4 PRECAUTIONS.....	38
4.1 PRECAUTIONS FOR THIS DEVICE*	38
4.1.1 PRECAUTIONS FOR USE OF THE COUCH	38
4.1.2 PRECAUTION FOR SERVICING AND MAINTENANCE OF THE EQUIPMENT.....	38
4.2 NON STANDARD USES	39

4.3	TRANSPORT AND STORAGE	40
4.3.1	TRANSPORT	40
4.3.2	STORAGE	40
4.4	WARRANTY	41
4.5	AFTER SALES SERVICE	41
4.6	ENVIRONNEMENT.....	41
4.6.1	USED EQUIPMENT.....	41

1 INTRODUCTION SPECIFICATIONS

1.1 PRODUCT PRESENTATION

1.1.1 DIVIA HEIGHT ADJUSTABLE EXAMINATION COUCH

TECHNICAL CHARACTERISTICS	
REFERENCE	DIVIA
DESIGNATION	HEIGHT ADJUSTABLE EXAMINATION COUCH
SAFE WORKING LOAD	250 kg (maximum allowable load including the total weight of the patient, mattress and accessories)
WEIGHT	100 kg
OVERALL DIMENSIONS	65 x 200 cm or 75 x 200 cm.
TRENDELENBURG / REVERSE TRENDelenburg	Electric (+/- 15°).
BACKREST	0 to 80°
EQUIPMENTS	2-sectional mattress base : 65x 200 cm or 75 x 200 cm Flat upholstery, high density foam 35 kg/m3. Width 6,5 cm Jersey with PVC coating, antibacterial and anti-fire (M1) treatment Adjustable headrest (0 - 90°) with 2 gas springs Electric height adjustment : 46 - 86 cm Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°) Steel structure, epoxy coated (light grey) Frame with 2 electric columns providing great stability Frame protection cover Remote control 4 lockable wheels Ø 125 mm Examination towel dispenser head end
AVAILABLE OPTIONS	Antistatic system Casters Ø 125 mm with central brake Double running castors Ø 150 mm with central brake Electric adjustable backrest Wireless remote control (in addition of the wired remote control) Foot pedal control (in addition of the wired remote control) Drop-down side rails (pair) Holder for clamps (pair) Fixation clamps for heels, armrest for blood sampling or leg support (pair) (DVSE required) Heels (pair) (DVETU required) Anatomic shaped armrest (DVETU necessary) Leg support (pair) (DVETU required) Protection cover for leg rest (pair) Clamps for IV pole holder (pair) Removable 2-stage IV pole, 2 hooks (DVTSR required) Comfort pillow (grey) Crystal protection cover for examination couch width 65cm Crystal protection cover for examination couch width 75cm Battery Examination lamp (DVETU required) Equipotential connection DIN rails

1.1.2 DIVIA MIXTE HEIGHT ADJUSTABLE EXAMINATION COUCH

TECHNICAL CHARACTERISTICS	
REFERENCE	DIVIA MIXTE
DESIGNATION	HEIGHT ADJUSTABLE COMBINED EXAMINATION COUCH
SAFE WORKING LOAD	250 kg (maximum allowable load including the total weight of the patient, mattress and accessories)
WEIGHT	90 kg
OVERALL DIMENSIONS	65 x 200 cm or 75 x 200 cm.
TRENDELENBURG / REVERSE TRENDelenburg	Electric (+/- 15°).
BACKREST	0 to 80°
EQUIPMENTS	<p>3-sectional mattress base: 65x 200 cm or 75 x 200 cm Flat upholstery, high density foam 35 kg/m3. Width 6,5 cm Jersey with PVC coating, antibacterial and anti-fire (M1) treatment</p> <p>Electric adjustable backrest (0 - 80°) Electric height adjustment: 46 - 86 cm Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°) Steel structure, epoxy coated (light grey) Frame with 2 electric columns providing great stability Frame protection cover Remote control 4 lockable double running wheels Ø 125 mm Examination towel dispenser head end</p>
AVAILABLE OPTIONS	<p>Antistatic system Casters Ø 125 mm with central brake Double running castors Ø 150 mm with central brake Wireless remote control (in addition of the wired remote control) Foot pedal control (in addition of the wired remote control) Drop-down side rails (pair) Holder for clamps (pair) Fixation clamps for heels, armrest for blood sampling or leg support (pair) (DVSE required) Heels (pair) (DVETU required) Anatomic shaped armrest (DVETU necessary) Leg support (pair) (DVETU required) Protection cover for leg rest (pair) Clamps for IV pole holder (pair) Removable 2-stage IV pole, 2 hooks (DVTSR required) Comfort pillow (grey) Battery Examination lamp (DVETU required) Recovery Tank (grey) Equipotential connection DIN rails</p>

OVERVIEW

1.2 AREA OF APPLICATION

This couch is an exam couch. It has equipment that makes it a perfect tool for gynaecological examinations.

1.3 PURPOSE OF THE USER MANUAL

Health professionals, this user manual has been created for the purpose of supplying you with all the necessary information and usage precautions to be followed to ensure that you can safely operate the product on your own.

The instructions for maintenance and product monitoring are intended to ensure its service life.

It is therefore essential for all staff to have access to them and consequently to be informed about the handling procedures and the risks associated with use of this couch. The couch should be handled by a person who is qualified and aware of the risks inherent in improper use of the product.

1.4 SAFEKEEPING

This user manual should be kept within close proximity to the product in a location designed for this purpose. It should be kept away from any item or substance that could obscure its legibility.

1.5 SYMBOLS AND STANDARDS USED

	Type B device
	Indicates that failure to follow the recommendation can lead to endangerment of the patient and the user, or damage to the equipment.
	Compliant with the essential requirements of Directive 93/42/CEE concerning medical devices.
	Protection against splashed water and sprayed water.
	Indicates the maximum weight of the patient and the safe working load.

This couch is a Class I medical device, non-invasive and active.

2 OPERATING INSTRUCTIONS

2.1 BATTERY

The battery housing has its own charger. It is advisable to charge the battery in a cool place for 24 hours using an electrical outlet as soon as it is received and prior to activating it.

For routine use, we strongly recommend keeping the couch plugged into the mains outlet for sedentary use. Reduced battery autonomy ensures that the couch can be operated in the event of a power failure or during transfer.

The battery ensures that the couch can be operated in the event of a power failure or enables the couch to be operated when not connected to an electric power system while being transferred.

However, once operation of the couch has ended, it is very important to remember to reconnect it (for a minimum of 12 hours a day). Never allow the battery to go completely flat (risk of diminished battery performance).

The life of the battery when charged is 25 cycles (powering up + powering down a function). When the sound signal indicates that the battery level is low, immediately reconnect the couch to the electrical outlet.

Conditions for connection to the electric power system and electromagnetic compatibility:

- Connection to the electric power system should be carried out in compliance with standards NF C 15-100 and NF C 15-211 (230 V electrical outlet). The power system must be protected by a circuit breaker with a maximum differential current of 30 mA.
- Operations performed on electrical parts must be carried out by authorized qualified staff.
- Electromagnetic interference from the couch may occur. It is therefore advisable to inhibit the couch functions and keep devices affected by the disturbances away from the couch.

2.2 ADJUSTABLE HEIGHT AND LOWERING


Before using these functions: Make sure that there are no objects or persons underneath the plinth area.

The height is adjusted by remote control

:

or

Raised position

Lowered position

2.3 BACKREST

Prior to activating the torso-lift, check that there is nothing hampering its movement (accessories, cables...)

The backrest is operated by remote control:

or

Or by manually activating the controls located under the torso-lift: press (with both hands for greater ease) to free the springs and adjust the torso-lift to the desired height.

Do not sit on the backrest

2.4 INTERMEDIATE BACKREST (DIVIA MIXTE)

The intermediate backrest is operated by remote control :


The intermediate backrest cannot be activated if the primary backrest is inclined to more than 60°, (see photo above).

2.5 TRENDelenburg POSITION

This is managed by remote control :

Flat position

Trendelenburg position

Do not activate the Trendelenburg position when the couch is in the lowered position

2.6 INDEPENDANT WHEELS BRAKES

Make sure to apply the couch brakes at all times, except during transport. Check that the brake is active by attempting to move the couch.

Do not park the couch on a sloping floor.

Press the pedal on the wheel to apply the brake and raise the pedal to release it.

Wheel brake released


Wheel brake applied

2.7 OPTIONS

2.7.1 HEAD REST

The head rest is easily adjusted via a counterbalancing system located at the rear of the backrest.

2.7.2 CENTRALLY CONTROLLED BRAKE AND STEERING PEDAL


Make sure to apply the couch brakes at all times, except during transport. Check that the brake is active by attempting to move the couch.

Do not park the couch on a sloping floor.

Lateral brake and steering pedals:

- **Brakes applied to couch** – red side in lowered position ①.
- **Directional steering of couch** – green side in lowered position ②.
- **Freehand steering of couch** - pedal at the middle position ③.

2.7.3 SIDERAILS


The bars are intended to reduce the risk of the patient suffering accidental falls. They are in no way designed to restrain the patient or to restrict the patient with strap-like immobilizing devices.

Before adjusting the bars, make sure that there is nothing hampering their movement.

Do not move the couch around by its bars.


Lowering the bars :

Unlock by pressing and holding in the button ①, and pushing the bar ②.

Raising the bars :


Pull up the upper bar until the button locks into place.

2.7.4 CLAMPS FOR FASTENING STIRRUPS, CRADLES, LEG SUPPORTS OR EXAMINATION LAMP

2.7.4.1 TO INSTALL THE CLAMPS

1. Slide the clamp onto the stand provided for this purpose (make sure beforehand that the thumbwheel does not get in the way of the clamp; if it does, unscrew it).


2. Tighten the clamp by turning the thumbwheel.

2.7.4.2 TO INSTALL THE ACCESSORIES

1. Slide the accessory onto the stand intended for this purpose (make sure beforehand that the thumbwheel does not get in the way of the accessory; if it does, unscrew it).

2. Tighten the accessory by turning the thumbwheel

2.7.4.3 TO INSTALL THE IV POLE

The IV holder is not inserted into the clamps provided for the stirrups, cradles or leg supports. Special clamps with an identical function to the one described above are provided.

Maximum safe working load: 10 kg(1)


When conveying the couch with the IV holder attached, make sure that the IV holder is not able to come into contact with any obstacle (e.g. : doorway, lighting fixture, etc.) or individual..

- (1) This indicative value is in keeping with safe use of the product. However, it is not representative of normal use with a much lower load.


2.7.5 LATERAL TECHNICAL RAILS

Slide the clamp on which your accessory is attached into place and secure the thumbscrew to tighten it.

Maximum Safe Working Load : 15 kg

2.7.6 RECOVERY TANKS

3 CLEANING AND DISINFECTION

3.1 INTRODUCTION

Regular proper cleaning ensures the service life of the device (cf. established cleaning schedule established when the device is put into service, specifying the resource reference staff).

Failure to carry out cleaning operations may lead to the risk of infection due to the presence of secretions and/or residue on the couch.

The following cleaning recommendations are not a substitute for your establishment's in-house protocol.

3.2 FREQUENCY

Cleaning operations should be performed:

- per patient, i.e. between patients
- after transporting an infected patient
- thoroughly on a weekly basis.

3.3 RECOMMENDED MATERIAL AND PRODUCTS

- Only use authorized, standardized cleaning products or disinfectants.
- Take note of products that are strictly prohibited by the manufacturer (Cf. list below).

List of strictly prohibited products:

Solvents :	Petrol / higher aliphatic and aromatic alcohols /ketones-esters / trichlorethylene-perchlorethylene
Base :	Soda 10% / Ammonia 10%
Acids :	Nitric 10% / Hydrochloric 30% / Citric 10%
Misc.:	Distilled water / Bleach / Abrasive powders / Any abrasive product

- Do not use abrasive utensils or mildly diluted bleach.
- Do not mix products.
- Follow the prescribe proportions and dilutions.
- Use clean disposable cloths.
- When cleaning and disinfecting, wear protective clothing and gloves.

3.4 DIRECTIONS FOR CLEANING AND DISINFECTION

3.4.1 CLEANING OPERATION BASED ON THE COMPONENTS

Before any cleaning operation is carried out on the couch, disconnect it from the electric power system and ensure that all internal parts are correctly connected.

COMPONENTS	OPERATIONS REQUIRED
ELECTRICAL PARTS	Wash with (very little) running water in cases where the springs are watertight.
	After cleaning, dry the couch
	Check that it is functioning properly
METAL PARTS (90% EPOXY COVERED)	Wash with soap and running water or with diluted alcohol
	After cleaning, dry the couch
	Check that it is functioning properly
UPHOLSTERY	Wash with soap and water
	After cleaning, dry the couch
	Check that it is functioning properly

3.4.2 DIRECTIONS TO FOLLOW

- All of the surfaces of the furniture must be cleaned.
- During cleaning or disinfection, it is advisable to:
 - o Apply the product by spraying or using a cloth.
 - o Start with the areas nearest to the patient, paying particular attention to the parts coming into contact with users.
 - o Work from the top down
- Allow the equipment to dry.

Products that stain (e.g. "blue product" or eosin) should be cleaned as quickly as possible in order to avoid it permeating the material.

4 PRECAUTIONS

4.1 PRECAUTIONS FOR THIS DEVICE*

4.1.1 PRECAUTIONS FOR USE OF THE COUCH

For effective use of the equipment, proceed only if the weight is properly distributed on the device.

In the case of a difficult or unsupervised patient, it is advisable to leave the splint area in the lowest position in order to reduce the risk of a potential fall.

Similarly, it is advisable to use the “adjustable height” function whenever handling the patient, in order to work at an ergonomic height.

The patient should be placed in a sitting or standing position on the base side, by lowering the bars. The patient should not be seated on the torso-lift. Risk of falling.

Under no circumstances can either of the bars be considered an area inaccessible to patients. The appropriate measures must be taken depending on the patient.

During conveyance, it is necessary to make sure that the power cord is disconnected and wrapped around the reel provided for this purpose before moving the couch.

When the couch does not need to be moved, it should be immobilized using the brakes, and then making sure that the wheels are locked by attempting to move it. Failure to observe this handling method can lead to a fall by someone attempting to lean against the couch.

All of the handling methods relating to the precautions for use mentioned above are described in section II– Operating instructions.

It is inadvisable to handle without appropriate training.

4.1.2 PRECAUTION FOR SERVICING AND MAINTENANCE OF THE EQUIPMENT

No modifications must be made to the device without the prior consent of the manufacturer: the modification may cause it to function in an unpredictable manner and thereby cause injury to the patient and the nursing staff.

The device requires regular inspection.

In order to maintain the device in good condition and ensure that it can be safely operated over the long term, it is necessary to service it in the prescribed manner.

Should abnormalities or damage which could compromise the operation and safety of the device be noted, the equipment must be immediately removed.

The device must be used exclusively in accordance with the instructions in this instruction manual.

The device must not be exposed to, much less come into contact with, combustion heat sources or flammable agents.

BMB MEDICAL shall not be held liable in the event of incidents or accidents resulting from improper use or modification of the product by the buyer or user.

4.2 NON STANDARD USES

 **This couch is designed and suitable for medical practices and procedures and therefore for specific usage. Non-standard use of this couch may result in accidents suffered by the user or even irreversible deterioration of the couch.**

The following list of non-standard uses of the couch is not an exhaustive one. BMB MEDICAL leaves it to users to apply common sense.

- Use by any person who has not read this user manual.
- Use by any person outside of the nursing staff who has not been authorized or advised by a qualified person.
- Use as a bed for a short person (under 1.30 m) or under the age of 12.
- Handling by a child.
- Use of the electric controls by more than one person.
- Use of the couch with a heavier weight than the safe working load (refer to the technical features).
- Use of the chassis to hold or store objects.
- Use of oxygen equipment under the plinth frame.
- Use in the presence of a flammable anesthetic mixture.
- Use outside of the establishment or on an unfinished floor surface.
- Use on a floor surface with a gradient greater than 10°.
- Washing by sluicing, high-pressure water or steam jet or passage through a tunnel washer.
- Rigorous use or jolting of the motors.
- Moving the couch while it is still connected to the electric power source.
- Connecting an electrical device to the couch by any means whatsoever.
- Any attempt to unscrew the ground terminal screw from the power supply box.
- Connecting the couch to a non-compliant power source.
- Any other use that does not qualify as a function for a hospital couch.

Information: The electrical springs are equipped with a thermal protection feature that prevents them from functioning if the motor overheats from rigorous use. They return to their original working condition within 2 hours of the couch being disconnected from the electric power source.

4.3 TRANSPORT AND STORAGE

4.3.1 TRANSPORT

When transporting the couch to its storage location, in order to avoid any risk of collision or falls it is advisable to:

- strap it to a support on the transport medium
- cover it
- lock the brakes
- adjust it in the lowered position

Any damage occasioned during transport of the couch are not covered by the warranty. Repair and replacement of the damaged parts are the responsibility of the client.

4.3.2 STORAGE

When the couch is stored, make sure that :

- the brakes are locked
- the couch is covered
- the couch is located in a cool, dry room (temperature: between +10°C and +40°C).

While in storage, it should not have any heavy weight placed on top of it.

The couch must not be treated as a support base for any type of equipment.

The device must not be exposed to, much less come into contact with sources of combustion heat or flammable agents.

In the event of prolonged storage, it is advisable to recharge the battery every 6 months as follows:

- Connect the power cord.
- Briefly press the control function to initiate charging.
- Leave the battery charging for 8 hours.

4.4 WARRANTY

Two-year warranty on parts and labor for the equipment.

However, the warranty is rendered null and void by:

- non-standard use ([Section 4.2](#))
- interference without prior notification to our technical department or servicing of:
 - o a mechanical part of the couch
 - o a spring or a motor
- failure to comply with any other condition of the warranty
- any inappropriate cleaning
- the use of accessories not sourced from the manufacturer
- failure to comply with the preventive servicing operations recommended on the servicing check-list

4.5 AFTER SALES SERVICE

For any information on installation, use or servicing of the product and any further questions, contact the After-Sales Service department at the following email address: sav@bmb-medical.com or on +33 (0)4 74 08 71 71 or your local distributor for International Clients.

You can access resource information in our “Customer Services Forum” and view maintenance tips and request our services directly through this secure interface.

The operator should carefully read the information contained in this user manual, particularly the information pertaining to precautions for safe use and methods of installation and use.

If you are unsure about correctly interpreting the instructions, kindly contact the After-Sales Service department which is always at your service.

For easy ordering of spare parts, you can request an exploded view of your products from our After-Sales Service department to obtain the references of the desired parts.

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4.6.1 USED EQUIPMENT

Recycling electrical equipment helps preserve natural resources and avoid the risk of pollution. To this end, BMB MEDICAL is fulfilling its obligations with respect to the expired life cycle of the couches it markets by funding the recycling subsidiary Recylum dedicated to WEEE products which it takes back free of cost (Find out more at www.recylum.com).

Tél: +33 (0)4 74 08 71 71

Fax: +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001
01604 TREVOUX Cedex
FRANCE