CAREXIA

MANUEL D'UTILISATION
USER MANUAL

 **BMB
MEDICAL**

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FRANÇAIS

BMB MEDICAL

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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www.bmb-medical.com

The purpose of this manual is to familiarize you with the use, maintenance and safety of your new medical device. We encourage you to read this document carefully and to follow the guidelines for optimum and safe use.

BMB MEDICAL offers a broad variety of options for the products covered by this user manual. Consequently, some features described in this manual may not apply to your chair.

Some illustrations are provided for information purposes and do not accurately represent the device in your possession.

Stick the product reference
label here

The information in this document are subject to change without notice and in no way engage the responsibility of BMB MEDICAL. This document may not be reproduced, in whole or in part, without prior written consent by BMB MEDICAL. This illicit representation or reproduction, by any means, constitutes counterfeiting. The illustrations in this document are not binding.

TABLE OF CONTENTS

1	INTRODUCTION SPECIFICATIONS	45
1.1	PRODUCT PRESENTATION	45
1.1.1	CAREXIA FP VERSATILE FIXED-HEIGHT CARE CHAIR	46
1.1.2	CAREXIA FPE VERSATILE FIXED-HEIGHT CARE CHAIR WITH ELECTRIC BACKREST AND LEGREST.....	48
1.1.3	CAREXIA FPV VERSATILE ELECTRICAL VARIABLE-HEIGHT CARE CHAIR	50
1.1.4	CAREXIA FPVE VERSATILE VARIABLE-HEIGHT CARE CHAIR WITH ELECTRIC BACKREST AND LEGREST.....	52
1.2	PURPOSE OF THE USER MANUAL.....	54
1.3	USERS	54
1.4	STORAGE.....	54
1.5	APPLIED PARTS	54
1.6	SYMBOLS AND STANDARDS USED	55
2	INSTRUCTIONS FOR USE	57
2.1	OPERATING CONDITIONS.....	57
2.2	INSTALLATION PROCEDURE	57
2.3	BATTERY (OPTIONAL).....	57
2.4	REMOTE CONTROL.....	59
2.4.1	CAREXIA FPV CHAIR.....	59
2.4.2	CAREXIA FPE CHAIR	59
2.4.3	CAREXIA FPVE CHAIR.....	60
2.4.4	TRANSPORT POSITION.....	60
2.4.5	RESET	61
2.5	VARIABLE HEIGHT AND LOWERING.....	61
2.6	BACKREST	61
2.6.1	CAREXIA GAS CYLINDER-ASSISTED BACKREST VERSION (FP AND FPV VERSIONS) 61	
2.6.2	CAREXIA ELECTRIC BACKREST VERSION (FPE AND FPVE VERSIONS)	62
2.7	LEGREST	62
2.7.1	CAREXIA GAS CYLINDER-ASSISTED BACKREST VERSION (FP AND FPV VERSIONS) 62	
2.8	FOOTREST.....	63
2.9	ARMRESTS.....	63
2.9.1	ERGONOMICALLY SHAPED REAR-FOLDING POLYURETHANE ARMRESTS.....	63
2.9.2	ANATOMICALLY SHAPED, FOLDING, SWIVEL AND HEIGHT ADJUSTABLE ARMRESTS	64

2.9.3	COMFORT ARMRESTS.....	65
2.10	BRAKE AND GUIDE PEDALS.....	66
2.10.1	CENTRAL BRAKING AND SIDE GUIDES.....	66
2.10.2	INDEPENDENT BRAKING.....	66
2.11	OPTIONS.....	67
2.11.1	IV POLE.....	67
2.11.2	FLEXIBLE TABLET HOLDER.....	68
2.11.3	PAPER ROLL HOLDER.....	69
2.11.4	PUSH BARS.....	69
2.11.5	READING LAMP.....	69
2.11.6	COMFORT CUSHION.....	70
2.11.7	USB CHARGING PORT.....	70
2.11.8	REMOVABLE CRYSTAL PROTECTION.....	70
3	CLEANING AND DISINFECTION.....	71
3.1	INTRODUCTION.....	71
3.2	FREQUENCY.....	71
3.3	CLEANING AND DISINFECTION ADVICE.....	71
3.3.1	CLEANING OPERATION.....	71
3.3.2	RECOMMENDED EQUIPMENT AND PRODUCTS.....	72
4	PRECAUTIONS.....	74
4.1	DEVICE-RELATED PRECAUTIONS.....	74
4.1.1	WARNINGS.....	74
4.1.2	CONTRAINDICATIONS.....	75
4.1.3	PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE.....	76
4.1.4	EMC-RELATED INFORMATION.....	77
4.2	TRANSPORT AND STORAGE.....	79
4.2.1	TRANSPORT.....	79
4.2.2	STORAGE.....	79
4.3	WARRANTY.....	80
4.4	AFTER-SALES SERVICE.....	81
4.5	ENVIRONMENT.....	81
4.5.1	SPENT EQUIPMENT.....	81

1 INTRODUCTION SPECIFICATIONS

1.1 PRODUCT PRESENTATION

This patient support device consists of a backrest and seat fitted onto a chassis with wheels. The device is equipped with armrests and, optionally, allows the temporary fitting of IV poles.

The device allows an operator to convey the patient within a healthcare establishment. The chair is braked by a mechanical system actuated by pressing a pedal.

The device is equipped with an electrically assisted actuator system actuated by a remote control to manage the movements of the following elements:

- backrest,
- legrest,
- footrest,
- electrical column used to adjust the product's height.

The products are designed to ensure that there are no unacceptable risks (essential performance) associated with patient support or caused by accidental movements.

1.1.1 CAREXIA FP VERSATILE FIXED-HEIGHT CARE CHAIR

TECHNICAL DATA	
REFERENCE	CAREXIA FP
DESIGNATION	VERSATILE FIXED-HEIGHT CARE CHAIR
SAFE WORKING LOAD	200 kg (indicates the maximum permissible load, including the total weight of the patient, mattress and accessories)
WEIGHT	65 to 85 kg depending on model
OVERALL DIMENSIONS	Width: 69.5 cm with armrests, or 74 to 90 cm with gutters. Seat width: 57 cm. Overall length: 76.4 cm in sitting position to 170 cm in lying position with long legrests. Overall height: 130 cm.
EQUIPMENT	Heavy-duty design mechanically welded steel frame. Oven-baked epoxy anti-corrosion coating. Fixed seat height: 55 cm. Ergonomically shaped seat and backrest. M1-rated non-flammable PVC coating. High-resilience foam, HR35 kg/m³ on rigid frame. Integrated headrest. Gas spring-assisted adjustable backrest from 80 to 60°. Mounted on 4 self-locking double-roller casters. ABS chassis cover with basket. Ergonomically shaped polyurethane armrests, can be folded back.
AVAILABLE OPTIONS	4 castors Ø 125 mm with central lock 4 double runing wheel castors Ø 150 mm with central lock Push bar Short legrest synchronized with the backrest. Trend position. 2 holders for IV pole (necessary for TSR2/TSR4) Removable IV pole, 2-stage, 2 hooks Removable IV pole, 4-stage, 4 hooks Comfort pillow (grey) Armrest comfort (pair) Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) Paper roll holder High comfort upholstery with 2 cm thickness memory foam Fixed I-Pad support with antitheft device

OVERVIEW

**BMB
MEDICAL**
dedicated solutions for your
patients' transfer & wellbeing

⚠ Safe working load (uniformly distributed): 200 kg

Integrated headrest

Gas spring-assisted backrest adjustable from 80° to 60° (or adjustable to recovery position with the legrest option)

ABS chassis cover with basket.

IV pole holder (optional)

Self-locking casters with double running wheels (optional central lock)

Removable and adjustable IV pole (optional)

Ergonomic seat and backrest with M1-rated non-flammable PVC coating

High-resilience foam, HR40 kg/m3 on rigid frame.

Armrests can be folded back

⚠ never apply undue force to the armrests

Adjustment control lever

Short legrest synchronized with the backrest (optional)

1.1.2 CAREXIA FPE VERSATILE FIXED-HEIGHT CARE CHAIR WITH ELECTRIC BACKREST AND LEGREST

TECHNICAL DATA	
REFERENCE	CAREXIA FPE
DESIGNATION	VERSATILE FIXED-HEIGHT CARE CHAIR WITH ELECTRIC BACKREST AND LEGREST
SAFE WORKING LOAD	200 kg (indicates the maximum permissible load, including the total weight of the patient, mattress and accessories)
WEIGHT	65 to 85 kg depending on model.
OVERALL DIMENSIONS	Width: 69.5 cm with armrests, or 74 to 90 cm with gutters. Seat width: 57 cm. Overall length: 76.4 cm in sitting position to 170 cm in lying position with long legrests. Overall height: 130 cm.
EQUIPMENT	Heavy-duty design mechanically welded steel frame. Oven-baked epoxy anti-corrosion coating. Fixed seat height: 55 cm. Ergonomically shaped seat and backrest. M1-rated non-flammable PVC coating. High-resilience foam, HR35 kg/m³ on rigid frame. Integrated headrest. Independent and synchronized electric backrest and legrest. Adjustable to emergency position - 1 key on the remote control. Mounted on 4 self-locking double-roller casters. ABS chassis cover with basket. Ergonomically shaped polyurethane armrests, can be folded back.
AVAILABLE OPTIONS	4 castors Ø 125 mm with central lock 4 double running wheel castors Ø 150 mm with central lock Push bar Long legrest Electric adjustable footrest (only for long legrest) Foot pedal control for relax and trend position 2 holders for IV pole (necessary for TSR2/TSR4) Removable IV pole, 2-stage, 2 hooks Removable IV pole, 4-stage, 4 hooks Comfort pillow (grey) Armrest comfort (pair) Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) Paper roll holder High comfort upholstery with 2 cm thickness memory foam Reading lamp DEL Battery 2 Ah USB Port Removable Crystal protection for footrest Removable Crystal protection for footrest Fixed I-Pad support with antitheft device

OVERVIEW

**BMB
MEDICAL**

dedicated solutions for your
patients' transfer & wellbeing

△ Safe working load (uniformly distributed): 200 kg

Paper roll holder
(optional)


Integrated headrest

Independent electric backrest:
adjustable to recovery position

ABS chassis cover with
basket.

Side brake and guide
pedal
(optional)

- green adjusted down:
directional guiding by the
2 front wheels.

- horizontal:
free guiding.

- red adjusted down:
locked position.

Self-locking casters with double running wheels
(optional central lock)

Ergonomic seat and backrest
with M1-rated non-flammable
PVC coating

High-resilience foam, HR40
kg/m3 on rigid frame.

Rear-folding armrests

(comfort armrest option on the photo)
△ never apply undue force to the
armrests

Independent short electric
legrest. Position adjustment by
remote control

1.1.3 CAREXIA FPV VERSATILE ELECTRICAL VARIABLE-HEIGHT CARE CHAIR

TECHNICAL DATA	
REFERENCE	CAREXIA FPV
DESIGNATION	VERSATILE ELECTRICAL VARIABLE-HEIGHT CARE CHAIR
SAFE WORKING LOAD	200 kg (maximum permissible load, including the total weight of the patient, mattress and accessories)
WEIGHT	65 to 85 kg depending on model.
OVERALL DIMENSIONS	Width: 73.2 cm with armrests, or 74 to 90 cm with gutters. Seat width: 57 cm. Overall length: 76.4 cm in sitting position to 170 cm in lying position with long legrests. Overall height: 130 cm.
EQUIPMENT	Heavy-duty design mechanically welded steel frame. Oven-baked epoxy anti-corrosion coating. Electric variable height from 54 to 78 cm (floor/seat height). Foot pedal variable height control on both sides. Foot pedal control on both sides. Ergonomically shaped seat and backrest. M1-rated non-flammable PVC coating. High-resilience foam, HR35 kg/m³ on rigid frame. Integrated headrest. Gas spring-assisted adjustable backrest from 80 to 60°. Mounted on 4 self-locking double-roller casters. ABS chassis cover with basket. Ergonomically shaped polyurethane armrests, can be folded back.
AVAILABLE OPTIONS	4 castors Ø 125 mm with central lock 4 double runing wheel castors Ø 150 mm with central lock Push bar Short legrest synchronized with the backrest. Trend position. 2 holders for IV pole (necessary for TSR2/TSR4) Removable IV pole, 2-stage, 2 hooks Removable IV pole, 4-stage, 4 hooks Comfort pillow (grey) Armrest comfort (pair) Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) Paper roll holder High comfort upholstery with 2 cm thickness memory foam Battery 2 Ah USB Port Fixed I-Pad support with antitheft device

OVERVIEW

**BMB
MEDICAL**

dedicated solutions for your
patients' transfer & wellbeing

⚠ Safe working load (uniformly distributed): 200 kg

Integrated headrest

Gas spring-assisted backrest
adjustable from 80° to 60°
(or adjustable to recovery position with
the legrest option)

Rear-folding armrests
(comfort armrest option on the photo)
⚠ never apply undue force to the
armrests

ABS chassis cover
with basket.

Foot pedal control on both sides.
Electric variable height from 54 to
78 cm (floor/seat height)

Ergonomic seat and backrest
with M1-rated non-flammable
PVC coating

High-resilience foam, HR 60
kg/m³ on rigid frame.

Adjustment control lever

Self-locking casters with
double running wheels
(optional central lock)


1.1.4 CAREXIA FPVE VERSATILE VARIABLE-HEIGHT CARE CHAIR WITH ELECTRIC BACKREST AND LEGREST

TECHNICAL DATA	
REFERENCE	CAREXIA FPVE
DESIGNATION	VERSATILE VARIABLE-HEIGHT CARE CHAIR WITH ELECTRIC BACKREST AND LEGREST
SAFE WORKING LOAD	200 kg (maximum permissible load, including the total weight of the patient, mattress and accessories)
WEIGHT	65 to 85 kg depending on model.
OVERALL DIMENSIONS	Width: 73.2 cm with armrests, or 74 to 90 cm with gutters. Seat width: 57 cm. Overall length: 76.4 cm in sitting position to 177 cm in lying position with long legrests. Overall height: 130 cm.
EQUIPMENT	Heavy-duty design mechanically welded steel frame. Oven-baked epoxy anti-corrosion coating. Electric variable height from 54 to 78 cm (floor/seat height). Ergonomically shaped seat and backrest. M1-rated non-flammable PVC coating. High-resilience foam, HR35 kg/m ³ on rigid frame. Integrated headrest. Independent and synchronized electric backrest and legrest. Adjustable to emergency position - 1 key on the remote control. Mounted on 4 self-locking double-roller casters. ABS chassis cover with basket. Ergonomically shaped polyurethane armrests, can be folded back.
AVAILABLE OPTIONS	4 castors Ø 125 mm with central lock 4 double runing wheel castors Ø 150 mm with central lock Push bar Long legrest Electric adjustable footrest (only for long legrest) Foot pedal control for relax and trend position Foot pedal control for variable height on both sides 2 holders for IV pole (necessary for TSR2/TSR4) Removable IV pole, 2-stage, 2 hooks Removable IV pole, 4-stage, 4 hooks Comfort pillow (grey) Armrest comfort (pair) Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) Paper roll holder High comfort upholstery with 2 cm thickness memory foam Reading lamp DEL Battery 2 Ah USB Port Removable Crystal protection for footrest Removable Crystal protection for footrest Fixed I-Pad support with antitheft device

OVERVIEW

BMB MEDICAL
dedicated solutions for your patients' transfer & wellbeing

⚠ Safe working load (uniformly distributed): 200 kg

1.2 PURPOSE OF THE USER MANUAL

Healthcare professionals, this user manual was created for you to provide you with the necessary information and precautions to observe for autonomous and safe use.

The product's service life is dependent upon observance of the maintenance and tracking instructions.

It is thus essential that all staff has access to and be informed of the handling procedures and risks associated with the use of this chair. The latter should be handled by a qualified person who is aware of the risks inherent to improper use. Additional training is not required as we consider that knowledge of the user manual is sufficient.

1.3 USERS

Healthcare professionals: physicians, caregivers and agents approved by Technical Services can operate all functions.

In use cases where the patient is the operator (**use of the remote control and armrests only**), the caregivers must first ascertain the patient's ability to use the product. The caregivers must also teach the patient how to safely use the product: adjustment to low position to step off the product.

1.4 STORAGE

This user manual should be kept near the product, in a place provided for this purpose. It must be shielded from all elements and substance that could compromise its perfect legibility.

1.5 APPLIED PARTS

The applied parts, according to the EN 60601-1 standard, are as follows:

- upholstery
- armrests
- remote controls
- reading lamp
- footrest

1.6 SYMBOLS AND STANDARDS USED

	Type B appliance
	Electrical class II equipment
	Electrical hazard, do not open
	Duty cycle (motor cycle) 2 min. operation / 18 min. stop
IPX6	Splash-resistant
	See the manual / instruction brochure
	Indicates the maximum permissible load, including the total weight of the patient, mattress (where applicable) and accessories.
	Conforms to the essential requirements of directive 93/42/EEC on medical devices.
	Serial number
	Manufacture date
	Manufacturer
	Alternating current
	Maximum patient weight
	Electronic waste: this symbol indicates that the product must not be discarded with unsorted municipal waste, but rather to a selective waste collection channel. Contact the local distributor for information concerning the return and/or collection systems available in the country.
	Protéger de l'humidité / keep dry
	Haut, tenir debout / stand up
	Fragile / fragile
	Recyclable / recycling
	Temperature range
	Humidity limit

General safety sign:

- Information highlighted by this symbol must be read and scrupulously observed.
- Indicates that failure to observe the guideline could endanger the patient or user, or damage the equipment.

	Sitting on the backrest is prohibited		Using the footrest as a footboard is prohibited
	Sitting on the legrest is prohibited		Leaning against the chair is prohibited
	Sitting on the armrests is prohibited		The chair must be set to its low position before being moved

“WARNING: This appliance must not be altered without the manufacturer’s written consent.”

Caution: only purpose-built accessories provided for this MD by BMB MEDICAL may be used.

The FP version of this chair is a Class I non-invasive and non-active medical device.

The FPE, FPV and FPVE versions of this chair are Class I non-invasive and active medical devices.

2 INSTRUCTIONS FOR USE

2.1 OPERATING CONDITIONS

AMBIENT CONDITIONS	OPERATION
Ambient temperature	+5 °C to +40 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

2.2 INSTALLATION PROCEDURE

Before commissioning the chair and depending on the model, check that the following are operational:

- Remote control
- Variable height
- Backrest
- Legrest
- Footrest
- Armrests
- Brake and guide pedals

Clean the product according to your establishment's internal protocol.

2.3 BATTERY (OPTIONAL)

The battery unit is fitted with its own charger. Before initial use, leave the product to stand for 12 h in a temperate environment (to allow the electrical and electronic components to reach normal operating temperature), then charge the battery for 12 hours from a mains outlet.

In routine use, we strongly recommend keeping the chair plugged into the mains during sedentary use.

The battery, with its limited operating life, can be used to operate the chair in the event of a power outage, or during transfer.

Once the manipulations are complete, remember to re-connect the chair to the mains (at least 12 hours per day). Never allow the battery to drain fully (risk of reduced battery performance).

“Warning: Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact.”

Battery characteristics		Battery safety symbols	
Battery type:	Lead	IPX6	Protection rating
Capacity:	1.2 Ah		CE marking
Input supply voltage:	29~45 VDC		Indoor use
Input supply voltage:	24 VDC, 1.2 A		Electronic waste

Power grid connection conditions and electromagnetic compatibility:

This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.

Position the chair such that it can be easily unplugged.

Unplug the product before moving it.

Do not leave the cables in contact with the floor or the casters,

Do not knot the cables,

Periodically check the condition of the cables.

Mains connection:	120 ~ 230 V depending on the country.
Frequency:	50/60 Hz
Rated current:	2 A
Protection rating / degree of electric shock protection:	Class 2 appliance
Power input:	525 W
Intermittent operating mode:	10%, i.e. 2 minutes of continuous operation followed by 18 minutes of rest.

The mains plug is the electrical grid insulation device.

Work on electrical parts should only be performed by qualified and authorized staff.

Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the EMC information [here](#).


2.4 REMOTE CONTROL

When using the remote control, all movements can be interrupted by simply releasing the corresponding key.


2.4.1 CAREXIA FPV CHAIR

Foot remote control


Variable height

2.4.2 CAREXIA FPE CHAIR

Hand-held remote control


2.4.3 CAREXIA FPVE CHAIR

4. Chair in patient installation position - press the remote control until the chair stops in the position
5. Chair in comfort position - press the remote control until the chair stops in the position
6. Emergency position - press the remote control and hold the button after arrival in comfort position until the chair reaches the emergency position

Foot remote control

Variable height

Backrest + legrest

Synchronized

Only activate the various positions after having ensured that no-one is touching any mechanical parts and that there are no objects under the seat, or between any of the moving structures.

2.4.4 TRANSPORT POSITION

When transporting the patient, the chair should be set to comfort position, with the seat in the low position (see Variable height [here](#))

2.4.5 RESET

This function should be used when the chair is not responding to the commands issued by the remote control. To reset the commands:

- Connect the chair to the mains
- Simultaneously press the keys indicated on the remote control until the LEDs start flashing.

2.5 VARIABLE HEIGHT AND LOWERING

Before using these functions, ensure that there are no objects or persons under the patient surface.

The height is adjusted using the remote control:

2.6 BACKREST

Before actuating the backrest, ensure that there are no obstacles to its movement (accessories, cables, etc.)

2.6.1 CAREXIA GAS CYLINDER-ASSISTED BACKREST VERSION (FP AND FPV VERSIONS)

The backrest is adjusted using the handle on the side of the seat and manually tilting it to the desired position. Backrest adjustable to emergency position.

The backrest movements are synchronized with those of the legrest, the adjustments are interdependent.

2.6.2 CAREXIA ELECTRIC BACKREST VERSION (FPE AND FPVE VERSIONS)

Before actuating the backrest, ensure that there are no obstacles to its movement (accessories, cables, etc.)

Long legrest
Independent and synchronized

Short legrest
Independent and synchronized

The short legrest can be fully folded under the seat, contrary to the long legrest.

The legrest is controlled using the remote control: or

If a long legrest is used, the patient should be seated from the side.

2.7 LEGREST

Before actuating the legrest, ensure that there are no obstacles to its movement (accessories, cables, etc.)

2.7.1 CAREXIA GAS CYLINDER-ASSISTED BACKREST VERSION (FP AND FPV VERSIONS)

The legrest movements are synchronized with those of the backrest, the adjustments are interdependent.


The short legrest can be fully folded under the seat when the chair is in patient installation position.

Short legrest
assisted and
synchronized

2.8 FOOTREST

The footrest is adjusted using the remote control:

Footrest in low position

Footrest in raised position

2.9 ARMRESTS

2.9.1 ERGONOMICALLY SHAPED REAR-FOLDING POLYURETHANE ARMRESTS

To raise or lower the armrests, simply lift or lower them manually:

Adjusting the distance between armrests:

To improve patient comfort, the distance between the armrests can be adjusted. Simply pull them out or push them in to the desired distance.

Pull or push the lower part of the armrest.

2.9.2 ANATOMICALLY SHAPED, FOLDING, SWIVEL AND HEIGHT ADJUSTABLE ARMRESTS

Armrest height adjustment - the armrests have 9 adjustment positions:

Lift the armrests to the desired position.

To lower the armrest, lift it to its high position,

then lower it.

Swiveling the armrest:

Press the paddle under the armrests. This paddle allows you to position the armrest as desired. The armrest can be rotated by up to 180°.

Adjusting the distance between armrests:

To improve patient comfort, the distance between the armrests can be adjusted. Simply pull them out or push them in to the desired distance.

Pull or push the lower part of the armrest.

2.9.3 COMFORT ARMRESTS

Adjusting the distance between armrests:

To improve patient comfort, the distance between the armrests can be adjusted. Simply pull them out or push them in to the desired distance.

2.10 BRAKE AND GUIDE PEDALS

Except during transport, always ensure that the chair's brake is engaged. Check whether the brake is engaged by attempting to move the chair.

Do not park the chair on a sloping floor.

2.10.1 CENTRAL BRAKING AND SIDE GUIDES

Brake and side guide pedals:

- Brake engaged - red side lowered ❶.
- Chair in directional guidance mode - green side lowered ❷.
- Chair in free guidance mode - pedal in median position ❸.

❶

❷

❸

2.10.2 INDEPENDENT BRAKING

Lock all 4 wheels by pressing the brake pedals.

Brake released

Brake engaged

2.11 OPTIONS

2.11.1 IV POLE

The IV pole is removable, it inserts into the holder on the armrests:

To deploy the IV pole:

- Turn the thumb wheel **①** counter-clockwise and extend the telescopic part of the IV pole **②** to the desired height.
- To lock the telescopic part, turn the thumb wheel **①** clockwise. Ensure that the IV pole is securely locked.

Safe working load: 6 kg.


This indicative value corresponds to safe product use, but is not representative of normal use for which the load is significantly lower.

When moving the chair fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.


Do not move the chair via the IV pole.

2.11.2 FLEXIBLE TABLET HOLDER

The flexible tablet holder is fitted onto the product's armrest:


Instructions for use:

Insert the tablet into the lower slot, then pull the handle until the tablet is locked.

The tablet holder is fitted with an anti-theft locking system. Use the key to lock the holder as shown in the image opposite.

2.11.3 PAPER ROLL HOLDER


Pull the tube to the left and raise it slightly to slide the paper roll on, then pull the tube to the left once more to reposition it in its housing.



Do not push or pull the product via the paper roll holder.

2.11.4 PUSH BARS

2.11.5 READING LAMP

To turn the lamp on or off, press:

2.11.6 COMFORT CUSHION

The comfort cushion can be easily adjusted via its counterweight system positioned behind the backrest.

2.11.7 USB CHARGING PORT

To use this option, connect the cable to the USB port under the chair.

2.11.8 REMOVABLE CRYSTAL PROTECTION

Secured with Velcro
at the back of the
legrest or on the
footrest

3 CLEANING AND DISINFECTION

3.1 INTRODUCTION

Regular and thorough cleaning will guarantee the device's life (see cleaning program drawn up at the time of device commissioning, designating the reference staff).

Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.

The following cleaning advice is not intended to replace your establishment's in-house protocol.

This chair has been designed to facilitate its cleaning thanks to its removable cover.

3.2 FREQUENCY

The chair should be cleaned:

- Between each patient.
- After transporting an infected patient.
- Completely each week

3.3 CLEANING AND DISINFECTION ADVICE

3.3.1 CLEANING OPERATION

Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that cause stains (e.g.: “blue product” or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

In the case of a product fitted with a USB charging port, ensure that no liquids enter the port.

3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS

- Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.
- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
 - o PVC, ABS, Polyamide, Polyurethane and Polypropylene
 - o Epoxy paint-coated metal surfaces
 - o Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

List of strictly prohibited products:

Solvents:	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
Bases:	10% Sodium hydroxide / 10% Ammonia
Acids:	10% Nitric / 30% Hydrochloric / 10% Citric
Miscellaneous:	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

High-pressure and steam cleaning are prohibited.

Any damage resulting from failure to observe these instructions for use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.

4 PRECAUTIONS

4.1 DEVICE-RELATED PRECAUTIONS

4.1.1 WARNINGS

This chair has been designed for and adapted to medical practices and procedures, and hence a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the common sense of the chair's users.

- The product must be connected to a mains power outlet compliant with standards in force.
- Once the brakes have been engaged, patient sitting and standing should be performed, with the seat in the low position:
 - o From the side, after raising the armrests (placing the chair in patient installation position)
 - o or from the front, if comfort armrests are used (placing the chair in patient installation position) with the help of a caregiver.
- When getting up, the patient must not rest on the articulated anatomically shaped armrests.
- The patient must not place his/her feet on the footrest.
- The patient must not sit on the legrest or backrest.
- In the case of an unsupervised patient, the seat should be left in its lowest position to reduce the risk of falling.
- In the case of an unassisted patient, the caregivers must inform the patient of the safety measures that must be implemented when using the product.
- Tell the patient not to sit on the ends (foot or head) of the chair. An excess weight applied to one of the ends could tip the patient surface and injure the patient.
- Tell the patient and any accompanying persons not to sit on the armrests.
- Before moving the product:
 - o Ensure that the power cord has been disconnected and placed on the product cover.
 - o Ensure that the brakes are released. It may be hazardous for the user and patient when attempting to move the product with the brakes engaged.
 - o Ensure that the patient's limbs do not extend beyond the armrests and legrest.
 - o The product must be placed in the lowered and comfort position.

- To ensure correct use of the product, only proceed if the weight is correctly distributed over the device.
- The “variable height” function should be used when handling the patient in order to work at an ergonomic height.
- Once movement is no longer the required, the chair should be immobilized with the brakes. Next, check that the wheels are locked by attempting to move the chair. Failure to observe this manipulation could lead to the fall of a person leaning against the chair.
- Do not use the MD on a floor with a slope of more than 10%.
- To avoid the risk of patient or user injury or IV pole damage when moving the chair, ensure that the IV pole is firmly secured.
- To avoid damaging the equipment during transport, ensure that the IV pole passes safely through all door frames and under all light fittings.
- Never use the chair in the presence of an anesthesia mixture flammable with air, with oxygen, or with nitrous oxide as this could cause a fire and/or explosion.
- To prevent failure, electromedical equipment requires special EMC-related precautions, described [here](#).
- Do not manually raise the legrest.
- Work on electrical parts should only be performed by qualified and authorized staff.
- Used to accommodate a person of small stature (less than 1.30 m tall), or under the age of 12 years.
- Do not alter the chair. Any changes made to the system could lead to unpredictable operation that could in turn jeopardize the safety of the user or patient. This would also void the chair's warranty.
- The chassis cover load capacity is of 35 kg. Do not sit or stand on the cover.
- To avoid material damage or injury, do not use the IV pole to pull or push the product.
- To avoid damaging the equipment, the maximum permissible load of the IV pole is 15 kg.
- Use by a person who has not first familiarized themselves with this user manual.
- Use by any non-caregiver who has not received the authorization and explanations from a qualified individual.
- Use of the chair with a load that exceeds the safe working load (see technical characteristics).
- Use outside of the establishment, or on unprepared ground.
- Washing with large amounts of water, with a high-pressure jet or with steam and passage through a wash tunnel.
- Connection, by any means, of an electric device to the chair.
- Any use that does not fall within the operating context of a hospital chair.

All manipulations relating to the aforementioned precautions for use are described in paragraph II - Instructions for use [HERE](#).

4.1.2 CONTRAINDICATIONS

There are no noteworthy contraindications for this device.

4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE

The device must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The device must be regularly inspected, at least once per year.

To maintain the device in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed. See the maintenance check-list that can be downloaded from your customer area at <https://www.bmb-medical.com> or from your distributor.

In the event of observed anomalies or damage that could jeopardize device operation and safety, the equipment must be immediately withdrawn.

The device should only be used according to the instructions in this user manual.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

Work on electrical parts should only be performed by qualified and authorized staff.

* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

4.1.4 EMC-RELATED INFORMATION

The CAREXIA care chair is suitable for use in all healthcare establishments, with the exception of facilities in the vicinity of HF ELECTROSURGICAL APPLIANCES, or near an electromagnetically shielded room for an EM magnetic resonance imaging system, in which the intensity of EM INTERFERENCES is high.

The CAREXIA care chair is not suitable for use in domestic homes or establishments or on any premises connected to the public low-voltage network supplying domestic buildings.

EM interference could cause the product to make unexpected movements.

This appliance should not be used next to, or stacked on other appliances as this could lead to incorrect operation. If this mode of operation is necessary, this and the other appliances should be observed to ensure their correct operation.

The use of accessories, transducers or cables other than those specified or supplied by the manufacturer can lead to an increase in electromagnetic emissions, or a decrease in this appliance's immunity, causing inappropriate operation.

Do not use portable RF communication devices (including peripherals such as antenna cables and external antennas) less than 30 cm (12 inches) from any part of the APPLIANCE, including the cables specified by the manufacturer. Otherwise, the performance of these devices could be impaired.

List of cables used:

4. Power cord, maximum length 150 cm.
5. Variable height column, maximum length 27 cm
6. Y cable for remote controls/motors, maximum length 25 cm

EMISSIONS	
Conducted and radiated RF emissions:	CISPR11 Group 1 Class A
Harmonic distortions IEC 61000-3-2:	Class A
Voltage fluctuations and Flicker IEC 61000-3-3:	Compliant

Immunity test	Test level IEC 60601	Compliance level																																	
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV in contact ± 2/4/8/15 kV in air	± 8 kV in contact ± 2/4/8/15 kV in air																																	
RF radiated electromagnetic fields IEC 61000-4-3	3 V/m (professional healthcare establishment) 80 MHz-2.7 GHz 80% AM at 1 kHz	10 V/m 80MHz-2.7GHz 80% AM at 1 kHz																																	
Proximity fields emitted by wireless RF communication devices IEC 61000-4-3	<table><tr><th>Frequency (MHz)</th><th>Modulation</th><th>Required level (V/m)</th><th>Compliance level (V/m)</th></tr><tr><td>385</td><td>Pulsed modulation: 18 Hz</td><td>27</td><td>27</td></tr><tr><td>450</td><td>Pulsed modulation: 18 Hz</td><td>28</td><td>28</td></tr><tr><td>710-745-780</td><td>Pulsed modulation: 217 Hz</td><td>9</td><td>9</td></tr><tr><td>810-871-930</td><td>Pulsed modulation: 18 Hz</td><td>28</td><td>28</td></tr><tr><td>1720-1845-1970</td><td>Pulsed modulation: 217 Hz</td><td>28</td><td>28</td></tr><tr><td>2450</td><td>Pulsed modulation: 217 Hz</td><td>28</td><td>28</td></tr><tr><td>5240-5500-5785</td><td>Pulsed modulation: 217 Hz</td><td>9</td><td>9</td></tr></table>			Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)	385	Pulsed modulation: 18 Hz	27	27	450	Pulsed modulation: 18 Hz	28	28	710-745-780	Pulsed modulation: 217 Hz	9	9	810-871-930	Pulsed modulation: 18 Hz	28	28	1720-1845-1970	Pulsed modulation: 217 Hz	28	28	2450	Pulsed modulation: 217 Hz	28	28	5240-5500-5785	Pulsed modulation: 217 Hz	9	9
Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)																																
385	Pulsed modulation: 18 Hz	27	27																																
450	Pulsed modulation: 18 Hz	28	28																																
710-745-780	Pulsed modulation: 217 Hz	9	9																																
810-871-930	Pulsed modulation: 18 Hz	28	28																																
1720-1845-1970	Pulsed modulation: 217 Hz	28	28																																
2450	Pulsed modulation: 217 Hz	28	28																																
5240-5500-5785	Pulsed modulation: 217 Hz	9	9																																
Rapid / burst electrical spikes IEC 61000-4-4	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz																																	
Shock waves IEC 61000-4-5	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV																																	
Conducted interference, induced by RF fields IEC 61000-4-6	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz																																	
Power frequency magnetic fields IEC 61000-4-8	30A/m	30A/m																																	
Voltage Dips and Interruptions: IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles																																	

NB: There is no deviation from the requirements of IEC 60601-1-2 ed. 4

NB: There are no other known measures for maintaining basic safety with respect to EMC phenomena.

4.2 TRANSPORT AND STORAGE

Storage and transport conditions:

Ambient conditions	Operation
Ambient temperature	-10 to +50 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

4.2.1 TRANSPORT

During chair transport with a view to its storage, the following is recommended to prevent any impacts or falls:

- Strap it to a surface of the means of transport
- Cover it
- Engage the brakes
- Set it to its low position

Damage caused during transport is not covered by the warranty. The costs of repair and replacement of damaged parts will be borne by the customer.

4.2.2 STORAGE

When the chair is in storage, ensure that:

- Its brakes are engaged
- It is covered

The product should be stored in a clean, ventilated, temperate and dry place, protected from the elements and from direct sunlight.

During storage, no heavy weights should be placed on the chair.

The chair must not be considered as a support surface for any other type of equipment.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

In the event of prolonged storage, the battery should be recharged every 6 months, proceeding as follows:

- Connect the mains power cord.
- Briefly press one of the remote control functions to initiate charging.
- Allow the battery to charge for 8 hours.

4.3 WARRANTY

Warranty: two years parts and labor according to the following terms:

BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.

This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL stretchers are designed for a projected service life of 10 years under normal conditions of use, with appropriate periodic maintenance as described in the maintenance check-list of each stretcher. BMB MEDICAL provides the original buyer with the guarantee that the welds on the stretchers are free from structural defects for the projected stretcher service life of 10 years, as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use (Warning paragraph [here](#))
- Work without prior notice to our technical services, or maintenance on:
 - A mechanical part of the chair
 - A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to observe the preventive maintenance operations recommended in the maintenance check-list.

4.4 AFTER-SALES SERVICE

For any information concerning product installation, use and maintenance, the technical documentation (wiring diagrams, list of components, maintenance procedures, etc.), or any other questions, please contact the After-Sales Service at the following address: **sav@bmb-medical.com**, or by phone at **+33 (0)4 74 08 71 71**, or your local distributor for overseas purchases.

From your “Customer Area” on the website www.bmb-medical.com, you can access your equipment, submit work requests and track these requests directly via this secure interface.

The operator must carefully read the information contained in this user manual, in particular the information concerning the precautions relating to safety and installation/use methods.

In the event of doubt concerning the correct interpretation of the instructions, do not hesitate to contact the After-Sales Service whose staff are always available.

To facilitate ordering spare parts, you can request the exploded diagrams of your products from the After-Sales Service. You can thus obtain the references of the desired parts.

4.5 ENVIRONMENT

4.5.1 SPENT EQUIPMENT

In Europe:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution.

You can obtain the appropriate disposal location from the authorities in charge of environmental protection.

France specificities: BMB MEDICAL fulfills its obligations concerning the end of life of CAREXIA range products that it markets by funding the Recylum recycling channel dedicated to Pro WEEE that collects them free of charge (more information at www.recylum.com).

Outside of Europe:

Dispose of the product or its components in accordance with local laws and regulations.

Contact the local distributor for information concerning the return and/or collection systems available in the country.

Tél: +33 (0)4 74 08 71 71

Fax: +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001
01604 TREVOUX Cedex
FRANCE