CARVI

MANUEL D'UTILISATION
USER MANUAL

 **BMB
MEDICAL**

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FRANÇAIS

BMB MEDICAL

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

The purpose of this manual is to inform you how to use, maintain and ensure the safety of your new medical device. We greatly advise that you read this document carefully and that you follow the recommendations for an optimum use in safe conditions.

BMB MEDICAL offers a wide variety of options for the stretcher-trolley referred to in this Operating Manual. As a result, some of the functionalities described herein may not be applicable to your stretcher-trolley.

Some of the presented illustrations are for indication purposes and are not an exact representation of the device that you have in your possession.

The information contained in this document can be modified without warning and BMB MEDICAL shall not be held responsible in any way. Any representation or reproduction, whether total or partial, carried out without BMB MEDICAL's authorisation is illicit. This illicit representation or reproduction, by any process, would be considered as counterfeiting. The illustrations included in this document are not contractual.

CONTENTS

1	SPECIFICATION INTRODUCTION	42
1.1	PRODUCT DESCRIPTION.....	42
1.1.1	CARVI STRETCHER TROLLEY.....	43
1.1.2	CARVI KIDS STRETCHER TROLLEY.....	45
1.2	SCOPE OF APPLICATION	47
1.3	PUPPOSE OF THE OPERATING MANUAL.....	47
1.4	USERS	47
1.5	STORAGE.....	47
1.6	APPLIED PARTS	47
1.7	SYMBOLS AND STANDARDS EMPLOYED.....	47
2	OPERATING INSTRUCTIONS.....	50
2.1	OPERATING CONDITIONS.....	50
2.2	INSTALLATION PROCEDURE	50
2.3	BATTERY (optional electric backrest and/or kneegatch)	50
2.4	PRODUCT POWER SUPPLY.....	51
2.5	ADJUSTABLE HEIGHT.....	52
2.5.1	RAISING THE BED.....	52
2.5.2	LOWERING THE BED.....	52
2.6	BACKREST LIFTING SYSTEM.....	53
2.7	SIDERAILS.....	53
2.8	BRAKING AND GUIDING PEDAL	54
2.9	BACKWARD / FORWARD TILTING SYSTEM	55
2.10	EMERGENCY POSITION.....	56
2.11	OPTIONS	56
2.11.1	RETRACTABLE PUSH BARS.....	56
2.11.2	IV POLE	57
2.11.3	PEDIATRIC SIDERAILS.....	58
2.11.4	HEAD SIDE PAEDIATRIC KIT.....	59
2.11.5	VERTICAL OXYGEN BOTTLE SUPPORT SYSTEM.....	60
2.11.6	SCOPE SUPPORT TRAY.....	61
2.11.7	X-RAY CASSETTE BELOW THE PATIENT SURFACE	62
2.11.8	CHART HOLDER BELOW THE HEADREST.....	62
2.11.9	ATTACHMENT HOOKS FOR BAGS ON BOTH SIDES BELOW THE PATIENT SURFACE.....	63

2.11.10	PAPER ROLL DISPENSER	63
2.11.11	ANTISTATIC SYSTEM	63
2.11.12	KNEE GATCH	64
.....		64
3	CLEANING AND DISINFECTION	65
3.1	INTRODUCTION	65
3.2	FREQUENCY	65
3.3	CLEANING AND DISINFECTION ADVICE.....	65
3.3.1	CLEANING OPERATION	65
3.3.2	RECOMMENDED EQUIPMENT AND PRODUCTS.....	66
3.4	DEVICE-RELATED PRECAUTIONS	68
3.4.1	WARNINGS	68
3.4.2	CONTRAINDICATIONS	70
3.4.3	PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE	71
3.4.4	EMC-RELATED INFORMATION.....	71
3.5	TRANSPORT AND STORAGE	73
3.5.1	TRANSPORT	73
3.5.2	STORAGE	74
3.6	WARRANTY	74
3.7	AFTER SALES SERVICE	75
3.8	ENVIRONMENT	76
3.8.1	SPENT EQUIPMENT.....	76

1 SPECIFICATION INTRODUCTION

1.1 PRODUCT DESCRIPTION

Variable-height patient support device on hydraulic columns and equipped with an electric-assisted actuator system that allows, depending on the options, the movements of the backrest, leg lifter and knee gatch. The equipment is mounted on casters.

The device allows an operator to convey the patient within a healthcare establishment. The 5th wheel option allows you to switch from straight in-line guidance during transport to the ability to pivot the product around a corner. The chair is braked by a mechanical system actuated by pressing a pedal.

The device is equipped with a hydraulic or electrically assisted (optional) actuator system operated by a remote control to ensure the movements of the following elements:

- backrest
- legrest
- knee gatch (optional)

The products are designed to ensure that there are no unacceptable risks (essential performance) associated with patient support or caused by accidental movements.

List of CARVI references :

- **CARVI 60** HYDRAULIC HEIGHT-ADJUSTABLE STRETCHER WITH 60 CM WIDE MATTRESS
- **CARVI 70** HYDRAULIC HEIGHT-ADJUSTABLE STRETCHER WITH 70 CM WIDE MATTRES
- **CARVI 80** HYDRAULIC HEIGHT-ADJUSTABLE STRETCHER WITH 80 CM WIDE MATTRES
- **CARVI KID** PEDIATRIC STRETCHER WITH HYDRAULIC HEIGHT ADJUSTMENT

1.1.1 CARVI STRETCHER TROLLEY

TECHNICAL CHARACTERISTICS	
REFERENCE	CARVI
DESIGNATION	HEIGHT-ADJUSTABLE HYDRAULIC STRETCHER-TROLLEY
SAFE OPERATING LOAD	250 kg (maximum acceptable load including the patient's total weight, mattress and fittings)
WEIGHT	90 kg (CARVI 60), 94 kg (CARVI 70), 98 kg (CARVI 80)
OVERALL DIMENSIONS	68 x 216 cm or 78 x 216 or 88 x 216 cm
FORWARD /BACKWARD	18 °
BACKREST	0 to 90 °
EQUIPMENT	2-sectional resin mattress base: 60 x 197 cm (CARVI 60), 70 x 197 cm (CARVI 70) or 80 x 197 cm (CARVI 80). Radiolucent material. Height adjustment range: from 62 to 94 cm (including 8cm thick mattress). Trendelenburg - Anti-Trendelenburg +/- 18°. Two-sided height adjustment and Trendelenburg / Anti-Trendelenburg adjustment. Removable mattress made of high quality and resilience foam (35kg/m³), thickness : 8 cm. Removable cover. PVC coated Jersey. Colour : grey. 6 holders for IV pole. Bumpers on all 4 corners. Backrest adjustable by 2 gas springs (0 - 90°). Castors with central braking system Ø 200 mm at foot-end. 2 directional wheels at the head-end (without 5th directional wheel). Drop-down side rails, 2 horizontal bars. Anodized aluminium. Foldable push bars at the foot-end 2 parts. Frame cover including utility tray and O₂ tank holder. Maximum patient weight : 185kg.
AVAILABLE OPTIONS	5th directional wheel. Reinforced castors. Anti-static system. Side rails 3 horizontal bars. Bed panel instead of push bars. Permanent, foldable 3-stage "comfort" IV pole, 2 hooks. Permanent, foldable 3-stage "comfort" IV pole, 4 hooks. Removable 2-stage IV pole, 2 hooks. Removable IV pole, 2-stage, 4 hooks. Foldable push bars head-end. Foldable push bar foot-end 1 part. X-ray cassette below the patient surface (double hard surface). Electric backrest with manual release. Electric backrest and knee gatch, manual release. Battery. Patient's remote control. Knee gatch. Oxygen bottle holder. Paper roll dispenser at the head-end.

	Chart holder (for paperwork) below the headrest. Attachment hooks for bags on both sides below the patient surface. Monitor tray at the foot-end. 2 locations available (1 tray in width 60cm / 1 or 2 trays in width 70 and 80 cm). Comfort mattress, high quality foam (35 kg/m³), height : 10 cm. High Comfort mattress, high quality foam (35 kg/m³) + 2 cm thick viscoelastic layer (total height : 12 cm). Equipotential connection. Rail DIN. Fixation clamp for armrest. Adjustable fixation clamp for armrest. Anatomic shaped armrest. Side rail padded cushions (pair). "Universal Meal tray. Extensible meal tray for side rails: from 65 to 88 cm. Fall protection side.
--	--

1.1.2 CARVI KIDS STRETCHER TROLLEY

TECHNICAL CHARACTERISTICS	
REFERENCE	CARVI KIDS
DESIGNATION	HEIGHT-ADJUSTABLE PAEDIATRIC HYDRAULIC STRETCHER-TROLLEY
SAFE OPERATING LOAD	250 kg (maximum acceptable load including the patient's total weight, mattress and fittings)
WEIGHT	85 kg
OVERALL DIMENSIONS	69 x 200 cm
FORWARD /BACKWARD	18 °
BACKREST	0 à 90 °
EQUIPMENT	2-sectional resin mattress base: 60 x 185 cm. Radiolucent material. Overall dimensions: 69 x 200 cm. Height adjustment: from 62 to 94 cm (including 8 cm thick mattress). Trendelenburg / Anti-Trendelenburg +/- 18°. Two-sided height adjustment and Trendelenburg / Anti-Trendelenburg adjustment. Removable mattress made of high quality and resilience foam (35kg/m³), height : 8 cm. Removable cover. PVC coated Jersey. Colour : grey. 6 holders for IV pole. Bumpers on all 4 corners. Backrest adjustable by 2 gas springs (0 - 90°). Castors with central braking system Ø 200 mm at foot-end. 2 directional wheels at the head-end (without 5th directional wheel). Drop-down side rails, 2 horizontal bars. Anodized aluminium. Foldable push bars at the foot-end 2 parts. Frame cover including utility tray and O₂ tank holder. Maximum patient weight : 185kg
AVAILABLE OPTIONS	5th directional wheel. Anti-static system. Side rails 3 horizontal bars. Bed panel instead of push bars. Polycarbonate paediatric side rails. Paediatric kit head-end. Permanent, foldable 3-stage "comfort" IV pole, 2 hooks. Permanent, foldable 3-stage "comfort" IV pole, 4 hooks. Removable 2-stage IV pole, 2 hooks. Removable IV pole, 2-stage, 4 hooks. Foldable push bars head-end. Foldable push bar foot-end 1 part. Paper roll dispenser at the head-end. Chart holder (for paperwork) below the headrest. Attachment hooks for bags on both sides below the patient surface. Monitor tray at the foot-end.

	Comfort mattress, high quality foam (35 kg/m³), height : 10 cm. Equipotential connection. Lateral protection cushion. Head protection cushion. Universal Meal tray. Extensible meal tray for side rails: from 65 to 88 cm. Fall protection side.
--	--

1.2 SCOPE OF APPLICATION

This stretcher-trolley is intended to be used, according to its specifications, options and fittings, for carrying casualties and performing the day care and emergency services provided in a hospital.

It is used for the transport, welcoming, transfer, rest and recovery of a patient.

1.3 PURPOSE OF THE OPERATING MANUAL

Healthcare professionals, this Operating Manual has been prepared for you in aim of providing all of the necessary information and all of the precautions to be implemented for a safe and independent use.

The instructions concerning the maintenance and monitoring of the product ensure its service life.

It is therefore essential that all of the personnel read this manual and be informed of the handling procedures and risks related to the use of the stretcher-trolley. It shall be handled by a qualified person, who is aware of the risks related to a bad use.

1.4 USERS

Healthcare professionals: physicians, caregivers and agents approved by Technical Services who are accustomed to working with similar equipment can operate all functions.

In use cases where the patient is the operator (**use of the remote control only**), the caregivers must first ascertain the patient's ability to use the product. The caregivers must also teach the patient how to safely use the product: adjustment to low position to step off the product.

1.5 STORAGE

This Operating Manual shall be kept near the product in a location provided for this purpose. It shall be kept away from any element and substance compromising its legibility.

1.6 APPLIED PARTS

The applied parts are as follows:

- mattress
- Side rails
- Remote controls

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	Type B appliance
	Electrical class II equipment
	Electrical hazard, do not open
	Duty cycle (motor cycle) 2 min. operation / 18 min. stop
IPX4	Splash-resistant
	See the manual / instruction brochure
	Indicates the maximum patient weight and safe working load
	Conforms to the essential requirements of directive 93/42/EEC on medical devices.
	Serial number
	Manufacture date
	Manufacturer
	Alternating current
	Maximum patient weight
	Electronic waste: this symbol indicates that the product must not be discarded with unsorted municipal waste, but rather to a selective waste collection channel. Contact the local distributor for information concerning the return and/or collection systems available in the country.
	Protéger de l'humidité / keep dry
	Haut, tenir debout / stand up
	Fragile / fragile
	Recyclable / recycling
	Temperature range
	Humidity limit

General safety sign:

- Information highlighted by this symbol must be read and scrupulously observed.
- Indicates that failure to observe the guideline could endanger the patient or user or damage the equipment.

	Sitting on the backrest is prohibited
	Sitting on the legrest is prohibited

“WARNING: This appliance must not be altered without the manufacturer’s written consent.”

“Caution: only purpose-built accessories provided for this MD by BMB MEDICAL may be used.”

2 OPERATING INSTRUCTIONS

2.1 OPERATING CONDITIONS

AMBIENT CONDITIONS	OPERATION
Ambient temperature	+5 °C to 30°C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

2.2 INSTALLATION PROCEDURE

Before commissioning the chair and depending on the model, check that the following are operational:

- [The product power supply](#)
- [The height adjustment](#)
- [Backrest](#)
- [Side rails](#)
- [Brake pedal / steer pedal](#)
- [Trendelenburg / reverse Trendelenburg](#)
- [CPR release mechanism](#)

Clean the product according to your establishment's internal protocol.

2.3 BATTERY (OPTIONAL ELECTRIC BACKREST AND/OR KNEEGATCH)

The battery unit is fitted with its own charger. Before initial use, leave the product to stand for 12 h in a temperate environment (to allow the electrical and electronic components to reach normal operating temperature), then charge the battery for 12 hours from a mains outlet.

In routine use, we strongly recommend keeping the chair plugged into the mains during sedentary use.

The battery, with its limited operating life, can be used to operate the stretcher-chair in the event of a power outage, or during transfer.

Once the manipulations are complete, remember to re-connect the stretcher-chair to the mains (at least 12 hours per day). Never allow the battery to drain fully (risk of reduced battery performance).

“Warning: Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact. »

Battery characteristics		Battery safety symbols	
Battery type:	Lead	IPX6	Protection rating
Capacity:	1,3 Ah		CE marking
Input supply voltage:	xxxx		Indoor use
Input supply voltage:	24 Vcc, 10A		Electronic waste

2.4 PRODUCT POWER SUPPLY

Power grid connection conditions and electromagnetic compatibility:

This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.

Position the product such that it can be easily unplugged.

Unplug the product before moving it.

Do not leave the cables in contact with the floor or the casters,

Do not knot the cables,

Periodically check the condition of the cables.

Mains connection:	120 ~ 240V depending on the country.
Frequency:	50 /60 Hz
Rated current:	3.3 A
Protection rating / degree of electric shock protection:	Class 2 appliance
Power input:	525 W
Intermittent operating mode:	10% i.e. 1 minutes of continuous operation followed by 19 minutes of rest.

The mains plug is the electrical grid insulation device.

Work on electrical parts should only be performed by qualified and authorized staff.

Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the [EMC information here](#).

2.5 ADJUSTABLE HEIGHT


Before using these functions:

- Make sure that no object nor any person is placed under the lying surface.
- Make sure that nobody has placed his/her foot under the different pedals (adjustable height, backward / forward tilt)

2.5.1 RAISING THE BED

The controls (grey pedals) are located on either side of the lying surface. The bed is raised by pumping the pedal:

2.5.2 LOWERING THE BED

The bed is lowered by pressing **simultaneously** on the two backward (yellow pedal) and forward (blue pedal) tilting pedals.

Make sure that your foot is placed so as to apply a homogeneous pressure on the 2 pedals.


2.6 BACKREST LIFTING SYSTEM

Before operating the backrest lifting system, check that nothing is blocking it (fittings, cables, etc.).

The backrest lifting system, assisted by 2 gas springs, is controlled by operating the bar located under the lying surface: apply a pressure (with 2 hands to make it easier) to unblock the cylinders and adjust the backrest lifting system at the chosen height.

Case of an electric trolley:

The backrest lifting system is adjusted with the remote control(s):

- Foot side remote control
- And patient remote control (optional)

2.7 SIDERAILS


The siderails are provided to reduce the risk of the patient accidentally falling from the stretcher. They are in no way designed to immobilise or maintain the patient with strap-type immobilisation devices.

Before handling the bars, make sure that nothing is blocking their movement. Do not use the bars to move the stretcher-trolley.


Lowering the bars:

To unlock, maintain the button pressed down ❶, and push the bar ❷.

Lifting the bars :

Lift the top bar until the button is locked.

2.8 BRAKING AND GUIDING PEDAL


Except during transport, make sure that you always apply the stretcher-trolley's brake. Check that the brake is active by trying to move the stretcher-trolley.

Do not park the stretcher-trolley on a slope.

The braking and guiding pedal is used for 3 functions:


- Locked stretcher-trolley function- pedal in low position ❶.
- Unlocked stretcher-trolley function, which is used to move it in several directions on 4 wheels (5th wheel is deactivated) - pedal in median position ❷.
- Stretcher-trolley in directional guiding mode which makes it easier to move the trolley in a straight line - side green pedal tilted downwards ❸.

2.9 BACKWARD / FORWARD TILTING SYSTEM


Backward tilt position:

In order to place the stretcher-trolley in the backward tilt position (stretcher head in low position), press on the yellow pedal.

Forward tilt position:

In order to place the stretcher-trolley in the forward tilt position (stretcher foot in low position), press on the blue pedal.

2.10 EMERGENCY POSITION

If it is necessary to rapidly lower the backrest, pull the red handle and follow the backrest by pulling it down.


If you use the emergency manual unlocking system, lower the backrest by using the remote control until it reaches its lowest position in order to reset its operation.

2.11 OPTIONS

2.11.1 RETRACTABLE PUSH BARS

In order to retract the handles, lift the locking ring up to the stop ②, the handles are then free to rotate downwards, away from the lying surface ③.

In order to reposition the handles, lift them and replace the ring in its low position. Check that they are positioned correctly.

2.11.2 IV POLE

Safe working load: 6 kg.

This indicative value corresponds to safe product use, but is not representative of normal use for which the load is significantly lower.

Do not affix anything to the telescopic parts of the IV pole (e.g. tape) that could hinder its movement.




When moving the stretcher fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.

Do not move the stretcher via the IV pole.

2.11.2.1 PERMANENTLY ATTACHED IV POLE

Installation :


1. Remove the IV pole from its storage location ❶ and raise into the vertical position until it clicks into place ❷.
2. To raise the height of the pole, (3 telescopic stages), pull up to either of the two upper stages until it locks into place ❸.

Storing the IV pole:


In order to fold the telescopic pole back into place, pull up the red locking bolts located below each stage.

Do not tilt the IV pole in its highest position. Lower to each of the various stages in succession until you reach the lowest position (one level) before inserting the pole back into its storage location.

2.11.2.2 REMOVABLE IV POLE

Insert the detachable IV pole into one of the corner locations provided for this purpose.


2.11.3 PEDIATRIC SIDERAILS


The pediatric bars are operated such as the bars described [here](#).

2.11.4 HEAD SIDE PAEDIATRIC KIT


In order to open the lateral panel, lift the red indexing finger ① and push the panel outward ②

In order to remove the head panel, pull the 2 red unlocking plates up to the stop ①. You can then free the head panel by removing it from its location ②


In order to replace the head panel in its position, make sure that the plates are at the stop by pulling them ①. Then insert the head panel into the slots provided for this purpose ②. And lock the system by pushing the plates to the stop ③.


Check that the head panel is secure by attempting to remove it.

2.11.5 VERTICAL OXYGEN BOTTLE SUPPORT SYSTEM

Rotate the oxygen bottle support system away from the lying surface.

Retracted bottle support system

Extended bottle support system

2.11.6 SCOPE SUPPORT TRAY

1 Insert the tray(s) into the slot(s) provided (right and/or left at the foot of the stretcher).	
2 Press the shutter release below the tray to unlock it	
3 Fold down the tray	
4 To use the tray, pull up it until it locks.	
5 Lift the tray to remove it from the holder	


When using Trendelenburg and reverse-Trendelenburg positioning, it is recommended that devices positioned on the fold-out monitor tray are attached using the strap provided for this purpose.


2.11.7 X-RAY CASSETTE BELOW THE PATIENT SURFACE


2.11.8 CHART HOLDER BELOW THE HEADREST

2.11.9 ATTACHMENT HOOKS FOR BAGS ON BOTH SIDES BELOW THE PATIENT SURFACE

2.11.10 PAPER ROLL DISPENSER


2.11.11 ANTISTATIC SYSTEM

Make sure the earthing tape touches the floor, is not bended and is not deformed.


2.11.12 KNEE GATCH

In order to place the stretcher in a knee bending position, pull on the yellow actuator ❶ while lifting the handle ❷ to the desired position.


In order to make the lying surface flat, pull on the yellow actuator ❶ while pressing the handle ❷.

Electric trolley option :

On the electric versions of the chair-stretchers, the knee bending position is adjusted with the remote control(s):

- Foot-side remote control
- and patient remote control (optional)

3 CLEANING AND DISINFECTION

3.1 INTRODUCTION

Regular and correct cleaning ensure the product's life service (see cleaning program prepared upon commissioning of the device by appointing the reference personnel).

If the cleaning operations are not carried out, this may cause a risk of infection due to the presence of secretions and/or residues on the stretcher-trolley.

The following cleaning recommendations do not replace your establishment's internal protocol.

3.2 FREQUENCY

The cleaning shall be performed:

- After each patient i.e. between two patients.
- After transporting an infected patient.
- Entirely on a weekly basis

3.3 CLEANING AND DISINFECTION ADVICE

3.3.1 CLEANING OPERATION

Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that cause stains (e.g.: "blue product" or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS

Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.

- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
 - PVC, ABS, Polyamide, Polyurethane and Polypropylene
 - Epoxy paint-coated metal surfaces
 - Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

List of strictly prohibited products:

Solvents:	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
Bases:	10% Sodium hydroxide / 10% Ammonia
Acids:	10% Nitric / 30% Hydrochloric / 10% Citric
Miscellaneous:	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

High-pressure and steam cleaning are prohibited.

Any damage resulting from failure to observe this instruction for use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.

3.4 DEVICE-RELATED PRECAUTIONS

3.4.1 WARNINGS

This MD has been designed for and adapted to medical practices and procedures, and hence a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the common sense of the stretcher's users.

- Do not use the product without first having familiarized yourself with this user manual.
- Do not alter the stretcher. Any changes made to the system could lead to unpredictable operation that could in turn jeopardize the safety of the user or patient. This would also void the stretcher's warranty.
- Used to accommodate a person of small stature (less than 1.30 m tall), or under the age of 12 years without adapted equipment (pediatric kit)
- To ensure correct use of the product, only proceed if the weight is correctly distributed over the device.
- The chassis cover load capacity is of 35 kg. Do not sit or stand on the cover.
- Only purpose-built accessories provided for this MD by BMB MEDICAL may be used
- Connection, by any means, of an electric device to the stretcher.
- The product must be connected to a mains power outlet compliant with standards in force.
- Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact.
- In case of prolonged interruption of use of your equipment, it is recommended to recharge the battery once a month. After more than 50 days of storage without recharging, battery performance may be impaired or even irreversibly damaged.
- This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.
- Position the product such that it can be easily unplugged.
- Unplug the product before moving it.
- Do not leave the cables in contact with the floor or the casters.
- Do not knot the cables.
- Periodically check the condition of the cables.
- Work on electrical parts should only be performed by qualified and authorized staff.
- Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the [EMC information here](#).
- When using the remote control, all movements can be interrupted by simply releasing the corresponding key.

- Only activate the various positions after having ensured that no-one is touching any mechanical parts and that there are no objects under the seat or patient surface, or between any of the moving structures.
- The side rails are intended to minimize the risk of the patient accidentally falling. They are not, in any way, designed to immobilize or maintain the patient.
- Before manipulating the side rails, ensure that nothing could hinder their movement.
- When lowering the side rail, take care not to trap the hands or feet of the patient or staff
- If the side rail locking operation is silent (no “click”), please revise the side rails.
- Except during transport, always ensure that the stretcher's brake is engaged. Check whether the brake is engaged by attempting to move the stretcher.
- The “variable height” function should be used when handling the patient in order to work at an ergonomic height.
- Make sure that no one has their feet under the different pedals (variable height, Trendelenburg, Reverse Trendelenburg).
- It is recommended to attach the devices to the monitor support plate with the strap provided for this purpose when positioning the stretcher in the Trendelenburg or Reverse Trendelenburg position.
- When using the manual emergency release, do not put the arm in the frame of the lying surface.
- When activating the knee gatch, to avoid any risk of pinching, do not place your arm or hand between the lying surface and the metal structure.
- Do not park the chair on a sloping floor.
- Use outside of the establishment, or on unprepared ground.
- Do not use the MD on a floor with a slope of more than 10%.
- Use of the chair with a load that exceeds the safe working load (see technical characteristics).
- Always ensure that the headrest is correctly locked in position.
- The product should be set to sitting position before changing the headrest.
- The weight of the patient's head rests on the headrest; it must thus be supported by the operator during adjustment. Failure to observe this recommendation could lead to patient injury.
- To avoid any risk of pinching during headrest adjustment, keep your fingers away from all articulated areas.
- Do not place your hand between the backrest and headrest to pull the unlock handle as you may injure your fingers.
- Always ensure that the headrest articulation is correctly locked.
- IV pole safe working load: 6 kg.
- When moving the chair fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.
- Do not move the chair via the IV pole
- In the case of an unsupervised patient, the product should be left in its lowest position to reduce the risk of falling.
- Tell the patient not to sit on the ends (foot or head) of the chair/stretcher. An excess weight applied to one of the ends could tip the patient surface and injure the patient.

- 
- Tell the patient and any accompanying persons not to sit on the armrests.
 - The equipotential connector must not, under any circumstances, be used as a protective earth connector
 - Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.
 - Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.
 - High-pressure and steam cleaning are prohibited.
 - Any use that does not fall within the operating context of a hospital chair.
 - To prevent failure, electromedical equipment requires special [EMC-related precautions](#).
 - Never use the chair in the presence of an anesthesia mixture flammable with air, with oxygen, or with nitrous oxide as this could cause a fire and/or explosion.
 - Use by any non-caregiver who has not received the authorization and explanations from a qualified individual.

All manipulations relating to the aforementioned precautions for use are described in [paragraph II - Instructions for use](#).

3.4.2 CONTRAINDICATIONS

There are no noteworthy contraindications for this device.

3.4.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE

The device must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The device must be regularly inspected, at least once per year. See the maintenance check-list. To maintain the device in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed

In the event of observed anomalies or damage that could jeopardize device operation and safety, the equipment must be immediately withdrawn.

The device should only be used according to the instructions in this user manual.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

Work on electrical parts should only be performed by qualified and authorized staff.

* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

3.4.4 EMC-RELATED INFORMATION

The products in the CARVI range are suitable for use in all healthcare establishments, with the exception of facilities in the vicinity of HF ELECTROSURGICAL APPLIANCES, or near an electromagnetically shielded room for an EM magnetic resonance imaging system, in which the intensity of EM INTERFERENCES is high.

The products in the CARVI range are not suitable for use in domestic homes or establishments or on any premises connected to the public low-voltage network supplying domestic buildings.

EM interference could cause the product to make unexpected movements.

This appliance should not be used next to, or stacked on other appliances as this could lead to incorrect operation. If this mode of operation is necessary, this and the other appliances should be observed to ensure their correct operation.

The use of accessories, transducers or cables other than those specified or supplied by the manufacturer can lead to an increase in electromagnetic emissions, or a decrease in this appliance's immunity, causing inappropriate operation.

Do not use portable RF communication devices (including peripherals such as antenna cables and external antennas) less than 30 cm (12 inches) from any part of products in the CARVI range, including the cables specified by the manufacturer. Otherwise, the performance of these devices could be impaired.

The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11, class A). If used in a residential environment (for which compliance

with the CISPR 11 class B standard is generally required), this equipment may not provide appropriate protection to radio frequency communication services. The user may need to take attenuation steps, such as moving or re-orienting the equipment.

List of cables used:

1. Power cord, maximum length 1000 cm.
2. Variable height column, maximum length 750 mm
3. Backrest cable: maximum length 1250 mm
4. Legrest cable: maximum length 750 mm
5. Footrest cable: maximum length 1500 mm

EMISSIONS	
Conducted and radiated RF emissions:	Conducted and radiated RF emissions:
Harmonic distortions IEC 61000-3-2:	Harmonic distortions IEC 61000-3-2:
Voltage fluctuations and Flicker IEC 61000-3-3:	Voltage fluctuations and Flicker IEC 61000-3-3:

Immunity test	Test level IEC 60601	Compliance level																																	
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV in contact ± 2/4/8/15 kV in air	± 8 kV in contact ± 2/4/8/15 kV in air																																	
RF radiated electromagnetic fields IEC 61000-4-3	3 V/m (professional healthcare establishment) 80MHz-2.7GHz 80% AM at 1 kHz	3V/m 80MHz-2.7GHz 80% AM at 1 kHz																																	
Proximity fields emitted by wireless RF communication devices IEC 61000-4-3	<table><tr><th>Frequency (MHz)</th><th>Modulation</th><th>Required level (V/m)</th><th>Compliance level (V/m)</th></tr><tr><td>385</td><td>Pulsed modulation: 18 Hz</td><td>27</td><td>27</td></tr><tr><td>450</td><td>Pulsed modulation: 18 Hz</td><td>28</td><td>28</td></tr><tr><td>710-745-780</td><td>Pulsed modulation: 217 Hz</td><td>9</td><td>9</td></tr><tr><td>810-870-930</td><td>Pulsed modulation: 18 Hz</td><td>28</td><td>28</td></tr><tr><td>1720-1845-1970</td><td>Pulsed modulation: 217 Hz</td><td>28</td><td>28</td></tr><tr><td>2450</td><td>Pulsed modulation: 217 Hz</td><td>28</td><td>28</td></tr><tr><td>5240-5500-5785</td><td>Pulsed modulation: 217 Hz</td><td>9</td><td>9</td></tr></table>			Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)	385	Pulsed modulation: 18 Hz	27	27	450	Pulsed modulation: 18 Hz	28	28	710-745-780	Pulsed modulation: 217 Hz	9	9	810-870-930	Pulsed modulation: 18 Hz	28	28	1720-1845-1970	Pulsed modulation: 217 Hz	28	28	2450	Pulsed modulation: 217 Hz	28	28	5240-5500-5785	Pulsed modulation: 217 Hz	9	9
Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)																																
385	Pulsed modulation: 18 Hz	27	27																																
450	Pulsed modulation: 18 Hz	28	28																																
710-745-780	Pulsed modulation: 217 Hz	9	9																																
810-870-930	Pulsed modulation: 18 Hz	28	28																																
1720-1845-1970	Pulsed modulation: 217 Hz	28	28																																
2450	Pulsed modulation: 217 Hz	28	28																																
5240-5500-5785	Pulsed modulation: 217 Hz	9	9																																
Rapid / burst electrical spikes IEC 61000-4-4	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz																																	
Shock waves IEC 61000-4-5	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV																																	
Conducted interference, induced by RF fields IEC 61000-4-6	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz																																	
Power frequency magnetic fields IEC 61000-4-8	30A/m	30A/m																																	
Voltage Dips and Interruptions: IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles																																	

NB: There is no deviation from the requirements of IEC 60601-1-2 ed. 4

NB: There are no other known measures for maintaining basic safety with respect to EMC phenomena.

3.5 TRANSPORT AND STORAGE

Storage and transport conditions:

Ambient conditions	Operation
Ambient temperature	-10 to +30 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

3.5.1 TRANSPORT

Before transporting the stretcher with a view to storing to ensure that there are no risks of impact or accidental falls during transport, we recommend that you:

- Fasten stretcher to the frame
- Place the cover over it
- Lock the brakes
- Place it in the lowest position

Damage to the appliance caused during transport is not covered by the guarantee. Repairs or replacement of damaged parts are the responsibility of the client.

3.5.2 STORAGE

When the chair is in storage, ensure that:

- Its brakes are engaged
- It is covered

The product should be stored in a clean, ventilated, temperate and dry place, protected from the elements and from direct sunlight.

During storage, no heavy weights should be placed on the chair.

The chair must not be considered as a support surface for any other type of equipment.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

In the event of prolonged storage, the battery should be recharged every 6 months, proceeding as follows:

- Connect the mains power cord.
- Briefly press one of the remote control functions to initiate charging.
- Allow the battery to charge for 8 hours.

3.6 WARRANTY

Warranty: two years parts and labor according to the following terms:

BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.

This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL stretchers are designed for a projected service life of 10 years under normal conditions of use, with appropriate periodic maintenance as described in the maintenance check-list of each stretcher. BMB MEDICAL provides the original buyer with the guarantee that the welds on the stretchers are free from structural defects for the projected stretcher service life of 10 years, as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use ([Warning paragraph](#))
- Work without prior notice to our technical services, or maintenance on:
 - A mechanical part of the chair
 - A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to observe the preventive maintenance operations recommended in the maintenance check-list.

3.7 AFTER SALES SERVICE

For any information relating to the installation, use and maintenance of the product and for any other questions, please contact the After-Sales Service at the following address: sav@bmb-medical.com or on **+33 (0)4 74 08 71 71** or **your local distributor if outside of France**. Thanks to our “Customer Space”, you can access your equipment items on-line, view the maintenance points to be attended to and request our services directly via this secure interface.

The operator must carefully read the instructions contained in this manual, paying particular attention to the safety precautions and the procedures to be followed for installation and correct use.

In the case of any doubts as to the correct interpretation of the instructions, please contact the After-Sales Service for any necessary clarifications.

To facilitate the ordering of spare parts, exploded parts drawings can be requested from the After-Sales Service to provide you with the references of the parts required.

3.8 ENVIRONMENT

3.8.1 SPENT EQUIPMENT

In France:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution. For this purpose, BMB MEDICAL fulfills its obligations concerning the end of life of CARVI range products that it markets by funding the Recylum recycling channel dedicated to Pro WEEE that collects them free of charge (more information on www.recylum.com).

In Europe:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution.

Outside Europe :

Dispose of the Product or its components in accordance with local laws and regulations.

Contact the local distributor to inquire about return and/or collection systems available in the country.

Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001
01604 TREVOUX Cedex
FRANCE

www.bmb-medical.com