
CLAUJA

MANUEL D'UTILISATION
USER MANUAL

 **BMB
MEDICAL**

FRANÇAIS

ENGLISH

FRANCAIS

BMB MEDICAL

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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Fax: +33(0)4 74 08 71 72

www.bmb-medical.com

The purpose of this manual is to familiarize you with the use, maintenance and safety of your new medical device. We encourage you to read this document carefully and to follow the guidelines for optimum and safe use.

BMB MEDICAL offers a broad variety of options for the products covered by this user manual. Consequently, some features described in this manual may not apply to your chair.

Some illustrations are provided for information purposes and do not accurately represent the device in your possession.

The information in this document are subject to change without notice and in no way engage the responsibility of BMB MEDICAL. This document may not be reproduced, in whole or in part, without prior written consent by BMB MEDICAL. This illicit representation or reproduction, by any means, constitutes counterfeiting. The illustrations in this document are not binding.

CONTENTS

ENGLISH	53
1 INTRODUCTION SPECIFICATIONS	56
1.1 PRODUCT PRESENTATION	56
1.1.1 CLAVIA LSA DAY SURGERY CHAIR	57
1.1.2 CLAVIA LSA CHIR STRETCHER-BED, DAY SURGERY VERSION	59
1.1.3 CLAVIA FCD CHEMOTHERAPY AND HEMODIALYSIS STRETCHER CHAIR	61
1.1.4 CLAVIA FSP MULTIPURPOSE CARE CHAIR	63
1.2 PURPOSE OF THE USER MANUAL	64
1.3 USERS	64
1.4 STORAGE	64
1.5 APPLIED PARTS	64
1.6 SYMBOLS AND STANDARDS USED	65
2 INSTRUCTIONS FOR USE	67
2.1 OPERATING CONDITIONS	67
2.2 INSTALLATION PROCEDURE	67
2.3 BATTERY (optional for the fcd version)	67
2.4 PRODUCT POWER SUPPLY	68
2.5 REMOTE CONTROL	69
2.5.1 PROGRAMMED POSITIONS	70
2.5.2 TRANSPORT POSITION	70
2.5.3 RESET	71
2.6 TRENDLENBURG / REVERSE TRENDLENBURG	73
2.7 EMERGENCY POSITION	73
2.8 VARIABLE HEIGHT AND LOWERING	73
2.9 BACKREST	74
2.10 LEGREST	74
2.11 STANDARDS AND PEDIATRICS SIDE RAILS	75
2.12 FOLDABLE PUSH BARS	76
2.13 ARMRESTS	77
2.13.1 SWIVELING THE ARMRESTS	77
2.13.2 ADJUSTING ARMREST HEIGHT	77
2.14 BRAKE AND GUIDE PEDALS	78
2.15 HEADRESTS	79
2.15.1 FITTING THE HEADRESTS	79
OPTIONS	80
2.15.2 OPHTHALMOLOGY HEADREST (for the CLAVIA LSA CHIR version)	80
2.15.3 REMOVABLE SURGERY HEADREST FOR STOMATOLOGY AND ENT (for the CLAVIA LSA CHIR version)	81
2.15.4 TECHNICAL SIDE RAILS	81
2.15.5 TECHNICAL RAIL-MOUNTED HEADREST SUPPORT	82
2.15.6 COMFORT CUSHION	82

2.15.7	ANTI-STATIC SYSTEM.....	83
2.15.8	FOOT PEDAL-CONTROLLED HEIGHT ADJUSTMENT (CLAVIA LSA SURGERY).....	83
2.15.9	DIRECTIONAL 5th WHEEL.....	84
2.15.10	IV POLE.....	85
2.15.11	60 CM WIDE BACKREST (CLAVIA LSA CHIR VERSION).....	86
2.15.12	FOOTREST (for the FCD and FSP versions)	86
2.15.13	IPAD HOLDER (CLAVIA FCD)	87
2.15.14	PAPER ROLL HOLDER	88
2.15.15	CRYSTAL PROTECTION.....	88
2.15.16	FIXED PUSH BAR (CLAVIA FCD & FSP versions).....	89
2.15.17	CABLE REEL.....	89
2.15.18	OXYGEN CYLINDER HOLDER UNDER BACKREST	89
2.15.19	EQUIPOTENTIAL CONNECTOR.....	90
2.15.20	SURGICAL FIELD HOLDER (CLAVIA LSA SURGERY VERSION)	90
2.15.21	PEDIATRIC KIT (for the CLAVIA LSA) version)	91
2.15.22	PROTECTIVE CUSHIONS	92
2.15.23	MEAL TRAY(CLAVIA FCD).....	92
3	CLEANING AND DISINFECTION	93
3.1	INTRODUCTION	93
3.2	FREQUENCY	93
3.3	CLEANING AND DISINFECTION ADVICE.....	93
3.3.1	CLEANING OPERATION	93
3.3.2	RECOMMENDED EQUIPMENT AND PRODUCTS.....	94
4	PRECAUTIONS.....	95
4.1	DEVICE-RELATED PRECAUTIONS	95
4.1.1	WARNINGS	95
4.1.2	CONTRAINDICATIONS	97
4.1.3	PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE	97
4.1.4	EMC-RELATED INFORMATION.....	98
4.2	TRANSPORT AND STORAGE	99
4.2.1	TRANSPORT	100
4.2.2	STORAGE	100
4.3	WARRANTY	101
4.4	AFTER-SALES SERVICE	102
4.5	ENVIRONMENT	102
4.5.1	SPENT EQUIPMENT.....	102

1 INTRODUCTION SPECIFICATIONS

1.1 PRODUCT PRESENTATION

This patient support device consists of a patient surface fitted onto a chassis with wheels. The device is equipped with side bars or armrests and, optionally, allows the temporary fitting of IV poles.

The device allows an operator to convey the patient within a healthcare establishment. The 5th wheel option allows you to switch from straight in-line guidance during transport to the ability to pivot the product around a corner. The chair is braked by a mechanical system actuated by pressing a pedal.

The device is equipped with an electrically assisted actuator system actuated by a remote control to manage the movements of the following elements:

- backrest,
- legrest,
- footrest,
- of 2 electrical columns used to adjust the product's height.

The products are designed to ensure that there are no unacceptable risks (essential performance) associated with patient support or caused by accidental movements.

- List of CLAVIA 5 (5th generation) references:
 - **CLAVIA 5 FCD 65** = Chemo Dialysis chair with a 65 cm wide mattress
 - **CLAVIA 5 FCD 75** = Chemo Dialysis chair with a 75cm wide mattress
 - **CLAVIA 5 FSP 65** = Sampling Care chair with a 65 cm wide mattress
 - **CLAVIA 5 FSP 75** = Sampling Care chair with a 75cm wide mattress
 - **CLAVIA 5 LSA 65** = Day Surgery Stretcher chair with a 65 cm wide mattress
 - **CLAVIA 5 LSA 75** = Day Surgery Stretcher chair with a 75cm wide mattress
 - **CLAVIA 5 LSA CHIR 65** = Day Surgery Stretcher chair with a 65 cm wide mattress
 - **CLAVIA 5 LSA CHIR 75** = Day Surgery Stretcher chair with a 75cm wide mattress

1.1.1 CLAVIA LSA STRETCHER BED

TECHNICAL DATA	
REFERENCE	CLAVIA LSA
DESIGNATION	Stretcher-bed for day surgery with electrical motion
SAFE WORKING LOAD	250 kg (maximum permissible load, including the total weight of the patient of 185 kg, mattress and accessories)
WEIGHT	84 kg (CLAVIA 65), 86 kg (CLAVIA 75)
OVERALL DIMENSIONS	219 x 80, 90 cm (reclined) / 155 x 80, 90 cm (seated)
DROP-DOWN SIDE RAILS	110 x 28 cm (height above the mattress)
TREND	Electrically assisted 11°
BACKREST	0 to 90° - Electrically assisted, manual release
EQUIPMENT	3-sectional mattress base : 65 x 200 cm, 75 x 200 cm Removable mattress made of high resilience foam (35 kg/m³), thickness 10 cm Removable cover of PVC coated Jersey, colour: slate grey Electrically assisted height adjustment on 2 LINAK columns: 59 - 99 cm Electrically (LINAK) assisted Trendelenburg (15°) / Reverse Trendelenburg (6°) Electrically (LINAK) assisted backrest, manual release Electrically (LINAK) assisted leg rest 0 - 60° Digital actuators allowing 6 programmed positions: Welcome, Relax, Reverse Trend, Secured Trendelenburg position, CPR Machine-welded steel frame, epoxy coated Frame protection cover Double running wheels with central braking system Ø 150 mm including 5th directional wheel Drop-down laminated side guards. Length 110 x Height 27cm above mattress 3 positions side rails. High security position. Lox position for easy access. Surgery position (ENT, Stomatology, Hand surgery, colonoscopy) Bumpers on all 4 corners fitted with holder for accessories 2 holders on siderails for IV pole Foldable push bars below backrest Battery
AVAILABLE OPTIONS	Anti-static system Comfort pillow (blue) Removable IV pole, 2-stage, 2 hooks Holder for oxygen cylinder (incompatible with the option CLOEC) Cable reel (incompatible with the option CLOSO2) Foldable push bars head-end Foldable push bars foot-end Mattress protection at foot end High Comfort mattress, high quality foam (35kg/m³), 2cm thick viscoelastic layer (total height : 12cm) Paper roll dispenser at the head-end (incompatible with BPSRB) Attachment hooks for bags on both sides below the patient surface Equipotential connection Side rail padded cushions (pair) CLAVIA side rails cushion/armrest, removable and right/left reversible (single)

AVAILABLE OPTIONS

Polycarbonate pediatric side rails (instead of the standard side rails)
Pediatric kit head-end
Pre-equipment for pediatric kit
Lateral protection cushion for pediatric side rails
Head protection cushion (CLOKCL obligatory)

SCOPE

This chair is intended, depending on its specifications, options and accessories, for a hospital's day surgery services.

1.1.2 CLAVIA LSA CHIR STRETCHER-BED, DAY SURGERY VERSION

TECHNICAL DATA	
REFERENCE	CLAVIA LSA CHIR
DESIGNATION	Day surgery stretcher-bed
SAFE WORKING LOAD	200 kg (indicates the maximum permissible load, including the total weight of the patient, mattress and accessories)
WEIGHT	113 kg (CLAVIA 65), 116 kg (CLAVIA 75)
OVERALL DIMENSIONS	213 x 80, 90 cm (reclined) / 165 x 80, 90 cm (seated)
DROP-DOWN SIDE RAILS	100 x 28 cm (height above the mattress)
TREND	Electrically assisted 11°
BACKREST	0 to 90° - Electrically assisted, manual release
EQUIPMENT	3-sectional mattress base : 65 x 200 cm, 75 x 200 cm Removable mattress made of high resilience foam (35 kg/m³), thickness 10 cm Removable cover of PVC coated Jersey, colour: slate grey 2 section backrest with a removable head section to fasten a surgery head rest for ophthalmology, stomatology and ENT Electrically assisted height adjustment on 2 LINAK columns : 59 - 99 cm Electrically (LINAK) assisted Trendelenburg (15°) / Reverse Trendelenburg (6°) Electrically (LINAK) assisted backrest, manual release Electrically (LINAK) assisted leg rest 0 - 60° Digital actuators allowing 6 programmed positions: Welcome, Relax, Reverse Trend, Secured Trendelenburg position, CPR Heavy duty design. Steel frame with powder coating ABS bonnet Double running wheels with central braking system Ø 150 mm including 5th directional wheel Drop-down laminated side guards. Length 110 x Height 27cm above mattress 3 positions side rails. High security position. Lox position for easy access. Surgery position (ENT, Stomatology, Hand surgery, colonoscopy) Bumpers on 2 corners at foot end fitted with holder for accessories 2 holders on siderails for IV pole Technical rails on backrest at both sides Foldable push bars below backrest Battery
AVAILABLE OPTIONS	Anti-static system Comfort pillow (blue) Without the standard headrest (depreciation) Removable IV pole, 2-stage, 2 hooks Holder for oxygen cylinder (incompatible with the option CLOEC) Cable reel (incompatible with the option CLOSO2) Foldable push bars foot-end Articulated removable surgery head rest for ophthalmology + hollow cushion Articulated removable surgery head rest for stomatology + flat cushion Removable surgery head rest for stomatology and ENT Additional hollow cushion for articulated removable surgery headrest Additional flat cushion for articulated removable surgery headrest

	Head rest support. Fixation on technical rails Width of the backrest 60 cm for CLAVIA LSA CHIR 70 only. Equipped with DIN rails Holder for operating field, with connection for oxygen Mattress protection at foot end Attachment hooks for bags on both sides below the patient surface Equipotential connection Foot pedal control for variable height (in addition of the wired remote control) CLAVIA side rails cushion/armrest, removable and right/left reversible (single)
--	--

SCOPE

This stretcher is intended, depending on its specifications, options and accessories, for a hospital's day surgery services.

1.1.3 CLAVIA FCD CHEMOTHERAPY AND HEMODIALYSIS STRETCHER CHAIR

TECHNICAL DATA	
REFERENCE	CLAVIA FCD
DESIGNATION	Chemotherapy / dialysis chair with electrically assisted movements
SAFE WORKING LOAD	250 kg (indicates the maximum permissible load, including the total weight of the patient 185 kg, mattress and accessories)
WEIGHT	90 kg (CLAVIA 65), 92 kg (CLAVIA 75).
OVERALL DIMENSIONS	202 x 98 or 108 cm (reclined) / 148 x 98 or 108 cm (seated)
TREND	Electrically assisted 11°
BACKREST	0 to 90° - Electrically assisted, manual release
EQUIPMENT	3-sectional mattress base : 65 x 200 cm or 75 x 200 cm Fixed mattress made of high resilience foam (35 kg/m³), width. 10 cm PVC coated Jersey Adjustable headrest Electrically assisted height adjustment on 2 LINAK columns : 59 - 99 cm Electrically (LINAK) assisted Trendelenburg (15°)/ Reverse Trendelenburg (6°) Electrically (LINAK) assisted backrest, manual release Electrically (LINAK) assisted leg rest 0 - 60° Digital actuators allowing 6 programmed positions: Welcome, Relax, Reverse Trend, Securedb Trendelenburg position, CPR Machine-welded steel frame, epoxy coated Frame protection cover 4 lockable castors Ø 125 mm Anatomic shaped armrest following the backrest inclination. Adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient 2 holders on siderails for IV pole
AVAILABLE OPTIONS	Anti-static system Castors Ø 125 mm with central brake Double running castors Ø 150 mm with central brake 5th wheel Battery (50 cycles) Removable IV pole, 2-stage, 2 hooks Fixed push bars head-end Cable reel with remote holder Paper roll dispenser at the head-end Adjustable electric footrest Removable Crystal protection for leg rest Removable Crystal protection for footrest High Comfort mattress, high quality foam (35 kg/m³) + 2 cm thick viscoelastic layer (total height : 12 cm) Fixed I-Pad support with antitheft device Attachment hooks for bags on both sides below the patient surface Equipotential connection Meal tray : 42 x 32 cm (CLOSPR obligatory) Meal tray support

SCOPE

This chair is intended, depending on its specifications, options and accessories, for chemotherapy and hemodialysis services.

It is used for patient chemotherapy and hemodialysis care.

1.1.4 CLAVIA FSP MULTIPURPOSE CARE CHAIR

TECHNICAL DATA	
REFERENCE	CLAVIA FSP
DESIGNATION	Multipurpose care chair
SAFE WORKING LOAD	250 kg (indicates the maximum permissible load, including the total weight of the patient, mattress and accessories)
WEIGHT	93 kg (CLAVIA 65), 95 kg (CLAVIA 75).
OVERALL DIMENSIONS	202 x 80 or 90 cm (reclined) / 148 x 80 or 90 (seated)
TREND	Electrically assisted 11°
BACKREST	0 to 90° - Electrically assisted, manual release
EQUIPMENT	3-sectional mattress base : 65 x 200 cm or 75 x 200 cm Fixed mattress made of high resilience foam (35 kg/m³), width. 10 cm PVC coated Jersey Adjustable headrest Electrically assisted height adjustment on 2 LINAK columns : 55 - 99 cm Electrically (LINAK) assisted Trendelenburg (15°) / Reverse Trendelenburg (6°) Electrically (LINAK) assisted backrest, manual release Electrically (LINAK) assisted leg rest 0 - 60° Digital actuators allowing 6 programmed positions: Welcome, Relax, Reverse Trend, Secured Trendelenburg position, CPR Push bars Machine-welded steel frame, epoxy coated Frame protection cover Castors with central braking system Ø 125 mm, including 2 directional ones Drop-down laminated side guards. Length 110 x Height 27cm above mattress 2 holders on siderails for IV pole Battery
AVAILABLE OPTIONS	Anti-static system Double running castors Ø 150 mm with central brake 5th wheel Removable IV pole, 2-stage, 2 hooks Cable reel with remote holder Paper roll dispenser at the head-end Adjustable electric footrest Removable Crystal protection for leg rest Removable Crystal protection for footrest High Comfort mattress, high quality foam (35 kg/m³) + 2cm thick viscoelastic layer (total height : 12 cm) Attachment hooks for bags on both sides below the patient surface Equipotential connection Side rail padded cushions (pair) CLAVIA side rails cushion/armrest, removable and right/left reversible (single)

1.2 PURPOSE OF THE USER MANUAL

Healthcare professionals, this user manual was created for you to provide you with the necessary information and precautions to observe for autonomous and safe use.

The product's service life is dependent upon observance of the maintenance and tracking instructions.

It is thus essential that all staff has access to and be informed of the handling procedures and risks associated with the use of this chair. This latter should be handled by a qualified person who is aware of the risks inherent to improper use. Additional training is not required as we consider that knowledge of the user manual is sufficient.

1.3 USERS

Healthcare professionals: physicians, caregivers and agents approved by Technical Services who are accustomed to working with similar equipment can operate all functions.

In use cases where the patient is the operator (**use of the remote control and armrests only**), the caregivers must first ascertain the patient's ability to use the product. The caregivers must also teach the patient how to safely use the product: adjustment to low position to step off the product.

1.4 STORAGE

This user manual should be kept near the product, in a place provided for this purpose. It must be shielded from all elements and substance that could compromise its perfect legibility.


1.5 APPLIED PARTS

The applied parts are as follows:

- Upholstery / mattress
- Armrests / side rails
- Remote controls
- Footrest

1.6 SYMBOLS AND STANDARDS USED

This chair is a Class I non-invasive and active medical device.

	Type B appliance
	Electrical class II equipment
	Electrical hazard, do not open
	Duty cycle (motor cycle) 1 min. operation / 19 min. stop
IPX4	Splash-resistant
	See the manual / instruction brochure
	Indicates the maximum patient weight and safe working load
	Conforms to the essential requirements of directive 93/42/EEC on medical devices.
	Serial number
	Manufacture date
	Manufacturer
	Alternating current
	Maximum patient weight
	Electronic waste: this symbol indicates that the product must not be discarded with unsorted municipal waste, but rather to a selective waste collection channel. Contact the local distributor for information concerning the return and/or collection systems available in the country.
	Protéger de l'humidité / keep dry
	Haut, tenir debout / stand up
	Fragile / fragile
	Recyclable / recycling
	Temperature range
	Humidity limit

General safety sign:

- Information highlighted by this symbol must be read and scrupulously observed.
- Indicates that failure to observe the guideline could endanger the patient or user, or damage the equipment.

Sitting on the backrest is prohibited

Using the footrest as a footboard is prohibited

Sitting on the legrest is prohibited

Sitting on the armrests is prohibited

Leaning against the chair is prohibited

The chair must be set to its low comfort (or lying) position before being moved.

“WARNING: This appliance must not be altered without the manufacturer’s written consent.”

“Caution: only purpose-built accessories provided for this MD by BMB MEDICAL may be used.”

2 INSTRUCTIONS FOR USE

2.1 OPERATING CONDITIONS

AMBIENT CONDITIONS	OPERATION
Ambient temperature	+5 °C to 30°C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

2.2 INSTALLATION PROCEDURE

Before commissioning the chair and depending on the model, check that the following are operational:

- [Product power supply](#)
- [Remote control](#)
- [Trendelenburg / Reverse Trendelenburg](#)
- [Emergency position](#)
- [Variable height and lowering](#)
- [Backrest](#)
- [Legrest](#)
- [Side rails](#)
- [Foldable push bars](#)
- [Armrests](#)
- [Brake and guide pedals](#)

Clean the product according to your establishment's internal protocol.

2.3 BATTERY (OPTIONAL FOR THE FCD VERSION)

The battery unit is fitted with its own charger. Before initial use, leave the product to stand for 12 h in a temperate environment (to allow the electrical and electronic components to reach normal operating temperature), then charge the battery for 12 hours from a mains outlet.

In routine use, we strongly recommend keeping the chair plugged into the mains during sedentary use.

The battery, with its limited operating life, can be used to operate the stretcher-chair in the event of a power outage, or during transfer.

Once the manipulations are complete, remember to re-connect the stretcher-chair to the mains (at least 12 hours per day). Never allow the battery to drain fully (risk of reduced battery performance).

“Warning: Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact. »

Battery characteristics		Battery safety symbols	
Battery type:	Lithium-Ion	IPX6	Protection rating
Capacity:	2.25 Ah	CE	CE marking
Input supply voltage:	29~45 VDC	House icon	Indoor use
Input supply voltage:	25.9 VDC, 10 A	Waste icon	Electronic waste

2.4 PRODUCT POWER SUPPLY

Power grid connection conditions and electromagnetic compatibility:

This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.

Position the product such that it can be easily unplugged.

Unplug the product before moving it.

Do not leave the cables in contact with the floor or the casters,

Do not knot the cables,

Periodically check the condition of the cables.

Mains connection:	100 ~ 240 V depending on the country.
Frequency:	50 /60 Hz
Rated current:	3.3 A
Protection rating / degree of electric shock protection:	Class 2 appliance
Power input:	525 W
Intermittent operating mode:	5 %, i.e. 1 minutes of continuous operation followed by 19 minutes of rest.

The mains plug is the electrical grid insulation device.

Work on electrical parts should only be performed by qualified and authorized staff.

Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the EMC information [here](#).

Power cable storage:

on CLAVIA LSA and LSA CHIR:

Optional on all versions: optional power cord and remote control stowage.

CLAVIA LSA


CLAVIA FCD

2.5 REMOTE CONTROL

When using the remote control, all movements can be interrupted by simply releasing the corresponding key.

2.5.1 PROGRAMMED POSITIONS

Only activate the various positions after having ensured that no-one is touching any mechanical parts and that there are no objects under the seat, or between any of the moving structures.

2.5.2 TRANSPORT POSITION

CLAVIA FCD transport position

When transporting the patient, the chair should be set to comfort position, with the seat in the low position (see [2.7 - VARIABLE HEIGHT](#)). The chair is moved via the head-end push bars

CLAVIA LSA / LSA CHIR transport position


When transporting the patient, the chair is moved via the foot-end push bar, with the backrest slightly raised.

2.5.3 RESET


2.5.3.1 GENERAL RESET

This function should be used when the chair is not responding to the commands issued by the remote control. After connecting the chair to the mains, perform the following 3 operations:


- Press the “height adjust - raise” key until the audible signal (repeated beeps) stops

- Press the “raise backrest” key until the audible signal (repeated beeps) stops

- Press the “raise legrest” key until the audible signal (repeated beeps) stops


The general reset procedure is complete.

2.5.3.2 BACKREST POSITION RESET


Perform only if the backrest emergency manual unlock (red handle) has been used.

Press and hold the “backrest” key and raise the backrest fully.

2.5.3.3 RESET FUNCTION

To restart the product following switching to safety mode after a power outage, over-current or surge, press and hold the “**RAISE LEGREST**” and “**LOWER LEGREST**” keys until the audible signal (repeated beeps) stops.

If the procedure fails, check the following connections:

- Power cord
- Battery cable
- Fastening the 5 cables to the unit under the chair
- Actuator cable under hard surface (less accessible)
- Ensure that the module is powered on when connected to the mains.


Repeat the reset procedure from the start.

If the chair moves to an unwanted or extreme position, stop immediately to avoid damaging the equipment. Contact the BMB MEDICAL Customer service for a repair intervention

2.6 TRENDLENBURG / REVERSE TRENDLENBURG

To set the product to the Trendelenburg position, press the Trendelenburg and Reverse Trendelenburg keys simultaneously.

This is a safety mechanism to ensure that the patient cannot set the product to the Trendelenburg position.

Secure Trendelenburg =
press both keys
simultaneously

To set the product to the reverse Trendelenburg position, press the reverse Trendelenburg key.

2.7 EMERGENCY POSITION

The emergency position is set using the remote control:

Simultaneously press the 2 height adjust keys

2.8 VARIABLE HEIGHT AND LOWERING


Before using these functions, ensure that there are no objects or persons under the lying surface.

The height is adjusted using the remote control:

2.9 BACKREST

Before actuating the backrest, ensure that there are no obstacles to its movement (accessories, cables, etc.)

The backrest is controlled using the remote control:

EMERGENCY: If it is necessary to rapidly lower the backrest, pull the red handle and accompany the backrest by pulling it down.

If the emergency manual unlock has been used, lower the backrest with the remote control to its lowest position to reset operation. See [2.5.3.2 BACKREST POSITION RESET](#).

2.10 LEGREST

Function activated by the remote control:

Legrest
lowered position

Legrest
intermediate position

Legrest
raised position

2.11 STANDARDS AND PEDIATRICS SIDE RAILS

The side rails are intended to minimize the risk of the patient accidentally falling. They are not, in any way, designed to immobilize or maintain the patient.


Before manipulating the side rails, ensure that nothing could hinder their movement.

When lowering the side rail, take care not to trap the hands or feet of the patient or staff

If the side rail locking operation is silent (no “click”), please revise the side rails.


Do not move the product via the side rails

Lowering the left- or right-hand side rail:

Side rail lowered at head-end

OR


Side rail lowered at foot-end

The side rail is locked in position by squeezing the handle **1** and pushing the side rail **2** down


Raising the side rails:

Squeeze the handle **1** and lift the side rail **2**. The side rail automatically locks in position


2.12 FOLDABLE PUSH BARS

To fold the push bars, lift the lock ring fully **②**, the handles can then swivel freely downwards, away from the patient surface **③**.

To reposition them, raise them and return the ring to its low position. Ensure that they are correctly positioned.

The version under the backrest for the CLAVIA LSA and CLAVIA LSA CHIR works in the same way:

2.13 ARMRESTS

2.13.1 SWIVELING THE ARMRESTS

To lower or raise the armrests, press the paddle under each armrest. This same paddle, thanks to its articulation, allows you to position the armrest as desired.

2.13.2 ADJUSTING ARMREST HEIGHT

You can also set the height of the armrest to one of 2 positions (photo. ① & photo. ③).

To lower the armrest, pull the index pin while simultaneously lifting the armrest (photo. ②)

To raise the armrest, lift as shown in photo ② without pulling the index pin; the former will then return to its initial position (photo ①).

2.14 BRAKE AND GUIDE PEDALS


Except during transport, always ensure that the chair's brake is engaged. Check whether the brake is engaged by attempting to move the chair.

Do not park the chair on a sloping floor.

The brake and guide pedals are located at the head of the product and on the sides.


The head-end brake and guide pedal offers 3 function:

- **Brake engaged** - pedal lowered❶.
- **Chair in directional guidance mode** (easier movement in a straight line) - pedal raised❷. (Excluding 5th wheel option)
- **Chair in free guidance mode** for multidirectional movement - pedal in median position❸.

Brake and side guide pedals:

- **Brake engaged** - red side lowered❶.
- **Chair in directional guidance mode** - green side lowered❷. (Excluding 5th wheel option)
- **Chair in free guidance mode** - pedal in median position❸.

2.15 HEADRESTS

Removable comfort headrest


Dual articulated removable surgery headrest for ophthalmology (optional)

Removable surgery headrest for stomatology and ENT (optional)

2.15.1 FITTING THE HEADRESTS

4. Lift the handle
5. Slide the headrest into the insert
6. Lock the headrest in position by lowering the handle

Always ensure that the headrest is correctly locked in position.
The product should be set to sitting position before changing the headrest.

The weight of the patient's head rests on the headrest; it must thus be supported by the operator during adjustment. Failure to observe this recommendation could lead to patient injury.

To avoid any risk of pinching during headrest adjustment, keep your fingers away from all articulated areas.

Do not place your hand between the backrest and headrest to pull the unlock handle as you may injure your fingers.

OPTIONS


2.15.2 OPHTHALMOLOGY HEADREST (for the CLAVIA LSA CHIR version)

To ensure ergonomic positioning, the headrest is fitted with a dual articulation.

The 1st articulation is adjusted by pushing the handle on the right under the headrest outwards.

Release the handles to lock (automatic lock) the headrest in the desired position.


The 2nd articulation is adjusted using the handle on the left. Release the handle to lock (automatic lock) the headrest in the desired position.

The cushion is fitted with a magnet. Simple, rapid, hygienic, and adjustable positioning.

2.15.3 REMOVABLE SURGERY HEADREST FOR STOMATOLOGY AND ENT (for the CLAVIA LSA CHIR version)

To adjust the headrest:


4. Unlock the headrest by loosening the handle.
5. Swivel the headrest to position it.
6. Lock the position by tightening the handle.

Always ensure that the headrest articulation is correctly locked.


2.15.4 TECHNICAL SIDE RAILS

Slide the clamp on which your accessory is fastened and tighten the thumb wheel to lock it in position.


- IV pole
- Remote control holder
- Surgical field holder with an oxygen inlet holder that can be fastened to the technical rail with this accessory

Clamp for fastening accessories to the technical rail (for service and day surgery versions).


2.15.5 TECHNICAL RAIL-MOUNTED HEADREST SUPPORT


2.15.6 COMFORT CUSHION

The headrest can be easily adjusted via its counterweight system positioned behind the backrest.

2.15.7 COUSSIN/ACCOUDOIR

The cushion/armrest is removable and reversible, adaptable to the CLAVIA LSA / LSA CHIR/ FSP solid barriers.


2.15.9 ANTI-STATIC SYSTEM

Regularly check that the bonding strip is in contact with the floor and that it is neither twisted nor deformed.

2.15.10 FOOT PEDAL-CONTROLLED HEIGHT ADJUSTMENT (CLAVIA LSA SURGERY)

In addition to the hand-held remote control. The pedal is positioned at the head-end.


To use the pedal:

- Press + to raise the chair
- Press - to lower the chair

Before using these functions, ensure that there are no objects or persons under the lying surface.

2.15.11 DIRECTIONAL 5th WHEEL

When the product is equipped with the directional 5th wheel:

Brake pedal:


Stretcher with brake engaged

Guide pedal:

- Stretcher in directional guidance mode: 5th wheel


- Stretcher in free guidance mode:

2.15.12 IV POLE

The removable IV POLE is placed in one of the 2 holders placed on either side of the barrier (LSA/ LSA SURGERY/FSP version) or the armrest (FCD version).

In each holder, three locations are available in order to adapt the position of the IV POLE to the positioning of the chair.

To deploy the IV pole:

- 3- Turn the thumb wheel ❶ counter-clockwise and extend the telescopic part of the IV pole ❷ to the desired height.
- 4- To lock the telescopic part, turn the thumb wheel ❶ clockwise.
Ensure that the IV pole is securely locked.

Safe working load: 6 kg.


This indicative value corresponds to safe product use, but is not representative of normal use for which the load is significantly lower.

When moving the chair fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.

Do not move the chair via the IV pole.

For CLAVIA version LSA and LSA CHIR products, the removable IV pole fits into one of the corner holders provided.

2.15.13 60 CM WIDE BACKREST (CLAVIA LSA CHIR VERSION)

Compatible only with the 75 cm wide CLAVIA LSA SURGERY. This narrow backrest is adapted for day surgery. It is equipped with a technical rail for mounting a hand surgery table.

2.15.14 FOOTREST (for the FCD and FSP versions)


The footrest is adjusted using the remote control:



The footrest must not, under any circumstances, be used by the patient to climb onto the product, it is a patient comfort element.

2.15.15 IPAD HOLDER (CLAVIA FCD)


The flexible tablet holder is fitted onto the product's armrest:

Instructions for use:

Insert the tablet into the lower slot, then pull the handle until the tablet is locked.

The tablet holder is fitted with an anti-theft locking system. Use the key to lock the holder as shown in the image opposite.

2.15.16 PAPER ROLL HOLDER


2.15.17 CRYSTAL PROTECTION

2.15.17.1 CLAVIA LSA VERSION

The protection can be easily positioned without the need for tools or pre-equipment. It prevents the mattress from moving when the product is in seat position.

The footrest protection is attached by a strap and clip mechanism and is adjustable.

Clip-strap for adjustments

2.15.17.2 CLAVIA FCD / FSP VERSION

Pre-equipment required for the order. Snap fasteners on rear of legrest.

Protection on footrest, attached with Velcro

2.15.18 FIXED PUSH BAR (CLAVIA FCD & FSP versions)

2.15.19 CABLE REEL

2.15.20 OXYGEN CYLINDER HOLDER UNDER BACKREST


2.15.21 EQUIPOTENTIAL CONNECTOR


The equipotential connector must not, under any circumstances, be used as a protective earth connector

The various metal parts of the patient surface are interconnected by a ground wire to ensure they are all at the same potential.

An equipotential bonding conductor is provided, to which the user can connect.

2.15.22 SURGICAL FIELD HOLDER (CLAVIA LSA SURGERY VERSION)


The surgical field holder is fitted with an oxygen inlet.


This holder can thus raise the surgical field, allowing a flow of oxygen at the holder outlet.

2.15.23 PEDIATRIC KIT (for the CLAVIA LSA) version)


To lock the panel, insert it into the slot provided.


To open the panel, turn the index pin ① and push the panel outwards ②.

To remove the head-end panel, pull the 2 unlock plates ① fully out. You can then release the head-end panel from its slot ②.

To re-fit the head-end panel, make sure that the plates are fully out by pulling them **1**. Insert the head-end panel into the slots provided **2**. Lock the system by pushing the plates fully in **3**.

Check that the head-end panel is securely in place by attempting to remove it.

2.15.24 PROTECTIVE CUSHIONS

Attached with Velcro.

2.15.25 MEAL TRAY (CLAVIA FCD)

The meal tray is removable. It can be installed in the space provided for this purpose.

3 locking notches. To change notches, lift the meal tray slightly and rotate it until you reach the following notch.

Do not place a weight greater than 10kg on the meal tray.

3 CLEANING AND DISINFECTION

3.1 INTRODUCTION

Regular and thorough cleaning will guarantee the device's life (see cleaning program drawn up at the time of device commissioning, designating the reference staff).

Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.

The following cleaning advice is not intended to replace your establishment's in-house protocol.

This chair has been designed to facilitate its cleaning thanks to its removable cover.

3.2 FREQUENCY

The chair should be cleaned:

- Between each patient.
- After transporting an infected patient.
- Completely each week

3.3 CLEANING AND DISINFECTION ADVICE

3.3.1 CLEANING OPERATION

Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that cause stains (e.g.: "blue product" or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

In the case of a product fitted with a USB charging port, ensure that no liquids enter the port.

3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS

Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.

- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
 - o PVC, ABS, Polyamide, Polyurethane and Polypropylene
 - o Epoxy paint-coated metal surfaces
 - o Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

List of strictly prohibited products:

Solvents:	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
Bases:	10% Sodium hydroxide / 10% Ammonia
Acids:	10% Nitric / 30% Hydrochloric / 10% Citric
Miscellaneous:	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

High-pressure and steam cleaning are prohibited.

Any damage resulting from failure to observe this instruction for use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.

4 PRECAUTIONS

4.1 DEVICE-RELATED PRECAUTIONS

4.1.1 WARNINGS

This MD has been designed for and adapted to medical practices and procedures, and hence a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the common sense of the chair's users.

- Do not use the product without first having familiarized yourself with this user manual.
- Do not alter the chair. Any changes made to the system could lead to unpredictable operation that could in turn jeopardize the safety of the user or patient. This would also void the chair's warranty.
- Used to accommodate a person of small stature (less than 1.30 m tall), or under the age of 12 years without adapted equipment (pediatric kit)
- To ensure correct use of the product, only proceed if the weight is correctly distributed over the device.
- The chassis cover load capacity is of 35 kg. Do not sit or stand on the cover.
- Only purpose-built accessories provided for this MD by BMB MEDICAL may be used
- Connection, by any means, of an electric device to the chair.
- The product must be connected to a mains power outlet compliant with standards in force.
- Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact.
- This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.
- Position the product such that it can be easily unplugged.
- Unplug the product before moving it.
- Do not leave the cables in contact with the floor or the casters.
- Do not knot the cables.
- Periodically check the condition of the cables.
- Work on electrical parts should only be performed by qualified and authorized staff.
- Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the EMC information [here EMC-RELATED INFORMATION](#).
- When using the remote control, all movements can be interrupted by simply releasing the corresponding key.

- Only activate the various positions after having ensured that no-one is touching any mechanical parts and that there are no objects under the seat or patient surface, or between any of the moving structures.
- If the emergency manual unlock has been used, lower the backrest with the remote control to its lowest position to reset operation. See [2.5.3.2 BACKREST POSITION RESET](#).
- The side rails are intended to minimize the risk of the patient accidentally falling. They are not, in any way, designed to immobilize or maintain the patient.
- Before manipulating the side rails, ensure that nothing could hinder their movement.
- When lowering the side rail, take care not to trap the hands or feet of the patient or staff.
- If the side rail locking operation is silent (no “click”), please revise the side rails.
- Except during transport, always ensure that the chair's brake is engaged. Check whether the brake is engaged by attempting to move the chair.
- Patient sitting and standing operations should be performed, with the chair in the low position, from the side, after raising the armrests (and placing the chair in the patient installation position for the CLAVIA FCD), or after lowering the side rails.
- When getting up, the patient must not rest on the articulated anatomically shaped armrests.
- The “variable height” function should be used when handling the patient in order to work at an ergonomic height.
- Do not park the chair on a sloping floor.
- Use outside of the establishment, or on unprepared ground.
- Do not use the MD on a floor with a slope of more than 10%.
- Use of the chair with a load that exceeds the safe working load (see technical characteristics).
- Always ensure that the headrest is correctly locked in position.
- The product should be set to sitting position before changing the headrest.
- The weight of the patient's head rests on the headrest; it must thus be supported by the operator during adjustment. Failure to observe this recommendation could lead to patient injury.
- To avoid any risk of pinching during headrest adjustment, keep your fingers away from all articulated areas.
- Do not place your hand between the backrest and headrest to pull the unlock handle as you may injure your fingers.
- Always ensure that the headrest articulation is correctly locked.
- IV pole safe working load: **6 kg**.
- When moving the chair fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.
- Do not move the chair via the IV pole
- The footrest must not, under any circumstances, be used by the patient to climb onto the product, it is a patient comfort element.
- In the case of an unsupervised patient, the product should be left in its lowest position to reduce the risk of falling.
- Tell the patient not to sit on the ends (foot or head) of the chair/stretcher. An excess weight applied to one of the ends could tip the patient surface and injure the patient.

- Tell the patient and any accompanying persons not to sit on the armrests.
- The equipotential connector must not, under any circumstances, be used as a protective earth connector
- Check that the head-end panel of the pediatric kit is securely in place by attempting to remove it.
- Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.
- Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.
- High-pressure and steam cleaning are prohibited.
- Any use that does not fall within the operating context of a hospital chair.
- To prevent failure, electromedical equipment requires special EMC-related precautions, described [here](#).
- Never use the chair in the presence of an anesthesia mixture flammable with air, with oxygen, or with nitrous oxide as this could cause a fire and/or explosion.
- Do not manually raise the legrest.
- Use by any non-caregiver who has not received the authorization and explanations from a qualified individual.
- In the case of an unsupervised patient, the product should be left in its lowest position to reduce the risk of falling.

All manipulations relating to the aforementioned precautions for use are described in paragraph [II - Instructions for use](#).

4.1.2 CONTRAINDICATIONS

There are no noteworthy contraindications for this device.

4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE

The device must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The device must be regularly inspected, at least once per year. See the maintenance check-list. To maintain the device in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed

In the event of observed anomalies or damage that could jeopardize device operation and safety, the equipment must be immediately withdrawn.

The device should only be used according to the instructions in this user manual.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

Work on electrical parts should only be performed by qualified and authorized staff.

* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

4.1.4 EMC-RELATED INFORMATION

The products in the CLAVIA range are suitable for use in all healthcare establishments, with the exception of facilities in the vicinity of HF ELECTROSURGICAL APPLIANCES, or near an electromagnetically shielded room for an EM magnetic resonance imaging system, in which the intensity of EM INTERFERENCES is high.

The products in the CLAVIA range are not suitable for use in domestic homes or establishments or on any premises connected to the public low-voltage network supplying domestic buildings.

EM interference could cause the product to make unexpected movements.

This appliance should not be used next to, or stacked on other appliances as this could lead to incorrect operation. If this mode of operation is necessary, this and the other appliances should be observed to ensure their correct operation.

The use of accessories, transducers or cables other than those specified or supplied by the manufacturer can lead to an increase in electromagnetic emissions, or a decrease in this appliance's immunity, causing inappropriate operation.

Do not use portable RF communication devices (including peripherals such as antenna cables and external antennas) less than 30 cm (12 inches) from any part of products in the CLAVIA range, including the cables specified by the manufacturer. Otherwise, the performance of these devices could be impaired.

The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11, class A). If used in a residential environment (for which compliance with the CISPR 11 class B standard is generally required), this equipment may not provide appropriate protection to radio frequency communication services. The user may need to take attenuation steps, such as moving or re-orienting the equipment.

List of cables used:

1. Power cord, maximum length 1000 cm.
2. Variable height column, maximum length 750 mm
3. Backrest cable: maximum length 1250 mm
4. Legrest cable: maximum length 750 mm
5. Footrest cable: maximum length 1500 mm

EMISSIONS	
Conducted and radiated RF emissions:	CISPR11 Group 1 Class A
Harmonic distortions IEC 61000-3-2:	Class A
Voltage fluctuations and Flicker IEC 61000-3-3:	Compliant

Immunity test	Test level IEC 60601	Compliance level		
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV in contact ± 2/4/8/15 kV in air	± 8 kV in contact ± 2/4/8/15 kV in air		
RF radiated electromagnetic fields IEC 61000-4-3	3 V/m (professional healthcare establishment) 80MHz-2.7GHz 80% AM at 1 kHz	3V/m 80MHz-2.7GHz 80% AM at 1 kHz		
Proximity fields emitted by wireless RF communication devices IEC 61000-4-3				
	Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)
	385	Pulsed modulation: 18 Hz	27	27
	450	Pulsed modulation: 18 Hz	28	28
	710-745-780	Pulsed modulation: 217 Hz	9	9
	810-870-930	Pulsed modulation: 18 Hz	28	28
	1720-1845-1970	Pulsed modulation: 217 Hz	28	28
	2450	Pulsed modulation: 217 Hz	28	28
5240-5500-5785	Pulsed modulation: 217 Hz	9	9	
Rapid / burst electrical spikes IEC 61000-4-4	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz		
Shock waves IEC 61000-4-5	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV		
Conducted interference, induced by RF fields IEC 61000-4-6	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz		
Power frequency magnetic fields IEC 61000-4-8	30A/m	30A/m		
Voltage Dips and Interruptions: IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles		

NB: There is no deviation from the requirements of IEC 60601-1-2 ed. 4

NB: There are no other known measures for maintaining basic safety with respect to EMC phenomena.

4.2 TRANSPORT AND STORAGE

Storage and transport conditions:

Ambient conditions	Operation
Ambient temperature	-10 to +30 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

4.2.1 TRANSPORT

During chair transport with a view to its storage, the following is recommended to prevent any impacts or falls:

- Strap it to a surface of the means of transport
- Cover it
- Engage the brakes
- Set it to its low position

Damage caused during transport is not covered by the warranty. The costs of repair and replacement of damaged parts will be borne by the customer.

4.2.2 STORAGE

When the chair is in storage, ensure that:

- Its brakes are engaged
- It is covered

The product should be stored in a clean, ventilated, temperate and dry place, protected from the elements and from direct sunlight.

During storage, no heavy weights should be placed on the chair.

The chair must not be considered as a support surface for any other type of equipment.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

In the event of prolonged storage, the battery should be recharged every 6 months, proceeding as follows:

- Connect the mains power cord.
- Briefly press one of the remote control functions to initiate charging.
- Allow the battery to charge for 8 hours.

4.3 WARRANTY

Warranty: two years parts and labor according to the following terms:

BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.

This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL stretchers are designed for a projected service life of 10 years under normal conditions of use, with appropriate periodic maintenance as described in the maintenance check-list of each stretcher. BMB MEDICAL provides the original buyer with the guarantee that the welds on the stretchers are free from structural defects for the projected stretcher service life of 10 years, as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use (Warning paragraph [here](#))
- Work without prior notice to our technical services, or maintenance on:
 - A mechanical part of the chair
 - A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to observe the preventive maintenance operations recommended in the maintenance check-list.

4.4 AFTER-SALES SERVICE

For any information concerning product installation, use and maintenance, the technical documentation (wiring diagrams, list of components, maintenance procedures, etc.), or any other questions, please contact the After-Sales Service at the following address: sav@bmb-medical.com, or by phone at +33 (0)4 74 08 71 71, or your local distributor for overseas purchases.

The operator must carefully read the information contained in this user manual, in particular the information concerning the precautions relating to safety and installation/use methods.

In the event of doubt concerning the correct interpretation of the instructions, do not hesitate to contact the After-Sales Service whose staff are always available.

To facilitate ordering spare parts, you can request the exploded diagrams of your products from the After-Sales Service. You can thus obtain the references of the desired parts.

4.5 ENVIRONMENT

4.5.1 SPENT EQUIPMENT

In France:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution. For this purpose, BMB MEDICAL fulfills its obligations concerning the end of life of CLAVIA range products that it markets by funding the Recylum recycling channel dedicated to Pro WEEE that collects them free of charge (more information on www.recylum.com).

In Europe:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution.

Obtain the appropriate disposal location from the authorities in charge of environmental protection.

Tél: +33 (0)4 74 08 71 71

Fax: +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001
01604 TREVOUX Cedex
FRANCE