DIVIA

MANUEL D'UTILISATION
USER MANUAL

 **BMB**
MEDICAL

 [FRANÇAIS](#)
[ENGLISH](#)

FRANÇAIS

BMB MEDICAL

100 Avenue du Formans
CS 11001
01604 TRÉVOUX CEDEX
Tél : +33(0)4 74 08 71 71
Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

This manual is for the purpose of familiarization with the use, maintenance and safe handling of your new medical device. We strongly advise you to carefully read this document and to follow the recommendations for safe, optimal use.

BMB MEDICAL offers a wide variety of options for both models featured in this user manual. Consequently, some features described in this manual may not apply to your couch.

Some illustrations are included for reference only and do not accurately represent the device you own.

Stick the
product reference
label here

The information in this document may be changed without notice, for which BMB MEDICAL can in no way be held responsible. Any representation or reproduction, whether in whole or in part, without permission from BMB MEDICAL, is illegal. Such illegal representation or reproduction, by any process whatsoever, shall constitute an infringement. The illustrations contained in this document are not contractually binding.

CONTENTS

ENGLISH	28
1 INTRODUCTION SPECIFICATIONS	31
1.1 PRODUCT PRESENTATION	31
1.1.1 DIVIA HEIGHT ADJUSTABLE EXAMINATION COUCH	31
1.1.2 DIVIA MIXTE HEIGHT ADJUSTABLE EXAMINATION COUCH	32
1.2 AREA OF APPLICATION	33
1.3 PURPOSE OF THE USER MANUAL	33
1.4 STORAGE	33
1.5 SYMBOLS AND STANDARDS USED	33
2 INSTRUCTIONS FOR USE	34
2.1 OPERATING CONDITIONS	34
2.2 INSTALLATION PROCEDURE	34
2.3 BATTERY	34
2.4 PRODUCT POWER SUPPLY	35
2.5 REMOTE CONTROL	36
2.5.1 RESET	36
2.6 ADJUSTABLE HEIGHT AND LOWERING	38
2.7 BACKREST	38
2.8 INTERMEDIATE BACKREST (DIVIA MIXTE)	39
2.9 TRENDLENBURG POSITION	39
2.10 INDEPENDANT WHEELS BRAKES	39
2.11 OPTIONS	40
2.11.1 HEAD REST	40
2.11.2 CENTRALLY CONTROLLED BRAKE AND STEERING PEDAL	40
2.11.3 SIDERAILS	40
2.11.4 CLAMPS FOR FASTENING STIRRUPS, CRADLES, LEG SUPPORTS OR EXAMINATION LAMP	41
2.11.5 LATERAL TECHNICAL RAILS	42
2.11.6 RECOVERY TANKS	43
3 CLEANING AND DISINFECTION	44
3.1 INTRODUCTION	44
3.2 FREQUENCY	44
3.3 CLEANING AND DISINFECTION ADVICE	44
3.3.1 CLEANING OPERATION	44

3.3.2	RECOMMENDED EQUIPMENT AND PRODUCTS	44
4	PRECAUTIONS	46
4.1	DEVICE RELATED PRECAUTIONS	46
4.1.1	WARNINGS	46
4.1.2	CONTRAINDICATIONS	46
4.1.3	PRECAUTIONS FOR USE OF THE COUCH	47
4.1.4	PRECAUTION FOR SERVICING AND MAINTENANCE OF THE EQUIPMENT	47
4.1.5	EMC-RELATED INFORMATION	48
4.2	TRANSPORT AND STORAGE	49
4.2.1	TRANSPORT	49
4.2.2	STORAGE	50
4.3	WARRANTY	50
4.4	AFTER SALES SERVICE	51
4.5	ENVIRONMENT	52
4.5.1	SPENT EQUIPMENT	52

1 INTRODUCTION SPECIFICATIONS

1.1 PRODUCT PRESENTATION

1.1.1 DIVIA HEIGHT ADJUSTABLE EXAMINATION COUCH

TECHNICAL CHARACTERISTICS	
REFERENCE	DIVIA
DESIGNATION	HEIGHT ADJUSTABLE EXAMINATION COUCH
SAFE WORKING LOAD	250 kg (maximum allowable load including the total weight of the patient, mattress and accessories)
WEIGHT	100 kg
OVERALL DIMENSIONS	65 x 200 cm or 75 x 200 cm.
TRENDELENBURG / REVERSE TRENDELENBURG BACKREST	Electric (+/- 15°).
EQUIPMENTS	<p>0 to 80°</p> <p>2-sectional mattress base: 65x 200 cm or 75 x 200 cm Flat upholstery, high density foam 35 kg/m³. Width 6,5 cm Jersey with PVC coating, antibacterial and anti-fire (M1) treatment Adjustable headrest (0 - 80°) with 2 gas springs Electric height adjustment: 46 - 86 cm Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°) Steel structure, epoxy coated (light grey) Frame with 2 electric columns providing great stability Frame protection cover Remote control 4 lockable wheels Ø 125 mm Examination towel dispenser head end</p>
AVAILABLE OPTIONS	<p>Antistatic system Casters Ø 125 mm with central brake Double running castors Ø 150 mm with central brake Electric adjustable backrest Wireless remote control (in addition of the wired remote control) Foot pedal control (in addition of the wired remote control) Drop-down side rails (pair) Holder for clamps (pair) Fixation clamps for heels, armrest for blood sampling or leg support (pair) Heels (pair) Anatomic shaped armrest Leg support (pair) Protection cover for leg rest (pair) Clamps for IV pole holder (pair) Removable 2-stage IV pole, 2 hooks Removable 2-stage IV pole, 4 hooks Comfort pillow (grey) Crystal protection cover for examination couch width 65cm Crystal protection cover for examination couch width 75cm Battery Examination lamp Equipotential connection DIN rails</p>

1.1.2 DIVIA MIXTE HEIGHT ADJUSTABLE EXAMINATION COUCH

TECHNICAL CHARACTERISTICS	
REFERENCE	DIVIA MIXTE
DESIGNATION	HEIGHT ADJUSTABLE COMBINED EXAMINATION COUCH
SAFE WORKING LOAD	250 kg (maximum allowable load including the total weight of the patient, mattress and accessories)
WEIGHT	90 kg
OVERALL DIMENSIONS	65 x 178 cm or 75 x 178 cm.
TRENDELENBURG / REVERSE TRENDELENBURG	Electric (+/- 15°).
BACKREST	0 to 80°
EQUIPMENTS	<p>3-sectional mattress base: 65x 178 cm or 75 x 178 cm Flat upholstery, high density foam 35 kg/m³. Width 6,5 cm Jersey with PVC coating, antibacterial and anti-fire (M1) treatment Electric adjustable backrest (0 - 80°) Electric height adjustment: 46 - 86 cm Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°) Steel structure, epoxy coated (light grey) Frame with 2 electric columns providing great stability Frame protection cover Remote control 4 lockable double running wheels Ø 125 mm Examination towel dispenser head end</p>
AVAILABLE OPTIONS	<p>Antistatic system Casters Ø 125 mm with central brake Double running castors Ø 150 mm with central brake Wireless remote control (in addition of the wired remote control) Foot pedal control (in addition of the wired remote control) Drop-down side rails (pair) Holder for clamps (pair) Fixation clamps for heels, armrest for blood sampling or leg support (pair) Heels (pair) Anatomic shaped armrest Leg support (pair) Protection cover for leg rest (pair) Clamps for IV pole holder (pair) Removable 2-stage IV pole, 2 hooks Removable 2-stage IV pole, 4 hooks Comfort pillow (grey) Battery Examination lamp Recovery Tank (grey) Equipotential connection DIN rails</p>

1.2 AREA OF APPLICATION

This couch is an exam couch. It has equipment that makes it a perfect tool for gynecological examinations.

1.3 PURPOSE OF THE USER MANUAL

Health professionals, this user manual has been created for the purpose of supplying you with all the necessary information and usage precautions to be followed to ensure that you can safely operate the product on your own.

The instructions for maintenance and product monitoring are intended to ensure its service life.

It is therefore essential for all staff to have access to them and consequently to be informed about the handling procedures and the risks associated with use of this couch. The couch should be handled by a person who is qualified and aware of the risks inherent in improper use of the product.

1.4 STORAGE

This user manual should be kept within close proximity to the product in a location designed for this purpose. It should be kept away from any item or substance that could obscure its legibility.

1.5 SYMBOLS AND STANDARDS USED

	Type B device
	Indicates that failure to follow the recommendation can lead to endangerment of the patient and the user, or damage to the equipment.
	Compliant with the essential requirements of Directive 93/42/CEE concerning medical devices.
IPx4	Protection against splashed water and sprayed water.
	Indicates the maximum weight of the patient and the safe working load.

This couch is a Class I medical device, non-invasive and active.

2 INSTRUCTIONS FOR USE

2.1 OPERATING CONDITIONS

AMBIENT CONDITIONS	OPERATION
Ambient temperature	+ 5 C° to 30°C
Relative humidity	20 to 90% at 30° C – no condensation
Atmospheric pressure	700 to 1060 ha

2.2 INSTALLATION PROCEDURE

Before commissioning the chair and depending on the model, check that the following are operational:

- [Product power supply](#)
- [Remote control](#)
- [Variable height and lowering](#)
- [Backrest](#)
- [Trendelenburg / Reverse Trendelenburg](#)
- [Brake and pedals](#)

Clean the product according to your establishment's internal protocol.

2.3 BATTERY

The battery unit is fitted with its own charger. Before initial use, leave the product to stand for 12 h in a temperate environment (to allow the electrical and electronic components to reach normal operating temperature), then charge the battery for 12 hours from a mains outlet.

In routine use, we strongly recommend keeping the chair plugged into the mains during sedentary use.

The battery, with its limited operating life, can be used to operate the stretcher-chair in the event of a power outage, or during transfer.

Once the manipulations are complete, remember to re-connect the stretcher-chair to the mains (at least 12 hours per day). Never allow the battery to drain fully (risk of reduced battery performance).

“Warning: Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact. »

Battery characteristics		Battery safety symbols	
Battery type:	Lithium-Ion	IPX6	Protection rating
Capacity:	2.25 Ah		CE marking
Input supply voltage:	29~45 VDC		Indoor use
Input supply voltage:	25.9 VDC, 10 A		Electronic waste

2.4 PRODUCT POWER SUPPLY

Power grid connection conditions and electromagnetic compatibility:

This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.

Position the product such that it can be easily unplugged.

Unplug the product before moving it.

Do not leave the cables in contact with the floor or the casters,

Do not knot the cables,

Periodically check the condition of the cables.

Mains connection:	100 ~ 240 V depending on the country.
Frequency:	50 /60 Hz
Rated current:	3.3 A
Protection rating / degree of electric shock protection:	Class 2 appliance
Power input:	525 W
Intermittent operating mode:	5 %, i.e. 1 minutes of continuous operation followed by 19 minutes of rest.

The mains plug is the electrical grid insulation device.

Work on electrical parts should only be performed by qualified and authorized staff.

Electromagnetic interference caused by the couch cannot be excluded. In this case, the couch's functions should be inhibited and those instruments subject to interference should be moved away. You can view the EMC information [here](#).

2.5 REMOTE CONTROL

When using the remote control, all movements can be interrupted by simply releasing the corresponding key.

2.5.1 RESET


2.5.1.1 GENERAL RESET

This function should be used when the couch is not responding to the commands issued by the remote control. After connecting the couch to the mains, perform the following 3 operations :


- Press the “height adjust - raise” key until the audible signal (repeated beeps) stops

- Press the « raise backrest » key until the end of the movement and then press the « lower backrest » key until the audible signal (repeated beeps) stops.


- Press the « raise second backrest » key until the end of the movement and then press the « lower second backrest » key until the audible signal (repeated beeps) stops.

The general reset procedure is complete.

2.5.1.2 **RESET FUNCTION**

To restart the product following switching to safety mode after a power outage, over-current or surge, press and hold the **“HEIGHT ADJUST RAISE”** and **“HEIGHT ADJUST LOW”** keys until the audible signal (repeated beeps) stops.

If the procedure fails, check the following connections:

- Power cord
- Battery cable
- Fastening the 5 cables to the unit under the couch
- Actuator cable under hard surface (less accessible)
- Ensure that the module is powered on when connected to the mains.

Repeat the reset procedure from the start.

If the couch moves to an unwanted or extreme position, stop immediately to avoid damaging the equipment. Contact the BMB MEDICAL Customer service for a repair intervention.

2.6 ADJUSTABLE HEIGHT AND LOWERING

Before using these functions: Make sure that there are no objects or persons underneath the plinth area.

The height is adjusted by remote control :

or

Raised position

Lowered position

2.7 BACKREST

Prior to activating the torso-lift, check that there is nothing hampering its movement (accessories, cables...)

The backrest is operated by remote control:

or


Or by manually activating the controls located under the torso-lift: press (with both hands for greater ease) to free the springs and adjust the torso-lift to the desired height.

2.8 INTERMEDIATE BACKREST (DIVIA MIXTE)


The intermediate backrest is operated by remote control:

The intermediate backrest cannot be activated if the primary backrest is inclined to more than 60°, (see photo above).

2.9 TRENDLENBURG POSITION

This is managed by remote control :

Reverse Trendelenburg position

Trendelenburg position

2.10 INDEPENDANT WHEELS BRAKES

Make sure to apply the couch brakes at all times, except during transport. Check that the brake is active by attempting to move the couch.

Do not park the couch on a sloping floor.

Press the pedal on the wheel to apply the brake and raise the pedal to release it.

Wheel brake released


Wheel brake applied

2.11 OPTIONS

2.11.1 HEAD REST

The head rest is easily adjusted via a counterbalancing system located at the rear of the backrest.

2.11.2 CENTRALLY CONTROLLED BRAKE AND STEERING PEDAL


Make sure to apply the couch brakes at all times, except during transport. Check that the brake is active by attempting to move the couch.

Do not park the couch on a sloping floor.

Lateral brake and steering pedals:

- **Brakes applied to couch** – red side in lowered position ❶.
- **Directional steering of couch** – green side in lowered position ❷.
- **Freehand steering of couch** - pedal at the middle position ❸.

2.11.3 SIDERAILS


The bars are intended to reduce the risk of the patient suffering accidental falls. They are in no way designed to restrain the patient or to restrict the patient with strap-like immobilizing devices.

Before adjusting the bars, make sure that there is nothing hampering their movement.

Do not move the couch around by its bars.


Lowering the bars:

Unlock by pressing and holding in the button ❶, and pushing the bar ❷.

Raising the bars :

Pull up the upper bar until the button locks into place.

2.11.4 CLAMPS FOR FASTENING STIRRUPS, CRADLES, LEG SUPPORTS OR EXAMINATION LAMP

2.11.4.1 TO INSTALL THE CLAMPS

1. Slide the clamp onto the stand provided for this purpose (make sure beforehand that the thumbwheel does not get in the way of the clamp; if it does, unscrew it).

2. Tighten the clamp by turning the thumbwheel.

2.11.4.2 TO INSTALL THE ACCESSORIES

1. Slide the accessory onto the stand intended for this purpose (make sure beforehand that the thumbwheel does not get in the way of the accessory; if it does, unscrew it).

2. Tighten the accessory by turning the thumbwheel

2.11.4.3 TO INSTALL THE IV POLE

The IV holder is not inserted into the clamps provided for the stirrups, cradles or leg supports. Special clamps with an identical function to the one described above are provided.


Maximum safe working load: 10 kg²

When conveying the couch with the IV holder attached, make sure that the IV holder is not able to come into contact with any obstacle (e.g.: doorway, lighting fixture, etc.) or individual.

2.11.5 LATERAL TECHNICAL RAILS

Slide the clamp on which your accessory is attached into place and secure the thumbscrew to tighten it.


Maximum Safe Working Load : 15 kg

² This indicative value is in keeping with safe use of the product,. However, it is not representative of normal use with a much lower load.

2.11.6 RECOVERY TANKS

3 CLEANING AND DISINFECTION

3.1 INTRODUCTION

Regular and thorough cleaning will guarantee the device's life (see cleaning program drawn up at the time of device commissioning, designating the reference staff).

Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.

The following cleaning advice is not intended to replace your establishment's in-house protocol.

This chair has been designed to facilitate its cleaning thanks to its removable cover.

3.2 FREQUENCY

The couch should be cleaned:

- Between each patient.
- After transporting an infected patient.
- Completely each week

3.3 CLEANING AND DISINFECTION ADVICE

3.3.1 CLEANING OPERATION

Before cleaning the couch, unplug it from the mains and ensure that all internal connections are correctly established..

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that cause stains (e.g.: "blue product" or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS

- Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds, and bleach solutions, etc.
- Observe the prescribed dosages and dilutions.

- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
 - PVC, ABS, Polyamide, Polyurethane and Polypropylene
 - Epoxy paint-coated metal surfaces
 - Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

List of strictly prohibited products:

Solvents:	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
Bases:	10% Sodium hydroxide / 10% Ammonia
Acids:	10% Nitric / 30% Hydrochloric / 10% Citric
Miscellaneous:	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

High-pressure and steam cleaning are prohibited.

Any damage resulting from failure to observe this instruction for use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.

4 PRECAUTIONS

4.1 DEVICE RELATED PRECAUTIONS

4.1.1 WARNINGS

This MD has been designed for and adapted to medical practices and procedures, and hence a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the common sense of the chair's users.

- Use by any person who has not read this user manual.
- Use by any person outside of the nursing staff who has not been authorized or advised by a qualified person.
- Use as a bed for a short person (under 1.30 m) or under the age of 12.
- Handling by a child.
- Use of the electric controls by more than one person.
- Use of the couch with a heavier weight than the safe working load (refer to the technical features).
- Use of the chassis to hold or store objects.
- Use of oxygen equipment under the plinth frame.
- Use in the presence of a flammable anesthetic mixture.
- Use outside of the establishment or on an unfinished floor surface.
- Use on a floor surface with a gradient greater than 10°.
- Washing by sluicing, high-pressure water or steam jet or passage through a tunnel washer.
- Rigorous use or jolting of the motors.
- Moving the couch while it is still connected to the electric power source.
- Connecting an electrical device to the couch by any means whatsoever.
- Any attempt to unscrew the ground terminal screw from the power supply box.
- Connecting the couch to a non-compliant power source.
- Any other use that does not qualify as a function for a hospital couch.

Information: The electrical springs are equipped with a thermal protection feature that prevents them from functioning if the motor overheats from rigorous use. They return to their original working condition within 2 hours of the couch being disconnected from the electric power source.

4.1.2 CONTRAINDICATIONS

There are no noteworthy contraindications for this device.

4.1.3 PRECAUTIONS FOR USE OF THE COUCH

For effective use of the equipment, proceed only if the weight is properly distributed on the device.

In the case of a difficult or unsupervised patient, it is advisable to leave the splint area in the lowest position in order to reduce the risk of a potential fall.

Similarly, it is advisable to use the “adjustable height” function whenever handling the patient, in order to work at an ergonomic height.

The patient should be placed in a sitting or standing position on the base side, by lowering the bars. The patient should not be seated on the torso-lift. Risk of falling.

Under no circumstances can either of the bars be considered an area inaccessible to patients. The appropriate measures must be taken depending on the patient.

During conveyance, it is necessary to make sure that the power cord is disconnected and wrapped around the reel provided for this purpose before moving the couch.

When the couch does not need to be moved, it should be immobilized using the brakes, and then making sure that the wheels are locked by attempting to move it. Failure to observe this handling method can lead to a fall by someone attempting to lean against the couch.

All manipulations relating to the aforementioned precautions for use are described in paragraph II - [Instructions for use](#).

It is inadvisable to handle without appropriate training.

4.1.4 PRECAUTION FOR SERVICING AND MAINTENANCE OF THE EQUIPMENT

No modifications must be made to the device without the prior consent of the manufacturer: the modification may cause it to function in an unpredictable manner and thereby cause injury to the patient and the nursing staff.

The device requires regular inspection.

In order to maintain the device in good condition and ensure that it can be safely operated over the long term, it is necessary to service it in the prescribed manner.

Should abnormalities or damage which could compromise the operation and safety of the device be noted, the equipment must be immediately removed.

The device must be used exclusively in accordance with the instructions in this instruction manual.

The device must not be exposed to, much less come into contact with, combustion heat sources or flammable agents.

BMB MEDICAL shall not be held liable in the event of incidents or accidents resulting from improper use or modification of the product by the buyer or user.

4.1.5 EMC-RELATED INFORMATION

The products in the DIVIA range are suitable for use in all healthcare establishments, with the exception of facilities in the vicinity of HF ELECTROSURGICAL APPLIANCES, or near an electromagnetically shielded room for an EM magnetic resonance imaging system, in which the intensity of EM INTERFERENCES is high.

The products in the DIVIA range are not suitable for use in domestic homes or establishments or on any premises connected to the public low-voltage network supplying domestic buildings.

EM interference could cause the product to make unexpected movements.

This appliance should not be used next to, or stacked on other appliances as this could lead to incorrect operation. If this mode of operation is necessary, this and the other appliances should be observed to ensure their correct operation.

The use of accessories, transducers or cables other than those specified or supplied by the manufacturer can lead to an increase in electromagnetic emissions, or a decrease in this appliance's immunity, causing inappropriate operation.

Do not use portable RF communication devices (including peripherals such as antenna cables and external antennas) less than 30 cm (12 inches) from any part of products in the DIVIA range, including the cables specified by the manufacturer. Otherwise, the performance of these devices could be impaired.

The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11, class A). If used in a residential environment (for which compliance with the CISPR 11 class B standard is generally required), this equipment may not provide appropriate protection to radio frequency communication services. The user may need to take attenuation steps, such as moving or re-orienting the equipment.

List of cables used:

1. Power cord, maximum length 5000 mm.
2. Variable height column head end cable, maximum length 460 mm
3. Variable height column foot end cable, maximum length 1700 mm
4. Backrest cable 1 : maximum length 850 mm
5. Backrest cable 2 : maximum length 1700 mm

EMISSION	
Conducted and radiated RF emissions:	CISPR11 Group 1 Class A
Harmonic distortions IEC 61000-3-2:	Class A
Voltage fluctuations and Flicker IEC 61000-3-3:	Compliant

Immunity test	Test level IEC 60601	Compliance level		
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV in contact ± 2/4/8/15 kV in air	± 8 kV in contact ± 2/4/8/15 kV in air		
RF radiated electromagnetic fields IEC 61000-4-3	3 V/m (professional healthcare establishment) 80MHz-2.7GHz 80% AM at 1 kHz	3V/m 80MHz-2.7GHz 80% AM at 1 kHz		
Proximity fields emitted by wireless RF communication devices IEC 61000-4-3	Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)
	385	Pulsed modulation: 18 Hz	27	27
	450	Pulsed modulation: 18 Hz	28	28
	710-745-780	Pulsed modulation: 217 Hz	9	9
	810-870-930	Pulsed modulation: 18 Hz	28	28
	1720-1845-1970	Pulsed modulation: 217 Hz	28	28
	2450	Pulsed modulation: 217 Hz	28	28
	5240-5500-5785	Pulsed modulation: 217 Hz	9	9
Rapid / burst electrical spikes IEC 61000-4-4	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz		
Shock waves IEC 61000-4-5	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV		
Conducted interference, induced by RF fields IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz		
Power frequency magnetic fields IEC 61000-4-8	30A/m	30A/m		
Voltage Dips and Interruptions: IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles		

NB: There is no deviation from the requirements of IEC 60601-1-2 ed. 4

NB: There are no other known measures for maintaining basic safety with respect to EMC phenomena.

4.2 TRANSPORT AND STORAGE

Storage and transport conditions:

Ambient conditions	Operation
Ambient temperature	-10 to +30 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

4.2.1 TRANSPORT

When transporting the couch to its storage location, in order to avoid any risk of collision or falls it is advisable to:

- strap it to a support on the transport medium
- cover it
- lock the brakes
- adjust it in the lowered position

Any damage occasioned during transport of the couch are not covered by the warranty. Repair and replacement of the damaged parts are the responsibility of the client.

4.2.2 STORAGE

When the couch is stored, make sure that :

- the brakes are locked
- the couch is covered

While in storage, it should not have any heavy weight placed on top of it.

The couch must not be treated as a support base for any type of equipment.

The device must not be exposed to, much less come into contact with sources of combustion heat or flammable agents.

In the event of prolonged storage, it is advisable to recharge the battery every 6 months as follows:

- Connect the power cord.
- Briefly press the control function to initiate charging.
- Leave the battery charging for 8 hours.

4.3 WARRANTY

Warranty: two years parts and labor according to the following terms:

BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.

This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL stretchers are designed for a projected service life of 10 years under normal conditions of use, with appropriate periodic maintenance as described in the maintenance check-list of each stretcher. BMB MEDICAL provides the original buyer with the guarantee that the welds on the stretchers are free from structural defects for the projected stretcher service life of 10 years, as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- non-standard use ([Section 4.2](#))
- interference without prior notification to our technical department or servicing of:
 - a mechanical part of the couch
 - a spring or a motor
- failure to comply with any other condition of the warranty
- any inappropriate cleaning
- the use of accessories not sourced from the manufacturer
- failure to comply with the preventive servicing operations recommended on the servicing check-list

4.4 AFTER SALES SERVICE

For any information concerning product installation, use and maintenance, the technical documentation (wiring diagrams, list of components, maintenance procedures, etc.), or any other questions, please contact the After-Sales Service at the following address: sav@bmb-medical.com, or by phone at +33 (0)4 74 08 71 71, or your local distributor for overseas purchases.

The operator must carefully read the information contained in this user manual, in particular the information concerning the precautions relating to safety and installation/use methods.

In the event of doubt concerning the correct interpretation of the instructions, do not hesitate to contact the After-Sales Service whose staff are always available.

To facilitate ordering spare parts, you can request the exploded diagrams of your products from the After-Sales Service. You can thus obtain the references of the desired parts.

4.5 ENVIRONMENT

4.5.1 SPENT EQUIPMENT

In France:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution. For this purpose, BMB MEDICAL fulfills its obligations concerning the end of life of CLAVIA range products that it markets by funding the Recylum recycling channel dedicated to Pro WEEE that collects them free of charge (more information on www.recylum.com).

In Europe:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution.

Obtain the appropriate disposal location from the authorities in charge of environmental protection.

Outside of Europe:

Dispose of the product or its components in accordance with local laws and regulations.

Contact the local distributor for information concerning the return and/or collection systems available in the country.

Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001

01604 TREVOUX Cedex

FRANCE