
I-CARE

MANUEL D'UTILISATION
USER MANUAL

 **BMB
MEDICAL**

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FRANÇAIS

BMB MEDICAL

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

This manual contains important information regarding the proper use, maintenance and safety of your new medical appliance. We recommend that you carefully review and comply with the recommendations contained in this manual to ensure optimum safe use.

BMB MEDICAL offers a wide range of accessories for use with the 3 models covered in this operating manual. As a consequence, some of the features described in this manual may not be applicable to your particular stretcher trolley.

Some illustrations are provided for information purposes and do not accurately represent the device in your possession.

The information contained in this document may be modified without any warning and is not intended as a commitment on behalf of BMB MEDICAL.

No part of this document can be photocopied or reproduced without the approval of BMB MEDICAL. Such reproduction or representation, by any process whatsoever, would constitute an infringement. All images in this document are for representational purposes only.

CONTENTS

ENGLISH	53
1 SPECIFICATION OVERVIEW	56
1.1 I-CARE PRESENTATION	56
1.1.1 I-CARE EMERGENCY	57
1.1.2 I-CARE AMBULATORY	59
1.1.3 I-CARE OPHTHALMOLOGY	61
1.2 I-CARE EZ-GO PRESENTATION	63
1.2.1 MOTORIZED STRETCHER I-CARE EZ-GO EMERGENCY	64
1.2.2 MOTORIZED STRETCHER I-CARE EZ-GO AMBULATORY	65
OVERVIEW	66
APPLICATION FIELD	67
1.3 OBJECTIVE OF THIS USER MANUAL	67
1.4 USERS	67
1.5 CONSERVATION	67
1.6 APPLIED PARTS	67
1.7 SYMBOLS AND STANDARDS USED	68
2 INSTRUCTION FOR USE	70
2.1 OPERATING CONDITIONS	70
2.2 INSTALLATION PROCEDURE	70
2.3 BATTERY (optional electric backrest and/or kneegatch)	70
2.4 PRODUCT POWER SUPPLY	71
2.5 HEIGHT ADJUSTMENT AND LOWERING	72
2.5.1 RAISING THE STRETCHER	72
2.5.2 LOWERING THE STRETCHER	72
2.6 BACKREST	73
2.6.1 MANUAL BACKREST	73
2.6.2 ELECTRIC BACKREST	73
2.7 SIDERAILS	73
2.7.1 I-CARE EMERGENCY AND OPHTHALMOLOGY	74
2.7.2 I-CARE AMBULATORY	75
2.8 BREAK PEDAL / STEER PEDAL (identical method of operation for all 3 models)	75
2.9 SELF POWERED MOTORISED (I-CARE EZ-GO)	76
2.9.1 INITIAL OPERATION	76
2.9.2 BATTERY	76
2.9.3 WHEELING THE STRETCHER	79
2.9.4 BRAKE SYSTEM	80
2.10 TRENDLENBURG / REVERSE TRENDLENBURG	81
2.11 CPR RELEASE MECHANISM	81
2.12 TILTING HEADREST (ONLY FOR THE ICARE OPHTHALMOLOGY MODEL)	82

2.12.1	INSTALLATION OF THE HEADREST	83
2.12.2	SURGERY HEADREST FOR OPHTHALMOLOGY.....	83
2.12.3	REMOVABLE SURGERY HEADREST FOR STOMATOLOGY AND ENT	84
2.13	OPTIONS	85
2.13.1	FOLDING PUSH BARS.....	85
2.13.2	FOLD OUT MONITOR TRAY	85
2.13.3	IV POLE	86
2.13.4	OXYGEN BOTTLE HOLDER BENEATH THE BACKREST	88
2.13.5	LATERAL TECHNICAL RAILS AND CLAMP	88
2.13.6	HEADREST SUPPORT, FIXATION ON TECHNICAL RAILS	88
2.13.7	PAPER ROLL DISPENSER	89
2.13.8	X RAY CASSETTE BELOW THE PATIENT SURFACE	89
2.13.9	ANTISTATIC SYSTEM.....	89
2.13.10	HOLDER FOR OPERATING FIELD WITH CONNECTION FOR OXYGEN (OPHTHALMOLOGY VERSION) 89	89
2.13.11	CHART HOLDER BELOW THE HEADREST.....	90
2.13.12	BAG HANGERS	90
2.13.13	MANUAL KNEE GATCH.....	90
2.13.14	ELECTRIC KNEE GATCH.....	91
3	CLEANING AND DISINFECTING	92
3.1	INTRODUCTION	92
3.2	FREQUENCY	92
3.3	CLEANING AND DISINFECTION ADVICE.....	93
3.3.1	CLEANING OPERATION	93
3.3.2	RECOMMENDED EQUIPMENT AND PRODUCTS.....	93
4	PRECAUTIONS.....	95
4.1	DEVICE-RELATED PRECAUTIONS	95
4.1.1	WARNINGS	95
4.1.2	CONTRAINDICATIONS	97
4.1.3	PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE	98
4.1.4	EMC-RELATED INFORMATION	98
4.2	TRANSPORT AND STORAGE.....	100
4.2.1	TRANSPORT	100
4.2.2	STORAGE.....	101
4.3	WARRANTY	101
4.4	AFTER SALES SERVICE	102
4.5	ENVIRONMENT.....	103
4.5.1	SPENT EQUIPMENT.....	103

1 SPECIFICATION OVERVIEW

1.1 I-CARE PRESENTATION

Variable-height patient support device on hydraulic columns and equipped with an electric-assisted actuator system that allows, depending on the options, the movements of the backrest, leg lifter and knee gatch. The equipment is mounted on casters.

The device allows an operator to convey the patient within a healthcare establishment. The 5th wheel option allows you to switch from straight in-line guidance during transport to the ability to pivot the product around a corner. The chair is braked by a mechanical system actuated by pressing a pedal.

The device is equipped with a hydraulic or electrically assisted (optional) actuator system operated by a remote control to ensure the movements of the following elements:

- backrest
- legrest
- knee gatch (optional)

The products are designed to ensure that there are no unacceptable risks (essential performance) associated with patient support or caused by accidental movements.

List of I-CARE references:

- **I-CARE 60 EMERGENCY** - HYDRAULIC HEIGHT-ADJUSTABLE I-CARE EMERGENCY STRETCHER WITH 60 CM WIDE MATTRESS
- **I-CARE 70 EMERGENCY** - HYDRAULIC HEIGHT-ADJUSTABLE I-CARE EMERGENCY WITH STRETCHER 70 CM WIDE MATTRESS
- **I-CARE 60 AMBULATORY** - HYDRAULICALLY HEIGHT-ADJUSTABLE I-CARE AMBULATORY STRETCHER WITH 60 CM WIDE MATTRESS
- **I-CARE 70 AMBULATORY** - HYDRAULICALLY HEIGHT-ADJUSTABLE I-CARE AMBULATORY STRETCHER WITH 70 CM WIDE MATTRESS
- **I-CARE 80 AMBULATORY** - HYDRAULIC HEIGHT-ADJUSTABLE I-CARE AMBULATORY STRETCHER WITH 80 CM WIDE MATTRESS
- **I-CARE 70 OPHTHALMOLOGY** - HYDRAULIC HEIGHT-ADJUSTABLE I-CARE OPHTHALMOLOGY STRETCHER WITH 70 CM WIDE MATTRESS

1.1.1 I-CARE EMERGENCY

TECHNICAL DATA	
REFERENCE	I-CARE EMERGENCY
DESIGNATION	STRETCHER TROLLEY WITH HEIGHT ADJUSTMENT ON HYDRAULIC COLUMNS
SAFE WORKING LOAD	320 kg (maximum load limit including the total weight of the patient, the mattress and accessories)
WEIGHT	80 kg
OVERALL DIMENSIONS	216 x 84 cm or 216 x 94 cm
FOLDABLE RAILS	115 x 26.5 cm (height above mattress)
TREND/REVERSE TREND	18°
BACKREST	0 to 90 °
EQUIPMENTS	2-sectional resin mattress base: 60 x 197 cm or 70 x 197 cm Radiolucent material Height adjustment range : from 64 to 96 cm (with 10 cm thick mattress) Trendelenburg / Anti-Trendelenburg 18° Two-sided + foot or head-end height adjustment and Trendelenburg / Anti-Trendelenburg adjustment Removable mattress made of high quality and resilience foam (35kg/m3), height : 10 cm Removable cover. PVC coated Jersey. Colour : grey 6 holders for IV pole Bumpers on all 4 corners. Backrest adjustable by 2 gas springs (0 - 90°) Castors with central braking system Ø 200 mm at foot- and head-end 5th retractable directional wheel Drop-down side rails. Anodized aluminium Foldable push bars at the foot-end 2 parts Frame cover (removable for easy cleaning and maintenance) including utility tray and O2 tank holder
OPTIONAL FEATURES	Anti-static system Permanent, foldable 3-stage «comfort» IV pole, 2 hooks Permanent, foldable 3-stage «comfort» IV pole, 4 hooks Removable IV pole, 2-stage, 2 hooks Foldable push bars head-end Foldable push bar foot-end 1 part X-ray cassette below the patient surface CORBE Electric backrest with manual release Electric backrest and knee gatch, manual release Battery Patient's remote control Knee gatch DIN rails Fixation clamp for armrest (on technical rails, CRORT necessary) Adjustable fixation clamp for armrest (on technical rails, CRORT necessary) Anatomic shaped armrest (ETU-R necessary) Vertical oxygen bottle holder Paper roll dispenser at the head-end Chart holder (for paperwork) below the headrest Attachment hooks for bags on both sides below the patient surface Monitor tray at the foot-end High Comfort mattress, high quality foam (35kg/m3), 2cm thick viscoelastic layer (total height : 12cm) Equipotential connection Side rail padded cushions (pair) Universal Meal tray / Extensible meal tray for side rails: from 65 to 88 cm / Fall protection side

OVERVIEW

BMB MEDICAL
dedicated solutions for your patients' transfer & wellbeing

⚠ Safe working load (SWL) : 320 kg

1.1.2 I-CARE AMBULATORY


TECHNICAL DATA	
REFERENCE	I-CARE AMBULATORY
DESCRIPTION	STRETCHER TROLLEY WITH HEIGHT ADJUSTMENT ON COLUMNS
SAFE WORKING LOAD	320 kg (maximum load limit including the total weight of the patient, the mattress and accessories)
WEIGHT	85 kg (I- CARE Ambulatory 60, 70), 90 kg (I-CARE 80)
OVERALL DIMENSIONS	69 x 216 cm, 79 x 216 cm or 89 x 216 cm.
FOLDABLE RAILS	126 x 30 cm (height above mattress)
TREND/REVERSE TREND	18°
BACKREST	0 to 90 °
EQUIPMENT	2-sectional resin mattress base: 60 x 197 cm, 70 x 197 cm or 80 x 197 cm Radiolucent material Overall dimensions: 69 x 216 cm, 79 x 216 cm or 89 x 216 cm Height adjustment range : from 64 to 96 cm (with 10 cm thick mattress) Trendelenburg / Anti-Trendelenburg 18° Two-sided + foot or head-end height adjustment and Trendelenburg/Anti-Trendelenburg adjustment Removable mattress made of high quality and resilience foam (35kg/m3), height : 10 cm Removable cover. PVC coated Jersey. Colour : grey 6 holders for IV pole Bumpers on all 4 corners Backrest adjustable by 2 gas springs (0 - 90°) Castors with central braking system Ø 200 mm at foot- and head-end 5th retractable directional wheel Drop-down side rails. Anodized aluminium Foldable push bars at the foot-end 2 parts Frame cover (removable for easy cleaning and maintenance) including utility tray and O2 tank holder
OPTIONAL FEATURES	Anti-static system Permanent, foldable 3-stage «comfort» IV pole, 2 hooks Permanent, foldable 3-stage «comfort» IV pole, 4 hooks Removable IV pole, 2-stage, 2 hooks Foldable push bars head-end Foldable push bar foot-end 1 part X-ray cassette below the patient surface Electric backrest with manual release Electric backrest and knee gatch, manual release Battery Patient's remote control Knee gatch Vertical oxygen bottle holder Paper roll dispenser at the head-end Chart holder (for paperwork) below the headrest Attachment hooks for bags on both sides below the patient surface Monitor tray at the foot-end High Comfort mattress, high quality foam (35kg/m3), 2cm thick viscoelastic layer (total height : 12cm) Equipotential connection DIN rails Fixation clamp for armrest (on technical rails, CRORT necessary) Adjustable fixation clamp for armrest Anatomic shaped armrest (ETU-R necessary) Side rail padded cushions (pair) Universal Meal tray / Extensible meal tray for side rails: from 65 to 88 cm / Fall protection side

OVERVIEW

BMB MEDICAL

dedicated solutions for your patients' transfer & wellbeing

⚠ Safe working load (SWL) : : 320 kg

1.1.3 I-CARE OPHTHALMOLOGY


TECHNICAL DATA	
REFERENCE	I-CARE OPHTHALMOLOGY
DESCRIPTION	STRETCHER TROLLEY WITH HEIGHT ADJUSTMENT ON COLUMNS
SAFE WORKING LOAD	200 kg (maximum load limit including the total weight of the patient, the mattress and accessories)
WEIGHT	90 kg
OVERALL DIMENSIONS	216 x 84 cm
FOLDABLE RAILS	115 x 26.5 cm (height above mattress)
TREND/REVERSE TREND	18 °
BACKREST	0 à 90 °
EQUIPMENT	2-sectional resin mattress base: 70 x 197 cm. Radiolucent material Height adjustment range: from 64 to 96 cm (with 8 cm thick mattress) Trendelenburg / Anti-Trendelenburg 18° Bilateral foot-operated variable height control / Trendelenburg / brakes Removable mattress made of high quality and resilience foam (35kg/m³), thickness : 8 cm Removable cover. PVC coated Jersey. Colour : grey Backrest adjustable by 2 gas springs (0 - 90°) Castors with central braking system Ø 200 mm at foot-end and head-end 5th retractable directional wheel Drop-down side rails. Anodized aluminium Foldable push bars at the foot-end 2 parts DIN rails Frame cover including utility tray and O2 tank holder
OPTIONAL FEATURES	Anti-static system Articulated removable surgery head rest for ophthalmology + hollow cushion Articulated removable surgery head rest for stomatology + flat cushion Removable surgery head rest for stomatology and ENT Additional hollow cushion for articulated removable surgery headrest Additional flat cushion for articulated removable surgery headrest Head rest support. Fixation on technical rails Removable IV pole, 2-stage, 2 hooks Fixation clamp for armrest Adjustable fixation clamp for armrest GT Anatomic shaped armrest (ETU-R necessary) Holder for operating field, with connection for oxygen Foldable push bar foot-end 1 part Knee gatch Equipotential connection Side rail padded cushions (pair) Universal Meal tray / Extensible meal tray for side rails: from 65 to 88 cm / Fall protection side

OVERVIEW

⚠ Safe working load (SWL) : 200 kg

BMB MEDICAL

dedicated solutions for your patients' transfer & wellbeing

1.2 I-CARE EZ-GO PRESENTATION

Variable-height patient support device on hydraulic columns and equipped with an electric-assisted actuator system that allows, depending on the options, the movements of the backrest, leg lifter and knee gatch.

The device allows the operator to transport the patient inside a health care facility. The equipment is mounted on wheels and is equipped with an electric wheel to assist in movement.

The braking system is a mechanical system operated by the support of a pedal.

The products are designed to ensure the absence of incapacitating risks (essential performance) related to patient support as well as the risk of unintentional movements.

List of I-CARE EZ-GO references) :

- **I-CARE EZ-GO 60 EMERGENCY –**
I-CARE EMERGENCY hydraulic variable height motorized stretcher with a 60 cm wide mattress.
- **I-CARE EZ-GO 70 EMERGENCY –**
I-CARE EMERGENCY hydraulic variable height motorized stretcher with a 70 cm wide mattress.
- **I-CARE EZ-GO 60 AMBULATORY –**
I-CARE AMBULATORY hydraulic variable height motorized stretcher with a 60 cm wide mattress.
- **I-CARE EZ-GO 70 AMBULATORY –**
I-CARE AMBULATORY hydraulic variable height motorized stretcher with a 70 cm wide mattress.
- **I-CARE EZ-GO 80 AMBULATORY –**
I-CARE AMBULATORY hydraulic variable height motorized stretcher with a 80 cm wide mattress.

1.2.1 MOTORIZED STRETCHER I-CARE EZ-GO EMERGENCY

TECHNICAL DATA	
REFERENCE	I-CARE EZ-GO EMERGENCY
DESCRIPTION	MOTORIZED STRETCHER
SAFE WORKING LOAD	320 kg (maximum load limit including the total weight of the patient, the mattress and accessories)
WEIGHT	112 to 117 kg
OVERALL DIMENSIONS	74 x 216 cm or 84 x 216 cm
FOLDABLE RAILS	115 x 26.5 cm (height above mattress)
TREND/REVERSE TREND	18 °
BACKREST	0 to 90 °
EQUIPMENT	Motorized drive system with autonomous wheel. Integrated battery. 15 km autonomy. Charge time: 5 hours. Lithium battery without memory effect. Forward/Backward twist grip with progressive speed acceleration (foot-end). Push bars foot-end. 2-sectional resin mattress base: 60 x 197 cm or 70 x 197 cm. Radiolucent material. Overall dimensions: 74 x 216 cm or 84 x 216 cm. Height adjustment range: from 64 to 96 cm (with 10 cm thick mattress). Trendelenburg / Anti-Trendelenburg 18°. Two-sided + foot or head-end height adjustment and Trendelenburg / Anti-Trendelenburg adjustment. Removable mattress made of high quality and resilience foam (35kg/m3), thickness: 10 cm. Removable cover. PVC coated Jersey. Colour : grey. 6 holders for IV pole. Bumpers on all 4 corners. Backrest adjustable by 2 gas springs (0 - 90°). Castors with central braking system Ø 200 mm at foot- and head-end. 5th retractable and directional wheel (motorized wheel). Antistatic system. Drop-down side rails. Anodized aluminium. Frame cover (removable for easy cleaning and maintenance) including utility tray and O2 tank holder.
OPTIONAL FEATURES	Permanent, foldable 3-stage «comfort» IV pole, 2 hooks Permanent, foldable 3-stage «comfort» IV pole, 4 hooks Removable IV pole, 2-stage, 2 hooks Foldable push bars head-end Electric backrest with manual release Electric backrest and knee gatch, manual release Battery Patient's remote control Knee gatch Vertical oxygen bottle holder Paper roll dispenser at the head-end Chart holder (for paperwork) below the headrest Attachment hooks for bags on both sides below the patient surface Monitor tray at the foot-end High Comfort mattress, high quality foam (35kg/m3), 2cm thick viscoelastic layer (total thickness : 12cm) Equipotential connection DIN rails Fixation clamp for armrest (on technical rails, CROORT necessary) Adjustable fixation clamp for armrest Anatomic shaped armrest (ETU-R necessary) Side rail padded cushions (pair)

1.2.2 MOTORIZED STRETCHER I-CARE EZ-GO AMBULATORY

TECHNICAL DATA	
REFERENCE	I-CARE EZ-GO AMBULATORY
DESCRIPTION	MOTORIZED STRETCHER
SAFE WORKING LOAD	320 kg (maximum load limit including the total weight of the patient, the mattress and accessories)
WEIGHT	112kg (I-CARE Ambulatory), 117 kg (I- CARE Ambulatory 70), 122 kg (I-CARE 80)
OVERALL DIMENSIONS	216 x 84 cm or 216 x 94 cm
FOLDABLE RAILS	126 x 30 cm (above the mattress)
TREND/REVERSE TREND	18°
BACKREST	0 to 90 °
EQUIPMENT	Motorized drive system with autonomous wheel Integrated battery. 15 km autonomy. Charge time: 5 hours. Lithium battery without memory effect. Forward/Backward twist grip with progressive speed acceleration (foot-end). 2-sectional resin mattress base: 60 x 197 cm, 70 x 197 cm or 80 x 197 cm. Radiolucent material. Overall dimensions: 69 x 216 cm, 79 x 216 cm or 89 x 216 cm. Height adjustment range : from 64 to 96 cm (with 10 cm thick mattress). Trendelenburg / Anti-Trendelenburg 18°. Two-sided + foot-end height adjustment and Trendelenburg / Anti-Trendelenburg adjustment. Removable mattress made of high quality and resilience foam (35kg/m3), thickness : 10 cm. Removable cover. PVC coated Jersey. Colour : grey. 6 holders for IV pole. Bumpers on all 4 corners. Backrest adjustable by 2 gas springs (0 - 90°). Castors with central braking system Ø 200 mm at foot- and head-end. 5th retractable directional wheel (motorized wheel). Antistatic system. Drop-down side rails. Anodized aluminium. Foldable push bars at the foot-end (1 or 2 parts). Frame cover (removable for easy cleaning and maintenance) including utility tray and O2 tank holder.
OPTIONAL FEATURES	Permanent, foldable 3-stage «comfort» IV pole, 2 hooks Permanent, foldable 3-stage «comfort» IV pole, 4 hooks Removable IV pole, 2-stage, 2 hooks Foldable push bars head-end Electric backrest with manual release Electric backrest and knee gatch, manual release Battery Patient's remote control Knee gatch Vertical oxygen bottle holder Paper roll dispenser at the head-end Chart holder (for paperwork) below the headrest Attachment hooks for bags on both sides below the patient surface Monitor tray at the foot-end High Comfort mattress, high quality foam (35kg/m3), 2cm thick viscoelastic layer (total thickness: 12cm) Equipotential connection DIN rails Fixation clamp for armrest (on technical rails, CRORT necessary) Adjustable fixation clamp for armrest Anatomic shaped armrest (ETU-R necessary) Side rail padded cushions (pair)

OVERVIEW

dedicated solutions for your patients' transfer & wellbeing

⚠ Safe working load (SWL) : 320 kg

VUE D'ENSEMBLE EZ-GO AMBULATORY

Forward/Backward twist grip with progressive speed acceleration

Backrest tilt (0 - 90°) with 2 gas springs

Retractable push handles (1 or 2 parts) head end (optional)

Mattress stop

Bumpers (4-corners)

Release tuck-away siderails

Hydraulic raise pedal (two-sided + foot end)

Brake/steer pedal foot end

Up : steering with 5th motor driven wheel

Mid-position : free steering

Down : braking position

Reverse Trend. pedal (two-sided + foot end)

Trend. pedal (two-sided + foot end)

Reverse Trend. tilt

Trendelenburg tilt

Low foot pedal

Oxygen tank holders

Brake/steer pedal head end

Self-sufficient motor driven wheel
Battery charger

APPLICATION FIELD

This stretcher trolley is intended according to its specifications, optional features and accessories for use in stretcher services, ambulatory care and hospital emergency units.

It is designed for transporting, treating and transferring patients and for their rest and recovery.

1.3 OBJECTIVE OF THIS USER MANUAL

This operating manual is intended for health professionals in order to furnish you with all the necessary information and all of precautions to be followed to ensure effective and safe use of this appliance.

The instructions relating to the maintenance and inspection of the product should be followed to ensure satisfactory equipment performance.

It is therefore essential that all staff have access to these instructions and are aware of the handling procedures and risks associated with the use of the stretcher trolley. The device must be handled by a qualified person who has been made aware of the risks associated with using the device incorrectly.

1.4 USERS

Healthcare professionals: physicians, caregivers and agents approved by Technical Services who are accustomed to working with similar equipment can operate all functions.

In use cases where the patient is the operator (**use of the remote control and armrests only**), the caregivers must first ascertain the patient's ability to use the product. The caregivers must also teach the patient how to safely use the product: adjustment to low position to step off the product.

1.5 CONSERVATION

This User Manual shall be kept near the product in a location provided for this purpose. It shall be kept away from any element and substance compromising its legibility.

1.6 APPLIED PARTS

The applied parts are as follows:

- Upholstery / mattress
- Armrests / side rails
- Remote controls

1.7 SYMBOLS AND STANDARDS USED

Electric I-CARE AMBULATORY and EMERGENCY models are Class I non-invasive and active medical devices.

I-CARE AMBULATORY, EMERGENCY and OPHTALMOLOGY models are Class I non-invasive and non-active medical devices.

Motorized I-CARE AMBULATORY EZ-GO and EMERGENCY EZ-GO models are Class I non-invasive and active medical devices.

	Type B appliance
	Electrical class II equipment
	Electrical hazard, do not open
	Duty cycle (motor cycle) 1 min. operation / 19 min. stop
IPX4	Splash-resistant
	See the manual / instruction brochure
	Indicates the maximum patient weight and safe working load
	Conforms to the essential requirements of directive 93/42/EEC on medical devices.
	Serial number
	Manufacture date
	Manufacturer
	Alternating current
	Maximum patient weight
	Electronic waste: this symbol indicates that the product must not be discarded with unsorted municipal waste, but rather to a selective waste collection channel. Contact the local distributor for information concerning the return and/or collection systems available in the country.
	Protéger de l'humidité / keep dry
	Haut, tenir debout / stand up
	Fragile / fragile

	Recyclable / recycling
	Temperature range
	Humidity limit

General safety sign:

- Information highlighted by this symbol must be read and scrupulously observed.
- Indicates that failure to observe the guideline could endanger the patient or user or damage the equipment.

Sitting on the backrest is prohibited

Sitting on the legrest is prohibited

Sitting on the armrests is prohibited

Leaning against the stretcher is prohibited

The stretcher must be set to its low comfort (or lying) position before being moved.

“WARNING: This appliance must not be altered without the manufacturer’s written consent.”

“Caution: only purpose-built accessories provided for this MD by BMB MEDICAL may be used.”

2 INSTRUCTION FOR USE

2.1 OPERATING CONDITIONS

AMBIENT CONDITIONS	OPERATION
Ambient temperature	+5 °C to 30°C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

2.2 INSTALLATION PROCEDURE

Before commissioning the chair and depending on the model, check that the following are operational:

- [The battery](#)
- [The product power supply](#)
- [The height adjustment and lowering](#)
- [Backrest](#)
- [Siderails](#)
- [Brake pedal / steer pedal](#)
- [Trendelenburg / reverse Trendelenburg](#)
- [CPR release mechanism](#)
- [Tilting headrest](#)

Clean the product according to your establishment's internal protocol.

2.3 BATTERY (OPTIONAL ELECTRIC BACKREST AND/OR KNEEGATCH)

The battery unit is fitted with its own charger. Before initial use, leave the product to stand for 12 h in a temperate environment (to allow the electrical and electronic components to reach normal operating temperature), then charge the battery for 12 hours from a mains outlet.

In routine use, we strongly recommend keeping the chair plugged into the mains during sedentary use.

The battery, with its limited operating life, can be used to operate the stretcher-chair in the event of a power outage, or during transfer.

Once the manipulations are complete, remember to re-connect the stretcher-chair to the mains (at least 12 hours per day). Never allow the battery to drain fully (risk of reduced battery performance).

“Warning: Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact. »

Battery characteristics		Battery safety symbols	
Battery type:	Lithium-Ion	IPX6	Protection rating
Capacity:	2.25 Ah		CE marking
Input supply voltage:	29~45 VDC		Indoor use
Input supply voltage:	25.9 VDC, 10 A		Electronic waste

2.4 PRODUCT POWER SUPPLY

Power grid connection conditions and electromagnetic compatibility:

This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.

Position the product such that it can be easily unplugged.

Unplug the product before moving it.

Do not leave the cables in contact with the floor or the casters,

Do not knot the cables,

Periodically check the condition of the cables.

Mains connection:	100 ~ 240 V depending on the country.
Frequency:	50 /60 Hz
Rated current:	3.3 A
Protection rating / degree of electric shock protection:	Class 2 appliance
Power input:	525 W
Intermittent operating mode:	5 %, i.e. 1 minutes of continuous operation followed by 19 minutes of rest.

The mains plug is the electrical grid insulation device.

Work on electrical parts should only be performed by qualified and authorized staff.

Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the [EMC information here](#).

2.5 HEIGHT ADJUSTMENT AND LOWERING

Before using these functions:


Make sure all persons and objects are clear of the area under the stretcher top,

Make sure that nobody's feet are placed under any of the pedals (height adjustment, trendelenburg/reverse trendelenburg).

2.5.1 RAISING THE STRETCHER

The controls (grey pedals), are located on both sides of the stretcher and there is a 3rd control at the foot- and head-ends.


To raise the stretcher, pump pedal:

2.5.2 LOWERING THE STRETCHER

To lower the stretcher, depress the Trendelenburg pedal (yellow pedal) and the reverse Trendelenburg pedal (blue pedal) together.

Make sure that feet are correctly positioned to ensure equal pressure is applied to both pedals.


2.6 BACKREST

Before adjusting the backrest, remove any equipment that may be in the way (accessories, cables etc.)

2.6.1 MANUAL BACKREST

The backrest adjustable by 2 gas springs is operated by pushing down on the bar located under the stretcher top: exert adequate pressure (using both hands) to activate the lift mechanism and adjust the backrest to the desired height.

2.6.2 ELECTRIC BACKREST

Backrest adjustment is operated by remote control(s):

- the foot-end remote control
- and the patient's remote controls (optional)

2.7 SIDERAILS

The siderails are intended to reduce the risk of patients from accidentally falling off the trolley. The siderails are not intended for use as patient restraints with strap-type securing devices.

When operating the siderails, remove any equipment out of the way.


Do not attempt to move the trolley when the siderails are raised

2.7.1 I-CARE EMERGENCY AND OPHTHALMOLOGY

Lowering the siderails :

To release the siderail, pull the lever downwards ❶, and push the siderail to the full down position ❷.

For easier handling, use one hand to operate the lever and place the other hand on the grab rail to guide the positioning.


Do not use excessive force when using the handle. Release the rail by moving it slightly to activate the lever.

Raising the siderails :

- Pull up the grab rail and raise it to the full up position until the latch engages.
- Locking is automatic.


Keep your hands away from the lever to avoid entrapment of fingers.

2.7.2 I-CARE AMBULATORY


Lowering the siderails :

To release the siderail, press down the button ❶, and push the siderail to the full down position ❷.

Raising the siderails:

Lift up the top rail until the button lock feature is activated.

2.8 BREAK PEDAL / STEER PEDAL (IDENTICAL METHOD OF OPERATION FOR ALL 3 MODELS)




When the stretcher trolley is not in motion, always make sure that the parking brake is engaged. Check that the brake is not activated when trying to move the stretcher trolley.

Do not park the stretcher trolley on sloping surfaces.

The brake-steer pedal has three positions:

- Stretcher trolley locked - pedal in the full down position ❶.
- Stretcher trolley engaged enabling multidirectional movement on 4 wheels (5th wheel is deactivated) - pedal in the centre position ❷.
- Stretcher trolley steered by the 5th wheel (in order to facilitate steering in a straight course) - pedal in the full up position ❸.

2.9 SELF POWERED MOTORISED (I-CARE EZ-GO)

2.9.1 INITIAL OPERATION

The storage and/or transport of the bed-chairs at very low temperatures (below 5°C) may impair the product's performance. This is why we recommend leaving the products for 24 hours in a warm environment before operating it for the first time in accordance with the conditions as described above.

2.9.2 BATTERY

The battery compartment is equipped with its own charger. We recommend, before operating the product for the 1st time, that you leave the product in a warm environment for 24 hours, before charging the battery for 5 hours from a mains outlet.

The battery allows the stretcher trolley to be operated in the event of a power cut or when AC power supply is not available during transfer.

In order to maintain an adequate charge level, we recommend that you connect the motorized stretchers to the mains supply whenever possible. You do not need to wait for the battery to fully discharge before recharging it.

“Warning: Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact. »

the appliance is removed from service for an extended period, we recommended that you recharge the battery once per month. After a period of 50 days of storage without recharging, the battery performance may be impaired or damaged irreversibly.

Battery characteristics		Battery safety symbols	
Battery type:	Lead	IPX6	Protection rating
Capacity:	1.3 Ah	CE	CE marking
Input supply voltage:	29~45 VDC		Indoor use
Input supply voltage:	24 Vcc, 10 A		Electronic waste

LED INDICATOR DISPLAY

Wheel display – *Image 1* as follows

		Green	The wheel is operating.
		Red	Error.
		Green	The wheel is charged.
		Flashing green	The wheel is charged but is still connected to the charger. It cannot be operated.
		Orange	The wheel is half charged.
		Flashing Orange	The wheel is charging
		Red	The wheel is charged at below one third of the full charge capacity.
		Flashing red	The wheel is fully discharged. The system cannot be operated. Recharge the wheel.
		Green	The wheel is communicating with the control unit
		Red	The control unit is not connected.
		Flashing red	The wheel is fully discharged. The communication between the wheel and control unit is disabled.


Control unit display – *Image 2* as follows

		Green	The control unit is operating
			The batteries of the wireless control module are discharged.
		Flashing red	There is an error with one of the wheels connected to the control unit.
		Green	The wheels are fully charged.
		Flashing green	The wheels are charged, but one of the system's wheels is connected to the charger. It cannot be operated.
		Orange	One of the system's wheels is half charged.
		Flashing Orange	One of the system's wheels is charging.
		Red	One of the system's wheels is charged at below one third of the full charge capacity.
		Flashing red	One of the system's wheels is fully discharged. The system cannot be operated. Recharge the wheel.
		Green	The control unit is communicating with all of the system's wheels.
		Red	The control unit is not connected.

RECHARGING THE BATTERY

The LED  indicates the charge status of the wheel.

The LED  indicates the charge status of the wheel.


The LED  indicates the charge status of the control unit.

LEDs   orange or   red → **RECHARGE THE BATTERY**

1. Connect the charger to the mains supply.


2. Unplug the connector from the charging station by pressing the safety lock.

3. Plug the connector into the wheel. Slide the latch into the groove on the socket.

LEDS  and  ● flashing green → **END OF CHARGING**

AN AUDIBLE ALARM SOUNDS AND THE ACTUATOR LOCK IS ACTIVATED IF YOU FORGET TO:

- **UNPLUG THE CONNECTOR FROM THE WHEEL**
- **PLUG THE CONNECTOR TO THE CHARGING STATION**

THE ACTUATOR LOCK IS ACTIVATED IF YOU FORGET TO UNPLUG THE DEVICE FROM THE MAINS SUPPLY

DO NOT ATTEMPT TO MOVE THE STRETCHER MANUALLY AND VERIFY THE ABOVE 3 POINTS

2.9.3 WHEELING THE STRETCHER

The electrical power traction drive using a self-propelled wheel facilitates the manoeuvrability of the stretcher.

1. Set the brake-steer pedal in the top position so that the motorised wheel is in contact with the floor.

2. Press the button located below the corner profile to activate the motorised wheel, the LED lights will then change to green.

3. Turn the knob to activate the stretcher.

Forward to move forward

Reverse to go backwards

Check regularly that there are no objects on the floor that may impair the rotation of the wheel.


2.9.4 BRAKE SYSTEM

Turn the knob back to stop the forward motion of the trolley, after releasing the knob the equipment is in brake mode. This will prevent the product from moving backwards.

2.10 TRENDLENBURG / REVERSE TRENDLENBURG


Trendelenburg position:

To tilt the stretcher trolley to the Trendelenburg position (head down), depress the yellow pedal.

Reverse Trendelenburg position:

To tilt the stretcher trolley to Reverse Trendelenburg position (foot down), depress the blue pedal.

2.11 CPR RELEASE MECHANISM

To lower the backrest quickly in an emergency situation, pull the red lever and lower the backrest to the full down position.


When using the manual CPR release, lower using the backrest remote control into the full down position before reinitiating operation.

Example of an electric trolley with knee gatch option :

To lower the patient to a horizontal position for emergency treatment, pull the red lever

When operating the manual CPR release do not place arms into the stretcher frame.

2.12 TILTING HEADREST (ONLY FOR THE ICARE OPHTHALMOLOGY MODEL)


Removable surgery head rest
for ophthalmology

Removable surgery head rest
for stomatology and ENT

2.12.1 INSTALLATION OF THE HEADREST

4. Raise the handle
5. Slide the headrest in the insert
6. Lock the headrest by putting the handle down

**Always ensure that the headrest is correctly locked in position.
The product should be set to sitting position before changing the headrest.**


The weight of the patient's head rests on the headrest; it must thus be supported by the operator during adjustment. Failure to observe this recommendation could lead to patient injury.

To avoid any risk of pinching during headrest adjustment, keep your fingers away from all articulated areas.



Do not place your hand between the backrest and headrest to pull the unlock handle as you may injure your fingers.

2.12.2 SURGERY HEADREST FOR OPHTHALMOLOGY


For ergonomic positioning, the headrest is equipped with two joints.

To adjust the first joint, push the right handle under the headrest outwards. Release the handles to lock (automatic locking) in the desired position.

To adjust the second joint, use the left-hand handle. Release the handles to lock (automatic locking) in the desired position.

The cushion is equipped with a magnet. Simple, quick, hygienic and adjustable positioning.

2.12.3 REMOVABLE SURGERY HEADREST FOR STOMATOLOGY AND ENT

In order to adjust the headrest:

1. To unlock the headrest, release the handle.
2. Spin the headrest to place it
3. Release the handle to lock the position


Be sure that the headrest is well locked.

2.13 OPTIONS

2.13.1 FOLDING PUSH BARS

To fold away the handles, lift up and slide the locking ring as far as it will go, the handles can then be swivelled downwards, away from the stretcher top.


To reposition the handles, move them back into place and place the locking ring back into the down position. Check that they are correctly positioned.

2.13.2 FOLD OUT MONITOR TRAY


The fold-out monitor tray at the foot-end of the stretcher trolley unfolds and folds by means of a pivoting mechanism.

To articulate the tray into its horizontal position ① to ④


①


②

3

4

When using Trendelenburg and reverse-Trendelenburg positioning, it is recommended that devices positioned on the fold-out monitor tray are attached using the strap provided for this purpose.

2.13.3 IV POLE

Safe working load: 6 kg.

This indicative value corresponds to safe product use, but is not representative of normal use for which the load is significantly lower.

Do not affix anything to the telescopic parts of the IV pole (e.g. tape) that could hinder its movement.




When moving the stretcher fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.

Do not move the stretcher via the IV pole.

2.13.3.1 PERMANENTLY ATTACHED IV POLE

Installation :


1. Remove the IV pole from its storage location ❶ and raise into the vertical position until it clicks into place ❷.
2. To raise the height of the pole, (3 telescopic stages), pull up to either of the two upper stages until it locks into place ❸.

Storing the IV pole:


In order to fold the telescopic pole back into place, pull up the red locking bolts located below each stage.

Do not tilt the IV pole in its highest position. Lower to each of the various stages in succession until you reach the lowest position (one level) before inserting the pole back into its storage location.


2.13.3.2 REMOVABLE IV POLE

Insert the detachable IV pole into one of the corner locations provided for this purpose.


2.13.4 OXYGEN BOTTLE HOLDER BENEATH THE BACKREST

Swing the oxygen bottle holder away from the stretcher top.


Bottle holder in fold-in position

Bottle holder in fold-out position

2.13.5 LATERAL TECHNICAL RAILS AND CLAMP


Slide the clamp onto which your accessory is attached to position and tighten the knob to lock in place.


Maximum safe working load: 15 kg.

2.13.6 HEADREST SUPPORT, FIXATION ON TECHNICAL RAILS

Insert the headrest in the support

2.13.7 PAPER ROLL DISPENSER

2.13.8 X RAY CASSETTE BELOW THE PATIENT SURFACE

2.13.9 ANTISTATIC SYSTEM


Make sure the earthing tape touches the floor, is not bended and is not deformed.


2.13.10 HOLDER FOR OPERATING FIELD WITH CONNECTION FOR OXYGEN (OPHTHALMOLOGY VERSION)

2.13.11 CHART HOLDER BELOW THE HEADREST

2.13.12 BAG HANGERS


2.13.13 MANUAL KNEE GATCH


To avoid the risk of entrapment, do not insert arms or hands between the stretcher top and the steel frame.

To raise the knee gatch, pull the yellow jack control ❶ while lifting the handle ❷ to the desired position.

To lower the knee gatch, pull the yellow jack control ❶ while pressing down on handle ❷.


2.13.14 ELECTRIC KNEE GATCH

To avoid the risk of entrapment, do not insert arms or hands between the stretcher top and the steel frame.

On the electric stretcher trolley versions, adjusting the knee gatch position is performed by remote control(s):

- the foot-end remote control
- and the patient's remote controls (optional))

3 CLEANING AND DISINFECTING

3.1 INTRODUCTION


Regular and correct cleaning is essential to ensure the satisfactory performance of the equipment (cf. a cleaning schedule should be established upon initial use of the device identifying the staff concerned).

If cleaning is not carried out regularly, there is a risk of infection caused by the presence of secretions and/or residues on the stretcher trolley.

The following cleaning instructions are not a substitute for the establishment's internal procedures.

The stretcher trolley is equipped with a removable cover cap to facilitate cleaning.

3.2 FREQUENCY

Cleaning must be carried out

- after each patient, that is, between two patients
- after transporting an infected patient
- in full on a weekly basis

3.3 CLEANING AND DISINFECTION ADVICE

3.3.1 CLEANING OPERATION

Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that cause stains (e.g.: “blue product” or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS

Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.

- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
 - o PVC, ABS, Polyamide, Polyurethane and Polypropylene
 - o Epoxy paint-coated metal surfaces
 - o Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

List of strictly prohibited products:

Solvents:	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
Bases:	10% Sodium hydroxide / 10% Ammonia
Acids:	10% Nitric / 30% Hydrochloric / 10% Citric
Miscellaneous:	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

High-pressure and steam cleaning are prohibited.

Any damage resulting from failure to observe this instruction for use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.

4 PRECAUTIONS

4.1 DEVICE-RELATED PRECAUTIONS

4.1.1 WARNINGS

This MD has been designed for and adapted to medical practices and procedures, and hence a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the common sense of the stretcher's users.

- Do not use the product without first having familiarized yourself with this user manual.
- Do not alter the stretcher. Any changes made to the system could lead to unpredictable operation that could in turn jeopardize the safety of the user or patient. This would also void the stretcher's warranty.
- Used to accommodate a person of small stature (less than 1.30 m tall), or under the age of 12 years without adapted equipment (pediatric kit)
- To ensure correct use of the product, only proceed if the weight is correctly distributed over the device.
- The chassis cover load capacity is of 35 kg. Do not sit or stand on the cover.
- Only purpose-built accessories provided for this MD by BMB MEDICAL may be used
- Connection, by any means, of an electric device to the stretcher.
- The product must be connected to a mains power outlet compliant with standards in force.
- Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact.
- In case of prolonged interruption of use of your equipment, it is recommended to recharge the battery once a month. After more than 50 days of storage without recharging, battery performance may be impaired or even irreversibly damaged.
- This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.
- Position the product such that it can be easily unplugged.
- Unplug the product before moving it.
- Do not leave the cables in contact with the floor or the casters.
- Do not knot the cables.
- Periodically check the condition of the cables.
- Work on electrical parts should only be performed by qualified and authorized staff.

- Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the [EMC information here](#).
- When using the remote control, all movements can be interrupted by simply releasing the corresponding key.
- Only activate the various positions after having ensured that no-one is touching any mechanical parts and that there are no objects under the seat or patient surface, or between any of the moving structures.
- The side rails are intended to minimize the risk of the patient accidentally falling. They are not, in any way, designed to immobilize or maintain the patient.
- Before manipulating the side rails, ensure that nothing could hinder their movement.
- When lowering the side rail, take care not to trap the hands or feet of the patient or staff
- If the side rail locking operation is silent (no "click"), please revise the side rails.
- Except during transport, always ensure that the stretcher's brake is engaged. Check whether the brake is engaged by attempting to move the stretcher.
- The "variable height" function should be used when handling the patient in order to work at an ergonomic height.
- Make sure that no one has their feet under the different pedals (variable height, Trendelenburg, Reverse Trendelenburg).
- It is recommended to attach the devices to the monitor support plate with the strap provided for this purpose when positioning the stretcher in the Trendelenburg or Reverse Trendelenburg position.
- When using the manual emergency release, do not put the arm in the frame of the lying surface.
- When activating the knee gatch, to avoid any risk of pinching, do not place your arm or hand between the lying surface and the metal structure.
- Do not park the chair on a sloping floor.
- Use outside of the establishment, or on unprepared ground.
- Do not use the MD on a floor with a slope of more than 10%.
- Use of the chair with a load that exceeds the safe working load (see technical characteristics).
- Always ensure that the headrest is correctly locked in position.
- The product should be set to sitting position before changing the headrest.
- The weight of the patient's head rests on the headrest; it must thus be supported by the operator during adjustment. Failure to observe this recommendation could lead to patient injury.
- To avoid any risk of pinching during headrest adjustment, keep your fingers away from all articulated areas.
- Do not place your hand between the backrest and headrest to pull the unlock handle as you may injure your fingers.
- Always ensure that the headrest articulation is correctly locked.
- IV pole safe working load: 6 kg.
- When moving the chair fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.

- 
- Do not move the chair via the IV pole
 - In the case of an unsupervised patient, the product should be left in its lowest position to reduce the risk of falling.
 - Tell the patient not to sit on the ends (foot or head) of the chair/stretcher. An excess weight applied to one of the ends could tip the patient surface and injure the patient.
 - Tell the patient and any accompanying persons not to sit on the armrests.
 - The equipotential connector must not, under any circumstances, be used as a protective earth connector
 - Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.
 - Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.
 - High-pressure and steam cleaning are prohibited.
 - Any use that does not fall within the operating context of a hospital chair.
 - To prevent failure, electromedical equipment requires special [EMC-related precautions](#).
 - Never use the chair in the presence of an anesthesia mixture flammable with air, with oxygen, or with nitrous oxide as this could cause a fire and/or explosion.
 - Use by any non-caregiver who has not received the authorization and explanations from a qualified individual.

All manipulations relating to the aforementioned precautions for use are described in [paragraph II - Instructions for use](#).

4.1.2 CONTRAINDICATIONS

There are no noteworthy contraindications for this device.

4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE

The device must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The device must be regularly inspected, at least once per year. See the maintenance check-list. To maintain the device in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed

In the event of observed anomalies or damage that could jeopardize device operation and safety, the equipment must be immediately withdrawn.

The device should only be used according to the instructions in this user manual.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

Work on electrical parts should only be performed by qualified and authorized staff.

* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

4.1.4 EMC-RELATED INFORMATION

The products in the I-CARE range are suitable for use in all healthcare establishments, with the exception of facilities in the vicinity of HF ELECTROSURGICAL APPLIANCES, or near an electromagnetically shielded room for an EM magnetic resonance imaging system, in which the intensity of EM INTERFERENCES is high.

The products in the I-CARE range are not suitable for use in domestic homes or establishments or on any premises connected to the public low-voltage network supplying domestic buildings.

EM interference could cause the product to make unexpected movements.

This appliance should not be used next to, or stacked on other appliances as this could lead to incorrect operation. If this mode of operation is necessary, this and the other appliances should be observed to ensure their correct operation.

The use of accessories, transducers or cables other than those specified or supplied by the manufacturer can lead to an increase in electromagnetic emissions, or a decrease in this appliance's immunity, causing inappropriate operation.

Do not use portable RF communication devices (including peripherals such as antenna cables and external antennas) less than 30 cm (12 inches) from any part of products in the I-CARE range, including the cables specified by the manufacturer. Otherwise, the performance of these devices could be impaired.

The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11, class A). If used in a residential environment (for which compliance

with the CISPR 11 class B standard is generally required), this equipment may not provide appropriate protection to radio frequency communication services. The user may need to take attenuation steps, such as moving or re-orienting the equipment.

List of cables used:

1. Power cord, maximum length 1000 cm.
2. Variable height column, maximum length 750 mm
3. Backrest cable: maximum length 1250 mm
4. Legrest cable: maximum length 750 mm
5. Footrest cable: maximum length 1500 mm

EMISSIONS	
Conducted and radiated RF emissions:	Conducted and radiated RF emissions:
Harmonic distortions IEC 61000-3-2:	Harmonic distortions IEC 61000-3-2:
Voltage fluctuations and Flicker IEC 61000-3-3:	Voltage fluctuations and Flicker IEC 61000-3-3:

Immunity test	Test level IEC 60601	Compliance level		
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV in contact ± 2/4/8/15 kV in air	± 8 kV in contact ± 2/4/8/15 kV in air		
RF radiated electromagnetic fields IEC 61000-4-3	3 V/m (professional healthcare establishment) 80MHz-2.7GHz 80% AM at 1 kHz	3V/m 80MHz-2.7GHz 80% AM at 1 kHz		
Proximity fields emitted by wireless RF communication devices IEC 61000-4-3				
	Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)
	385	Pulsed modulation: 18 Hz	27	27
	450	Pulsed modulation: 18 Hz	28	28
	710-745-780	Pulsed modulation: 217 Hz	9	9
	810-870-930	Pulsed modulation: 18 Hz	28	28
	1720-1845-1970	Pulsed modulation: 217 Hz	28	28
	2450	Pulsed modulation: 217 Hz	28	28
5240-5500-5785	Pulsed modulation: 217 Hz	9	9	
Rapid / burst electrical spikes IEC 61000-4-4	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz		
Shock waves IEC 61000-4-5	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV		
Conducted interference, induced by RF fields IEC 61000-4-6	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz		
Power frequency magnetic fields IEC 61000-4-8	30A/m	30A/m		
Voltage Dips and Interruptions: IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles		

NB: There is no deviation from the requirements of IEC 60601-1-2 ed. 4

NB: There are no other known measures for maintaining basic safety with respect to EMC phenomena.

4.2 TRANSPORT AND STORAGE

Storage and transport conditions:

Ambient conditions	Operation
Ambient temperature	-10 to +30 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

4.2.1 TRANSPORT

Before transporting the stretcher with a view to storing to ensure that there are no risks of impact or accidental falls during transport, we recommend that you:

- Fasten stretcher to the frame
- Place the cover over it
- Lock the brakes
- Place it in the lowest position

Damage to the appliance caused during transport is not covered by the guarantee. Repairs or replacement of damaged parts are the responsibility of the client.

4.2.2 STORAGE

When the chair is in storage, ensure that:

- Its brakes are engaged
- It is covered

The product should be stored in a clean, ventilated, temperate and dry place, protected from the elements and from direct sunlight.

During storage, no heavy weights should be placed on the chair.

The chair must not be considered as a support surface for any other type of equipment.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

In the event of prolonged storage, the battery should be recharged every 6 months, proceeding as follows:

- Connect the mains power cord.
- Briefly press one of the remote control functions to initiate charging.
- Allow the battery to charge for 8 hours.

4.3 WARRANTY

Warranty: two years parts and labor according to the following terms:

BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.

This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL stretchers are designed for a projected service life of 10 years under normal conditions of use, with appropriate periodic maintenance as described in the maintenance check-list of each stretcher. BMB MEDICAL provides the original buyer with the guarantee that the welds on the stretchers are free from structural defects for the projected stretcher service life of 10 years, as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use ([Warning paragraph](#))
- Work without prior notice to our technical services, or maintenance on:
 - A mechanical part of the chair
 - A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to observe the preventive maintenance operations recommended in the maintenance check-list.

4.4 AFTER SALES SERVICE

For any information relating to the installation, use and maintenance of the product and for any other questions, please contact the After-Sales Service at the following address: sav@bmb-medical.com or on **+33 (0)4 74 08 71 71** or **your local distributor if outside of France**. Thanks to our “Customer Space”, you can access your equipment items on-line, view the maintenance points to be attended to and request our services directly via this secure interface.

The operator must carefully read the instructions contained in this manual, paying particular attention to the safety precautions and the procedures to be followed for installation and correct use.

In the case of any doubts as to the correct interpretation of the instructions, please contact the After-Sales Service for any necessary clarifications.

To facilitate the ordering of spare parts, exploded parts drawings can be requested from the After-Sales Service to provide you with the references of the parts required.

4.5 ENVIRONMENT

4.5.1 SPENT EQUIPMENT

In France:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution. For this purpose, BMB MEDICAL fulfills its obligations concerning the end of life of I-care range products that it markets by funding the Recylum recycling channel dedicated to Pro WEEE that collects them free of charge (more information on www.recylum.com).

In Europe:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution.

Outside Europe :

Dispose of the Product or its components in accordance with local laws and regulations.

Contact the local distributor to inquire about return and/or collection systems available in the country.

Tél: +33 (0)4 74 08 71 71

Fax: +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001
01604 TREVOUX Cedex
FRANCE