MATERNA

MANUEL D'UTILISATION
USER MANUAL

 **BMB
MEDICAL**

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FRANÇAIS

BMB MEDICAL

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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www.bmb-medical.com

The purpose of this manual is to familiarize you with the use, maintenance and safety of your new medical device. We encourage you to read this document carefully and to follow the guidelines for optimum and safe use.

BMB MEDICAL offers a broad variety of options for the products covered by this user manual. Consequently, some features described in this manual may not apply to your product.

Some illustrations are provided for information purposes and do not accurately represent the device in your possession.

The information in this document are subject to change without notice and in no way engage the responsibility of BMB MEDICAL. This document may not be reproduced, in whole or in part, without prior written consent by BMB MEDICAL. This illicit representation or reproduction, by any means, constitutes counterfeiting. The illustrations in this document are not binding.

TABLE OF CONTENTS

| | |
|---|----|
| ENGLISH | 17 |
| 1 INTRODUCTION SPECIFICATIONS | 19 |
| 1.1 PRODUCT PRESENTATION | 19 |
| 1.1.1 CRADLE | 19 |
| 1.2 PURPOSE OF THE USER MANUAL..... | 19 |
| 1.3 USERS | 20 |
| 1.4 STORAGE..... | 20 |
| 1.5 SYMBOLS AND STANDARDS USED | 20 |
| 2 INSTRUCTIONS FOR USE | 22 |
| 2.1 INSTALLATION PROCEDURE | 22 |
| 2.2 BRAKING AND MOVING | 22 |
| 2.3 HEIGHT ADJUSTMENT | 22 |
| 2.4 REVERSE TRENDELENBURG..... | 23 |
| 2.5 REMOVABLE CARRIER..... | 24 |
| 2.5.1 U-shaped stand | 24 |
| 3 CLEANING AND DISINFECTION | 25 |
| 3.1 INTRODUCTION | 25 |
| 3.2 FREQUENCY | 25 |
| 3.3 CLEANING AND DISINFECTION ADVICE..... | 25 |
| 3.3.1 CLEANING OPERATION | 25 |
| 3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS..... | 26 |
| 4 PRECAUTIONS..... | 27 |
| 4.1 DEVICE-RELATED PRECAUTIONS..... | 27 |
| 4.1.1 WARNINGS | 27 |
| 4.1.2 CONTRAINDICATIONS | 27 |
| 4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE | 28 |
| 4.2 TRANSPORT AND STORAGE | 29 |
| 4.2.1 TRANSPORT | 29 |
| 4.2.2 STORAGE | 29 |
| 4.3 WARRANTY | 29 |
| 4.4 AFTER-SALES SERVICE | 31 |
| 4.5 ENVIRONMENT | 31 |
| 4.5.1 SPENT EQUIPMENT..... | 31 |

1 INTRODUCTION SPECIFICATIONS

1.1 PRODUCT PRESENTATION

Infant support device, composed of a shell mounted on a wheel frame. The device is equipped with an independent braking wheel system.

The device is fitted with a lever-actuated cylinder-powered height adjustment system. Reverse Trendelenburg is actuated by an indexer.

The product is designed to ensure that there are no unacceptable risks (essential performance) related to patient support and the risk of unexpected movements.

1.1.1 CRADLE

| TECHNICAL DATA | |
|---------------------------|--|
| REFERENCE | MATERNA |
| DESIGNATION | Cradle |
| SAFE WORKING LOAD | 9 kg (indicates the maximum admissible load including the total weight of the patient and accessories) |
| WEIGHT | 20 kg depending on the model |
| OVERALL DIMENSIONS | 105 x 58 Height 70 cm to 100 cm |
| EQUIPMENT | <ul style="list-style-type: none">* Overall dimensions : 105 cm x 58 cm* Removable cradle made of transparent thermoplastic polymer 75 cm x 35 cm x 29 cm* Pneumatic variable height from 70 cm to 100cm* Trendelenburg 10 °* Structure made of steel epoxy coated with* Drilled bottom for better ventilation* 4 self lockable double running castors (Ø 80 mm)* «U» base* M1 fire-resistant mattress 5 cm, removable cover |

1.2 PURPOSE OF THE USER MANUAL

Healthcare professionals, this user manual was created for you to provide you with the necessary information and precautions to observe for autonomous and safe use.

The product's service life is dependent upon observance of the maintenance and tracking instructions.

It is thus essential that all staff has access to and be informed of the handling procedures and risks associated with the use of this device. The latter should be handled by a qualified person who is aware of the risks inherent to improper use. Additional training is not required as we consider that knowledge of the user manual is sufficient.

1.3 USERS

Healthcare professionals: physicians, caregivers and agents approved by Technical Services can operate all functions.

In use cases where the patient is the operator, the nursing staff must judge the patient's ability to use the product. The nursing staff must also train the patient in the safe use of the product.

1.4 STORAGE

This user manual should be kept near the product, in a place provided for this purpose. It must be shielded from all elements and substance that could compromise its perfect legibility.

1.5 SYMBOLS AND STANDARDS USED

| | |
|---|---|
|  | See the manual / instruction brochure |
|  | Indicates the maximum permissible load, including the total weight of the patient, mattress (where applicable) and accessories. |
|  | Conforms to the essential requirements of directive 93/42/EEC on medical devices. |
|  | Serial number |
|  | Product reference |
|  | Manufacture date |
|  | Manufacturer |
|  | Maximum patient weight |
|  | Keep dry |
|  | Top, stand up |
|  | Fragile |
|  | Recyclable |
|  | Temperature range |

| | |
|---|----------------|
|  | Humidity limit |
|---|----------------|

General safety sign:

- Information highlighted by this symbol must be read and scrupulously observed.
- Indicates that failure to observe the guideline could endanger the patient or user, or damage the equipment.

“WARNING: This appliance must not be altered without the manufacturer’s written consent.”

Caution: only purpose-built accessories provided for this MD by BMB MEDICAL may be used.

The MATERNA is a Class I non-invasive and non-active medical device.

2 INSTRUCTIONS FOR USE

2.1 INSTALLATION PROCEDURE

Before commissioning the product, ensure, depending on the model, that the following are operational:

- **Casters with independent brakes**
- **Variable height**
- **Carrier reverse Trendelenburg**

Clean the product according to your establishment's internal protocol.

2.2 BRAKING AND MOVING

Except during transport, always ensure that the cradle's brakes are engaged. Ensure that the brake is engaged by attempting to move the cradle.

Engage the wheel brakes as soon as the product is immobilized.

Do not park the product on a sloping floor.


Wheel in free position

Wheel in locked position

2.3 HEIGHT ADJUSTMENT

To adjust the height, activate the cylinder using the lever (travel: 29 cm)

2.4 REVERSE TRENDELENBURG

The cradle is placed in reverse Trendelenburg position using the bouton at the top of the product's column

Bed in horizontal position

Bed in reverse Trendelenburg position max. 10°

Placing the bed in reverse Trendelenburg

Ensure that the carrier is locked by shaking it gently.

2.5 REMOVABLE CARRIER

The carrier can be removed from the frame; to do so, pull on the carrier

Do not remove or re-fit the carrier with the infant in it.

Ensure that the hooks are securely clipped by lifting the carrier.

2.5.1 U-shaped stand

Casters with independent brakes

3 CLEANING AND DISINFECTION

3.1 INTRODUCTION

Regular and thorough cleaning will guarantee the device's life (see cleaning program drawn up at the time of device commissioning, designating the reference staff).

Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the device.

The following cleaning advice is not intended to replace your establishment's in-house protocol.

This device has been designed to facilitate cleaning.

3.2 FREQUENCY

The chair should be cleaned:

- Between each patient.
- Completely each week

3.3 CLEANING AND DISINFECTION ADVICE

3.3.1 CLEANING OPERATION

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that cause stains (e.g.: “blue product” or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS

- Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.
- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
 - o PVC, ABS, Polyamide, Polyurethane and Polypropylene
 - o Epoxy paint-coated metal surfaces
 - o Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

List of strictly prohibited products:

| | |
|-----------------------|--|
| Solvents: | Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene |
| Bases: | 10% Sodium hydroxide / 10% Ammonia |
| Acids: | 10% Nitric / 30% Hydrochloric / 10% Citric |
| Miscellaneous: | Pure or weakly diluted bleach / Abrasive powders / All abrasive products |

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

High-pressure and steam cleaning are prohibited.

Any damage resulting from failure to observe this Instructions for Use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.

4 PRECAUTIONS

4.1 DEVICE-RELATED PRECAUTIONS

4.1.1 WARNINGS

This device is designed for and adapted to medical practices and procedures, therefore to a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.

The following list of abnormal device uses is not exhaustive and BMB MEDICAL counts on the users' common sense.

- Not for use with children over the age of 6 months (or over 10 kg).
- Use of the device with a load that exceeds the safe working load (see technical characteristics).
- This product should no longer be used for a child that can sit, kneel or stand up on their own.
- The product should not be used if any part is broken, torn or missing.
- When the product is used, the casters must always be locked.
- All adjustment, fixing and positioning of the product must be carried out exclusively by an adult.
- Do not place any objects inside the product that could reduce its depth.
- Use only one mattress in the cradle.
- Do not adjust the product when the child is inside.
- Do not move the product when the child is inside.
- Do not use the product with more than one child at a time.
- The cradle should be placed on a flat, level floor.
- Do not use the MD on a floor with a slope of more than 5°.
- Use only the mattress supplied with the product.
- To avoid any risk of fire, do not place the cradle near an intense heat source.
- Do not alter the product. Any changes made to the system could lead to unpredictable operation that could in turn jeopardize the safety of the user or patient. This would void the device's warranty.
- Use by a person who has not first familiarized themselves with this user manual.
- Washing with large amounts of water, with a high-pressure jet or with steam and passage through a wash tunnel.
- Any use that does not fall within the operating context of a hospital cradle.

All manipulations relating to the aforementioned precautions for use are described in paragraph II - Instructions for Use [HERE](#).

4.1.2 CONTRAINDICATIONS

There are no other contraindications for this device.

4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE

The device must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn mattress upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The device must be regularly inspected, at least twice per year.

To maintain the device in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed:

- Ensure that all wheels are securely attached.
- Visually check all fastenings.
- Check the cylinder fastenings.
- Ensure that the variable height, reverse Trendelenburg and wheel braking all function correctly.
- Lubricate all the moving mechanisms.

To keep the system operating optimally, repeat this check regularly.

If you have any doubts about the functioning and integrity of the product, take it out of service and inform your technical service or BMB MEDICAL after-sales service.

In the event of observed anomalies or damage that could jeopardize device operation and safety, the equipment must be immediately withdrawn.

The device should only be used according to the instructions in this user manual.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

Work on electrical parts should only be performed by qualified and authorized staff.

* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

4.2 TRANSPORT AND STORAGE

Storage and transport conditions:

| Ambient conditions | Operation |
|----------------------|--------------------------------------|
| Ambient temperature | -10 to +50 °C |
| Relative humidity | 20 to 90% at 30 °C - no condensation |
| Atmospheric pressure | 700 to 1060 ha |

4.2.1 TRANSPORT

During device transport with a view to its storage, the following is recommended to prevent any impacts or falls:

- Strap it to a surface of the means of transport
- Cover it
- Engage the brakes
- Set it to its low position

Damage caused during transport is not covered by the warranty. The costs of repair and replacement of damaged parts will be borne by the customer.

4.2.2 STORAGE

When the product is in storage, ensure that:

- Its brakes are engaged
- It is covered

The product should be stored in a clean, ventilated, temperate and dry place, protected from the elements and from direct sunlight.

During storage, no heavy weights should be placed on the chair.

The product must not be considered as a support surface for any other type of equipment.

The product must not be exposed to and even less be in contact with, sources of heat, combustion and flammable agents.

4.3 WARRANTY

Warranty: two years parts and labor according to the following terms:

BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.

This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL products are designed for an expected service life of 10 years under normal conditions of use and with appropriate periodic maintenance as described above. BMB MEDICAL guarantees the original buyer that the welds on its products will be free from structural defects during the expected 10-year service life of the product as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use (*Warning paragraph [here](#)*)
- Work without prior notice to our technical services, or maintenance on:
 - o A mechanical part of the product
 - o A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to comply with the recommended preventive maintenance operations.

4.4 AFTER-SALES SERVICE

For any information concerning product installation, use and maintenance, the technical documentation (wiring diagrams, list of components, maintenance procedures, etc.), or any other questions, please contact the After-Sales Service at the following address: **sav@bmb-medical.com**, or by phone at **+33 (0)4 74 08 71 71**, or your local distributor for overseas purchases.

From your “Customer Area” on the website www.bmb-medical.com, you can access your equipment, submit work requests and track these requests directly via this secure interface.

The operator must carefully read the information contained in this user manual, in particular the information concerning the precautions relating to safety and installation/use methods.

In the event of doubt concerning the correct interpretation of the instructions, do not hesitate to contact the After-Sales Service whose staff are always available.

To facilitate ordering spare parts, you can request the exploded diagrams of your products from the After-Sales Service. You can thus obtain the references of the desired parts.

4.5 ENVIRONMENT

4.5.1 SPENT EQUIPMENT

To be disposed of separately from household waste.

Dispose of the product or its components in accordance with local and national laws and regulations on environmental protection and recycling of raw materials

You can obtain the appropriate disposal location from the authorities in charge of environmental protection.

France specificities: BMB MEDICAL fulfills its obligations concerning the end of life of the products in the MATERNA range, which it puts on the market by financing the recycling sector of Valdelia which takes them back free of charge (More information on www.valdelia.org).

Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001

01604 TREVOUX Cedex

FRANCE

www.bmb-medical.com