# SYREETA ARETHA

MANUEL D'UTILISATION  
USER MANUAL

 **BMB  
MEDICAL**

FRANÇAIS

ENGLISH

# FRANÇAIS

---

## **BMB MEDICAL**

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

[www.bmb-medical.com](http://www.bmb-medical.com)

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[www.bmb-medical.com](http://www.bmb-medical.com)

The purpose of this manual is to familiarize you with the use, maintenance and safety of your new medical device. We encourage you to read this document carefully and to follow the guidelines for optimum and safe use.

**BMB MEDICAL** offers a broad variety of options for the products covered by this user manual. Consequently, some features described in this manual may not apply to your chair.

Some illustrations are provided for information purposes and do not accurately represent the device in your possession.

The information in this document are subject to change without notice and in no way engage the responsibility of BMB MEDICAL. This document may not be reproduced, in whole or in part, without prior written consent by BMB MEDICAL. This illicit representation or reproduction, by any means, constitutes counterfeiting. The illustrations in this document are not binding.

## TABLE OF CONTENTS

ENGLISH .....	17
1 INTRODUCTION SPECIFICATIONS .....	19
1.1 PRODUCT PRESENTATION .....	19
1.1.1 SYREETA SOFA BED .....	19
1.1.2 ARETHA SOFA BED .....	20
1.2 PURPOSE OF THE USER MANUAL .....	20
1.3 STORAGE .....	20
1.4 SYMBOLS AND STANDARDS USED .....	21
2 INSTRUCTIONS FOR USE .....	22
2.1 INSTALLATION PROCEDURE .....	22
2.2 OPEN/CLOSE CHAIR .....	22
2.3 BRAKING & MOVEMENT .....	23
2.4 OPTIONS .....	24
2.4.1 POWER-ASSISTED ADJUSTABLE BACKREST .....	24
2.4.2 READING LAMP AND USB PORT .....	24
2.4.3 HEADREST .....	25
2.4.4 ARMREST CUSHIONS .....	25
2.4.5 STORAGE BOX .....	25
3 CLEANING AND DISINFECTION .....	26
3.1 INTRODUCTION .....	26
3.2 FREQUENCY .....	26
3.3 CLEANING AND DISINFECTION ADVICE .....	26
3.3.1 CLEANING OPERATION .....	26
3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS .....	27
4 PRECAUTIONS .....	28
4.1 DEVICE-RELATED PRECAUTIONS .....	28
4.1.1 WARNINGS .....	28
4.1.2 CONTRAINDICATIONS .....	29
4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE .....	29
4.2 TRANSPORT AND STORAGE .....	29
4.2.1 TRANSPORT .....	29
4.2.2 STORAGE .....	30
4.3 WARRANTY .....	30
4.4 AFTER-SALES SERVICE .....	31
4.5 ENVIRONMENT .....	31
4.5.1 SPENT EQUIPMENT .....	31

# 1 INTRODUCTION SPECIFICATIONS

## 1.1 PRODUCT PRESENTATION

The SYREETA and FERGUS sofa beds are as comfortable as a lounge chair. Ideal for fathers accompanying their wives during maternity stays, this multi-functional convertible armchair offers pleasant seating and the comfort of a real bed.

### 1.1.1 SYREETA SOFA BED

TECHNICAL DATA	
<b>REFERENCE</b>	SYREETA
<b>DESIGNATION</b>	COMPANION CHAIR-BED
<b>SAFE WORKING LOAD</b>	160 kg (indicates the maximum admissible load including the total weight of the patient and accessories)
<b>WEIGHT</b>	65 kg depending on the model
<b>OVERALL DIMENSIONS</b>	<b>Width:</b> 70 cm <b>Overall length:</b> 90 cm (folded) 230 cm (open)
<b>EQUIPMENT</b>	<ul style="list-style-type: none"><li>* Fold-out sofa bed - steel frame</li><li>* Thin armrests</li><li>* Retractable panels and cushions for cleaning</li><li>* Mattress dimensions : 70 cm x 188 cm</li><li>* High resiliency complex foam seat 35 Kg/m<sup>3</sup> dacron + HR</li><li>* 4 self lockable casters ø 80 mm</li><li>* M2 fire-resistant coating Valencia, 10 colours</li><li>* Removable mattress made of articulated foam (35 kg/m<sup>3</sup>), thickness 10 cm.</li><li>* Bed base: 13 resin slats for better hygiene and easy cleaning</li></ul>
<b>AVAILABLE OPTIONS</b>	Reading lamp + USB port Light beams Fixed headrest 2 soft cushions Power backrest with adjustable angle Storage box

## 1.1.2 ARETHA SOFA BED

TECHNICAL DATA	
<b>REFERENCE</b>	ARETHA
<b>DESIGNATION</b>	COMPANION CHAIR-BED
<b>SAFE WORKING LOAD</b>	160 kg (indicates the maximum admissible load including the total weight of the patient and accessories)
<b>WEIGHT</b>	65 kg depending on the model
<b>OVERALL DIMENSIONS</b>	<b>Width:</b> 89 cm <b>Overall length:</b> 100 cm (folded) 240 cm (open)
<b>EQUIPMENT</b>	<ul style="list-style-type: none"><li>* Armchair convertible into bed with wooden armrests and steel frame covered with epoxy.</li><li>* All panels and cushions are retractable for cleaning or repairs.</li><li>* Mattress dimensions : 70 cm x 188 cm</li><li>* Independent seat</li><li>* Full armrests</li><li>* Seat made from 35 kg/m<sup>3</sup> high resiliency complex foams (Dacron and HR)</li><li>* 4 wheels (ø 80 mm) with independent brakes</li><li>* M2 flame retardant cover, 10 color options</li><li>* HR foam mattress with a density of 35 kg/m<sup>3</sup>, 10 cm thick</li><li>* Bed base: 13 resin slats for better hygiene and easy cleaning</li></ul>
<b>AVAILABLE OPTIONS</b>	Reading lamp + USB port Fixed headrest Reclining backrest

## 1.2 PURPOSE OF THE USER MANUAL

We created this user manual for you in order to provide you with all the necessary information and all the usual precautions to be observed for fully autonomous and completely safe use.

The product's service life is dependent upon observance of the maintenance and tracking instructions.

## 1.3 STORAGE

This user manual should be kept near the product, in a place provided for this purpose. It must be shielded from all elements and substance that could compromise its perfect legibility.

## 1.4 SYMBOLS AND STANDARDS USED

	See the manual / instruction brochure		Protéger de l'humidité / keep dry
	Indicates the maximum permissible load, including the total weight of the patient, mattress (where applicable) and accessories.		Haut, tenir debout / stand up
	Serial number		Fragile / fragile
	Product reference		Recyclable / recycling
	Manufacture date		Temperature range
	Manufacturer		Humidity limit

### General safety sign:

- Information highlighted by this symbol must be read and scrupulously observed.
- Indicates that failure to observe the guideline could endanger the user, or damage the equipment.

**“WARNING: This appliance must not be altered without the manufacturer’s written consent.”**

## 2 INSTRUCTIONS FOR USE

---

### 2.1 INSTALLATION PROCEDURE

Before commissioning the chair and depending on the model, ensure that the product functions correctly.

Clean the product according to your establishment's internal protocol.

### 2.2 OPEN/CLOSE CHAIR

- 1 Tilt the armrest cushions outwards

- 2 Facing the chair, grab the tube located under the seat cushion. Pull the seat towards you.

- 3 Unfold the bed base

- 4 Repeat the operation in reverse to close the chair.

In order not to damage the cover, the mattress must be put in the axis when the chair is folded. This will also allow the base part to be reattached so as not to damage the upholstery.

## 2.3 BRAKING & MOVEMENT

Wheel in locked position

Wheel in free position

**Except during transport, always ensure that the chair's brake is engaged. Check whether the brake is engaged by attempting to move the chair.**

**Make sure to lock all 4 wheels.**

**Do not park the chair on a sloping floor.**

**Do not move the chair with the wheels locked.**

**Do not move the chair when the bed mechanism is unfolded.**

**Do not move the chair with the wheel brakes locked.**

**Do not move the chair with a person on it.**

**NOTE: the wheels are not intended for long distance movement, simply to allow the chair to be moved for cleaning.**

## 2.4 OPTIONS

### 2.4.1 POWER-ASSISTED ADJUSTABLE BACKREST

Before actuating the backrest, ensure that there are no obstacles to its movement (accessories, cables, etc.)

Do not sit on the backrest.

#### **SYREETA backrest adjustment:**

The backrest is adjusted with the trigger:


### 2.4.2 READING LAMP AND USB PORT

The lamp plugs into an AC outlet.

The lamp is turned on and off simply by touching the silver part.

Light intensity is controlled by prolonged contact with the silver part.

### 2.4.3 HEADREST

### 2.4.4 ARMREST CUSHIONS

The cushions can be removed and repositioned using snap buttons.

### 2.4.5 STORAGE BOX

## 3 CLEANING AND DISINFECTION

---

### 3.1 INTRODUCTION

Regular and thorough cleaning will guarantee the device's life (see cleaning program drawn up at the time of device commissioning, designating the reference staff).

**Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.**

The following cleaning advice is not intended to replace your establishment's in-house protocol.

This chair has been designed to facilitate cleaning.

### 3.2 FREQUENCY

The chair should be cleaned:

- For each patient/accompanying person, i.e. between patients.
- Completely each week

### 3.3 CLEANING AND DISINFECTION ADVICE

#### 3.3.1 CLEANING OPERATION

**Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.**

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that can cause stains (e.g.: "blue product" or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

In the case of a product fitted with a USB charging port, ensure that no liquids enter the port.

### 3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS

- Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.
- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
  - o PVC, ABS, Polyamide, Polyurethane and Polypropylene
  - o Epoxy paint-coated metal surfaces
  - o Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

#### List of strictly prohibited products:

<b>Solvents:</b>	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
<b>Bases:</b>	10% Sodium hydroxide / 10% Ammonia
<b>Acids:</b>	10% Nitric / 30% Hydrochloric / 10% Citric
<b>Miscellaneous:</b>	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

**High-pressure and steam cleaning are prohibited.**

**Any damage resulting from failure to observe this Instructions for Use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.**

## 4 PRECAUTIONS

---

### 4.1 DEVICE-RELATED PRECAUTIONS

#### 4.1.1 WARNINGS

**This chair has been designed for and adapted to medical practices and procedures, and hence a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.**

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the common sense of the chair's users.

- The chair's mechanisms could cause lesions; ensure that no-one is touching the mechanical elements before activating the chair.
- Do not use the product on a floor with a slope greater than 10%.
- Do not alter the chair. Any changes made to the system could lead to unpredictable operation that could in turn jeopardize the safety of the user or patient. This would void the product's warranty.
- Use by a person who has not first familiarized themselves with this user manual.
- Washing with large amounts of water, with a high-pressure jet or with steam and passage through a wash tunnel.
- Any use that does not fall within the operating context of a hospital chair.
- Using the chair with a load greater than that recommended.
- Do not overload the chair, even if unfolded.
- Moving the chair when the wheels are locked.
- Impacts of varying force.
- Sitting on the armrests.
- Storing the chair in a damp (risk of component oxidation) or extremely sunny environment (risk of deterioration of components such as upholstery and risk of burns to the patient by contact on the upholstery).
- Moving the chair when the bed mechanism is unfolded.
- Moving the chair with the wheel brakes engaged.
- The bed mechanism is intended for use by one person only.
- Moving the chair with a person on it.
- Using the chair as a transport trolley.

**All manipulations relating to the aforementioned precautions for use are described in [paragraph II - Instructions for Use HERE](#).**

## 4.1.2 CONTRAINDICATIONS

There are no other contraindications for this device.

## 4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE

The device must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The device must be regularly inspected, at least twice per year.

To maintain the device in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed:

- Ensure that all wheels are securely attached.
- Visually check all fastenings.
- Check the cylinder fastenings.
- Ensure that the backrest functions correctly and is securely attached.
- Lubricate all the moving mechanisms.

In the event of observed anomalies or damage that could jeopardize device operation and safety, the equipment must be immediately withdrawn.

The device should only be used according to the instructions in this user manual.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

Work on electrical parts should only be performed by qualified and authorized staff.

\* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

## 4.2 TRANSPORT AND STORAGE

### 4.2.1 TRANSPORT

During chair transport with a view to its storage, the following is recommended to prevent any impacts or falls:

- Strap it to a surface of the means of transport
- Cover it
- Engage the brakes
- Set it to its low position

Damage caused during transport is not covered by the warranty. The costs of repair and replacement of damaged parts will be borne by the customer.

## 4.2.2 STORAGE

When the chair is in storage, ensure that:

- Its brakes are engaged
- It is covered

The product should be stored in a clean, ventilated, temperate and dry place, protected from the elements and from direct sunlight.

During storage, no heavy weights should be placed on the chair.

The chair must not be considered as a support surface for any other type of equipment.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

## 4.3 WARRANTY

Warranty: two years parts and labor according to the following terms:

*BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.*

*This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.*

*Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.*

*Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.*

*BMB MEDICAL products are designed for a projected service life of 10 years under normal conditions of use and with appropriate periodic maintenance as described in the maintenance check-list of each product. BMB MEDICAL guarantees the original buyer that the welds on its products will be free from structural defects during the expected 10-year service life of the product as long as the original buyer is in possession of the product.*

*No guarantees are provided concerning appropriateness for a specific use.*

*This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.*

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use ([Warning paragraph here](#))
- Work without prior notice to our technical services, or maintenance on:
  - A mechanical part of the product
  - A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to observe the preventive maintenance operations recommended in the maintenance check-list.

## 4.4 AFTER-SALES SERVICE

For any information concerning product installation, use and maintenance, the technical documentation (wiring diagrams, list of components, maintenance procedures, etc.), or any other questions, please contact the After-Sales Service at the following address: **sav@bmb-medical.com**, or by phone at **+33 (0)4 74 08 71 71**, or your local distributor for overseas purchases.

From your “Customer Area” on the website [www.bmb-medical.com](http://www.bmb-medical.com), you can access your equipment, submit work requests and track these requests directly via this secure interface.

The operator must carefully read the information contained in this user manual, in particular the information concerning the precautions relating to safety and installation/use methods.

In the event of doubt concerning the correct interpretation of the instructions, do not hesitate to contact the After-Sales Service whose staff are always available.

To facilitate ordering spare parts, you can request the exploded diagrams of your products from the After-Sales Service. You can thus obtain the references of the desired parts.

## 4.5 ENVIRONMENT

### 4.5.1 SPENT EQUIPMENT

To be disposed of separately from household waste.

Dispose of the product or its components in accordance with local and national laws and regulations on environmental protection and recycling of raw materials

You can obtain the appropriate disposal location from the authorities in charge of environmental protection.

**France specificities:** BMB MEDICAL fulfills its obligations concerning the end of life of the products in the SYREETA & FERGUS range, which it puts on the market by financing the recycling sector of Valdelia which takes them back free of charge (More information on [www.valdelia.org](http://www.valdelia.org)).

Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

[info@bmb-medical.com](mailto:info@bmb-medical.com)

100 Av. du Formans – CS 11001

01604 TREVOUX Cedex

FRANCE