# CLEM

MANUEL D'UTILISATION  
USER MANUAL

 **BMB**  
**MEDICAL**

 [FRANÇAIS](#)  
[ENGLISH](#)

## **BMB MEDICAL**

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

[www.bmb-medical.com](http://www.bmb-medical.com)

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[www.bmb-medical.com](http://www.bmb-medical.com)

The purpose of this manual is to familiarize you with the use, maintenance and safety of your new product. We encourage you to read this document carefully and to follow the guidelines for optimum and safe use.

**BMB MEDICAL** offers a broad variety of options for the products covered by this user manual. Consequently, some features described in this manual may not apply to your chair.

Some illustrations are provided for information purposes and do not accurately represent the product in your possession.

The information in this document are subject to change without notice and in no way engage the responsibility of BMB MEDICAL. This document may not be reproduced, in whole or in part, without prior written consent by BMB MEDICAL. This illicit representation or reproduction, by any means, constitutes counterfeiting. The illustrations in this document are not binding.

## TABLE OF CONTENTS

---

ENGLISH .....	17
1 INTRODUCTION SPECIFICATIONS .....	19
1.1 PRODUCT PRESENTATION .....	19
1.1.1 CLEM FOLDING RECLINER ARMCHAIR .....	19
1.2 PURPOSE OF THE USER MANUAL.....	19
1.3 STORAGE.....	20
1.4 SYMBOLS AND STANDARDS USED .....	20
2 INSTRUCTIONS FOR USE .....	22
2.1 INSTALLATION PROCEDURE .....	22
2.2 PRODUCT FOLDING AND UNFOLDING.....	22
2.3 POSITION ADJUSTMENT .....	24
2.4 RETRACTABLE HEADREST .....	25
3 CLEANING AND DISINFECTION .....	26
3.1 INTRODUCTION .....	26
3.2 CLEANING AND DISINFECTION ADVICE.....	26
3.2.1 CLEANING OPERATION .....	26
3.2.2 RECOMMENDED EQUIPMENT AND PRODUCTS.....	26
4 PRECAUTIONS.....	28
4.1 DEVICE-RELATED PRECAUTIONS .....	28
4.1.1 WARNINGS .....	28
4.1.2 CONTRAINDICATIONS .....	28
4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE .....	28
4.1.4 STORAGE .....	29
4.2 WARRANTY .....	29
4.3 AFTER-SALES SERVICE .....	30
4.4 ENVIRONMENT .....	31

# 1 INTRODUCTION SPECIFICATIONS

## 1.1 PRODUCT PRESENTATION

Multi-position relax chair.

The CLEM is a relaxing chair for accompanying persons, maternity units or even for the staff.

Thanks to its generous dimensions, the CLEM offers ease of use and space savings.

### 1.1.1 CLEM FOLDING RECLINER ARMCHAIR

TECHNICAL DATA	
<b>REFERENCE</b>	CLEM
<b>DESIGNATION</b>	CLEM FOLDING RECLINER ARMCHAIR
<b>SAFE WORKING LOAD</b>	140 kg (indicates the maximum admissible load including the total weight of the patient and accessories)
<b>WEIGHT</b>	10 kg
<b>OVERALL DIMENSIONS</b>	<b>Seat width:</b> 52 cm. <b>Overall length:</b> from 90 cm in sitting position to 163 cm in reclined position. <b>Overall height:</b> 115 cm (chair position) 86 cm (lying position) 98 cm (folded)
<b>EQUIPMENT</b>	* Manual folding relaxation chair * Reclining, with ultra-comfort position : places the legs slightly above the heart to provide optimal well-being * Steel structure ø 20 mm * Foam HR 25 Kg - Thickness 7cm * Fixed black polypropylene armrests * Retractable headrest * High stability * Foldable * VALENCIA color range : M2 fireproof coating, 10 colors

## 1.2 PURPOSE OF THE USER MANUAL

We created this user manual for you in order to provide you with all the necessary information and all the usual precautions to be observed for fully autonomous and completely safe use.

The product's service life is dependent upon observance of the maintenance and tracking instructions.

It is thus essential that all staff has access to and be informed of the handling procedures and risks associated with the use of this chair. The latter should be handled by a qualified person who is aware of the risks inherent to improper use. Additional training is not required as we consider that knowledge of the user manual is sufficient.

### 1.3 STORAGE

This user manual should be kept near the product, in a place provided for this purpose. It must be shielded from all elements and substance that could compromise its perfect legibility.

### 1.4 SYMBOLS AND STANDARDS USED

	See the manual / instruction brochure
	Indicates the maximum permissible load, including the total weight of the patient, mattress (where applicable) and accessories.
	Conforms to the essential requirements of directive 93/42/EEC on medical devices.
	Serial number
	Product reference
	Manufacture date
	Manufacturer
	Keep dry
	Top, stand up
	Fragile
	Recyclable
	Temperature range
	Humidity limit

#### General safety sign:

- Information highlighted by this symbol must be read and scrupulously observed.
- Indicates that failure to observe the guideline could endanger the user, or damage the equipment.

	Sitting on the backrest is prohibited		Leaning against the chair is prohibited
	Sitting on the legrest is prohibited		Sitting on the armrests is prohibited

**“WARNING: This appliance must not be altered without the manufacturer’s written consent.”**

**Caution: only purpose-built accessories provided for this MD by BMB MEDICAL may be used.**

## 2 INSTRUCTIONS FOR USE

---



### 2.1 INSTALLATION PROCEDURE

Before commissioning the chair, check that it is operational.  
Clean the product according to your establishment's internal protocol.

### 2.2 PRODUCT FOLDING AND UNFOLDING

#### **Unfolding:**


- A. Place the CLEM on a flat surface with the base on the floor and the backrest up.
- B. Unfold the CLEM by leaning on the end of the backrest and the armrest.
- C. Pull the backrest backwards and the armrest downwards until the legrest is completely unfolded.
- D. Place the product in vertical position.

### **Folding:**

- E. Place the locking system in the neutral position
- F. Place the product in vertical position.
- G. Push the armrests towards the backrest until they are completely folded.
- H. Lock the position

**Beware of pinching and shear risks during installation.**

**To use this item safely, please make sure that it is properly unfolded. During unfolding and folding, take all the necessary precautions to avoid pinching.**

## 2.3 POSITION ADJUSTMENT

In order to use this item safely, make sure it is completely unfolded.

To switch to the relaxed position, push the back of the chair with your back.

To return to the sitting position, carry the weight forward by pressing with your feet on the frame.

To lock the position, turn the locking system to the horizontal position.

## 2.4 RETRACTABLE HEADREST

## 3 CLEANING AND DISINFECTION

---

### 3.1 INTRODUCTION

Regular and thorough cleaning will guarantee the product's life (see cleaning program drawn up at the time of product commissioning, designating the reference staff).

**Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.**

The following cleaning advice is not intended to replace your establishment's in-house protocol.

This chair has been designed to facilitate cleaning.

### 3.2 CLEANING AND DISINFECTION ADVICE

#### 3.2.1 CLEANING OPERATION

For routine cleaning of upholstery, painted stainless steel or aluminum surfaces, or plastic parts, a mild detergent such as soapy water can be used followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product with a clean disposable cloth.
- Insist on the parts in contact with the users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that cause stains (e.g.: "blue product" or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

#### 3.2.2 RECOMMENDED EQUIPMENT AND PRODUCTS

- Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.
- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
  - PVC, ABS, Polyamide, Polyurethane and Polypropylene
  - Epoxy paint-coated metal surfaces
  - Stainless steel and aluminum surfaces

- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

**List of strictly prohibited products:**

<b>Solvents:</b>	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
<b>Bases:</b>	10% Sodium hydroxide / 10% Ammonia
<b>Acids:</b>	10% Nitric / 30% Hydrochloric / 10% Citric
<b>Miscellaneous:</b>	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

**High-pressure and steam cleaning are prohibited.**

**Any damage resulting from failure to observe these Instructions for Use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.**

## 4 PRECAUTIONS

---

### 4.1 DEVICE-RELATED PRECAUTIONS

#### 4.1.1 WARNINGS

**Abnormal use of the chair can cause accidents for the user or even irreparable damage to the product.**

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the users' common sense.

- Failure to respect the load for which it was designed, i.e. up to 160 kg, which excludes, among other things, use by two people simultaneously, even with temporary or partial support.
- Sitting on the chair when the backrest is not in the vertical position.
- Letting oneself "drop" into the chair.
- Changing position too quickly.
- Use by children.
- Product storage in a damp place.
- Never use the chair in the presence of an anesthesia mixture flammable with air, with oxygen, or with nitrous oxide as this could cause a fire and/or explosion.

**All manipulations relating to the aforementioned precautions for use are described in [paragraph II - Instructions for Use](#).**

#### 4.1.2 CONTRAINDICATIONS

There are no other contraindications for this product.

#### 4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE

The product must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The product must be regularly inspected, at least once per year.

To maintain the product in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed:

- Visually check all fastenings.
- Check the cylinder fastenings.
- Check that the backrest, legrest and all fastenings are working correctly.
- Lubricate all the moving mechanisms.

To keep the system operating optimally, repeat this check regularly.

If you have any doubts about the functioning and integrity of the product, take it out of service and inform your technical service or BMB MEDICAL after-sales service.

In the event of observed anomalies or damage that could jeopardize product operation and safety, the equipment must be immediately withdrawn.

The product should only be used according to the instructions in this user manual.

The product must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

\* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

#### **4.1.4 STORAGE**

The product should be stored in a clean, ventilated, temperate and dry place, protected from the elements and from direct sunlight.

During storage, no heavy weights should be placed on the chair.

The chair must not be considered as a support surface for any other type of equipment.

The product must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

## **4.2 WARRANTY**

Warranty: two years parts and labor according to the following terms:

*BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.*

*This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.*

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL products are designed for an expected service life of 10 years under normal conditions of use and with appropriate periodic maintenance as described above. BMB MEDICAL guarantees the original buyer that the welds on its products will be free from structural defects during the expected 10-year service life of the product as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use ([Warning paragraph](#))
- Work without prior notice to our technical services, or maintenance on:
  - o A mechanical part of the product
  - o A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to comply with the recommended preventive maintenance operations.

### **4.3 AFTER-SALES SERVICE**

For any information concerning product installation, use and maintenance, the technical documentation (wiring diagrams, list of components, maintenance procedures, etc.), or any other questions, please contact the After-Sales Service at the following address: **sav@bmb-medical.com**, or by phone at **+33 (0)4 74 08 71 71**, or your local distributor for overseas purchases.

From your “Customer Area” on the website [www.bmb-medical.com](http://www.bmb-medical.com), you can access your equipment, submit work requests and track these requests directly via this secure interface.

The operator must carefully read the information contained in this user manual, in particular the information concerning the precautions relating to safety and installation/use methods.

In the event of doubt concerning the correct interpretation of the instructions, do not hesitate to contact the After-Sales Service whose staff are always available.

To facilitate ordering spare parts, you can request the exploded diagrams of your products from the After-Sales Service. You can thus obtain the references of the desired parts.

#### 4.4 ENVIRONMENT

To be disposed of separately from household waste.

Dispose of the product or its components in accordance with local and national laws and regulations on environmental protection and recycling of raw materials

You can obtain the appropriate disposal location from the authorities in charge of environmental protection.

**France specificities:** BMB MEDICAL fulfills its obligations concerning the end of life of the products in the CLEM range, which it puts on the market by financing the recycling sector of Valdelia which takes them back free of charge (More information on [www.valdelia.org](http://www.valdelia.org)).

Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

[info@bmb-medical.com](mailto:info@bmb-medical.com)

100 Av. du Formans – CS 11001

01604 TREVOUX Cedex

FRANCE

[www.bmb-medical.com](http://www.bmb-medical.com)