# I-MOVE

MANUEL D'UTILISATION  
USER MANUAL

BMB  
MEDICAL

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# FRANÇAIS

---

## **BMB MEDICAL**

100 Avenue du Formans  
CS 11001  
01604 TRÉVOUX CEDEX  
Tél : +33(0)4 74 08 71 71  
Fax : +33(0)4 74 08 71 72

[www.bmb-medical.com](http://www.bmb-medical.com)

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[www.bmb-medical.com](http://www.bmb-medical.com)

The purpose of this manual is to familiarize you with the use, maintenance, and safety procedures for your new medical equipment. We strongly advise you to read this document carefully and apply the recommendations it gives, for optimum and safe use.

**BMB MEDICAL** offers a wide range of options for the models covered by this handbook. Some of the features described in this manual might therefore not apply to your care chair.

Some illustrations are provided for information purposes and do not accurately represent the device in your possession.

The information in this document may be modified without notice and does not imply any responsibility of any sort of BMB MEDICAL. Any representation or copying, whether complete or partial, without authorization from BMB MEDICAL is illegal. This illegal representation or copying, by any method at all, would constitute an infringement of rights. The illustrations in this document are not contractual.

# CONTENTS

ENGLISH .....	38
1 INTRODUCTION .....	41
1.1 PRODUCT PRESENTATION .....	41
1.1.1 I-MOVE TRANSFER CHAIR .....	42
1.1.2 I-MOVE TRANSFER CHAIR FOR OVERWEIGHT PATIENTS .....	43
1.1.3 MOTORIZED I-MOVE EZ-GO TRANSFER CHAIR .....	44
1.1.4 1.1.4 MOTORIZED I-MOVE EZ-GO TRANSFER CHAIR FOR OVERWEIGHT PATIENTS .....	45
1.1.5 OUTDOOR I-MOVE .....	46
1.2 APPLICATION FIELD .....	47
1.3 PURPOSE OF THE USER MANUAL .....	47
1.4 USERS .....	47
1.5 STORAGE .....	47
1.6 APPLIED PARTS .....	47
1.7 SYMBOLS AND STANDARDS USED .....	48
2 INSTRUCTION FOR USE .....	50
2.1 OPERATING CONDITIONS .....	50
2.2 INSTALLATION PROCEDURE .....	50
2.3 MOTORIZED WHEEL .....	50
2.3.1 START-UP .....	50
2.3.2 BATTERY .....	50
2.4 PRODUCT POWER SUPPLY .....	52
2.4.1 MOTORISED DRIVE .....	53
2.4.2 STANDARD I-MOVE .....	54
2.5 ARMREST .....	54
2.6 FOOTREST .....	55
2.7 OPTIONS .....	55
2.7.1 LEGREST .....	55
2.7.2 ANTISTATIC SYSTEM .....	56
2.7.3 REMOVABLE IV POLE / ANTI-THEFT BAR .....	57
2.7.4 HOLDER FOR OXYGEN CYLINDER .....	57
2.7.5 ACCESSORY BASKET .....	58
2.7.6 COIN ACCEPTOR AND WALL FIXATION .....	59
2.7.7 HOOKS TO HOLD BAG .....	59

3	CLEANING AND DISINFECTION .....	60
3.1	INTRODUCTION .....	60
3.2	FREQUENCY .....	60
3.3	CLEANING AND DISINFECTION ADVICE.....	60
3.3.1	CLEANING OPERATION .....	60
3.4	RECOMMENDED EQUIPMENT AND PRODUCTS .....	61
3.4.1	CLEANING OF THE OUTDOOR CHAIR.....	62
4	PRECAUTIONS.....	63
4.1	PRECAUTIONS RELATIVES A CE DISPOSITIF .....	63
4.1.1	AVERTISSEMENTS.....	63
4.1.2	<b>CONTRAINdications.....</b>	64
4.1.3	PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE .....	65
4.1.4	MC-RELATED INFORMATION.....	65
4.2	TRANSPORT AND STORAGE .....	67
4.2.1	TRANSPORT .....	67
4.2.2	STORAGE .....	67
4.3	WARRANTY .....	67
4.4	AFTER SALES SERVICE .....	68
4.5	ENVIRONMENT.....	69
4.5.1	SPENT EQUIPMENT.....	69

# 1 INTRODUCTION

---

## 1.1 PRODUCT PRESENTATION

Patient support device. The frame of the chair is mounted on casters to allow patient transfer, and it can be equipped with an electric wheel to assist movement. The MD is equipped with a mechanical braking system.

The device allows an operator to convey the patient within a healthcare establishment.

The products are designed to ensure that there are no unacceptable risks (essential performance) associated with patient support or caused by accidental movements.

List of I-Move references:

- I-MOVE – Stackable transfer chair I-MOVE
- I-MOVE BARIATRIC – Bariatric transfer chair I-MOVE
- I-MOVE EZ-GO – motorized transfer chair I-MOVE
- I-MOVE EZ-GO BARIATRIC – Motorized bariatric transfer chair I-MOVE
- I-MOVE OUTDOOR – Outdoor transfer chair I-MOVE

### 1.1.1 I-MOVE TRANSFER CHAIR

TECHNICAL CHARACTERISTICS	
<b>REFERENCE</b>	I-MOVE
<b>DESCRIPTION</b>	PATIENT TRANSFER CHAIR
<b>SAFE WORKING LOAD</b>	200 kg (maximum acceptable weight comprising the weight of the patient, the mattress and the accessories)
<b>WEIGHT</b>	35 kg
<b>OVERALL DIMENSION</b>	105 x 75 x 100 cm
<b>CUSHION AND BACK SIZE</b>	46 x 46 cm
<b>BACK HEIGHT</b>	55 cm
<b>EQUIPMENT</b>	<ul style="list-style-type: none"> <li>* Overall dimensions : 105 (length) x 75 (width) x 100 (height) cm</li> <li>* Dimensions seat and back: 46 x 46 cm</li> <li>* Seat height : 55 cm</li> <li>* Steel frame epoxy coated, white RAL 9003</li> <li>* Patient file holder</li> <li>* Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person</li> <li>* Front swivel wheels Ø 125 mm</li> <li>* Sliding and self-blocking foot rest</li> <li>* Polyurethan armrests, can be folded up</li> <li>* Anti-shock protection</li> <li>* Cushion in foam - Thickness 5 cm</li> <li>* Seamless finish, vinyl coating</li> <li>* Colour range : fuchsia, green, blue</li> </ul>
<b>OPTIONS</b>	<ul style="list-style-type: none"> <li>Removable IV pole, 2-stage, 2 hooks.</li> <li>Fixed IV pole, 2-stage, 2 hooks.</li> <li>Removable IV pole, 2-stage, 4 hooks.</li> <li>Anti-static system.</li> <li>U-shape anti-theft, fitted with 2 hooks infusion holder.</li> <li>Accessory basket.</li> <li>Holder for oxygen cylinder.</li> <li>Equipment for leg rest (pair).</li> <li>Adjustable leg rest (1 piece).</li> <li>Coin acceptor for 1€ coin or equivalent.</li> <li>Coin acceptor for specific token.</li> <li>Additional specific token.</li> <li>Wall fixation.</li> <li>Hooks to hold bags.</li> </ul>

## 1.1.2 I-MOVE TRANSFER CHAIR FOR OVERWEIGHT PATIENTS

TECHNICAL CHARACTERISTICS	
<b>REFERENCE</b>	I-MOVE -Transfer chair for overweight patient
<b>DESCRIPTION</b>	OVERWEIGHT PATIENT TRANSFER CHAIR
<b>SAFE WORKING LOAD</b>	250 kg (shows the maximum acceptable weight comprising the weight of the patient, the mattress and the accessories)
<b>WEIGHT</b>	50 kg
<b>OVERALL DIMENSION</b>	105 x 90 x 102 cm
<b>CUSHION AND BACK SIZE</b>	46 x 65 cm
<b>BACK HEIGHT</b>	55 cm
<b>EQUIPMENT</b>	<ul style="list-style-type: none"> <li>* Overall dimensions : 105 (length) x 90 (width) x 102 (height) cm</li> <li>* Dimensions seat and back: 46 x 65 cm</li> <li>* Seat height : 55 cm</li> <li>* Steel frame epoxy coated, white RAL 9003</li> <li>* Patient file holder</li> <li>* Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person</li> <li>* Front swivel wheels Ø 125 mm</li> <li>* Sliding and self-blocking foot rest</li> <li>* Polyurethan armrests, can be folded up</li> <li>* Anti-shock protection</li> <li>* Cushion in foam - Thickness 8 cm</li> <li>* PVC coated Jersey</li> <li>* Colour range : lavender, green, squirrel</li> </ul>
<b>OPTIONS</b>	<ul style="list-style-type: none"> <li>Removable IV pole, 2-stage, 2 hooks.</li> <li>Fixed IV pole, 2-stage, 2 hooks.</li> <li>Removable IV pole, 2-stage, 4 hooks.</li> <li>Anti-static system.</li> <li>U-shape anti-theft, fitted with 2 hooks infusion holder.</li> <li>Accessory basket.</li> <li>Holder for oxygen cylinder.</li> <li>Equipment for leg rest (pair).</li> <li>Adjustable leg rest (1 piece).</li> <li>Coin acceptor for 1€ coin or equivalent.</li> <li>Coin acceptor for specific token.</li> <li>Additional specific token.</li> <li>Wall fixation.</li> <li>Hooks to hold bags.</li> </ul>

### 1.1.3 MOTORIZED I-MOVE EZ-GO TRANSFER CHAIR

TECHNICAL CHARACTERISTICS	
<b>REFERENCE</b>	I-MOVE EZ-GO - Motorized
<b>DESCRIPTION</b>	MOTORIZED PATIENT TRANSFER CHAIR
<b>SAFE WORKING LOAD</b>	200 kg (shows the maximum acceptable weight comprising the weight of the patient, the mattress and the accessories)
<b>WEIGHT</b>	52 kg
<b>OVERALL DIMENSION</b>	105 x 70 x 102 cm
<b>CUSHION AND BACK SIZE</b>	46 x 46 cm
<b>BACK HEIGHT</b>	55 cm
<b>EQUIPMENT</b>	<ul style="list-style-type: none"> <li>* Overall dimensions: 105 (length) x 70 (width) x 102 (height) cm</li> <li>* Dimensions seat and back: 46 x 46 cm</li> <li>* Seat height: 55 cm</li> <li>* Motorized drive system with autonomous wheel</li> <li>* Switchable battery - 5 km autonomy per battery - Charge time : 5 hours</li> <li>* Lithium battery without memory effect</li> <li>* Forward/Backward twist grip with progressive speed acceleration</li> <li>* Steel frame epoxy coated</li> <li>* Patient file holder</li> <li>* Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person</li> <li>* Front swivel wheels Ø 125 mm</li> <li>* Sliding and self-blocking footrest</li> <li>* Polyurethan armrests, can be folded up</li> <li>* Accessory basket</li> <li>* Anti-shock protection</li> <li>* Supplied with 2 batteries and external charger</li> <li>* Cushion in foam - Thickness 5 cm</li> <li>* Seamless finish, vinyl coating</li> <li>* Color range: fuchsia, green, blue</li> </ul>
<b>OPTIONS</b>	<ul style="list-style-type: none"> <li>Removable IV pole, 2-stage, 2 hooks.</li> <li>Removable IV pole, 2-stage, 4 hooks.</li> <li>U-shape anti-theft, fitted with 2 hooks infusion holder.</li> <li>Holder for oxygen cylinder.</li> <li>Hooks to hold bags.</li> <li>Equipment for leg rest (pair).</li> <li>Adjustable leg rest (1 piece).</li> <li>Coin acceptor for 1€ coin or equivalent.</li> <li>Coin acceptor for specific token.</li> <li>Additional specific token.</li> <li>Wall fixation.</li> <li>Without battery (depreciation).</li> <li>Additional battery.</li> <li>External wall charger</li> <li>External charger mobile version</li> <li>Charging station 1 charger</li> <li>Charging station 2 chargers (one multi-socket included)</li> <li>Charging station 3 chargers (one multi-socket included)</li> <li>Charging station 4 chargers (one multi-socket included)</li> </ul>

## 1.1.4 1.1.4 MOTORIZED I-MOVE EZ-GO TRANSFER CHAIR FOR OVERWEIGHT PATIENTS

CARACTERISTIQUES TECHNIQUES	
<b>REFERENCE</b>	I-MOVE EZ-GO BARIATRIC - Transfer chair for overweight patient
<b>DESCRIPTION</b>	OVERWEIGHT PATIENT TRANSFER CHAIR
<b>SAFE WORKING LOAD</b>	250 kg (shows the maximum acceptable weight comprising the weight of the patient, the mattress and the accessories)
<b>WEIGHT</b>	65 kg
<b>OVERALL DIMENSION</b>	105 x 90 x 102 cm
<b>CUSHION AND BACK SIZE</b>	46 x 65 cm
<b>BACK HEIGHT</b>	55 cm
<b>EQUIPMENT</b>	<ul style="list-style-type: none"> <li>* Overall dimensions: 105 (length) x 90 (width) x 102 (height) cm</li> <li>* Dimensions seat and back: 46 x 65 cm</li> <li>* Seat height : 55 cm</li> <li>* Motorized drive system with autonomous wheel</li> <li>* Switchable battery - 5 km autonomy per battery - Charge time : 5 hours</li> <li>* Lithium battery without memory effect.</li> <li>* Forward/Backward twist grip with progressive speed acceleration</li> <li>* Steel frame epoxy coated RAL 9010</li> <li>* Patient file holder</li> <li>* Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person</li> <li>* Front swivel wheels Ø 125 mm</li> <li>* Sliding and self-blocking foot rest</li> <li>* Polyurethan armrests, can be folded up</li> <li>* Accessory basket</li> <li>* Supplied with 2 batteries and external charger</li> <li>* Anti-shock protection</li> <li>* Cushion in foam HR 35kg/m3 - Thickness 8 cm</li> <li>* Seamless finish, vinyl coating</li> <li>* Colour range : lavender, green, squirrel</li> </ul>
<b>OPTIONS</b>	<p>Removable IV pole, 2-stage, 2 hooks.      Removable IV pole, 2-stage, 4 hooks.      U-shape anti-theft, fitted with 2 hooks infusion holder.      Holder for oxygen cylinder.      Hooks to hold bags.      Equipment for leg rest (pair).      Adjustable leg rest (1 piece).      Coin acceptor for 1€ coin or equivalent.      Coin acceptor for specific token.      Additional specific token.      Wall fixation.      Without battery (depreciation).      Additional battery.      External wall charger      External charger mobile version      Charging station 1 charger      Charging station 2 chargers (one multi-socket included)      Charging station 3 chargers (one multi-socket included)      Charging station 4 chargers (one multi-socket included)</p>

### 1.1.5 OUTDOOR I-MOVE

CARACTERISTIQUES TECHNIQUES	
<b>REFERENCE</b>	I-MOVE OUTDOOR – Outdoor transfer chair
<b>DESCRIPTION</b>	TRANSFER CHAIR
<b>SAFE WORKING LOAD</b>	200 kg (shows the maximum acceptable weight comprising the weight of the patient, the mattress and the accessories)
<b>WEIGHT</b>	35 kg
<b>OVERALL DIMENSION</b>	105 x 75 x 100 cm
<b>CUSHION AND BACK SIZE</b>	46 x 46 cm
<b>BACK HEIGHT</b>	55 cm
<b>EQUIPMENT</b>	<ul style="list-style-type: none"> <li>* Overall dimensions : 105 (length) x 75 (width) x 100 (height) cm</li> <li>* Dimensions seat and back : 46 x 46 cm</li> <li>* Seat height : 55 cm</li> <li>* Steel frame epoxy coated, white RAL 9003</li> <li>* Patient file holder</li> <li>* Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person</li> <li>* Front swivel wheels Ø 125 mm</li> <li>* Sliding and self-blocking foot rest</li> <li>* Polyurethan armrests, can be folded up</li> <li>* Anti-shock protection</li> <li>* Polypropylene plastic fitting</li> <li>* Color : blue</li> </ul>
<b>OPTIONS</b>	<ul style="list-style-type: none"> <li>Removable IV pole, 2-stage, 2 hooks.</li> <li>Fixed IV pole, 2-stage, 2 hooks.</li> <li>Removable IV pole, 2-stage, 4 hooks.</li> <li>Anti-static system.</li> <li>U-shape anti-theft, fitted with 2 hooks infusion holder.</li> <li>Accessory basket.</li> <li>Holder for oxygen cylinder.</li> <li>Coin acceptor for 1€ coin or equivalent.</li> <li>Coin acceptor for specific token.</li> <li>Additional specific token.</li> <li>Wall fixation.</li> <li>Hooks to hold bags.</li> </ul>

## **1.2 APPLICATION FIELD**

This chair is intended for the reception and transfer of the patient in sitting position. Its technical characteristics also allow an outdoor use.

## **1.3 PURPOSE OF THE USER MANUAL**

Healthcare professionals, this user manual was created for you to provide you with the necessary information and precautions to observe for autonomous and safe use.

The product's service life is dependent upon observance of the maintenance and tracking instructions.

It is thus essential that all staff has access to and be informed of the handling procedures and risks associated with the use of this chair. This latter should be handled by a qualified person who is aware of the risks inherent to improper use. Additional training is not required as we consider that knowledge of the user manual is sufficient.

## **1.4 USERS**

Healthcare professionals: physicians, caregivers and agents approved by Technical Services who are accustomed to working with similar equipment can operate all functions.

## **1.5 STORAGE**

This user manual should be kept near the product, in a place provided for this purpose. It must be shielded from all elements and substance that could compromise its perfect legibility.

## **1.6 APPLIED PARTS**

The applied parts are as follows:

- Upholstery
- Armrests
- Footrest

## 1.7 SYMBOLS AND STANDARDS USED

I-MOVE and I-MOVE transfer chairs for overweight patients are Class I, non-invasive and non-active medical devices.

The I-MOVE EZ-GO transfer chair is a Class I, non-invasive and active medical device.

	Type B appliance
	Electrical class II equipment
	Electrical hazard, do not open
	Duty cycle (motor cycle) 1 min. operation / 19 min. stop
<b>IPX5</b>	Splash-resistant
	See the manual / instruction brochure
	Indicates the maximum patient weight and safe working load
	Complies with the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council
<b>SN</b>	Serial number
	Manufacture date
	Manufacturer
	Alternating current
	Maximum patient weight
	Electronic waste: this symbol indicates that the product must not be discarded with unsorted municipal waste, but rather to a selective waste collection channel. Contact the local distributor for information concerning the return and/or collection systems available in the country.
	Protéger de l'humidité / keep dry
	Haut, tenir debout / stand up
	Fragile / fragile
	Recyclable / recycling
	Temperature range

	Humidity limit
---	----------------

**General safety sign:**

- Information highlighted by this symbol must be read and scrupulously observed.
- Indicates that failure to observe the guideline could endanger the patient or user, or damage the equipment.

Sitting on the armrests is prohibited

Leaning against the chair is prohibited

**“WARNING: This appliance must not be altered without the manufacturer’s written consent.”**

**“Caution: only purpose-built accessories provided for this MD by BMB MEDICAL may be used.”**

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### 2.2 INSTALLATION PROCEDURE

Before putting the chair into service, depending on the model, check the operation of:

- [Motorized drive](#)
- [Braking](#)
- [Armrests](#)
- [Footrest](#)

Clean the product according to your establishment's internal protocol.

### 2.3 MOTORIZED WHEEL

The motorized traction by the wheel facilitates the movement of the chair.

#### 2.3.1 START-UP

Storing and/or carrying motorized chairs in very low temperature conditions (below 5 °C) may alter the operation of the product. This is why it is recommended to leave the products for 24 hours in a temperate environment before putting them into service following the conditions described above.

#### 2.3.2 BATTERY

The battery box  comes with its own external charger. 

Before first use, leave the product in a warm environment for 24 hours and then charge the battery for 6 to 8 hours on a mains socket.

In order to maintain an acceptable level of charge, it is recommended to charge the motorized chairs as soon as possible. It is not necessary to wait until the battery is completely discharged.


**« Warning: Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact. »**

**In case of prolonged interruption of use of your equipment, it is recommended to recharge the battery once a month. Beyond 50 days of storage without recharging, the performance of the battery may be impaired or even irreversibly damaged.**

Battery specifications		Battery safety symbols	
Battery type:	Lithium-Ion	IPX4 and IPX5	Protection rating
Capacity:	4.50 Ah	CE	CE marking
Input supply voltage:	29~45V CC	House	Indoor use

Label with information on the operation of the batteries.

**Rechargeable battery  
(Option : EZSBA Without battery (depreciation))**

## 2.4 PRODUCT POWER SUPPLY

**Power grid connection conditions and electromagnetic compatibility:**

**This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.**

**Position the product such that it can be easily unplugged.**

**Unplug the product before moving it.**

**Do not leave the cables in contact with the floor or the casters,**

**Do not knot the cables,**

**Periodically check the condition of the cables.**

Mains connection:	100~240V depending on the country.
Frequency:	50 /60 Hz
Rated current:	3.3 A
Protection rating / degree of electric shock protection:	Class 2 appliance
Power input:	525 W
Intermittent operating mode:	5% i.e., 1 minute of continuous operation followed by 19 minutes of rest.

The mains plug is the electrical grid insulation device.


Work on electrical parts should only be performed by qualified and authorized staff.

Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the [EMC information here](#).

#### 2.4.1 MOTORISED DRIVE

In forward motion there are two speeds:  
"turtle" being slow and "rabbit" being fast.

In reverse there is only one speed.

##### DEMARRAGE

- Press the start button (F)
  - o The battery charge level indicator lights up (E)
- Fully press the brake control (A) and then release it.  
When a green or red arrow is displayed (D), the chair is started and ready for use.

##### SHIFTING

###### Forward:

Select the drive direction (B).  
Select speed (C).

Turn the control handle forward (G). The speed is controlled by this handle.

### **Reverse:**

Select the drive direction (B).

Turn the control handle forward (G). The speed is controlled by this handle. The selected speed (C) has no impact in reverse.

For safety reasons, the command reaction time is 0.5 seconds

**Check regularly that nothing on the ground is hindering the rotation of the wheel.**

### **BRAKING**

- Fully press the brake control (A) to activate the brake.

For safety reasons, the command reaction time is 0.5 seconds.

#### **2.4.2 STANDARD I-MOVE**

- Push with your foot the green pedal to unblock the chair.
- Push with your foot the red pedal to immobilize the chair.

### **2.5 ARMREST**

To raise or lower the armrests, simply raise them or lower them manually.


## 2.6 FOOTREST

The opening at the front of the tray allows a foot to be inserted to advance or push back the footrest.

To deploy the footrest, extend as far as the stop to lock it. To retract the footrest, lift it and push it.

- When moving a patient make sure that the footrest is deployed.
- Do not step on the foot rest to sit in the transfer chair.
- Do not stand up on the footrest.


## 2.7 OPTIONS

### 2.7.1 LEGREST

To install the leg rest, insert the bar in the jaws provided for the purpose under the chair.

It is possible to incline the leg rest by placing its bar in the 2nd notch.


Inclined position

Straight position

For storage, place the leg rest upside down as shown below

- Do not sit on the leg rest
- Do not apply a load greater than 20kg on the leg rest.

## 2.7.2 ANTISTATIC SYSTEM

The I-MOVE EZ-GO is equipped with antistatic wheel (optional for I-MOVE standard version)

In order to maintain the antistatic functions of the wheels, you must clean it regularly (no dust )

### 2.7.3 REMOVABLE IV POLE / ANTI-THEFT BAR

The back frame of the transfer chair has two places to insert the removable IV Pole and the anti-theft bar.

The removable IV Pole is placed in one of the 2 holders placed on both sides of the back of the chair.

To deploy the IV pole:

- 1- Turn the thumb wheel ① counter-clockwise and extend the telescopic part of the IV pole ② to the desired height.
- 2- To lock the telescopic part, turn the thumb wheel ① clockwise.  
Ensure that the IV pole is securely locked.

**Safe working load: 6 kg.**

This indicative value corresponds to safe product use, but is not representative of normal use for which the load is significantly lower.

When moving the chair fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.

**Do not move the chair via the IV pole.**


### 2.7.4 HOLDER FOR OXYGEN CYLINDER

Put the oxygen cylinder in the holder:

The holder is compatible with oxygen cylinder in diameter 10cm and 14cm.

In order to adjust the holder for oxygen cylinder in diameter 10cm or 14cm, you have to reposition the retaining plate to the good place unscrewing and retightening it in the hole. Check if the screws are tight. Ask your technical department or our after-sales department doing it.

### 2.7.5 ACCESSORY BASKET

Removable for easy cleaning.

## 2.7.6 COIN ACCEPTOR AND WALL FIXATION

Token or coin. Put the coin in the space provided for this purpose to unlock the chair.

Fix the wall fixation on the wall and connect the chain with the coin acceptor.

## 2.7.7 HOOKS TO HOLD BAG

## 3 CLEANING AND DISINFECTION

---

### 3.1 INTRODUCTION

Regular and proper cleaning ensures the service life of the appliance (cf. cleaning program drawn up at the time of entry into service of the appliance with the identities of the personnel responsible).

**Failure to perform cleaning operations may cause a risk of infection due to the presence of secretions and/or residues on the transfer chair.**

The following cleaning advice is not intended to replace your establishment's in-house protocol.

### 3.2 FREQUENCY

The chair should be cleaned:

- Between each patient.
- After transporting an infected patient.
- Completely each week

### 3.3 CLEANING AND DISINFECTION ADVICE

#### 3.3.1 CLEANING OPERATION

**Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.**

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that cause stains (e.g.: "blue product" or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

### 3.4 RECOMMENDED EQUIPMENT AND PRODUCTS

Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.

- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
  - o PVC, ABS, Polyamide, Polyurethane and Polypropylene
  - o Epoxy paint-coated metal surfaces
  - o Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

#### List of strictly prohibited products:

Solvents:	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
Bases:	10% Sodium hydroxide / 10% Ammonia
Acids:	10% Nitric / 30% Hydrochloric / 10% Citric
Miscellaneous:	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

N.B :

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.


**High-pressure and steam cleaning are prohibited.**

**Any damage resulting from failure to observe this instruction for use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.**

### 3.4.1 CLEANING OF THE OUTDOOR CHAIR

*Cleaning and disinfection recommendations for the I-MOVE  
OUTDOOR seat and backrest*

## 4 PRECAUTIONS

---

### 4.1 PRECAUTIONS RELATIVES A CE DISPOSITIF

#### 4.1.1 AVERTISSEMENTS

**This MD has been designed for and adapted to medical practices and procedures, and hence a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.**

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the common sens of the chair's users.

- Do not use the product without first having familiarized yourself with this user manual.
- Do not alter the chair. Any changes made to the system could lead to unpredictable operation that could in turn jeopardize the safety of the user or patient. This would also void the chair's warranty.
- The chassis cover load capacity is of 35 kg. Do not sit or stand on the cover.
- Only purpose-built accessories provided for this MD by BMB MEDICAL may be used
- Connection, by any means, of an electric device to the chair
- The product must be connected to a mains power outlet compliant with standards in force.
- Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact.
- This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.
- Position the product such that it can be easily unplugged.
- Unplug the product before moving it.
- Do not leave the cables in contact with the floor or the casters.
- Do not knot the cables.
- Periodically check the condition of the cables.
- Work on electrical parts should only be performed by qualified and authorized staff.
- Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view the EMC information [here](#).
- Except during transport, always ensure that the chair's brake is engaged. Check whether the brake is engaged by attempting to move the chair.
- Do not park the chair on a sloping floor.

- Use outside the facility on an undeveloped floor or Do not use the MD on a floor with a slope greater than 10% or When moving down a slope (maximum 10°), the patient must be positioned in the opposite direction of travel.
- When moving the patient, the footrest must be in the extended position and both armrests must be lowered to prevent a fall.
- When transferring the patient (moving the patient from bed to chair or from chair to bed), the chair should be positioned as close to the bed as possible to ensure maximum patient safety. Adjust the height of the bed to match the height of the transfer chair.
- Do not sit on the leg rest
- Do not apply a load greater than 20kg on the leg rest.
- Use of the chair with a load that exceeds the safe working load (see technical characteristics).
- IV pole safe working load: 6 kg.
- When moving the chair fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.
- Do not move the chair via the IV pole
- The footrest must not, under any circumstances, be used by the patient to climb onto the product, it is a patient comfort element.
- Tell the patient and any accompanying persons not to sit on the armrests.
- Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.
- Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.
- High-pressure and steam cleaning are prohibited.
- Any use that does not fall within the operating context of a hospital chair.
- To prevent failure, electromedical equipment requires special [EMC-related precautions, described here](#).
- Never use the chair in the presence of an anesthesia mixture flammable with air, with oxygen, or with nitrous oxide as this could cause a fire and/or explosion.
- Do not manually raise the legrest.
- Use by any non-caregiver who has not received the authorization and explanations from a qualified individual.
- In the case of an unsupervised patient, the product should be left in its lowest position to reduce the risk of falling.

**All manipulations relating to the aforementioned precautions for use are described in [paragraph II - Instructions for use](#).**

#### **4.1.2 CONTRAINDICATIONS**

There are no noteworthy contraindications for this device.

#### **4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE**

The device must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The device must be regularly inspected, at least once per year. See the maintenance check-list. To maintain the device in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed.

In the event of observed anomalies or damage that could jeopardize device operation and safety, the equipment must be immediately withdrawn.

The device should only be used according to the instructions in this user manual.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion, and flammable agents.

Work on electrical parts should only be performed by qualified and authorized staff.

\* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

#### **4.1.4 MC-RELATED INFORMATION**

The products in the I-MOVE range are suitable for use in all healthcare establishments, with the exception of facilities in the vicinity of HF ELECTROSURGICAL APPLIANCES, or near an electromagnetically shielded room for an EM magnetic resonance imaging system, in which the intensity of EM INTERFERENCES is high.

The products in the I-MOVE range are not suitable for use in domestic homes or establishments or on any premises connected to the public low-voltage network supplying domestic buildings.

**EM interference could cause the product to make unexpected movements.**

**This appliance should not be used next to, or stacked on other appliances as this could lead to incorrect operation. If this mode of operation is necessary, this and the other appliances should be observed to ensure their correct operation.**

**The use of accessories, transducers or cables other than those specified or supplied by the manufacturer can lead to an increase in electromagnetic emissions, or a decrease in this appliance's immunity, causing inappropriate operation.**

**Do not use portable RF communication devices (including peripherals such as antenna cables and external antennas) less than 30 cm (12 inches) from any part of products in the I-MOVE range, including the cables specified by the manufacturer. Otherwise, the performance of these devices could be impaired.**

**The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11, class A). If used in a residential environment (for which compliance with the CISPR 11 class B standard is generally required), this equipment may not provide appropriate protection to radio frequency communication services. The user may need to take attenuation steps, such as moving or re-orienting the equipment.**

List of cables used:

1. Power cord, maximum length 1000 cm.
2. Variable height column, maximum length 750 mm
3. Backrest cable: maximum length 1250 mm
4. Legrest cable: maximum length 750 mm
5. Footrest cable: maximum length 1500 mm

EMISSION	
Conducted and radiated RF emissions:	CISPR11 Group 1 Class A
Harmonic distortions IEC 61000-3-2:	Class A
Voltage fluctuations and Flicker IEC 61000-3-3:	Compliant

Immunity test	Test level IEC 60601		Compliance level	
Electrostatic discharges (ESD) IEC 61000-4-2	$\pm 8$ kV in contact $\pm 2/4/8/15$ kV in air		$\pm 8$ kV in contact $\pm 2/4/8/15$ kV in air	
RF radiated electromagnetic fields IEC 61000-4-3	3 V/m (professional healthcare establishment) 80MHz-2.7GHz 80% AM at 1 kHz		3V/m 80MHz-2.7GHz 80% AM at 1 kHz	
Proximity fields emitted by wireless RF communication devices IEC 61000-4-3				
Rapid / burst electrical spikes IEC 61000-4-4	Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)
	385	Pulsed modulation: 18 Hz	27	27
	450	Pulsed modulation: 18 Hz	28	28
	710-745-780	Pulsed modulation: 217 Hz	9	9
	810-870-930	Pulsed modulation: 18 Hz	28	28
	1720-1845-1970	Pulsed modulation: 217 Hz	28	28
	2450	Pulsed modulation: 217 Hz	28	28
Shock waves IEC 61000-4-5	5240-5500-5785		9	9
	Pulsed modulation: 217 Hz			
Rapid / burst electrical spikes IEC 61000-4-4	Power supply: $\pm 2$ kV Input/output lines: $\pm 1$ kV Repeat frequency: 100 kHz		Power supply: $\pm 2$ kV Input/output lines: $\pm 1$ kV Repeat frequency: 100 kHz	
Shock waves IEC 61000-4-5	Between phases: $\pm 0.5$ kV, $\pm 1$ kV Between earth and phases $\pm 0.5$ kV, $\pm 1$ kV, $\pm 2$ kV		Between phases: $\pm 0.5$ kV, $\pm 1$ kV Between earth and phases $\pm 0.5$ kV, $\pm 1$ kV, $\pm 2$ kV	
Conducted interference, induced by RF fields IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz		3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Power frequency magnetic fields IEC 61000-4-8	30A/m		30A/m	
Voltage Dips and Interruptions: IEC 61000-4-11	0% UT; 0.5 cycle At $0^\circ$ , $45^\circ$ , $90^\circ$ , $135^\circ$ , $180^\circ$ , $225^\circ$ , $270^\circ$ and $315^\circ$ 0% UT; 1 cycle at $0^\circ$ 70% UT; 25/30 cycles at $0^\circ$ 0% UT; 250/300 cycles		0% UT; 0.5 cycle At $0^\circ$ , $45^\circ$ , $90^\circ$ , $135^\circ$ , $180^\circ$ , $225^\circ$ , $270^\circ$ and $315^\circ$ 0% UT; 1 cycle at $0^\circ$ 70% UT; 25/30 cycles at $0^\circ$ 0% UT; 250/300 cycles	

NB: There is no deviation from the requirements of IEC 60601-1-2 ed. 4

NB: There are no other known measures for maintaining basic safety with respect to EMC phenomena.

## 4.2 TRANSPORT AND STORAGE

Storage and transport conditions :

Ambient conditions	Operation
Ambient temperature	-10 to +30 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

### 4.2.1 TRANSPORT

During chair transport with a view to its storage, the following is recommended to prevent any impacts or falls:

- Strap it to a surface of the means of transport
- Cover it
- Engage the brakes
- Set it to its low position

Damage caused during transport is not covered by the warranty. The costs of repair and replacement of damaged parts will be borne by the customer.

### 4.2.2 STORAGE

When the chair is in storage, ensure that:

- Its brakes are engaged
- It is covered
- In a temperate and not humid room (temperature: from +10°C to +40°C).

During storage, no heavy weights should be placed on the chair.

The chair must not be considered as a support surface for any other type of equipment.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

In case of prolonged storage, it is recommended to recharge the battery every 50 days.

## 4.3 WARRANTY

Warranty: two years parts and labor according to the following terms:

BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.

*This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.*

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL chairs are designed for a projected service life of 10 years under normal conditions of use, with appropriate periodic maintenance as described in the maintenance check-list of each chair. BMB MEDICAL provides the original buyer with the guarantee that the welds on the chairs are free from structural defects for the projected chair service life of 10 years, as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use (Warning paragraph [here](#))
- Work without prior notice to our technical services, or maintenance on:
  - o A mechanical part of the chair
  - o A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to observe the preventive maintenance operations recommended in the maintenance check-list.

#### **4.4 AFTER SALES SERVICE**

For any information concerning product installation, use and maintenance, the technical documentation (wiring diagrams, list of components, maintenance procedures, etc.), or any other questions, please contact the After-Sales Service at the following address: [sav@bmb-medical.com](mailto:sav@bmb-medical.com), or by phone at +33 (0)4 74 08 71 71, or your local distributor for overseas purchases.

From your "Customer Area" on the website [www.bmb-medical.com](http://www.bmb-medical.com), you can access your equipment, submit work requests and track these requests directly via this secure interface.

The operator must carefully read the information contained in this user manual, in particular the information concerning the precautions relating to safety and installation/use methods.

In the event of doubt concerning the correct interpretation of the instructions, do not hesitate to contact the After-Sales Service whose staff are always available.

To facilitate ordering spare parts, you can request the exploded diagrams of your products from the After-Sales Service. You can thus obtain the references of the desired parts.

## 4.5 ENVIRONMENT

### 4.5.1 SPENT EQUIPMENT

#### In France:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution. For this purpose, BMB MEDICAL fulfills its obligations concerning the end of life of I-MOVE range products that it markets by funding the Recylum recycling channel dedicated to Pro WEEE that collects them free of charge (more information on [www.recylum.com](http://www.recylum.com)).

#### In Europe:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution.

You can obtain the appropriate disposal location from the authorities in charge of environmental protection.

#### Outside of Europe:

Dispose of the product or its components in accordance with local laws and regulations.

Contact the local distributor for information concerning the return and/or collection systems available in the country.

Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

[info@bmb-medical.com](mailto:info@bmb-medical.com)

100 Av. du Formans – CS 11001  
01604 TREVOUX Cedex  
FRANCE