KYELLE

MANUEL D'UTILISATION
USER MANUAL

 **BMB
MEDICAL**

 [FRANÇAIS](#)
[ENGLISH](#)

FRANÇAIS

BMB MEDICAL

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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www.bmb-medical.com

The purpose of this manual is to give you an understanding of the use, maintenance and safety of your new medical device. We strongly advise that you read this document carefully and follow its recommendations for optimal, safe use.

BMB MEDICAL provides a wide range of options for the products covered in this user manual. As a result, some of the features described herein may not be applicable to your recliner.

Some illustrations are provided solely for information purposes and do not accurately represent the device in your possession.

The information in this document is subject to change without notice; the BMB MEDICAL Company accepts no responsibility whatsoever for said information. Any representation or reproduction, in whole or in part, without the consent of BMB MEDICAL is unlawful. Said unlawful representation or reproduction, by any means whatsoever, shall constitute an infringement. The illustrations in this document are not binding.

TABLE DES MATIERES

ENGLISH	24
1 INTRODUCTION SPECIFICATIONS	27
1.1 PRODUCT PRESENTATION	27
1.1.1 KYELLE RECLINER.....	28
1.2 SCOPE	29
1.3 PURPOSE OF THE USER MANUAL.....	29
1.4 USERS	29
1.5 STORAGE.....	29
1.6 SYMBOLS AND STANDARDS USED	30
2 INSTRUCTIONS FOR USE	31
2.1 INSTALLATION PROCEDURE	31
2.2 ARMRESTS.....	31
2.3 BACKREST	32
2.4 LEG REST.....	33
2.5 OPTIONS	33
2.5.1 GAS-SPRING ASSISTED, INDEPENDENT LEG REST	33
2.5.2 FIXED PUSH BAR.....	34
2.5.3 REVERSIBLE PUSH BAR	34
2.5.4 HEADREST CUSHION.....	35
2.5.5 FOOTREST	35
2.6 BRAKING.....	35
2.6.1 WHEELS WITH INDEPENDENT BRAKE.....	35
2.6.2 WHEEL WITH DIRECTIONAL BRAKE	36
2.6.3 CENTRALIZED BRAKING.....	36
3 CLEANING AND DISINFECTION	37
3.1 INTRODUCTION	37
3.2 FREQUENCY	37
3.3 CLEANING AND DISINFECTION RECOMMENDATIONS.....	37
3.3.1 CLEANING OPERATION	37
3.4 RECOMMENDED EQUIPMENT AND PRODUCTS	38
4 PRECAUTIONS.....	39
4.1 PRECAUTIONS WITH REGARD TO THIS DEVICE.....	39
4.1.1 WARNINGS	39
4.1.2 CONTRAINDICATIONS	40

4.1.3	PRECAUTION WHEN MAINTAINING AND KEEPING THE EQUIPMENT IN GOOD CONDITION.....	40
4.2	TRANSPORT AND STORAGE	41
4.2.1	TRANSPORT	41
4.2.2	STORAGE	41
4.3	WARRANTY	42
4.4	CUSTOMER SERVICE.....	43
4.5	ENVIRONMENT	43

1 INTRODUCTION SPECIFICATIONS

1.1 PRODUCT PRESENTATION

Patient support device, consisting of a backrest and a seat mounted on a chassis on wheels (depending on the options). The device is equipped with armrests and, as an option, supports the temporary installation of IV poles. It is used to care for and transfer patients. The MD is equipped with a mechanical braking system actuated by pressing on a pedal or brake bar (depending on the options).

The device allows the operator to carry the patient inside a healthcare facility.

The products are designed to prevent the risk of incapacitation associated with patient support as well as the risk of involuntary movement.

1.1.1 KYELLE RECLINER

TECHNICAL FEATURES	
REFERENCE	KYELLE
DESCRIPTION	EMERGENCY AND REST RECLINER
SAFE WORKING LOAD	200kg (indicates the maximum permissible load, including the total weight of the patient, mattress and accessories)
WEIGHT	27 to 34 kg depending on the model
OVERALL DIMENSION	Width: 73 cm Length: 76 cm in seated position, 161 cm when lying down Height: 120 cm
SEAT WIDTH	52 cm
SEAT HEIGHT	48 cm
EQUIPMENT	<ul style="list-style-type: none"> * Fixed backrest inclined at 85° * Epoxy coated steel frame * Assembled on 4 pads (2 back pads with integrated casters) * Fixed polyurethane armrests * Foam upholstery thickness 8cm * Unclippable seat - Ergonomically shaped backrest - M1 fireproof PVC covering * Range of colors URBAN
OPTIONS AVAILABLE	<p>Version assembled on 4 casters of 125mm diameter: 2 free casters at the front, 2 rear casters with independent brakes "</p> <p>"Version assembled on 4 casters : 2 free casters at the front of 125mm diameter, 2 centrally braked rear casters of 75mm diameter"</p> <p>Front directional wheel (KYRF or KYRS mandatory) Ergonomically shaped polyurethane retractable armrests. Backrest adjustable from 85 to 10° assisted by gas spring. (Gas cylinders) Leg rest synchronized with the back rest adjustable from 0 to 100°. (KYRB mandatory) (Gas cylinders) Independent leg rest adjustable from 0 to 100° assisted by gas spring. (Gas cylinders) Retractable footrest (KYRF or KYRS mandatory) Fixed push bar (KYRF or KYRS mandatory) Reversible push bar (sitting and lying position) - (KYRF or KYRS mandatory) Comfort headrest cushion Removable IV pole, 2-stage, 2 hooks (KYSTR mandatory) Removable IV pole, 2-stage, 4 hooks (KYSTR mandatory) IV Pole holder (on both sides of the product) Support for meal tray on the right side of the patient (compatible with the KYSPR-G) Support for meal tray on the left side of the patient (compatible with the KYSPR-D) Small Shelf : 42 x 32 cm (KYSPR-D or KYSPR-G mandatory) Large Shelf : 81 x 41 cm (KYSPR-D or KYSPR-G mandatory)</p>

1.2 SCOPE

This recliner is intended to care for and transfer patients.

1.3 PURPOSE OF THE USER MANUAL

This user manual is intended to provide healthcare professionals with all the information they need and all the precautions to be observed for use with complete autonomy and safety.

Product maintenance and tracking instructions guarantee its long service life.

It is therefore imperative that all personnel have access to this manual and be informed of the handling procedures and risks associated with the use of the recliner. Said recliner must be handled by a qualified person, aware of the risks inherent in its misuse. Further training is not required; we consider it sufficient to read the user manual.

1.4 USERS

Healthcare professionals: physician, caregivers and authorized Technical Officers used to working with similar equipment can handle all features.

1.5 STORAGE

This user manual should be kept near the product in a place intended for this purpose. It must be kept away from any element and any substance likely to affect its perfect legibility.

1.6 SYMBOLS AND STANDARDS USED

The KYELLE recliner is a non-invasive, non-active Class I medical device.

	Please refer to the instruction brochure/manual
	Indicates the maximum patient weight and the safe working load
	Compliant with the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council
	Serial number
	Date of manufacture
	Manufacturer
	Maximum patient weight
	Protect from moisture
	This end up
	Fragile
	Recyclable

General safety symbol:

- Information marked with this symbol must be read and strictly adhered to.
- Indicates that failure to comply with the recommendation may endanger patients, users or damage the equipment.

	Do not sit on the backrest		Do not sit on the armrest
	Do not sit on the leg rest		Do not lean on the recliner

“WARNING: Do not alter this device without the written consent of the manufacturer.”

“Note: only the use of accessories intended and provided for this MD by BMB MEDICAL is authorized”

2 INSTRUCTIONS FOR USE

2.1 INSTALLATION PROCEDURE

Before bringing the recliner into service, check the functioning of the following items, depending on the model:

- [Braking](#)
- [Armrests](#)
- [Backrest](#)
- [Leg rest](#)

Clean the product in accordance with the protocol specific to your facility.

2.2 ARMRESTS


Before handling the armrests, make sure nothing hinders their movement.

Pay attention to the risk of entrapment.

Lowering the armrests:

The armrests are released by pulling the indexing pin.

Raising the armrests:

Lift the armrest until the indexing pin locks and you hear a “click”

Before handling the armrests, make sure nothing hinders their movement.

Before raising the armrest, take care not to get the hands of the patient or personnel caught.

If you hear no “click” when your armrest locks into place, you should have your armrest serviced.


Do not use the armrests to move the product.

2.3 BACKREST

Before activating the backrest, make sure nothing hinders its movement (accessories, cables, etc.)

The backrest is adjusted using the lever to the right of the seat. To unlock the backrest, pull the lever up and adjust the backrest by manually tilting it to the desired position. Release the lever to secure the position. The backrest can be lowered into a flat position.

2.4 LEG REST

2.4.1.1 LEG REST SYNCHRONIZED WITH THE BACKREST

Before activating the leg rest, make sure nothing hinders its movement (accessories, cables, etc.)

The movements of the leg rest are synchronized with those of the backrest, as part of an interdependent adjustment system ([see paragraph 2.3 Backrest](#))

The short leg rest is fully retractable under the seat when in sitting position.

2.5 OPTIONS

2.5.1 GAS-SPRING ASSISTED, INDEPENDENT LEG REST

The leg rest is adjusted using the lever to the left of the seat. To unlock the leg rest, pull the lever up and adjust the leg rest until you reach the desired position.

Keep the lever in position to raise the leg rest. Release the lever to secure the position.


Keep the lever in position and push the leg rest to lower it.

2.5.2 FIXED PUSH BAR

2.5.3 REVERSIBLE PUSH BAR

Position 1

Position 2

When carrying a patient in supine position, take care to protect their feet.

2.5.4 HEADREST CUSHION

2.5.5 FOOTREST

Pull or push the yellow lever to raise or retract the footrest

2.6 BRAKING

2.6.1 WHEELS WITH INDEPENDENT BRAKE

Brake with the 2 back wheels by depressing the brake pedals.

Brake in released position

Brake in locked position

2.6.2 WHEEL WITH DIRECTIONAL BRAKE

Wheels with directional brake facilitate travel in a straight line

Free guidance

Directional guidance activated

2.6.3 CENTRALIZED BRAKING

Stop the recliner by pressing the brake bar down with your foot

Release the recliner by lifting the brake bar

3 CLEANING AND DISINFECTION

3.1 INTRODUCTION

Regular and proper cleaning guarantees the long service life of the device (see cleaning program established when bringing the device into service by appointing reference personnel).

Failure to perform cleaning operations may lead to a risk of infection due to the presence of secretions and/or residue on the recliner.

The cleaning instructions hereafter are not a substitute for your facility's internal protocol.

3.2 FREQUENCY

The device must be cleaned:

- For every patient, i.e. between two patients.
- After carrying an infected patient.
- Thoroughly once a week

3.3 CLEANING AND DISINFECTION RECOMMENDATIONS

3.3.1 CLEANING OPERATION

For the routine cleaning of upholstery, painted, stainless steel or aluminum surfaces, plastic parts and control elements, a mild detergent such as soapy water can be used, before rinsing with clear water. After cleaning, thoroughly dry the product.

All furniture surfaces must be cleaned.

When cleaning or disinfecting, it is advisable to:

- Spray or apply the cleaning product onto the MD using a clean, single-use mop.
- Start with areas closest to the patient, focusing on parts in contact with users.
- Work from the top down
- Allow the equipment to dry.

Stain-forming products (e.g.: "blue product" or eosin) must be cleaned as soon as possible to avoid any impregnation.

Airborne disinfection of surfaces (ADS) is authorized for this product.

3.4 RECOMMENDED EQUIPMENT AND PRODUCTS

- Only use authorized cleaners or disinfectants, standard quaternary, phenolic cleaners and bleach solutions...
- Comply with the recommended dosage and dilution.
- Comply with the conditions of application (contact time, deposited amount, rinsing, etc.)
- Do not mix products.
- Make sure the detergents and disinfectants used are compatible with:
 - o PVC, ABS, Polyamide, Polyurethane, Polypropylene
 - o Epoxy painted metal surfaces
 - o Metal surfaces made of stainless steel and aluminum
- Take note of the products strictly prohibited by the manufacturer (see list below).

List of strictly prohibited products:

Solvents:	Gasoline/higher aliphatic and aromatic alcohols/ketones-esters / trichlorethylene-perchlorethylene
Base:	10% soda/10% ammonia
Acids:	10% nitric/30% hydrochloric/10% citric
Miscellaneous:	Bleach, pure or diluted in a little water / Abrasive powders / Any abrasive product

Note:

- **Quaternary cleaners** are made up of ingredients containing sentence fragments such as "...yl" and "ammonium chloride"
- **Quat./isopropyl cleaners** are made up of a quaternary ingredient as mentioned above as well as isopropyl alcohol
- **Phenolic cleaners** are made up of ingredients containing the "-phenol" suffix
- **Bleach** is known under the generic name of "sodium hypochlorite"

Do not use abrasive utensils.

Use clean, single-use mops.

When cleaning and disinfecting, protective clothing and gloves must be worn.

High-pressure and steam cleaning is prohibited.

Damage caused by failure to comply with the instructions for use of a detergent/disinfectant shall under no circumstances be covered by BMB MEDICAL under the warranty.

4 PRECAUTIONS

4.1 PRECAUTIONS WITH REGARD TO THIS DEVICE

4.1.1 WARNINGS

This MD is designed for and tailored to medical practices and procedures, i.e. a specific use. Any abnormal use thereof may endanger users, patients or even cause irremediable deterioration.

The following list of abnormal uses of the recliner is not exhaustive, and the BMB MEDICAL Company counts on the users' common sense.

- Do not use the product without having read this user manual first.
- Used to place a person of short stature (under 4 ft 3') or under 12 years of age in supine position.
- Do not alter the recliner. Any system alteration may lead to unpredictable operation, which would compromise the safety of the user or patient. In this case, the recliner's warranty would be null and void.
- Only the use of accessories intended and provided for this MD by BMB MEDICAL is authorized
- When it no longer has to be moved, the recliner must be secured with the brakes. In this case, make sure the wheels are locked by attempting to move it. Failure to perform this operation may cause a person to fall as they try to lean on the recliner.
- Do not park the recliner on a slope.
- Do not use the MD on a slope steeper than 10%.
- When moving a patient, both armrests must be raised to prevent falls.
- When transferring a patient (patient moved from a bed to the recliner or from the recliner to a bed), the recliner must be positioned as close to the bed as possible to ensure maximum patient safety. Adjust the height of the bed so that it is identical to that of the recliner.
- When moving a patient in supine position, take care to protect their feet
- The patient must not sit on the leg rest or backrest.
- Tell the patient not to sit on the ends (leg or top) of the recliner. Excessive weight placed at either end may cause the lying surface to tilt over and injure the patient.
- Use of the recliner with a load exceeding the safe working load (please refer to the technical features).
- Tell the patient and his/her companions not to sit on the armrests.
- Before moving the product:
 - Make sure the brakes are disengaged. Trying to move the secured product may be dangerous for users as well as patients.
 - Make sure the patient's limbs are not outside the armrests and leg rest.

- 
- The recliner's mechanisms may cause lesions; make sure nobody touches mechanical parts before activating the recliner.
 - For the effective use of the equipment, only start moving if the weight is well distributed on the device.
 - Failure to perform cleaning operations may lead to a risk of infection due to the presence of secretions and/or residue on the recliner.
 - Any other use outside the recommended operation of the recliner.
 - Do not lift the leg rest manually.
 - Use by anyone other than caregivers who has not received clearance and explanations from a qualified person.
 - Use outside the facility or on unsuitable ground.
 - Washdown, high-pressure jet or steam washing and use of a tunnel washer.
 - Any other use outside the recommended operation of a hospital recliner.

All handling operations relating to the above precautions for use are described in paragraph II– Instructions for use.

4.1.2 CONTRAINDICATIONS

There is no contraindication for this device

4.1.3 PRECAUTION WHEN MAINTAINING AND KEEPING THE EQUIPMENT IN GOOD CONDITION

No alteration shall be made to the device without the prior consent of the manufacturer, as said alteration could lead to unpredictable operation, thus causing damage to the patient and caregivers.

Torn upholstery no longer provides an effective antibacterial barrier and must therefore be replaced without delay.

The device requires regular inspections, at least once a year. Please refer to the maintenance checklist. To keep the device in good condition while guaranteeing its safe and durable operation, the prescribed maintenance must be performed

If anomalies or damage are found, which are likely to affect the operation and safety of the device, the equipment must be immediately removed.

The device must be used exclusively in accordance with the instructions of this user manual.

The device must not be exposed to, let alone in contact with combustion sources and flammable materials.

* The BMB MEDICAL Company shall not be held liable for any incidents or accidents arising from the improper use or alteration of the product by the purchaser or the user.

4.2 TRANSPORT AND STORAGE

4.2.1 TRANSPORT

When transporting the recliner with a view to storage, to prevent any risk of impact or fall, it is advisable to:

- Strap it to a support of the means of transport
- Cover it
- Engage the brake
- Place it in the lowest position

Damage caused during transport is not covered by the warranty. The repair and replacement of damaged parts shall be borne by the customer.

4.2.2 STORAGE

When the recliner is in storage, make sure it is:

- Secure
- Covered
- In a warm, non-damp room (temperature: +10°C to +40°C).

No heavy weight must be placed on it during storage.

The recliner may not be used as a bearing plane for any type of equipment.

The device must not be exposed to, let alone in contact with combustion sources and flammable materials.

4.3 WARRANTY

2-year warranty, parts and labor, under the following terms:

BMB MEDICAL grants its customers a two-year commercial warranty from the delivery date within the scope defined hereafter.

This commercial warranty is not a substitute for the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer Code or that relating to the defects of the item sold, under the conditions provided for in articles 1641 to 1648 and 2232 of the French Civil Code.

Beyond the legal guarantee, the BMB MEDICAL Company guarantees the exchange and replacement of parts found to be unserviceable, due to a machining, assembly or material defect. Parts rendered unusable by normal wear and tear or by improper use or unauthorized repair are excluded from warranty coverage.

Where appropriate, at the request of BMB MEDICAL, any product or part for which a warranty claim has been made must be returned carriage paid to the BMB MEDICAL factory. Any improper use, alteration or repair made with parts that have not been supplied or authorized by BMB MEDICAL shall void this warranty.

BMB MEDICAL products are designed for an expected operating life of 10 years under normal conditions of use and with appropriate periodic maintenance, as described in the maintenance checklist of each product. BMB MEDICAL guarantees to the original purchaser that the welds on its products shall be free from structural defects for the product's 10-year expected operating life as long as the original purchaser is in possession of the product.

No warranty is made with respect to fitness for a particular purpose.

This warranty excludes disposable items, IV poles, mattresses, batteries and damage caused by misuse.

However, the following shall render the product warranty null and void:

- Abnormal use ([Warning paragraph here](#))
- Work without advance notification to our technical or maintenance departments on:
 - o A mechanical part
 - o A jack
- Failure to comply with any other warranty condition
- Any unsuitable cleaning
- Use of accessories that have not been supplied by the manufacturer
- Failure to adhere to the preventive maintenance operations recommended in the maintenance checklist.

4.4 CUSTOMER SERVICE

For any information on the installation, use and maintenance of the product or any other question, please contact the Customer Service department at: sav@bmb-medical.com or your local distributor outside France.

With your “Customer Space” on the www.bmb-medical.com website, you can access your equipment online, request and monitor interventions directly via this secure interface.

Operators should carefully read the information contained in this user manual, in particular that relating to precautions for use in terms of safety as well as installation and operating methods.

If you have any doubt about the correct interpretation of instructions, feel free to contact the Customer Service department who is always at your disposal.

To make it easier to order spare parts, you can ask the Customer Service department for the exploded views of your products, to obtain the references of the parts you need.

4.5 ENVIRONMENT

To be disposed of separately from household waste.

Dispose of the product or its components in accordance with the local and national laws and regulations governing environmental protection and raw materials recycling.

Ask the authorities responsible for environmental protection where to find the appropriate waste disposal site.

Specific to France: BMB MEDICAL fulfills its obligations relating to the end of life of the products in the TIM range it places on the market by funding the Valdelia recycling channel which takes them back free of charge (Further information at www.valdelia.org).

Tel.: +33 (0)4 74 08 71 71

Fax: +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001

01604 TREVOUX Cedex

FRANCE