TIM VELVET SWANN

MANUEL D'UTILISATION
USER MANUAL

 **BMB
MEDICAL**

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FRANÇAIS

BMB MEDICAL

100 Avenue du Formans

CS 11001

01604 TRÉVOUX CEDEX

Tél : +33(0)4 74 08 71 71

Fax : +33(0)4 74 08 71 72

www.bmb-medical.com

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Fax: +33(0)4 74 08 71 72

www.bmb-medical.com

The purpose of this manual is to familiarize you with the use, maintenance and safety of your new medical device. We encourage you to read this document carefully and to follow the guidelines for optimum and safe use.

BMB MEDICAL offers a broad variety of options for the products covered by this user manual. Consequently, some features described in this manual may not apply to your chair.

Some illustrations are provided for information purposes and do not accurately represent the device in your possession.

The information in this document are subject to change without notice and in no way engage the responsibility of BMB MEDICAL. This document may not be reproduced, in whole or in part, without prior written consent by BMB MEDICAL. This illicit representation or reproduction, by any means, constitutes counterfeiting. The illustrations in this document are not binding.

TABLE OF CONTENTS

ENGLISH	36
1 INTRODUCTION SPECIFICATIONS	39
1.1 PRODUCT PRESENTATION	39
1.1.1 TIM VERSATILE FIXED-HEIGHT CARE CHAIR	40
1.1.2 VELVET VERSATILE FIXED-HEIGHT CARE CHAIR	41
1.1.3 SWANN VERSATILE FIXED-HEIGHT CARE CHAIR	42
1.2 PURPOSE OF THE USER MANUAL	43
1.3 USERS	43
1.4 STORAGE	43
1.5 APPLIED PARTS	43
1.6 SYMBOLS AND STANDARDS USED	44
2 INSTRUCTIONS FOR USE	46
2.1 OPERATING CONDITIONS	46
2.2 INSTALLATION PROCEDURE	46
2.3 BATTERY (OPTIONAL) & CONNECTION TO MAINS POWER	46
2.4 REMOTE CONTROL	48
2.5 PROGRAMMED POSITIONS	48
2.5.1 TRENDLENBURG POSITION	48
2.5.2 VARIABLE HEIGHT (optional)	49
2.6 BACKREST & LEGREST	49
2.6.1 SYNCHRONIZED ADJUSTMENT OF BACKREST & LEGREST (TIM)	49
2.6.2 MANUAL BACKREST (TIM)	50
2.6.3 ELECTRIC BACKREST (VELVET & SWANN)	50
2.6.4 ELECTRIC LEGREST	50
2.7 PIVOTING HEADREST	51
2.8 BRAKE AND GUIDE PEDAL	51
2.8.1 CENTRAL BRAKING AND SIDE GUIDES (OPTIONAL)	52
2.9 ARMRESTS	52
2.9.1 REMOVABLE COMFORT ARMRESTS	52
2.9.2 FIXED ALUMINUM ARMRESTS	52
2.9.3 RETRACTABLE ARMRESTS	53
2.10 OPTIONS	53
2.10.1 IV POLE	53
2.10.2 PAPER ROLL HOLDER	54
2.10.3 READING LAMP + USB PORT	54
2.10.4 ANATOMICAL ARMREST FOR SAMPLING	54

2.10.5	CYLINDRICAL NECK CUSHION.....	55
2.10.6	MEAL TRAY	55
2.10.7	PUSH BARS.....	55
2.10.8	STAINLESS STEEL BASE.....	55
2.10.9	CRYSTAL PROTECTION.....	56
2.10.10	LIFTER OPTION	56
3	CLEANING AND DISINFECTION	58
3.1	INTRODUCTION	58
3.2	FREQUENCY	58
3.3	CLEANING AND DISINFECTION ADVICE.....	58
3.3.1	CLEANING OPERATION	58
3.3.2	RECOMMENDED EQUIPMENT AND PRODUCTS.....	59
4	PRECAUTIONS.....	60
4.1	DEVICE-RELATED PRECAUTIONS	60
4.1.1	WARNINGS	60
4.1.2	CONTRAINDICATIONS	61
4.1.3	PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE	61
4.1.4	EMC-RELATED INFORMATION.....	62
4.2	TRANSPORT AND STORAGE	64
4.2.1	TRANSPORT	64
4.2.2	STORAGE	64
4.3	WARRANTY	64
4.4	AFTER-SALES SERVICE	66
4.5	ENVIRONMENT	66
4.5.1	SPENT EQUIPMENT.....	66

1 INTRODUCTION SPECIFICATIONS

1.1 PRODUCT PRESENTATION

This patient support device consists of a backrest and seat fitted onto a chassis with wheels. The device is equipped with armrests and, optionally, allows the temporary fitting of IV poles.

The device allows an operator to convey the patient within a healthcare establishment. The chair is braked by a mechanical system actuated by pressing a pedal.

The device is equipped with an electrically assisted actuator system actuated by a remote control to manage the movements of the following elements:

- backrest,
- legrest,
- electrical column used to adjust the product's height.

The products are designed to ensure that there are no unacceptable risks (essential performance) associated with patient support or caused by accidental movements.

1.1.1 TIM VERSATILE FIXED-HEIGHT CARE CHAIR

TECHNICAL DATA	
REFERENCE	TIM
DESIGNATION	TIM MANUAL RELAXATION AND REST CHAIR
SAFE WORKING LOAD	130 kg (indicates the maximum admissible load including the total weight of the patient and accessories)
WEIGHT	35 kg depending on the model
OVERALL DIMENSIONS	Width: 73 cm Seat width: 57 cm Overall length: from 85 cm in sitting position to 153 cm in reclined position. Overall height: 115 cm.
EQUIPMENT	* Assisted synchronized backrest and legrest until the confort position (Gas assist cylinder) * The backrest goes independently of the legrest until the recovery position * Reinforced steel frame * Rack adjustment headrest * Seat on coil springs * Removable comfort armrests * Front exit * Retractable footrest under the seat * Tubular frame : - Front: 2 adjustable support pads - Back : 2 casters (Ø50 mm) * Seat and back cushions in HR 35 Kg density foam, fireproof, and covered * VALENCIA color range: M2 fireproof coating, 10 colors
AVAILABLE OPTIONS	Reading lamp + USB Port Fixed aluminium armrests 4 self lockable casters (Ø50 mm) Inox base Ø 690mm Star base : 3 back casters + 2 front fixed pads Neck cylindrical cushion Legrest crystal protection

1.1.2 VELVET VERSATILE FIXED-HEIGHT CARE CHAIR

TECHNICAL DATA	
REFERENCE	VELVET
DESIGNATION	VELVET ELECTRIC CHAIR
SAFE WORKING LOAD	160 kg (indicates the maximum admissible load including the total weight of the patient and accessories)
WEIGHT	40 kg depending on the model
OVERALL DIMENSIONS	Width (with ALU armrests) : 82 cm Width (with comfort armrests) : 74 cm Seat width: 57 cm. Overall length: from 100 cm in sitting position to 177 cm in reclined position. Overall height: 125 cm.
EQUIPMENT	Electric independant backrest and legrest (Electric cylinder) Reinforced steel frame, epoxy coated Rack adjustment headrest Seat on coil springs Foamed fireproof M1 vinyl, with silicone dacron Armrests with 7 cm HR foam padding and memory foam Front exit Retractable footrest under the seat Fixed height Star base with self lockable casters $\varnothing$ 75 mm Remote control holder Seat and back cushions in HR 35 Kg density foam, fireproof VALENCIA color range: M2 fireproof coating, 10 colors
AVAILABLE OPTIONS	* Emergency electric recovery position on remote control * Reading lamp + USB Port * Bilateral holder for IV Pole (incompatible FA022) * Paper roll dispenser * Removable aluminium armrests 2 positions * Removable comfort armrests - incompatible with FA018V-G/D * Neck cylindrical cushion * Small Shelf : 42 x 32 cm (FA018V-D or FA018V-G necessary) * Large Shelf : 81 x 41 cm (FA018V-D or FA018V-G necessary) * Support for meal tray on the right side of the patient (compatible with the FA018V-G) * Support for meal tray on the right side of the patient (compatible with the FA018V-D) * 4 self lockable casters ($\varnothing$50 mm) * 2 back casters + 2 front fixed pads * Inox base $\varnothing$ 690mm * Legrest crystal protection * Removable IV pole, 2-stage, 2 hooks (FA0088VE mandatory) * Removable IV pole, 2-stage, 4 hooks (FA0088VE mandatory)

1.1.3 SWANN VERSATILE FIXED-HEIGHT CARE CHAIR

TECHNICAL DATA	
REFERENCE	SWANN
DESIGNATION	SWANN ELECTRIC CHAIR
SAFE WORKING LOAD	160 kg (indicates the maximum admissible load including the total weight of the patient and accessories)
WEIGHT	80 kg depending on the model
OVERALL DIMENSIONS	Width: 80 cm Seat width: 57 cm. Overall length: from 104 cm in sitting position to 183 cm in reclined position. Overall height: 131 cm.
EQUIPMENT	Fixed height tubular frame 2 pivoting self lockable back casters ø 125 mm 2 fixed front casters ø 125 mm Emergency electric recovery position on the remote control Aluminium armrest Rack adjustment headrest 2 holders for IV pole Electric backrest and legrest (Electric cylinders) Push bar PVC coated jersey VALENCIA color range: M2 fireproof coating, 10 colors
AVAILABLE OPTIONS	* Reading lamp + USB Port * IV Pole * Battery * Paper roll dispenser * Phlebotomy armrest - Dimensions 400/200/40 mm (FA012S + FA015 or FA015B required) * Support for Phlebotomy armrest- The pair (FA015 or FA015B required) * Electric height adjustment (Electric cylinders) * Variable height with lift (Electric cylinders) * Two-sided foot pedals for Trendelenburg adjustment (Electric height adjustment necessary FA014) * Retractable armrests 3 positions (FA014 necessary) * Retractable aluminium armrests 2 positions * Small Shelf : 42 x 32 cm (FA018S-D or FA018S-G necessary) * Large Shelf : 81 x 41 cm (FA018S-D or FA018S-G necessary) * Support for meal tray on the right side of the patient (compatible with the FA018S-G) * Support for meal tray on the right side of the patient (compatible with the FA018S-D) * Back double running castors Ø 150 mm with central lock. * 2 pivoting self lockable front casters ø 100 mm * Beam of light * Neck cylindrical cushion * Removable legrest crystal protection

1.2 PURPOSE OF THE USER MANUAL

Healthcare professionals, this user manual was created for you to provide you with the necessary information and precautions to observe for autonomous and safe use.

The product's service life is dependent upon observance of the maintenance and tracking instructions.

It is thus essential that all staff has access to and be informed of the handling procedures and risks associated with the use of this chair. The latter should be handled by a qualified person who is aware of the risks inherent to improper use. Additional training is not required as we consider that knowledge of the user manual is sufficient.

1.3 USERS

Healthcare professionals: physicians, caregivers and agents approved by Technical Services can operate all functions.

In use cases where the patient is the operator, the caregivers must first ascertain the patient's ability to use the product. The caregivers must also teach the patient how to safely use the product: adjustment to low position to step off the product.

1.4 STORAGE

This user manual should be kept near the product, in a place provided for this purpose. It must be shielded from all elements and substance that could compromise its perfect legibility.

1.5 APPLIED PARTS

The applied parts, according to the EN 60601-1 standard, are as follows:

- upholstery
- armrests
- remote control
- reading lamp

1.6 SYMBOLS AND STANDARDS USED

	Type B appliance
	Electrical class II equipment
	Electrical hazard, do not open
	Duty cycle (motor cycle) 2 min. operation / 18 min. stop
IPX4	Splash-resistant
	See the manual / instruction brochure
	Indicates the maximum permissible load, including the total weight of the patient, mattress (where applicable) and accessories.
	Conforms to the essential requirements of directive 93/42/EEC on medical devices.
	Serial number
	Product reference
	Manufacture date
	Manufacturer
	Alternating current
	Electronic waste: this symbol indicates that the product must not be discarded with unsorted municipal waste, but rather to a selective waste collection channel. Contact the local distributor for information concerning the return and/or collection systems available in the country.
	Protéger de l'humidité / keep dry
	Haut, tenir debout / stand up
	Fragile / fragile
	Recyclable / recycling
	Temperature range
	Humidity limit

General safety sign:

- **Information highlighted by this symbol must be read and scrupulously observed.**
- **Indicates that failure to observe the guideline could endanger the patient or user, or damage the equipment.**

	Sitting on the backrest is prohibited		Leaning against the chair is prohibited
	Sitting on the legrest is prohibited		The chair must be set to its low position before being moved
	Sitting on the armrests is prohibited		

“WARNING: This appliance must not be altered without the manufacturer’s written consent.”

Caution: only purpose-built accessories provided for this MD by BMB MEDICAL may be used.

The TIM is a Class I non-invasive and non-active medical device.

The VELVET is a Class I non-invasive and active medical device.

The SWANN is a Class I non-invasive and active medical device.

2 INSTRUCTIONS FOR USE

2.1 OPERATING CONDITIONS

AMBIENT CONDITIONS	OPERATION
Ambient temperature	+5 °C to +40 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

2.2 INSTALLATION PROCEDURE

Before commissioning the chair and depending on the model, check that the following are operational:

- Remote control
- Variable height
- Backrest
- Legrest
- Armrests
- Brake and guide pedals

Clean the product according to your establishment's internal protocol.

2.3 BATTERY (OPTIONAL) & CONNECTION TO MAINS POWER

The electronic unit is equipped with its own charger (the battery is charged with the chair cable). Before initial use, leave the product to stand for 12 h in a temperate environment (to allow the electrical and electronic components to reach normal operating temperature), then charge the battery for 12 hours from a mains outlet.

In routine use, we strongly recommend keeping the chair plugged into the mains during sedentary use.

The battery, with its limited operating life, can be used to operate the chair in the event of a power outage, or during transfer.

Once the manipulations are complete, remember to re-connect the chair to the mains (at least 12 hours per day). Never allow the battery to drain fully (risk of reduced battery performance).

“Warning: Do not short-circuit the batteries and do not leave the system unused and without charging for extensive periods of time as this could cause self-discharge, leading to medium-term damage. This will also have an environmental impact. »

Battery characteristics		Battery safety symbols	
Battery type:	Lead gel	IPX6	Protection rating
Capacity:	1.2 Ah		CE marking
			Indoor use
Input supply voltage:	24 VDC, 1.2 A		Electronic waste

Power grid connection conditions and electromagnetic compatibility:

This medical device's power cord may represent an obstacle that could cause a fall. Be mindful of its presence when moving around the appliance, when the latter is connected to the mains.

Position the chair such that it can be easily unplugged.

Unplug the product before moving it.

Do not leave the cables in contact with the floor or the casters,

Do not knot the cables,

Periodically check the condition of the cables.

Mains connection:	220 V
Frequency:	50 Hz
Rated current:	3.5 A
Protection rating / degree of electric shock protection:	Class 2 appliance
Power input:	770 W
Intermittent operating mode:	10%, i.e. 2 minutes of continuous operation followed by 18 minutes of rest.

The mains plug is the electrical grid insulation device.

Work on electrical parts should only be performed by qualified and authorized staff.

Electromagnetic interference caused by the chair cannot be excluded. In this case, the chair's functions should be inhibited and those instruments subject to interference should be moved away. You can view [the EMC information here](#).

2.4 REMOTE CONTROL


When using the remote control, all movements can be interrupted by simply releasing the corresponding key.

Only activate the various positions after having ensured that no-one is touching any mechanical parts and that there are no objects under the seat, or between any of the moving structures.

VELVET

SWANN

SWANN with
Variable height with lift

2.5 PROGRAMMED POSITIONS

Sitting position

Reclining position

2.5.1 TRENDLENBURG POSITION

The Trendelenburg position is obtained by a long press on

2.5.1.1 BILATERAL TRENDLENBURG FOOT PEDAL (OPTIONAL)

Before activating Trendelenburg, check that nothing hinders the movement (accessories, cables, etc.)

Press the blue pedal to move to the Trendelenburg position.


2.5.2 VARIABLE HEIGHT (optional)


Before using these functions, ensure that there are no objects or persons under the patient surface.

Height is adjusted using the remote control:

2.6 BACKREST & LEGREST

Before actuating the legrest and backrest, check that there are no obstacles to their movement (accessories, cables, etc.)

2.6.1 SYNCHRONIZED ADJUSTMENT OF BACKREST & LEGREST (TIM)

❶ The transition from patient installation position to the comfort position is done by gravity-driven by leaning back and pushing on the armrests.

❷ Transition from the comfort position to the patient installation position is done by moving the legs to fold the leg rest under the seat.

❶

Patient installation

❷


Comfort position

2.6.2 MANUAL BACKREST (TIM)

The backrest is adjusted by tilting the torso back and pushing on the armrests.

2.6.3 ELECTRIC BACKREST (VELVET& SWANN)

Before actuating the backrest, ensure that there are no obstacles to its movement (accessories, cables, etc.)

The backrest is adjusted with the remote control:

2.6.4 ELECTRIC LEGREST

Before actuating the legrest, ensure that there are no obstacles to its movement (accessories, cables, etc.)

The legrest is controlled using the remote control:


2.7 PIVOTING HEADREST

The position of the headrest is adjusted manually:

Pull the headrest forwards to tilt it.

To tilt it backwards, move it forwards, ❶ then backwards ❷.

2.8 BRAKE AND GUIDE PEDAL

Except during transport, always ensure that the chair's brake is engaged. Check whether the brake is engaged by attempting to move the chair.

Do not park the chair on a sloping floor.

Star base on casters with independent brakes

Tubular base 2 casters with independent brakes + 2 front glides

Wheel in locked position

Wheel in free position

To move the TIM: manually unfold the legrest and disengage the rear wheel brakes. The wheels are free.

2.8.1 CENTRAL BRAKING AND SIDE GUIDES (OPTIONAL)



Except during transport, always ensure that the chair's brake is engaged. Check whether the brake is engaged by attempting to move the chair.

Do not park the chair on a sloping floor.

Brake and side guide pedals:

- Brake engaged - red side lowered **①**.
- Chair in directional guidance mode - green side lowered **②**.
- Chair in free guidance mode - pedal in median position **③**.

2.9 ARMRESTS

2.9.1 REMOVABLE COMFORT ARMRESTS

2.9.2 FIXED ALUMINUM ARMRESTS

2.9.3 RETRACTABLE ARMRESTS

To raise or lower the armrests, simply pull the index pin and raise or lower the armrest,

To ensure that the armrest is locked in place, check that the index pin has returned to its housing.

2.10 OPTIONS

2.10.1 IV POLE

The IV pole is removable; it fits into the holder at the back of the product. This holder is rotary.

To deploy the IV pole, pull on the telescopic part.

To lower the IV pole, lift the ring

Safe working load: 6 kg.

This indicative value corresponds to safe product use, but is not representative of normal use for which the load is significantly lower.

When moving the chair fitted with the IV pole, ensure that the latter does not strike any obstacles (e.g. door frame, lighting, etc.) or individuals.

Do not move the chair via the IV pole.

2.10.2 PAPER ROLL HOLDER


Do not push or pull the product via the paper roll holder.

2.10.3 READING LAMP + USB PORT

The lamp plugs into an AC outlet.

The lamp is turned on and off simply by touching the silver part.

Light intensity is controlled by prolonged contact with the silver part.

2.10.4 ANATOMICAL ARMREST FOR SAMPLING

The sampling armrest fits into the tube located on the side of the seat:

- slide the armrest rod into the tube,
- and lock it by tightening the knob.


The armrest is directional thanks to its ball joint system. To adjust its position:

- loosen the knob,
- place the armrest in the desired position,
- then tighten the knob to lock the armrest in position.

2.10.5 CYLINDRICAL NECK CUSHION

2.10.6 MEAL TRAY

The tray swivels on its axis. Maximum weight: 15 kg .

2.10.7 PUSH BARS

2.10.8 STAINLESS STEEL BASE


2.10.9 CRYSTAL PROTECTION

The protection can be easily positioned without the need for tools or pre-equipment. The footrest protection is attached by a strap and clip mechanism and is adjustable

2.10.10 LIFTER OPTION

The "lift" option on the remote control is **locked** automatically.

To unlock it, press the buttons marked ① and ② until you hear a long "beep".

Once the option is unlocked, you can activate the lifting aid by pressing the button marked ①.

To return the chair to the sitting position you can press the button marked ②

Effortless lifting assistance

Before using the lift assist feature, make sure no limbs are in the area or in contact with the lift mechanism. Also ensure that no objects are under the chair.

If it is necessary to remove objects from underneath the wheelchair, you must unplug the power supply.

Use the lift function only on a level surface to ensure the stability of the product.

Always brake the product before operating the lift.

Do not sit on the wheelchair in the lift position.

Be sure to leave the chair in the sitting position when not in use.

Be sure to unfold the dining tray before activating the lift position to avoid the risk of entrapment

3 CLEANING AND DISINFECTION

3.1 INTRODUCTION

Regular and thorough cleaning will guarantee the device's life (see cleaning program drawn up at the time of device commissioning, designating the reference staff).

Failure to perform these cleaning operations can lead to a risk of infection caused by the presence of secretions and/or residues on the chair.

The following cleaning advice is not intended to replace your establishment's in-house protocol. This chair has been designed to facilitate cleaning.

3.2 FREQUENCY

The chair should be cleaned:

- Between each patient.
- After transporting an infected patient.
- Completely each week.

3.3 CLEANING AND DISINFECTION ADVICE

3.3.1 CLEANING OPERATION

Before cleaning the chair, unplug it from the mains and ensure that all internal connections are correctly established.

When cleaning, check that the upholstery does not show any cuts, tears, punctures, cracks or pulled threads. Never use a mattress or upholstery with a damaged cover. Replace this latter before the internal part of the mattress or upholstery becomes contaminated. It should also be replaced if stained.

For routine cleaning of the upholstery, painted stainless steel or aluminum surfaces, plastic parts and control elements, a gentle detergent such as soapy water may be used, followed by rinsing with clear water. After cleaning, dry the product thoroughly.

All of the product's surfaces must be cleaned.

During cleaning or disinfection, it is recommended to:

- Spray or apply the cleaning product to the MD with a clean disposable cloth.
- Start with the areas closest to the patient, focusing on parts in contact with users.
- Proceed from top to bottom.
- Allow the equipment to dry.

Products that can cause stains (e.g.: "blue product" or eosin) must be cleaned as rapidly as possible to prevent them from soaking in.

Aerial disinfection of surfaces is allowed for this product.

In the case of a product fitted with a USB charging port, ensure that no liquids enter the port.

3.3.2 RECOMMENDED EQUIPMENT AND PRODUCTS

- Use only authorized and standardized cleaning products or disinfectants: quaternary cleaning products, phenolic compounds and bleach solutions, etc.
- Observe the prescribed dosages and dilutions.
- Observe the application conditions (contact time, quantity applied, rinsing, etc.)
- Do not mix the products.
- Ensure that all disinfectant detergents used are compatible with:
 - o PVC, ABS, Polyamide, Polyurethane and Polypropylene
 - o Epoxy paint-coated metal surfaces
 - o Stainless steel and aluminum surfaces
- Familiarize yourself with the products strictly prohibited by the manufacturer (see list below).

List of strictly prohibited products:

Solvents:	Gasoline / aliphatic and aromatic higher alcohols / ketones-esters / trichloroethylene-perchloroethylene
Bases:	10% Sodium hydroxide / 10% Ammonia
Acids:	10% Nitric / 30% Hydrochloric / 10% Citric
Miscellaneous:	Pure or weakly diluted bleach / Abrasive powders / All abrasive products

NB:

- **Quaternary cleaning products** consist of ingredients containing “-yl” suffix elements and ammonium chloride
- **Quat./isopropyl cleaning products** consist of a quaternary ingredient, as described above, along with isopropyl alcohol
- **Phenolic cleaning products** consist of ingredients containing the “-pheno” suffix
- **Bleach** is known under the generic name “sodium hypochlorite”

Never use abrasive utensils.

Use clean disposable cloths.

Wear protective clothing and gloves during cleaning and disinfection.

High-pressure and steam cleaning are prohibited.

Any damage resulting from failure to observe this Instructions for Use of a detergent-disinfectant will not be covered by BMB MEDICAL under the terms of the warranty.

4 PRECAUTIONS

4.1 DEVICE-RELATED PRECAUTIONS

4.1.1 WARNINGS

This chair has been designed for and adapted to medical practices and procedures, and hence a specific use. Abnormal use of the latter can lead to user or patient accidents, or even to irreversible damage.

The following list of abnormal chair uses is not exhaustive and BMB MEDICAL counts on the common sense of the chair's users.

- The product must be connected to a mains power outlet compliant with standards in force.
- Once the brakes have been engaged, patient sitting and standing should be performed, with the seat in the low position:
 - o From the side by lowering the armrests (and returning the chair to the patient installation position)
 - o or from the front by returning the chair to the patient installation position with the help of the nursing staff.
- When getting up, the patient must not rest on the articulated anatomical armrests.
- The chair's brakes must be engaged when the patient sits down and stands up
- The patient must not sit on the legrest or backrest.
- In the case of an unassisted patient, the caregivers must inform the patient of the safety measures that must be implemented when using the product.
- Tell the patient not to sit on the ends (foot or head) of the chair. An excess weight applied to one of the ends could tip the patient surface and injure the patient.
- Tell the patient and any accompanying persons not to sit on the armrests.
- Before moving the product:
 - o Ensure that the power cord is unplugged.
 - o Ensure that the brakes are released. It may be hazardous for the user and patient when attempting to move the product with the brakes engaged.
 - o Ensure that the patient's limbs do not extend beyond the armrests and legrest.
 - o The product must be placed in the lowered and comfort position.
- The chair's mechanisms could cause lesions; ensure that no-one is touching the mechanical elements before activating the chair.
- Before using the lift assist feature, make sure no limbs are in the area or in contact with the lift mechanism. Also ensure that no objects are under the chair. If it is necessary to remove objects from underneath the wheelchair, you must unplug the power supply.
- -Use the lift function only on a level surface to ensure the stability of the product.
- Always brake the product before operating the lift.
- Do not sit on the wheelchair in the lift position. Be sure to leave the chair in the seated position when not in use.
- Be sure to unfold the dining tray before activating the "lift" position in order to

- to avoid the risk of entrapment.
- To ensure correct use of the product, only proceed if the weight is correctly distributed over the device.
- The “variable height” function should be used when handling the patient in order to work at an ergonomic height.
- Once movement is no longer required, the chair should be immobilized with the brakes. Next, check that the wheels are locked by attempting to move the chair. Failure to observe this manipulation could lead to the fall of a person leaning against the chair.
- Do not use the MD on a floor with a slope of more than 10%.
- To avoid the risk of patient or user injury of IV pole damage when moving the chair, ensure that the IV pole is firmly secured.
- To avoid damaging the equipment during transport, ensure that the IV pole passes safely through all door frames and under all light fittings.
- Never use the chair in the presence of an anesthesia mixture flammable with air, with oxygen, or with nitrous oxide as this could cause a fire and/or explosion.
- To prevent failure, electromedical equipment requires special [EMC-related precautions, described here](#).
- Do not manually raise the legrest.
- Work on electrical parts should only be performed by qualified and authorized staff.
- Used to accommodate a person of small stature (less than 1.30 m tall), or under the age of 12 years.
- Do not alter the chair. Any changes made to the system could lead to unpredictable operation that could in turn jeopardize the safety of the user or patient. This would also void the chair’s warranty.
- To avoid material damage or injury, do not use the IV pole to pull or push the product.
- To avoid damaging the equipment, the maximum permissible load of the IV pole is 6 kg.
- Use by a person who has not first familiarized themselves with this user manual.
- Use by any non-caregiver who has not received the authorization and explanations from a qualified individual.
- Use of the chair with a load that exceeds the safe working load (see technical characteristics).
- Use outside of the establishment, or on unprepared ground.
- Washing with large amounts of water, with a high-pressure jet or with steam and passage through a wash tunnel.
- Connection, by any means, of an electric device to the chair.
- Any use that does not fall within the operating context of a hospital chair.

All manipulations relating to the aforementioned precautions for use are described in [paragraph II - Instructions for Use HERE](#).

4.1.2 CONTRAINDICATIONS

There are no other contraindications for this device

4.1.3 PRECAUTION ASSOCIATED WITH EQUIPMENT CARE AND MAINTENANCE

The device must not be altered in any way without the manufacturer's prior consent: the alteration could lead to unpredictable operation, potentially injuring the patient and caregivers.

Torn upholstery no longer presents an effective antibacterial barrier and must hence be replaced as soon as possible.

The device must be regularly inspected, at least twice per year.

To maintain the device in good condition and to guarantee its safe and durable operation, the prescribed maintenance must be performed:

- Ensure that all wheels are securely attached.
- Visually check all fastenings.
- Check the cylinder fastenings.
- Check that the backrest, legrest and all fastenings are working correctly.
- Lubricate all the moving mechanisms.

To keep the system operating optimally, repeat this check regularly.

If you have any doubts about the functioning and integrity of the product, take it out of service and inform your technical service or BMB MEDICAL after-sales service.

In the event of observed anomalies or damage that could jeopardize device operation and safety, the equipment must be immediately withdrawn.

The device should only be used according to the instructions in this user manual.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

Work on electrical parts should only be performed by qualified and authorized staff.

* BMB MEDICAL will accept no liability in the event of incidents or accidents resulting from improper use or alteration of the product by the buyer or user.

4.1.4 EMC-RELATED INFORMATION

VELVET & SWANN care chairs are suitable for use in all healthcare establishments, with the exception of facilities in the vicinity of HF ELECTROSURGICAL APPLIANCES, or near an electromagnetically shielded room for an EM magnetic resonance imaging system, in which the intensity of EM INTERFERENCES is high.

VELVET & SWANN care chairs are suitable for use in domestic homes or establishments or on any premises connected to the public low-voltage network supplying domestic buildings.

EM interference could cause the product to make unexpected movements.

This appliance should not be used next to, or stacked on other appliances as this could lead to incorrect operation. If this mode of operation is necessary, this and the other appliances should be observed to ensure their correct operation.

The use of accessories, transducers or cables other than those specified or supplied by the manufacturer can lead to an increase in electromagnetic emissions, or a decrease in this appliance's immunity, causing inappropriate operation.

Do not use portable RF communication devices (including peripherals such as antenna cables and external antennas) less than 30 cm (12 inches) from any part of the APPLIANCE, including the cables specified by the manufacturer. Otherwise, the performance of these devices could be impaired.

List of cables used:

1. Power cord, maximum length 270 cm.
2. Variable height cylinder cable, maximum length 70 cm
3. Y cable for remote controls/motors, maximum length 130 cm

EMISSION	
Conducted and radiated RF emissions:	CISPR11 Group 1 Class A
Harmonic distortions IEC 61000-3-2:	Class A
Voltage fluctuations and Flicker IEC 61000-3-3:	Compliant

Immunity test	Test level IEC 60601	Compliance level																																	
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV in contact ± 2/4/8/15 kV in air	± 8 kV in contact ± 2/4/8/15 kV in air																																	
RF radiated electromagnetic fields IEC 61000-4-3	3 V/m (professional healthcare establishment) 80 MHz-2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2.7 GHz 80% AM at 1 kHz																																	
Proximity fields emitted by wireless RF communication devices IEC 61000-4-3	<table><tr><th>Frequency (MHz)</th><th>Modulation</th><th>Required level (V/m)</th><th>Compliance level (V/m)</th></tr><tr><td>385</td><td>Pulsed modulation: 18 Hz</td><td>27</td><td>27</td></tr><tr><td>450</td><td>Pulsed modulation: 18 Hz</td><td>28</td><td>28</td></tr><tr><td>710-745-780</td><td>Pulsed modulation: 217 Hz</td><td>9</td><td>9</td></tr><tr><td>810-871-930</td><td>Pulsed modulation: 18 Hz</td><td>28</td><td>28</td></tr><tr><td>1720-1845-1970</td><td>Pulsed modulation: 217 Hz</td><td>28</td><td>28</td></tr><tr><td>2450</td><td>Pulsed modulation: 217 Hz</td><td>28</td><td>28</td></tr><tr><td>5240-5500-5785</td><td>Pulsed modulation: 217 Hz</td><td>9</td><td>9</td></tr></table>			Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)	385	Pulsed modulation: 18 Hz	27	27	450	Pulsed modulation: 18 Hz	28	28	710-745-780	Pulsed modulation: 217 Hz	9	9	810-871-930	Pulsed modulation: 18 Hz	28	28	1720-1845-1970	Pulsed modulation: 217 Hz	28	28	2450	Pulsed modulation: 217 Hz	28	28	5240-5500-5785	Pulsed modulation: 217 Hz	9	9
	Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)																															
	385	Pulsed modulation: 18 Hz	27	27																															
	450	Pulsed modulation: 18 Hz	28	28																															
	710-745-780	Pulsed modulation: 217 Hz	9	9																															
	810-871-930	Pulsed modulation: 18 Hz	28	28																															
	1720-1845-1970	Pulsed modulation: 217 Hz	28	28																															
	2450	Pulsed modulation: 217 Hz	28	28																															
5240-5500-5785	Pulsed modulation: 217 Hz	9	9																																
Rapid / burst electrical spikes IEC 61000-4-4	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz	Power supply: ± 2 kV Input/output lines: ± 1 kV Repeat frequency: 100 kHz																																	
Shock waves IEC 61000-4-5	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases ± 0.5 kV, ± 1 kV, ± 2 kV																																	
Conducted interference, induced by RF fields IEC 61000-4-6	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	3 V 0.15 MHz — 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz																																	
Power frequency magnetic fields IEC 61000-4-8	30A/m	30A/m																																	
Voltage Dips and Interruptions: IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle at 0° 70% UT; 25/30 cycles at 0° 0% UT; 250/300 cycles																																	

NB: There is no deviation from the requirements of IEC 60601-1-2 ed. 4

NB: There are no other known measures for maintaining basic safety with respect to EMC phenomena.

4.2 TRANSPORT AND STORAGE

Transport and storage conditions :

Ambient conditions	Operation
Ambient temperature	0 to +50 °C
Relative humidity	20 to 90% at 30 °C - no condensation
Atmospheric pressure	700 to 1060 ha

4.2.1 TRANSPORT

During chair transport with a view to its storage, the following is recommended to prevent any impacts or falls:

- Strap it to a surface of the means of transport
- Cover it
- Engage the brakes
- Set it to its low position

Damage caused during transport is not covered by the warranty. The costs of repair and replacement of damaged parts will be borne by the customer.

4.2.2 STORAGE

When the chair is in storage, ensure that:

- Its brakes are engaged
- It is covered

The product should be stored in a clean, ventilated, temperate and dry place, protected from the elements and from direct sunlight.

During storage, no heavy weights should be placed on the chair.

The chair must not be considered as a support surface for any other type of equipment.

The device must not be exposed to, or even less be in contact with, sources of heat, combustion and flammable agents.

In the event of prolonged storage, the battery should be recharged every 6 months, proceeding as follows:

- Connect the mains power cord.
- Briefly press one of the remote control functions to initiate charging.
- Allow the battery to charge for 8 hours.

4.3 WARRANTY

Warranty: two years parts and labor according to the following terms:

BMB MEDICAL grants its customers a commercial warranty of two (2) years from the date of delivery, within the limits defined below.

This commercial warranty is not intended to replace the legal guarantee of conformity mentioned in articles L.211-4 to L.211-13 of the French Consumer code, or that pertaining defects of the purchased object, covered by articles 1641 to 1648 and 2232 of the French Code of civil law.

Beyond the legal guarantee, BMB MEDICAL guarantees the exchange and replacement of parts accepted to be out of order, due to machining fault, assembly fault or material failure. The scope of the warranty does not include parts made unusable either by normal wear, or by improper use or non-compliant repair.

Where applicable, at the request of BMB MEDICAL, any product or part that is the subject of a warranty claim must be returned postage paid to the BMB MEDICAL plant. Any improper use, alteration or repair performed with parts not supplied or certified by BMB MEDICAL will void this warranty.

BMB MEDICAL products are designed for an expected service life of 10 years under normal conditions of use and with appropriate periodic maintenance as described above. BMB MEDICAL guarantees the original buyer that the welds on its products will be free from structural defects during the expected 10-year service life of the product as long as the original buyer is in possession of the product.

No guarantees are provided concerning appropriateness for a specific use.

This warranty excludes any disposable parts, IV poles, mattresses, batteries, along with damage incurred by abuse.

The product warranty shall be considered expired, however, by any of the following:

- Abnormal use ([Warning paragraph here](#))
- Work without prior notice to our technical services, or maintenance on:
 - A mechanical part of the product
 - A cylinder or motor
- Failure to observe any other term of the warranty
- Any inappropriate washing
- The use of accessories not provided by the manufacturer
- Failure to comply with the recommended preventive maintenance operations.

4.4 AFTER-SALES SERVICE

For any information concerning product installation, use and maintenance, the technical documentation (wiring diagrams, list of components, maintenance procedures, etc.), or any other questions, please contact the After-Sales Service at the following address: **sav@bmb-medical.com**, or by phone at **+33 (0)4 74 08 71 71**, or your local distributor for overseas purchases.

From your “Customer Area” on the website www.bmb-medical.com, you can access your equipment, submit work requests and track these requests directly via this secure interface.

The operator must carefully read the information contained in this user manual, in particular the information concerning the precautions relating to safety and installation/use methods.

In the event of doubt concerning the correct interpretation of the instructions, do not hesitate to contact the After-Sales Service whose staff are always available.

To facilitate ordering spare parts, you can request the exploded diagrams of your products from the After-Sales Service. You can thus obtain the references of the desired parts.

4.5 ENVIRONMENT

4.5.1 SPENT EQUIPMENT

4.5.1.1 TIM

To be disposed of separately from household waste.

Dispose of the product or its components in accordance with local and national laws and regulations on environmental protection and recycling of raw materials

You can obtain the appropriate disposal location from the authorities in charge of environmental protection.

France specificities: BMB MEDICAL fulfills its obligations concerning the end of life of the products in the TIM range, which it puts on the market by financing the recycling sector of Valdelia which takes them back free of charge (More information on www.valdelia.org).

4.5.1.2 VELVET & SWANN

In Europe:

To be disposed of separately from household waste.

In accordance with Directive 2012/19/EU of the European Parliament and Council of July 4, 2012 on waste electrical and electronic equipment (WEEE), this symbol indicates that the product must not be disposed of along with non-sorted municipal waste, but rather discarded to selective collection.

The recycling of electrical equipment serves to preserve natural resources and avoids any risks of pollution.

You can obtain the appropriate disposal location from the authorities in charge of environmental protection.

France specificities: BMB MEDICAL fulfills its obligations concerning the end of life of VELVET range products that it markets by funding the Recylum recycling channel dedicated to Pro WEEE that collects them free of charge (more information at www.recylum.com).

Outside of Europe:

Dispose of the product or its components in accordance with local laws and regulations.

Contact the local distributor for information concerning the return and/or collection systems available in the country.

Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans – CS 11001

01604 TREVoux Cedex

FRANCE