

CARUI

Patient trolley with height adjustment on columns
DEDICATED SOLUTIONS FOR YOUR PATIENTS' TRANSFER & WELLBEING





CARVI 60 - 70 - 80 Patient trolley with height adjustment on columns

DESCRIPTION

CARVI is a patient trolley with height adjustment on hydraulic columns.

Light, robust, and user-friendly, it is designed for patient transportation, emergency and PACU settings.

Various optional features are available to meet all your needs.

Overall width : 69 cm / 79 cm / 89 cm.

Overall length : 216 cm.

Overall height : 62 to 94 cm.

Weight : 90 kg / 94 kg / 98 kg/

2 year warranty

Technical data CARVI

2-sectional resin mattress base: 60 x 197 cm (CARVI 60), 70 x 197 cm (CARVI 70) or 80 x 197 cm (CARVI 80).

Radiolucent material.

Height adjustment range: from 62 to 94 cm (including 8cm thick mattress).

Trendelenburg / Anti-Trendelenburg 18°.

Two-sided height adjustment and Trendelenburg / Anti-Trendelenburg adjustment.

Removable mattress made of high quality and resilience foam (35kg/m³), height : 8 cm.

Removable cover. PVC coated Jersey. Colour : grey.

6 holders for IV pole.

Bumpers on all 4 corners.

Backrest adjustable by 2 gas springs (0 - 90°).

Castors with central braking system Ø 200 mm at foot-end.

2 directional wheels at the head-end (without 5th directional wheel).

Drop-down side rails, 2 horizontal bars. Anodized aluminium.

Foldable push bars at the foot-end (2 parts).

Frame cover including utility tray and O₂ tank holder.

Safe Working Load : 250 kg.

Trendelenburg/ Reverse
Trend hand-free control



Release button side rails



CARVI KIDS

Patient trolley

with height adjustment on columns

DESCRIPTION

CARVI is a patient trolley with height adjustment on hydraulic columns.

CARVI KIDS complies, with its specific side rails and paediatric kit at head-end, with the requirements from April, 26th 2010 of the French Medical Safety Agency ANSM with regard to children under 146 cm.

CARVI KIDS is designed to provide absolute safety while transporting children.

Light, robust, and user-friendly, it is designed for patient transportation, emergency and PACU settings. Various optional features are available to meet all your needs.

Overall width : 69 cm.

Overall length : 200 cm.

Overall height : 62 to 94 cm.

Weight : 85 kg.

2 year warranty

Technical data CARVI

2-sectional resin mattress base: 60 x 185 cm

Radiolucent material.

Height adjustment range: from 62 to 94 cm (including 8 cm thick mattress).

Trendelenburg / Anti-Trendelenburg 18°.

Two-sided height adjustment and Trendelenburg / Anti-Trendelenburg adjustment.

Removable mattress made of high quality and resilience foam (35kg/m³), height : 8 cm.

Removable cover. PVC coated Jersey. Colour : grey.

6 holders for IV pole.

Bumpers on all 4 corners.

Backrest adjustable by 2 gas springs (0 - 90°).

Castors with central braking system Ø 200 mm at foot-end.

2 directional wheels at the head-end (without 5th directional wheel).

Drop-down side rails, 2 horizontal bars. Anodized aluminium.

Foldable push bars at the foot-end (2 parts).

Frame cover including utility tray and O2 tank holder.

Safe Working Load : 250 kg.



Paediatric kit head-end (optional)

Polycarbonate paediatric side rails (optional)



OPTIONS: OPTIONAL FEATURES

Ref.	Description	CARVI	CARVI KIDS
CRO5E	5th directional wheel	•	•
CLORA	Anti-static system	•	•
CROB3B	Side rails 3 horizontal bars	•	•
MLOPL	Bed panel instead of push bars	•	•
CKOBP	Polycarbonate paediatric side rails		•
CKOKT	Paediatric kit head-end		•
TSE2C	Permanent, foldable 3-stage «comfort» IV pole, 2 hooks	•	
TSE4C	Permanent, foldable 3-stage «comfort» IV pole, 4 hooks	•	
TSR2	Removable IV pole, 2-stage, 2 hooks	•	•
CROBPT	Foldable push bars head-end	•	•
CROBPP1	Foldable push bar foot-end 1 part	•	•
CROK7	X-ray cassette below the patient surface (incompatible with CRORJ + CRODOC)	•	
CORBE	Electric backrest with manual release	•	
CORJE	Electric backrest and knee gatch, manual release	•	
COBBA	Battery	•	
COTLC	Patient's remote control	•	
CRORJ	Knee gatch	•	
CROSO2	Vertical oxygen bottle holder	•	•
CROSP	Paper roll dispenser at the head-end	•	•
CRODOC	Chart holder (for paperwork) below the headrest (incompatible with CROK7)	•	
PSU	Attachment hooks for bags on both sides below the patient surface	•	•



Ref TSE4C
Permanent IV pole, 4 hooks



Ref TSR2
Permanent IV pole, 2 hooks



Ref CROK7
X-Ray cassette below the patient surface



Ref CRO5E
5th directional wheel

OPTIONS (following) :


Réf.	Description	CARVI	CARVI KIDS
CROPSC	Monitor tray at the foot-end	●	
MC	Comfort mattress, high quality foam (35kg/m ³), height : 10cm	●	●
MCC	High Comfort mattress, high quality foam (35kg/m ³), 2cm thick viscoelastic layer (total height : 12cm)	●	
EQUI	Equipotential connection	●	●
CKOCL	Lateral protection cushion		●
CKOCT	Head protection cushion		●



Réf CROBPT
Foldable push bars head-end



Réf CROPSC
Monitor tray at foot-end



Ref CROSO2
Vertical oxygen bottle holder



Ref CRODOC
Chart holder (for paperwork) below the headrest

CARVI



100 Av. du Formans - CS 11001
01604 TREVOUX Cedex
FRANCE

Phone. : +33 (0)4 74 08 71 71
Fax : +33 (0)4 74 08 71 72
E-mail : info@bmb-medical.com
www.bmb-medical.com

© 2017 -  - Tél. 06 77 96 82 22 / www.krea-communication.fr

