

CLICK & MOVE

Motorized tractor for stretchers and chairs
DEDICATED SOLUTIONS FOR YOUR PATIENTS' TRANSFERT & WELLBEING





Click & Move
Motorized drive system
with autonomous wheel.

CLICK & MOVE

DESCRIPTION

CLICK & MOVE is a motorized traction unit with an autonomous wheel.

Perfect complement to a large park of stretcher, CLICK & MOVE requires a connexion on the stretcher (available on after sales).

It fits exclusively to BMB Medical products.

Overall width : 60 cm.

Overall length : 80 cm.

Overall height : 121 cm.

Weight : 50 kg.

2 year warranty

Technical data Click & Move

Motorized drive system with autonomous wheel.

Integrated battery. 15 km autonomy. Charge time : 5 hours.

Lithium battery without memory effect.

Integrated charger.

Forward/Backward twist grip with progressive speed acceleration.

Connection by pedal on the furniture.

Connection system requested on the furniture.

5th wheel obligatory.

Charge de fonctionnement en sécurité : 320 kg



Autonomous motorized wheel with integrated battery



Forward / backward twist grip with progressive speed acceleration.

Click & Move

OPTIONS: OPTIONAL FEATURES

Ref.	Description
BDM- BRA-CAR	Connection for Click & Move for CARVI
BDM- BRA-CLV	Connection for Click & Move for CLAVIA
BDM- BRA-ICA	Connection for Click & Move for I-CARE



Réf BDM-BRA-CAR/CLV/ICA
Connection for
Click & Move



Click & Move and I-Care



Click & Move and Clavia




Click & Move



100 Av. du Formans - CS 11001
01604 TREVoux Cedex
FRANCE

Phone : +33 (0)4 74 08 71 71
Fax : +33 (0)4 74 08 71 72
E-mail : info@bmb-medical.com
www.bmb-medical.com

© 2015 -  - Tél. 06 77 96 82 22 / www.krea-communication.fr

