

Accessories Range BMB Medical

DEDICATED SOLUTIONS FOR YOUR PATIENTS' TRANSFER & WELLBEING





Assistant
Height adjustment with gas spring
foldable Table-Top

Assistant
with integrated Drip Stand (optional)

DESCRIPTION

The ASSISTANT is a multipurpose product for your patients and your staff.

This is the perfect shelf in your one day care unit: use it for your patient effects, a snack or the small medical equipment. It is as little as a IV pole.

You can adapt it to every product with its star base.

It has a variable height and an optional removable 2-stage IV pole.

Overall width : from 57,5 to 63,3 cm

Overall length : from 57,5 to 99,7 cm

Overall height : 82 to 112 cm

Weight : 17Kg

Warranty : 2 years.

Technical data Assistant

2 sections folding table: 75 x 35 cm unfolded / 37,5 x 35 cm (folded)

Height adjustment range: from 82 to 112 cm.

No brake system.

5 rotary wheels Ø 40mm.

Storage space under the table: 27 x 25 x 10 cm.


Handbag / backpack holder (max 15kg).

IV pole holder.

Maximum load : 15 Kg with folded table / 5 Kg at the table end

OPTIONS:

Ref.	Description
TSR2	removable 2-stage IV pole



Storage space under the table

Réf TSR2
removable 2-stage IV pole



**Bedside cabinet / Locker
Included wardrobe**



**Bedside cabinet / Locker
Included shoes tray**

DESCRIPTION

The BEDSIDE CABINET / LOCKER is the perfect product to welcome your patient in a one day care unit.

It is equipped with a secured locker for the storage of your patients' personal effects.

Mobil, resistant and easy to clean, it joins perfectly a versatile care chair.

Overall width : 47,5 cm

Overall length : 59 cm

Overall height : 126 cm

Weight : 25 Kg

Warranty : 2 years.

Technical data Bedside cabinet / Locker

Stratified compact, acrylic resin finishing touch.

4 double running wheel castors Ø 100mm.

Removable shoes tray in white ABS.

Key lock system (2 keys included + optional universal key).

Push handle.

Drawer in ABS.

Wardrobe for clothes.

Identification system for patient / care chair / locker.

Maximum load : 15 kg.

Colors : White, black (optional), anthracite grey (optional)

OPTIONS:

Ref.	Description
VEST_C	Optional universal key
VEST_SD	Optional digital door lock



Identification system for patient / care chair / locker



Digital door and drawer lock.



Drawer for the storage of your patients' personal effects



Universal meal tray
on CLAVIA LSA SURGERY



Universal meal tray
on I-CARE trolley

DESCRIPTION

The universal meal tray is an accessory which can be fixed on the side rails.

Light and easy, it fits with all width of products

Overall width : 43,2 cm

Overall length : 72,5 cm

Overall height : 63 cm

Weight : 3,5 Kg

Warranty : 2 years.

Technical data universal meal tray

Extensible meal tray for side rails: from 65 to 88 cm.

Fall protection side.

Compatible with standardized tray meal :

- GastroNorm - GN1/1 : 53 x 32.5 cm
- EuroNorm : 53 x 37 cm

It fits with all width of products.

Maximum load : 10Kg.



Universal meal tray
on CLAVIA LSA
with laminated
side guards



Universal meal tray
fits with all width of products





100 Av. du Formans - CS 11001
01604 TREVOUX Cedex
FRANCE

Phone : +33 (0)4 74 08 71 71
Fax : +33 (0)4 74 08 71 72
E-mail : info@bmb-medical.com
www.bmb-medical.com

© 2017 -  - Tél. 06 77 96 82 22 / www.krea-communication.fr

