

# CAREXIA RANGE

## Versatile care chair

DEDICATED SOLUTIONS FOR YOUR PATIENTS' TRANSFER & WELLBEING





**Carexia FP**  
Versatile care chair  
synchronized legrest (optional)



**Carexia FP**  
Versatile care chair  
in trend position

Versatile care chair  
with fixed height

CAREXIA FP

## DESCRIPTION

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.

**Overall width :** 69,5 cm with armrests  
90 cm with gutters

**Seat width:** 57 cm.

**Overall length :** 76.4 cm in the sitting position at 170 cm lying down.

**Overall height :** 130 cm.

**Weight :** 65 to 85 kg depending on model.

**2 year warranty**

## Technical data Carexia FP

Heavy duty design. Steel frame with powder coating. ABS cover.

Fixed height : 55 cm.

AErgonomic shaped backrest and seat. PVC coated jersey.

High resilience foam 40 kg/m<sup>3</sup>.

Integrated headrest.

Adjustable back rest from 80° to 60° with 1 gaz spring

Trend position with optional legrest.

4 self lockable casters with double running wheels

ABS basket.

Polyurethan armrests, can be folded up.

Safe working load : 200 kg.



Carexia FP with comfort armrest.



**Carexia FPE**  
*Versatile care chair*  
*Electric backrest and legrest*



**Carexia FPE**  
*Versatile care chair*  
*with comfort armrest*

## DESCRIPTION

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.

**Overall width :** 69,5 cm with armrests  
 90 cm with gutters

**Seat width:** 57 cm.

**Overall length :** 76.4 cm in the sitting position at 177 cm lying down with long legrest.

**Overall height :** 130 cm.

**Weight :** 65 to 85 kg depending on model.

**2 year warranty**

## Technical data Carexia FPE

Heavy duty design. Steel frame with powder coating. ABS cover.

Fixed height : 55 cm.

Ergonomic shaped backrest and seat. PVC coated jersey.

High resilience foam 60 kg/m<sup>3</sup>.

Integrated headrest.

Electric backrest and legrest.

Trend position.

4 self lockable casters with double running wheels

ABS basket.

Polyurethan armrests, can be folded up.

Safe working load : 200 kg.



Carexia FPE in trend position.

Versatile care chair with fixed height  
 and electric backrest and leg rest





Carexia FPV  
Versatile care chair  
with electric variable height

CAREXIA FPV

## DESCRIPTION

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.

**Overall width :** 73,2 cm with armrests  
90 cm with gutters

**Seat width :** 57 cm.

**Overall length :** 76.4 cm in the sitting position at 170 cm lying down.

**Overall height :** 130 cm.

**Weight :** 65 to 85 kg depending on model.

**2 year warranty**

## Technical data Carexia FPV

Heavy duty design. Steel frame with powder coating. ABS cover.

Electric variable height from 54 to 78 cm.

Foot pedal control on both sides.

Ergonomic shaped backrest and seat. PVC coated jersey.

High resilience foam 60 kg/m<sup>3</sup>.

Integrated headrest.

Adjustable back rest from 80° to 60° with 1 gaz spring

Trend position with optional legrest.

4 self lockable casters with double running wheels

ABS basket.

Polyurethan armrests, can be folded up.

Safe working load : 200 kg.



Carexia FPV with legrest in trend position.



Foot pedal control.



**Carexia FPVE**  
*Versatile care chair with electric variable height,  
 with comfort armrest (optional)*



**Carexia FPVE**  
*Versatile care chair  
 with long legrest and footrest (optional)*

Versatile care chair  
 with electric movement

CAREXIA FPVE

## DESCRIPTION

CAREXIA is a versatile treatment chair. Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient. Pure and harmonious design is the vector of a positive image for your unit. CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.

**Overall width :** 73,2 cm with armrests  
 90 cm with gutters

**Seat width :** 57 cm.

**Overall length :** 76.4 cm in the sitting position at 177 cm lying down with long legrest.

**Overall height :** 130 cm.

**Weight :** 65 to 85 kg depending on model.

**2 year warranty**

## Technical data Carexia FPVE

Heavy duty design. Steel frame with powder coating. ABS cover.

Electric variable height from 54 to 78 cm.

Commande par pédale bilatérale.

Ergonomic shaped backrest and seat. PVC coated jersey.

High resilience foam 60 kg/m3.

Integrated headrest.

Electric backrest and legrest.

Trend position.

4 self lockable casters with double running wheels

ABS basket.

Polyurethan armrests, can be folded up.

Safe working load : 200 kg.



Carexia FPVE in trend position.



Foot pedal control (optional).

# CAREXIA RANGE

## OPTIONS: OPTIONAL FEATURES

Ref.	Description	CAREXIA FP	CAREXIA FPE	CAREXIA FPU	CAREXIA FPU€
FPORS1	4 castors Ø 125 mm with central lock	●	●	●	●
FPORS2	4 double runing wheel castors Ø 150 mm with central lock	●	●	●	●
FPOBP	Push bar	●	●	●	●
FPORJ	Short legrest synchronized with the backrest. Trend position.	●		●	
FPORJL	Long legrest		●		●
FPORP	Electric adjustable footrest (only for long legrest)		●		●
FPOCP1	Foot pedal control for relax and trend position		●		●
FPOCP2	Foot pedal control for variable height on both sides				●
FPOSA	2 holders for IV pole (necessary for TSR2)	●	●	●	●
TSR2	Removable IV pole, 2-stage, 2 hooks	●	●	●	●
FPOAC4	Armrest comfort (pair)	●	●	●	●
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)	●	●	●	●
FPOSR	Paper roll holder	●	●	●	●
FPOLAMP	Reading lamp DEL		●		●
FPOBA	Battery 2 Ah		●	●	●
FPOGC	High comfort upholstery with 3 cm thickness memory foam	●	●	●	●

## OPTIONS: OPTIONAL FEATURES

Ref.	Description	CAREXIA FP	CAREXIA FPE	CAREXIA FPU	CAREXIA FPUE
ACC-CTG	Comfort pillow (grey)	•	•	•	•
SUP-USB	USB Port		•	•	•
PROTEC RJL	Removable Crystal protection for footrest		•		•
PROTEC RP	Removable Crystal protection for footrest		•		•
SUPTAB	Fixed I-Pad support with antitheft device	•	•	•	•



Ref SUP-USB  
USB Port



Ref FPOSA  
2 holders for IV pole  
(necessary for TSR2)



Ref FPORS2  
4 double runing wheel  
castors Ø 150 mm



Ref FPORP  
Electric adjustable footrest  
(only for long legrest)



Ref FPOAC3  
Anatomic shaped armrest  
adjustable in height  
and in all directions.



Ref FPOSR  
Paper roll holder



Ref FPOLAMP  
Reading lamp DEL.



Ref TSRZ  
Removable IV pole, 2-stage,  
2 hooks

CAREXIA RANGE

