

# I-CARE Patient trolley

DEDICATED SOLUTIONS FOR YOUR PATIENTS' TRANSFER & WELLBEING





**I-CARE Emergency 60 / Emergency 70**  
Patient trolley  
With hydraulic height adjustment



**I-CARE Ambulatory 60 / Ambulatory 70 / Ambulatory 80**  
Patient trolley  
With hydraulic height adjustment

I-CARE

DESCRIPTION

I-CARE is a patient trolley with height adjustment on hydraulic columns.

Stable and robust, it is ideally designed to an intensive use.

Its features and its wide range of accessories make it especially suitable for Emergency, OR, PACU or ASC settings.

**Overall width :** Emergency 74 / 84 cm.  
Ambulatory 69 / 79 / 89 cm.

**Overall length :** 216 cm.

**Overall height :** 64 to 96 cm.

**Weight :** 85 kg / 90 kg

**2 year warranty.**

Technical data I-CARE

2-sectional resin mattress base : 60 x 197 cm / 70 x 197 cm / 80 x 197 cm (according to the model).

Radiolucent material.

Two-sided + foot-end (or head-end, free choice) height adjustment and Trendelenburg / Anti-Trendelenburg adjustment.

Height adjustment range: from 64 to 96 cm (including 10 cm thick mattress).

Hand-free Trendelenburg / Anti-Trendelenburg (18°)

Removable mattress made of high quality and resilience foam (35kg/m3), thickness : 10 cm.

Removable cover. PVC coated Jersey. Colour : grey.

6 holders for IV pole.

Bumpers on all 4 corners.

Backrest adjustable by 2 gas springs (0 - 90°).

Castors with central braking system Ø 200 mm at foot- and head-end.

5th retractable and directional wheel.

Drop-down side rails. Anodized aluminium.

Frame cover (removable for easy cleaning and maintenance) including utility tray and O2 tank holder.

Central braking at foot- and head-end.

Foldable push handles at foot- end (in 2 parts)

Max. load : 320 kg.



Detail view foldable push handles





## I-CARE OPHTHALMOLOGY

### Patient trolley

With hydraulic height adjustment

## DESCRIPTION

I-CARE is a patient trolley with height adjustment on hydraulic columns.

With its removable surgery headrests, the product is adapted to one day surgery.

I-CARE OPHT is for ophthalmic surgery, Stomatology, ENT, hand surgery...

**Overall width :** 74 cm.

**Overall length :** 200 cm (without headrest)  
235 cm (with headrest) .

**Overall height :** 64 to 96 cm (including mattress)

**Weight :** 90 kg  
**2 year warranty.**

## Technical data I-CARE

2-sectional resin mattress base : 70 x 197 cm

Radiolucent material.

Two-sided + foot-end height adjustment and Trendelenburg / Anti-Trendelenburg adjustment.

Height adjustment range: from 64 to 96 cm (including 8 cm thick mattress).

Hand-free Trendelenburg / Anti-Trendelenburg (18°)

Removable mattress made of high quality and resilience foam (35kg/m<sup>3</sup>), thickness : 8 cm.

Removable cover. PVC coated Jersey. Colour : grey.

Backrest adjustable by 2 gas springs (0 - 90°).

Castors with central braking system Ø 200 mm at foot- and head-end.

5th retractable and directional wheel.

Drop-down side rails. Anodized aluminium.

Frame cover (removable for easy cleaning and maintenance) including utility tray and O2 tank holder.

Foldable push handles at foot- end (in 2 parts)

Max. load : 200 kg.



Ophthalmologic headrest (optional)

Stomatologic and ENT headrest (optional)





## OPTIONS: OPTIONAL FEATURES

Réf.	Description	I-CARE AMBU	I-CARE EMER	I-CARE OPHT	I-CARE AMBU EZ-GO	I-CARE EMER EZ-GO
CLORA	Anti-static system	●	●	●	INCLUDED	INCLUDED
TSE2C	Permanent, foldable 3-stage «comfort» IV pole, 2 hooks	●	●		●	●
TSE4C	Permanent, foldable 3-stage «comfort» IV pole, 4 hooks	●	●		●	●
TSR2	Removable IV pole, 2-stage, 2 hooks	●	●	●	●	●
CROBPT	Foldable push bars head-end	●	●		●	●
CROBPP1	Foldable push bar foot-end 1 part	●	●	●		
CROK7	X-ray cassette below the patient surface (incompatible with CRORJ + CRODOC)	●	●			
CORBE	Electric backrest with manual release	●	●		●	●
CORJE	Electric backrest and knee gatch, manual release	●	●		●	●
COBBA	Battery	●	●		●	●
COTLC	Patient's remote control	●	●		●	●
CRORJ	Knee gatch	●	●		●	●
CRORT	Side technical rails		●	INCLUDED		●
ETU	Fixation clamp for armrest (on technical rails)		●	●		●
CROSO2	Vertical oxygen bottle holder	●	●		●	●
CROSP	Paper roll dispenser at the head-end	●	●		●	●
CRODOC	Chart holder (for paperwork) below the headrest (incompatible with CROK7)	●	●		●	●
PSU	Attachment hooks for bags on both sides below the patient surface	●	●		●	●
CROPSC	Monitor tray at the foot-end	●	●		●	●
MMM	High Comfort mattress, high quality foam (35kg/m <sup>3</sup> ), 2cm thick viscoelastic layer (total height : 12cm)	●	●		●	●
EQUI	Equipotential connection	●	●	●	●	●



## OPTIONS (following) : \_\_\_\_\_


Réf.	Description	I-CARE AMBU	I-CARE EMER	I-CARE OPHT	I-CARE AMBU EZ-GO	I-CARE EMER EZ-GO
LCA-OPH	Removable surgery head rest for ophthalmology			●		
LCA-STO	Removable surgery head rest for stomatology and ENT			●		
LCA-SUP	Head rest support. Fixation on technical rails.			●		
FSC	Holder for operating field, with connection for oxygen (ETU necessary)			●		

I-CARE



100 Av. du Formans - CS 11001  
01604 TREVoux Cedex  
FRANCE

Phone : +33 (0)4 74 08 71 71  
Fax : +33 (0)4 74 08 71 72  
E-mail : [info@bmb-medical.com](mailto:info@bmb-medical.com)  
[www.bmb-medical.com](http://www.bmb-medical.com)

© 2017 -  - Tél. 06 77 96 82 22 / [www.krea-communication.fr](http://www.krea-communication.fr)

