

I-MOVE

Patient transport chair

DEDICATED SOLUTIONS FOR YOUR PATIENTS' TRANSFER & WELLBEING





Patient transport chair I-MOVE
For the patients' reception



Stackable patient transport chair I-MOVE
for optimized tidying

I-MOVE

DESCRIPTION

I-MOVE is a transport chair for your patients' reception and sitting transport.

Manageable, robust and attractive, it is now stackable. It is the modern and colored object for your department.

With our big range of accessories, you can adapt your use.

Overall width: 75 cm.

Overall length: 105 cm.

Overall height: 100 cm.

Weight: 35 kg.

2 years warranty.

Technical data I-Move

Steel frame epoxy coated.

Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person.

Front swivel wheels Ø 125 mm.

Cushion in foam width 5 cm.

Seamless finish, vinyl coating.

Color range: fuchsia, green, blue.

Dimensions seat and back: 46 x 46 cm

Seat height : 55 cm.


Sliding and self-blocking foot rest.

Polyurethan armrests, can be folded up.

Safe Working load : 200 kg.



Two positions leg rest (optional)



Different colors



Sliding and self-blocking foot rest





**Bariatric patient transport chair I-MOVE
for the patients' reception**

I-MOVE

DESCRIPTION

I-MOVE is a bariatric transport chair for your patients' reception and sitting transport.

Manageable, robust and attractive, it is the modern and colored object for your department.

With our big range of accessories, you can adapt your use.

Overall width: 90 cm.

Overall length: 105 cm.

Overall height: 102 cm.

Weight: 50 kg.

2 years warranty.

Technical data I-Move

Steel frame epoxy coated.

Twin back wheels Ø 300 mm with central braking system, can be operated in stand position.

Front swivel wheels Ø 125 mm.

Cushion in foam width 8 cm.

PVC coated Jersey.

Color range: lavender, green, squirrel

Dimensions seat and back: 46 x 65 cm

Seat height : 55 cm.

Sliding and self-blocking foot rest.

Polyurethan armrests, can be folded up.

Safe Working load : 250 kg.



Sliding and self-blocking foot rest



Different colors



Zoom on the twin back wheel Ø 300 mm with central braking system





Motorized patient transport chair
I-MOVE EZ-GO

I-MOVE

DESCRIPTION

I-MOVE is a transport chair for your patients' reception and sitting transport.

I-MOVE EZ-GO is equipped with an autonomous motorized wheel (including its own battery). Thanks to I-MOVE EZ-GO, patients transfer is child's play!

With our big range of accessories, you can adapt your use.

Overall width: 70 cm or 90 cm.

Overall length: 105 cm.

Overall height: 102 cm.

Weight: 50 kg or 65 kg.

2 years warranty.

Technical data I-Move

Steel frame epoxy coated.

Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person.

Front swivel wheels Ø 125 mm.

Cushion in foam width 5 cm.

Seamless finish, vinyl coating.

Color range: fuchsia, green, blue (standard chair) / lavender, green, squirrel (bariatric chair).

Dimensions seat and back: 46 x 46 cm or 46 x 65 cm.

Seat height : 55 cm.

Sliding and self-blocking foot rest.

Polyurethan armrests, can be folded up.

Antistatic system.

Accessory basket.

Motorized drive system with autonomous wheel.

Forward/Backward twist grip with progressive speed acceleration.

Safe Working load : 200 kg.



Motorized wheel including battery.



OPTIONS: OPTIONAL FEATURES

Ref.	Description	I-MOVE	I-MOVE BARIATRIC	I-MOVE EZ-GO
TSR2	Removable IV pole, 2-stage, 2 hooks	●	●	●
CLORA	Anti-Static system	●	●	INCLUDED
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with TSR2)	●	●	●
TROPO	Accessory basket	●	●	INCLUDED
TROSO2	Holder for oxygen cylinder	●	●	●
TROPE	Equipment for leg rest (pair)	●	●	●
TRORJ	Adjustable leg rest (1 piece)	●	●	●
TROMO	Coin acceptor	●	●	●
TROMOP	Wall fixation	●	●	●
PSU	Hooks to hold bag	●	●	●



Ref TROMO
Coin acceptor



Ref TSR2
IV pole, 2 hooks



Ref TRORJ
Adjustable leg rest




Ref PSU
Hooks to hold bag

I-MOVE



100 Av. du Formans - CS 11001
01604 TREVoux Cedex
FRANCE

Tél : +33 (0)4 74 08 71 71
Fax : +33 (0)4 74 08 71 72
E-mail : info@bmb-medical.com
www.bmb-medical.com

© 2017 -  - Tél. 06 77 96 82 22 / www.krea-communication.fr

