

TRANSFER CHAIR I-MOVE





CONTENTS

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position.

Manageable, robust and attractive it will bring for a long time a touch of modernity and color to your department.

A wide range of accessories allows it to adapt to each use.



I - M O V E



I - M O V E B A R I A T R I C



I - M O V E E Z - G O



I - M O V E E Z - G O B A R I A T R I C



Sturdy

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Easy to handle & Stackable for optimal storage



TRANSFER CHAIR I-MOVE

I-MOVE is a transfer chair, it allows the reception of the patient and his transport in sitting position.





I - M O V E

I-MOVE DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	75 cm	100 cm	H. 55 cm D. 46 cm W. 46 cm	35 Kg

Safe Working Load : 200 Kg



Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

I-MOVE transfer chair for patient welcoming

- Steel frame epoxy coated
- Patient file holder
- Back wheels Ø 300 mm with central braking system
- Can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Cushion in foam width 5 cm
- Seamless finish, vinyl coating
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up

Back wheels Ø 300 mm with central braking system



Front swivel wheels Ø 125 mm






Polyurethan armrests, can be folded up



Sliding foot rest



Different colors




Composition : PVC 83 % / Coton 17 %
Thickness: 1,2 mm +/- 0,1 mm
Light resistance: >= 6



FIRE RATING
M1 / AM 18 / EN 1021.1-2-3.

ACCESSORIES

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR2-F	Fixed IV pole, 2-stage, 2 hooks
CLORA	Anti-Static system
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with TSR2)
TROPO	Accessory basket
TROSO2	Holder for oxygen cylinder
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO -EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
PSU	Hooks to hold bag



Removable IV pole,
2-stage, 2 hooks
TSR2



Accessory basket
TROPO



Holder for oxygen
cylinder
TROSO2



Adjustable leg rest
(1 piece)
TRORJ



Coin acceptor
TROMO



Hooks to hold bag
PSU

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit



I-MOVE BARIATRIC TRANSFER CHAIR

I-MOVE is a bariatric transport chair for your patients' reception and sitting transport.




I-MOVE BARIATRIC

I-MOVE BARIATRIC DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	90 cm	102 cm	H. 55 cm D. 46 cm W. 46 cm	50 Kg

Safe Working Load : 250 Kg



Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position.


I-MOVE transfer chair for patient welcoming

- Steel frame epoxy coated
- Patient file holder
- Back wheels Ø 300 mm with central braking system
- Can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Cushion in foam thickness 8 cm.
- PVC coated Jersey
- Sliding and self-blocking foot rest.
- Polyurethan armrests, can be folded up.

Back wheels Ø 300 mm with central braking system



Polyurethan armrests, can be folded up.



Front swivel wheels Ø 125 mm



Sliding foot rest.



Differents colors



ORANGE



GREEN



LAVANDER

Composition : PVC 83 % / Coton 17 %
Thickness: 1,2 mm +/- 0,1 mm
Light resistance: >= 6



FIRE RATING
M1 / AM 18 / EN 1021.1-2-3.

I-MOVE BARIATRIC

ACCESSORIES

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR2-F	Fixed IV pole, 2-stage, 2 hooks
CLORA	Anti-Static system
TROPO	Accessory basket
TROSO2	Holder for oxygen cylinder
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO -EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
PSU	Hooks to hold bag



IV pole, 2-stage,
2 hooks
TSR2



Accessory basket
TROPO




Holder for oxygen
cylinder
TROSO2



Adjustable leg rest
(1 piece)
TRORJ



Coin acceptor
TROMO



Hooks to hold bag
PSU

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Motorized Very easy to use, progressive speed



TRANSFER CHAIR I-MOVE EZ-GO

I-MOVE EZ-GO is a transport chair for your patients' reception and sitting transport. It is equipped with an autonomous motorized wheel (including its own battery). Thanks to I-MOVE EZ-GO, patients transfer is child's play!




I - MOVE EZ - GO

I-MOVE EZ-GO DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	70 cm	102 cm	H. 55 cm D. 46 cm W. 46 cm	50 Kg

Safe Working Load : 200 Kg




Motorized wheel including battery

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

Motorized transfer chair I-MOVE EZ-GO for patient welcoming

- Steel frame epoxy coated.
- Patient file holder
- Back wheels Ø 300 mm with central braking system
- Can be operated in stand position by the accompanying person.
- Front swivel wheels Ø 125 mm.
- Cushion in foam thickness 5 cm.
- Seamless finish, vinyl coating
- Sliding and self-blocking foot rest.
- Polyurethan armrests, can be folded up.
- Accessory basket.
- Antistatic system
- Motorized drive system with autonomous wheel.
- Forward/Backward twist grip with progressive speed acceleration

Polyurethan armrests, can be folded up



Back wheels Ø 300 mm with central braking system



Sliding footrest




Front swivel wheels Ø 125 mm



Accessory basket



Differents colors



FUSCHIA



GREEN



BLUE



Composition : PVC 83 % / Coton 17 %
Thickness: 1,2 mm +/- 0,1 mm
Light resistance: >= 6



FIRE RATING
M1 / AM 18 / EN 1021.1-2-3.

I-MOVE EZ-GO



ACCESSORIES


REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with TSR2)
TROSO2	Holder for oxygen cylinder
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROBA	Additional battery
PSU	Hooks to hold bag



Removable IV pole,
2-stage, 2 hooks
TSR2



Hooks to hold bag
PSU



Holder for oxygen
cylinder
TROSO2



Adjustable leg rest
(1 piece)
TRORJ

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Motorized Very easy to use, progressive speed



TRANSFER CHAIR I-MOVE EZ-GO BARIATRIC

I-MOVE EZ-GO is a transport chair for your patients' reception and sitting transport. It is equipped with an autonomous motorized wheel (including its own battery).

Thanks to I-MOVE EZ-GO, patients transfer is child's play! It is particularly suitable for the transfer of overweight patients.




I-MOVE EZ-GO BARIATRIC

I-MOVE EZ-GO BARIATRIC DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	90 cm	102 cm	H. 55 cm D. 46 cm W. 46 cm	65 Kg

Safe Working Load: 250 Kg




Motorized wheel including battery

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

Motorized transfer chair I-MOVE EZ-GO for patient welcoming

- Steel frame epoxy coated.
- Patient file holder
- Back wheels Ø 300 mm with central braking system
- Can be operated in stand position by the accompanying person.
- Front swivel wheels Ø 125 mm.
- Cushion in foam thickness 8 cm.
- Seamless finish, vinyl coating
- Sliding and self-blocking foot rest.
- Polyurethan armrests, can be folded up.
- Accessory basket.
- Antistatic system
- Motorized drive system with autonomous wheel.
- Forward/Backward twist grip with progressive speed acceleration

Polyurethan armrests, can be folded up



Back wheels Ø 300 mm with central braking system



Sliding footrest



Front swivel wheels Ø 125 mm



Accessory basket



Différents coloris



ORANGE



GREEN



LAVANDER

Composition : PVC 83 % / Coton 17 %
Thickness: 1,2 mm +/- 0,1 mm
Light resistance: >= 6



FIRE RATING
M1 / AM 18 / EN 1021.1-2-3.

I-MOVE EZ-GO BARIATRIC



ACCESSORIES


REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TROSO2	Holder for oxygen cylinder
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROBA	Additional battery
PSU	Hooks to hold bag



Removable IV pole,
2-stage, 2 hooks
TSR2



Hooks to hold bag
PSU



Holder for oxygen
cylinder
TROSO2



Adjustable leg rest
(1 piece)
TRORJ



Tél : +33 (0)4 74 08 71 71
Fax : +33 (0)4 74 08 71 72
info@bmb-medical.com

100 Av. du Formans - CS 11001
01604 TREVOUX Cedex
FRANCE

www.bmb-medical.com