

ACCESSORIES
MOBILES SUPPORTS
HOSPITAL SCREENS





CONTENTS

Our "MOBILE SUPPORTS & HOSPITAL SCREENS" range is suitable for all types of services in hospitals.

Combining comfort, ergonomics and style, our products are designed to accomodate and meet the needs of companions, parents, and caregivers.

- ▶ MOBILE SUPPORTS
 - ▶ TABOLYS
 - ▶ MATERNA
 - ▶ ASSISTANT
 - ▶ INFORMALYS
 - ▶ LOCKER
 - ▶ LARGE MODEL
 - ▶ SMALL MODEL
 - ▶ MODULYS
- ▶ HOSPITAL SCREENS
- ▶ ACCESSORIES

TABOLYS

Height-adjustable mobile table for patients and hospital staff. It accommodates patient effects, a snack or small medical equipment. Its H-shaped base is designed to fit under all types of beds, furniture and armchairs.

TABOLYS HPL

- HPL tray (white compact)
- «H» base
- 2 back bearing castors Ø 50 mm
- 2 front castors Ø 30 mm
- Paint Blue Grey Anod B711



TABOLYS ABS

- ABS top (Thermoformed grey)
- «H» base
- 2 back bearing castors Ø 50 mm
- 2 front castors Ø 30 mm
- Paint Blue Grey Anod B711



TABOLYS HPL DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	82.5 cm	53 cm	67 to 107 cm
BASE	64 cm	53 cm	4 cm
TRAY	82.5 cm	46.5 cm	0.8 cm

TABOLYS ABS DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	81 cm	53 cm	67 to 107 cm
BASE	64 cm	53 cm	4 cm
TRAY	81 cm	41 cm	4 cm

SAFE WORKING LOAD : 15 Kg

ACCESSORIES

REFERENCE	DESIGNATION
TAHPL	Tilting tray for TABOLYS HPL (compact)
TAABS	Tilting tray for TABOLYS ABS (Thermoformed)



Tilting tray
TAHPL



Tilting tray TAABS

Variable height cradle. The Materna is the ideal cradle to accompany parents during their stay in maternity/ neonatology. *Its base is designed to fit under all types of healthcare beds, allowing the mother to have her baby as close as possible to her, while making the baby safe.*

- Pneumatic variable height from 88 cm to 109 cm
- Trendelenburg 10 °
- Structure made of steel epoxy coated (Oxide ou gray)
- Removable cradle made of transparent plexiglass
- Drilled bottom for better ventilation
- M1 fire-resistant mattress 3 cm
- 4 self lockable double running castors (Ø 75 mm)

M A T E R N A DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	98 cm	61 cm	88 to 109 cm

SAFE WORKING LOAD : 10 Kg



REFERENCE	DESCRIPTION
MATERNAU	"U" base



"U" base




ASSISTANT

ACCESSORIES

Multifunctional day care assistant for both patients and hospital staff. *For patient effects, a snack or small medical equipment, it is an ideal support for your day surgery services.* Compact, it takes up no more space than an IV pole. Its star-shaped base fits perfectly with all our products.

- 2 sections folding table : 75 x 35 cm unfolded / 37,5 x 35 cm folded). Weight : 17kg
- Height adjustment range : from 82 to 112 cm
- Unbraked product
- 4 rotary wheels Ø 40mm
- Storage space under the table : 27 x 25 x 10 cm.
- Handbag / backpack holder (maximum 15kg)
- IV Pole holder

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks



Removable IV pole, 2-stage,
2 hooks
TSR2

ASSISTANT DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	57,5 to 99,7 cm	57,5 to 63,3 cm	82 to 112 cm
BASE			7.5 cm

Maximum load : 15 Kg folded / 5Kg at the table end



Storage space



INFORMALYS

Mobile computer table with variable height. Informalys is the ideal product to accompany the medical staff in their units. **Designed for laptops, it has a handle to facilitate its movement, as well as a mouse pad.** Mobile, resistant and easy to clean, it is the perfect companion for your computer equipment.

- Computer table Trespa matt white tray
- For laptops
- Flexibility - All input-output ports on the laptop are accessible
- Front handle
- 3 self lockable pivoting casters Ø 80 mm
- Retractable mouse pad

INFORMALYS DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	60 cm	60 cm	74 to 113 cm
TRAY	57 cm	54 cm	2.4 cm
BASE	60 cm	60 cm	16 cm

SAFE WORKING LOAD : 15 Kg



LOCKER LARGE MODEL

ACCESSORIES

The bedside cabinet locker is the ideal product for the reception of the patient in an day surgery unit. Intended for the storage of your patients' personal belongings, the product is equipped with a secure locking system. Mobile, resistant and easy to clean, it is the perfect companion to a care chair.

- Stratified compact, acrylic resin finishing touch, white color
- 4 double running wheel castors Ø 100mm
- Removable shoes tray in white ABS
- Key lock system (2 keys included)
- Push handle - Drawer in ABS - Wardrobe for clothes
- Identification system for patient / care chair / locker
- Anti-fall edge
- Door opening on the right (Left on request)
- Weight : 25 kg

REFERENCE	DESCRIPTION
VEST_C	Universal key for key lock system
VEST_SD	Digital door lock with code (1 universal key included)
VEST_CSD	Additional universal key for digital lock

LOCKER DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	59 cm	47.5 cm	126 cm

SAFE WORKING LOAD : 15 Kg



Universal key for key lock system
VEST_C



Digital door lock with code
VEST_SD



Additional universal key for digital lock
VEST_CSD



LOCKER SMALL MODEL

ACCESSORIES

The bedside cabinet locker is the ideal product for the reception of the patient in an day surgery unit. **Intended for the storage of your patients' personal belongings, the product is equipped with a secure locking system.** Mobile, resistant and easy to clean, it is the perfect companion to a care chair.

- Stratified compact, acrylic resin finishing touch, white color
- 4 double running wheel castors Ø 80mm
- Removable shoes tray in white ABS
- Key lock system (2 keys included)
- Push handle - Drawer in ABS - Wardrobe for clothes
- Anti-fall edge
- Door opening on the right (Left on request)
- Weight : 21 kg

LOCKER DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	59 cm	47.5 cm	92 cm

SAFE WORKING LOAD : 15 Kg



REFERENCE	DESCRIPTION
VEST_C	Universal key for key lock system
VEST_SD	Digital door lock with code (1 universal key included)
VEST_CSD	Additional universal key for digital lock
VEST_PB	Waste basket for 10 L bag (only available for VEST_P)




Waste basket for 10 L bag
VEST_C



Digital door lock with code
VEST_SD



Additional universal key for
digital lock
VEST_CSD



Waste basket for 10 L bag
VEST_PB

Locker large model
VEST



Locker small model
VEST_P



System allowing the installation of a cyclor and its consumables for daily or home peritoneal dialysis.

Allows the installation of 4 types of known devices : Sleep-Safe Harmony - HomeChoice Claria - S3 - NXstage

- Central part in Trespa and steel materials
- Top part in Trespa materials
- Structure made of steel. Epoxy coated.
- 3 casters Ø 75mm
- Manual variable height
- Top plate which can support until 80kg, room for a 60x40cm machine
- Easy adjustable height, even the machine on it, from 48cm to 78cm
- Weight : 15 KG - Can be easily taken to pieces

REFERENCE	DESCRIPTION
MODULSUP	Attachment hooks for bags
MODULTUBS	Tubes for the top part
MODU-DESS	Removable additional table for consumables. Retractable legs. 520 x 800 x 390 mm Weight : 5 kg
MODUSAC	Transport bag

MODULYS COMPOSITION

MODULCO


Top part :
600 x 410 mm



+ column
+ down part:
H 620 x 550 mm
weight : 15 kg
Down height: 490mm
Up height: 780 mm

MODU-DESS

520 x 800 x 390 mm
Weight : 5kg



MODULTUBS

Tubes for the top part




310 x 555 mm
Weight: 1,2 kg

MODULSUP


Attachment hooks for bags




310 x 65 x 40 mm
Weight : 400 g



Attachment hooks for bags
MODULSUP



Tubes for the top part
MODULTUBS



MODU-DESS





Versatile, adapted to different hospital units

Adaptable, many options available

Two models Wall-mounted and two mobile bases

Easy to use, thanks to its mobility



HOSPITAL SCREENS

ACCESSORIES

Rigid panel screens.

- M1 white PVC panels
- From 5 to 13 panels
- Wall mounting (including 10 cm spacer / on request 20 or 30 cm) or two mobile bases
- Finishing profiles available in different colours: Blue, Grey and Orange
- Strap for closing the screen

REFERENCE	DESCRIPTION
TENENT1	Strut of 10 cm (included with wall mounting)
TENENT2	Strut of 20 cm (On request with wall mounting)
TENENT3	Strut of 30 cm (On request with wall mounting)

SCREENS DIMENSIONS

HEIGHT	LENGTH (meter)	NUMBER OF PANELS								
		5	6	7	8	9	10	11	12	13
1,65 M		1,365	1,575	1,785	1,995	2,205	2,415	2,625	2,835	3,045



Strut of 10, 20 or 30 cm
TENENT1 2 3
(included with wall mounting)

FINISHING PROFILES BLUE - GREY - ORANGE



ACCESSORIES



Lightweight Arm & Handsurgery table in Carbon Fibre with twin clamps (hook on)
One (1) glassfiber supporting leg included for left or right, mounting depending of location of surgeon or C-arm/microscope.
Cushion included.
Length: 815 mm. Weight: 5 kgs.
Metalfree space.

CHIRMAIN



Safety belt for stretchers, chairs and armchairs

ACC-CTP



Protection tunel from rain for CARVI and I-CARE

PROPLUIE



Roll board / patient transfert device
Size: 177 X 55cm ($\pm 3\%$) (unfolded) / 88 x 55cm ($\pm 3\%$) (folded)
Weight : 2.4 kg ($\pm 3\%$)
100% Nylon
Aluminum Handle
SAFE WORKING LOAD : 200 kg

SWITCHEASY



Tél : +33 (0)4 74 08 71 71
Fax : +33 (0)4 74 08 71 72
info@bmb-medical.com

100 Av. du Formans - CS 11001
01604 TREVOUX Cedex
FRANCE

www.bmb-medical.com