

ACCESSORIES  
MOBILES SUPPORTS  
HOSPITAL SCREENS





## CONTENTS

Our "MOBILE SUPPORTS & HOSPITAL SCREENS" range is suitable for all types of services in hospitals.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

- ▶ MOBILE SUPPORTS
  - ▶ TABOLYS
  - ▶ MATERNA
  - ▶ ASSISTANT
  - ▶ INFORMALYS
  - ▶ LOCKER
    - ▶ LARGE MODEL
    - ▶ SMALL MODEL
  - ▶ MODULYS
- ▶ HOSPITAL SCREENS
- ▶ ACCESSORIES



# TABOLYS

Height-adjustable mobile table for patients and hospital staff. It accommodates patient effects, a snack or small medical equipment. Its H-shaped base is designed to fit under all types of beds, furniture and armchairs.

## TABOLYS HPL

- HPL tray (white compact)
- «H» base
- 2 back bearing castors Ø 50 mm
- 2 front castors Ø 30 mm
- Paint Blue Grey Anod B711



## TABOLYS ABS

- ABS top (Thermoformed grey)
- «H» base
- 2 back bearing castors Ø 50 mm
- 2 front castors Ø 30 mm
- Paint Blue Grey Anod B711



## TABOLYS HPL DIMENSIONS

	LENGHT	WIDTH	HEIGHT
<b>OVERALL</b>	82.5 cm	53 cm	67 to 107 cm
<b>BASE</b>	64 cm	53 cm	4 cm
<b>TRAY</b>	82.5 cm	46.5 cm	0.8 cm

## TABOLYS ABS DIMENSIONS

	LENGHT	WIDTH	HEIGHT
<b>OVERALL</b>	81 cm	53 cm	67 to 107 cm
<b>BASE</b>	64 cm	53 cm	4 cm
<b>TRAY</b>	81 cm	41 cm	4 cm

**SAFE WORKING LOAD : 15 Kg**

## ACCESSORIES

REFERENCE	DESIGNATION
<b>TAHPL</b>	Tilting tray for TABOLYS HPL (compact)
<b>TAABS</b>	Tilting tray for TABOLYS ABS (Thermoformed)



Tilting tray  
TAHPL



Tilting tray TAABS


Variable height cradle. The Materna is the ideal cradle to accompany parents during their stay in maternity/ neonatology. *Its base is designed to fit under all types of healthcare beds, allowing the mother to have her baby as close as possible to her, while making the baby safe.*

- Pneumatic variable height from 70cm to 100cm
- Trendelenburg 10 °
- Structure made of steel epoxy coated blue grey
- Removable cradle made of transparent plexiglass: 75 cm x 35 cm x 29 cm
- Drilled bottom for better ventilation
- M1 fire-resistant mattress 5 cm, removable cover
- 4 self lockable double running castors (Ø 80 mm)

## M A T E R N A DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	105 cm	58 cm	70 to 100 cm

**SAFE WORKING LOAD : 9Kg**



REFERENCE	DESCRIPTION
<b>MATDOC</b>	Patient file holder
<b>MATPA</b>	Storage basket



Storage basket MATPA



Patient file holder

# ASSISTANT

# ACCESSORIES

Multifunctional day care assistant for both patients and hospital staff. *For patient effects, a snack or small medical equipment, it is an ideal support for your day surgery services.* Compact, it takes up no more space than an IV pole. Its star-shaped base fits perfectly with all our products.

- 2 sections folding table : 75 x 35 cm unfolded / 37,5 x 35 cm folded). Weight : 17kg
- Height adjustment range : from 82 to 112 cm
- Unbraked product
- 4 rotary wheels Ø 40mm
- Storage space under the table : 27 x 25 x 10 cm.
- Handbag / backpack holder (maximum 15kg)
- IV Pole holder

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks



Removable IV pole, 2-stage, 2 hooks  
TSR2

## ASSISTANT DIMENSIONS

	LENGHT	WIDTH	HEIGHT
<b>OVERALL</b>	57,5 to 99,7 cm	57,5 to 63,3 cm	82 to 112 cm
<b>BASE</b>			7.5 cm

Maximum load : 15 Kg folded / 5Kg at the table end



Storage space

# INFORMALYS

Mobile computer table with variable height. Informalys is the ideal product to accompany the medical staff in their units. **Designed for laptops, it has a handle to facilitate its movement, as well as a mouse pad.** Mobile, resistant and easy to clean, it is the perfect companion for your computer equipment.

- Computer table Trespa matt white tray
- For laptops
- Flexibility - All input-output ports on the laptop are accessible
- Front handle
- 3 self lockable pivoting casters Ø 80 mm
- Retractable mouse pad

## INFORMALYS DIMENSIONS

	LENGHT	WIDTH	HEIGHT
OVERALL	60 cm	60 cm	74 to 113 cm
TRAY	57 cm	54 cm	2.4 cm
BASE	60 cm	60 cm	16 cm

**SAFE WORKING LOAD : 15 Kg**





# LOCKER LARGE MODEL

# ACCESSORIES

The bedside cabinet locker is the ideal product for the reception of the patient in an day surgery unit. Intended for the storage of your patients' personal belongings, the product is equipped with a secure locking system. Mobile, resistant and easy to clean, it is the perfect companion to a care chair.

- Stratified compact, acrylic resin finishing touch, white color
- 4 double running wheel castors Ø 100mm
- Removable shoes tray in white ABS
- Key lock system (2 keys included)
- Push handle - Drawer in ABS - Wardrobe for clothes
- Identification system for patient / care chair / locker
- Anti-fall edge
- Door opening on the right (Left on request)
- Weight : 25 kg

REFERENCE	DESCRIPTION
<b>VEST_C</b>	Universal key for key lock system
<b>VEST_SD</b>	Digital door lock with code (1 universal key included)
<b>VEST_CSD</b>	Additional universal key for digital lock

## LOCKER DIMENSIONS

	LENGHT	WIDTH	HEIGHT
<b>OVERALL</b>	59 cm	47.5 cm	126 cm

**SAFE WORKING LOAD : 15 Kg**



Universal key for key lock system  
VEST\_C



Digital door lock with code  
VEST\_SD



Additional universal key for digital lock  
VEST\_CSD



# LOCKER SMALL MODEL

# ACCESSORIES

The bedside cabinet locker is the ideal product for the reception of the patient in an day surgery unit. **Intended for the storage of your patients' personal belongings, the product is equipped with a secure locking system.** Mobile, resistant and easy to clean, it is the perfect companion to a care chair.

- Stratified compact, acrylic resin finishing touch, white color
- 4 double running wheel castors Ø 80mm
- Removable shoes tray in white ABS
- Key lock system (2 keys included)
- Push handle - Drawer in ABS - Wardrobe for clothes
- Anti-fall edge
- Door opening on the right (Left on request)
- Weight : 21 kg

## LOCKER DIMENSIONS

	LENGHT	WIDTH	HEIGHT
<b>OVERALL</b>	59 cm	47.5 cm	92 cm

**SAFE WORKING LOAD : 15 Kg**



REFERENCE	DESCRIPTION
<b>VEST_C</b>	Universal key for key lock system
<b>VEST_SD</b>	Digital door lock with code (1 universal key included)
<b>VEST_CSD</b>	Additional universal key for digital lock
<b>VEST_PB</b>	Waste basket for 10 L bag (only available for VEST_P)




Waste basket for 10 L bag  
VEST\_C



Digital door lock with code  
VEST\_SD



Additional universal key for  
digital lock  
VEST\_CSD



Waste basket for 10 L bag  
VEST\_PB



Locker large model  
VEST



Locker small model  
VEST\_P



System allowing the installation of a cyclor and its consumables for daily or home peritoneal dialysis.

Allows the installation of 4 types of known devices : Sleep-Safe Harmony - HomeChoice Claria - S3 - NXstage

- Central part in Trespa and steel materials
- Top part in Trespa materials
- Structure made of steel. Epoxy coated.
- 3 casters Ø 75mm
- Manual variable height
- Top plate which can support until 80kg, room for a 60x40cm machine
- Easy adjustable height, even the machine on it, from 48cm to 78cm
- Weight : 15 KG - Can be easily taken to pieces

REFERENCE	DESCRIPTION
<b>MODULSUP</b>	Attachment hooks for bags
<b>MODULTUBS</b>	Tubes for the top part
<b>MODU-DESS</b>	Removable additional table for consumables. Retractable legs. 520 x 800 x 390 mm Weight : 5 kg
<b>MODUSAC</b>	Transport bag

## MODULYS COMPOSITION

### MODULCO


Top part :  
600 x 410 mm



+ column  
+ down part:  
H 620 x 550 mm  
weight : 15 kg  
Down height: 490mm  
Up height: 780 mm

### MODU-DESS

520 x 800 x 390 mm  
Weight : 5kg



### MODULTUBS

Tubes for the top part




310 x 555 mm  
Weight: 1,2 kg

### MODULSUP


Attachment hooks for bags




310 x 65 x 40 mm  
Weight : 400 g



Attachment hooks for bags  
MODULSUP



Tubes for the top part  
MODULTUBS



MODU-DESS





**Versatile,** adapted to different hospital units

**Adaptable,** many options available

**Two models** Wall-mounted and two mobile bases

**Easy to use,** thanks to its mobility





# HOSPITAL SCREENS

# ACCESSORIES

## Rigid panel screens.

- M1 white PVC panels
- From 5 to 13 panels
- Wall mounting (including 10 cm spacer / on request 20 or 30 cm) or two mobile bases
- Finishing profiles available in different colours: Blue, Grey and Orange
- Strap for closing the screen

REFERENCE	DESCRIPTION
<b>TENENT1</b>	Strut of 10 cm (included with wall mounting)
<b>TENENT2</b>	Strut of 20 cm (On request with wall mounting)
<b>TENENT3</b>	Strut of 30 cm (On request with wall mounting)

## SCREENS DIMENSIONS

HEIGHT	LENGTH (meter)	NUMBER OF PANELS								
		5	6	7	8	9	10	11	12	13
<b>1,65 M</b>		1,365	1,575	1,785	1,995	2,205	2,415	2,625	2,835	3,045



Strut of 10, 20 or 30 cm  
TENENT1 2 3  
(included with wall mounting)

## FINISHING PROFILES BLUE - GREY - ORANGE





# ACCESSORIES



Lightweight Arm & Handsurgery table in Carbon Fibre with twin clamps (hook on)  
One (1) glassfibre supporting leg included for left or right, mounting depending of location of surgeon or C-arm/microscope.  
Cushion included.  
Length: 815 mm. Weight: 5 kgs.  
Metalfree space.

**CHIRMAIN**



Safety belt for stretchers, chairs and armchairs

**ACC-CTP**



Protection tunel from rain for CARVI and I-CARE

**PROPLUIE**



Roll board / patient transfert device  
Size: 177 X 55cm ( $\pm 3\%$ ) (unfolded) / 88 x 55cm( $\pm 3\%$ ) (folded)  
Weight : 2.4 kg ( $\pm 3\%$ )  
100% Nylon  
Aluminum Handle  
SAFE WORKING LOAD : 200 kg

**SWITCHEASY**



Tél : +33 (0)4 74 08 71 71  
Fax : +33 (0)4 74 08 71 72  
[info@bmb-medical.com](mailto:info@bmb-medical.com)

100 Av. du Formans - CS 11001  
01604 TREVOUX Cedex  
FRANCE

[www.bmb-medical.com](http://www.bmb-medical.com)