

EXAMINATION COUCH





CONTENTS

Our range of EXAMINATION COUCH is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and esthetics, our products are designed to accommodate and meet the needs of patients and caregivers.




DIVIA



DIVIA MIXTE



COLOR CHART



Stability

Intensive use


Easy to use electric height adjustment



DIVIA VARIABLE-HEIGHT COUCH

DIVIA is an examination couch with electric height adjustment on columns.

The couch is particularly suitable for emergency services and intensive use.



DIVIA DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	205 cm	75 to 85 cm	46 to 86 cm
SEAT	65 x 200 cm or 75 x 200 cm		
INCLINE FOLDER	0 to 90°		

PRODUCT WEIGHT : 80 Kg
SAFE WORKING LOAD: 250 kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of EXAMINATION COUCH is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and esthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Examination couch with electric height adjustment on columns.

- Electric height adjustment : 46 - 86 cm
- Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°)
- Adjustable headrest (0 - 90°) with 2 gas springs
- 2-sectional mattress base : 65x 200 cm or 75 x 200 cm
- Flat upholstery, high density foam 35 kg/m³ - Thickness 6,5 cm
- Jersey with PVC coating, antibacterial and anti-fire (M1) treatment
- Steel structure, epoxy coated (light grey)
- Frame with 2 electric columns providing great stability
- Frame protection cover
- Remote control
- 4 lockable wheels Ø 125 mm
- Examination towel dispenser head end



Remote control




Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°)

Electric height adjustment : 46 - 86 cm



Examination towel dispenser head end



Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 85 % / Coton 15 %
Thickness: 1,2 mm +/- 0,1 mm
Effective against bacteria (MRSA)
Phthalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



FIRE PROOF
M1

ACCESSORIES



Antistatic system
DVRA



Casters Ø 125 mm
with central brake
DVRS1



Double running castors Ø 150 mm
with central brake
DVRS2



Electric adjustable backrest
DVRBE



Wireless remote control
(in addition of the wired
remote control)
DVTSF




Foot pedal control (in addition of the wired
remote control) If the options DVTAP and
DVTSF are chosen, the wired remote control
is removed)
DVTAP



Drop-down side rails (pair)
DVBR




Holder for clamps (pair)
DVSE




Fixation clamps for heels, armrest for
blood sampling or leg support (pair)
(DVSE required)
DVETU



Heels (pair) (DVETU required)
DVET



Leg support (pair) (DVETU required)
DVSJ




Removable 2-stage IV pole, 2 hooks
(DVTSR required)
TSR2

REFERENCE	DESCRIPTION
DVRA	Antistatic system
DVRS1	Casters Ø 125 mm with central brake
DVRS2	Double running castors Ø 150 mm with central brake
DVRBE	Electric adjustable backrest
DVTSF	Wireless remote control (in addition of the wired remote control)
DVTAP	Foot pedal control (in addition of the wired remote control) If the options DVTAP and DVTSF are chosen, the wired remote control is removed)
DVBR	Drop-down side rails (pair)
DVSE	Holder for clamps (pair)
DVETU	Fixation clamps for heels, armrest for blood sampling or leg support (pair) (DVSE required)
DVET	Heels (pair) (DVETU required)
GT	Anatomic shaped armrest (DVETU necessary)
DVSJ	Leg support (pair) (DVETU required)
DVHSJ	Protection cover for leg rest (pair)
DVTSR	Clamps for IV pole holder (pair)
TSR2	Removable 2-stage IV pole, 2 hooks (DVTSR required)
TSR4	Removable IV pole, 2-stage, 4 hooks (DVTSR required)
ACC-CTG	Comfort pillow (grey)
DVPC7	Crystal protection cover for examination couch width 65cm
DVPC8	Crystal protection cover for examination couch width 75cm
DVBBA	Battery
DVLAMP	Examination lamp (DVETU and DVSE required)
EQUI	Equipotential connection
DVRT	DIN rails



Comfort pillow (grey)
ACC-CTG



Examination lamp
(DVETU required)
DVLAMP

Intensive use

Versatile suitable for gynecological examinations

Easy to use Electric backrest and height adjustment



GYNECOLOGICAL COUCH DIVIA MIXTE

DIVIA is a gynecological examination couch with electric height adjustment on columns.

It has a second electrically operated intermediate backrest for gynecological examinations.



DIVIA MIXTE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	178 cm	75 to 85 cm	46 to 86 cm
SEAT	65 x 178 cm or 75 x 178 cm		
INCLINE FOLDER	0 to 80°		

PRODUCT WEIGHT : 90 Kg
SAFE WORKING LOAD: 250 kg




NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of EXAMINATION COUCH is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and esthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Examination couch with variable height for gynaecological examinations.

- Electric height adjustment : 46 - 86 cm
- Electric adjustable backrest (0 - 80°)
- Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°)
- 3-sectional mattress base : 65x 178 cm or 75 x 178 cm
- Flat upholstery, high density foam 35 kg/m3. Width 6,5 cm
- Steel structure, epoxy coated (light grey)
- Frame with 2 electric columns providing great stability
- Frame protection cover
- Remote control
- 4 lockable double running wheels Ø 125 mm
- Examination towel dispenser head end
- Jersey with PVC coating, antibacterial and anti-fire (M1) treatment



Remote control




Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°)

Electric height adjustment : 46 - 86 cm



Examination towel dispenser head end



Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA



PLUM



DEEP OCEAN



PETROL



GREY

URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 85 % / Coton 15 %
Thickness: 1,2 mm +/- 0,1 mm
Effective against bacteria (MRSA)
Phthalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



FIRE PROOF
M1

ACCESSORIES



Antistatic system
DVRA



Casters Ø 125 mm
with central brake
DVRS1



Double running castors Ø 150 mm
with central brake
DVRS2



Wireless remote control
(in addition of the wired
remote control)
DVTSF



Foot pedal control (in addition of the wired
remote control) If the options DVTAP and
DVTSF are chosen, the wired remote control
is removed)
DVETAP



Drop-down side rails (pair)
DVBR



Holder for clamps (pair)
DVSE



Fixation clamps for heels, armrest for
blood sampling or leg support (pair)
(DVSE required)
DVETU



Heels (pair) (DVETU required)
DVET



Leg support (pair) (DVETU required)
DVSJ




Removable 2-stage IV pole, 2 hooks
(DVTSR required)
TSR2




Comfort pillow (grey)
ACC-CTG

REFERENCE	DESCRIPTION
DVRA	Antistatic system
DVRS1	Casters Ø 125 mm with central brake
DVRS2	Double running castors Ø 150 mm with central brake
DVTSF	Wireless remote control (in addition of the wired remote control)
DVTAP	Foot pedal control (in addition of the wired remote control) If the options DVTAP and DVTSF are chosen, the wired remote control is removed)
DVBR	Drop-down side rails (pair)
DVSE	Holder for clamps (pair)
DVETU	Fixation clamps for heels, armrest for blood sampling or leg support (pair) (DVSE required)
DVET	Heels (pair) (DVETU required)
GT	Anatomic shaped armrest (DVETU necessary)
DVSJ	Leg support (pair) (DVETU required)
DVHSJ	Protection cover for leg rest (pair)
DVTSR	Clamps for IV pole holder (pair)
TSR2	Removable 2-stage IV pole, 2 hooks (DVTSR required)
TSR4	Removable IV pole, 2-stage, 4 hooks (DVTSR required)
ACC-CTG	Comfort pillow (grey)
DVBBA	Battery
DVLAMP	Examination lamp (DVETU and DVSE required)
DVMB	Recovery Tank (grey) (can hold up to 0,4 Liter)
EQUI	Equipotential connection
DVRT	DIN rails



Examination lamp
(DVETU required)
DVLAMP



Recovery Tank (grey)
DVMB



URBAN TEXTILE M1

DIVIA RANGE



ABSINTHE



TOMATO



PANSY



FUCHSIA



PLUM



DEEP OCEAN



PETROL



GREY

COLOR CHART




CLASSIFICATION
M1
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 85 % / Coton 15 %
Thickness: 1,2 mm +/- 0,1 mm



Effective
against bacteria
(MRSA)



Solvent-free



Phtalate-free



Martindale
abrasion 100000
cycles



Waterproof

A medical professional must always rely on his or her own professional clinical judgment when deciding whether to use a product(s) when treating a particular patient.
A healthcare professional anyway, refers to the package insert, product label and/or instructions for use before using any BMB MEDICAL product.
The products depicted are CE marked in accordance with EU Regulations and Directives.
The products are Class 1 Medical Devices.



Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans - CS 11001

01604 TREVOUX Cedex

FRANCE

www.bmb-medical.com