

TRANSFER CHAIR I-MOVE





CONTENTS

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position.

Manageable, robust and attractive it will bring for a long time a touch of modernity and color to your department.

A wide range of accessories allows it to adapt to each use.



I - MOVE



I - MOVE OUTDOOR



I - MOVE BARIATRIC



I - MOVE EZ-GO



I - MOVE EZ-GO BARIATRIC





TRANSFER CHAIR I-MOVE

Sturdy

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Easy to handle & Stackable for optimal storage





I-MOVE is a transfer chair, it allows the reception of the patient and his transport in sitting position.



I-MOVE DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	75 cm	100 cm	H. 55 cm D. 46 cm W. 46 cm	35 Kg

PRODUCT WEIGHT : 35 Kg
SAFE WORKING LOAD : 200 Kg



Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

I-MOVE transfer chair for patient welcoming

- Steel frame epoxy coated
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Anti-shock protection
- Cushion in foam - Thickness 5 cm
- Seamless finish, vinyl coating



Back wheels Ø 300 mm with central braking system



Front swivel wheels Ø 125 mm
Sliding foot rest



Patient file holder



Polyurethan armrests, can be folded up


Different colors



FUSCHIA



GREEN




BLUE

Composition : PVC 83 % / Coton 17 %
Thickness: 1,2 mm +/- 0,1 mm
Light resistance: >= 6



FIRE RATING
M1 / AM 18 / EN 1021.1-2-3.

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
TSR2-F	Fixed IV pole, 2-stage, 2 hooks
CLORA	Anti-Static system
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with TSR2)
TROPO	Accessory basket
TROSO2	Holder for oxygen cylinder
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO -EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
PSU	Hooks to hold bag



Removable IV pole,
2-stage, 2 hooks
TSR2



Accessory basket
TROPO



Holder for oxygen
cylinder
TROSO2



Adjustable leg rest
(1 piece)
TRORJ



Coin acceptor
TROMO



Hooks to hold bag
PSU





TRANSFER CHAIR I-MOVE OUTDOOR

Sturdy It is designed to withstand weather conditions

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Easy to handle & Stackable For optimal storage

I-MOVE is a transfer chair, it allows the reception of the patient and his transport in sitting position, outdoor.

Its coating is weather resistant, so it can be stored outside.




I - M O V E O U T D O O R

I-MOVE OUTDOOR DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	75 cm	100 cm	H. 55 cm D. 46 cm W. 46 cm	35 Kg

PRODUCT WEIGHT : 35 Kg
SAFE WORKING LOAD : 200 Kg



Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position indoor and outdoor.

I-MOVE OUTDOOR transfer chair for patient welcoming

- Steel frame epoxy coated
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Anti-shock protection
- Polypropylene plastic fitting
- Colour : blue



Back wheels Ø 300 mm with central braking system



Polyurethan armrests, can be folded up



Patient file holder




Front swivel wheels Ø 125 mm
Sliding and self-blocking foot rest



Composition : Polypropylene

I-MOVE OUTDOOR



Removable IV pole,
2-stage, 2 hooks
TSR2



U-shape anti-theft, fitted with 2 hooks
infusion holder
TROAA



Accessory basket
TROPO



Holder for oxygen
cylinder
TROSO2



Wall fixation
TROMOP



Hooks to hold bags
PSU

ACCESSORIES

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
TSR2-F	Fixed IV pole, 2-stage, 2 hooks
CLORA	Anti-Static system
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (incompatible with TSR2)
TROPO	Accessory basket
TROSO2	Holder for oxygen cylinder
TROMO -EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
PSU	Hooks to hold bag

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit



TRANSFER CHAIR I-MOVE BARIATRIC

I-MOVE is a bariatric transport chair for your patients' reception and sitting transport.




I-MOVE BARIATRIC

I-MOVE BARIATRIC DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	90 cm	102 cm	H. 55 cm D. 46 cm W. 65 cm	50 Kg

PRODUCT WEIGHT : 50Kg
SAFE WORKING LOAD : 250 Kg



Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position.

I-MOVE transfer chair for patient welcoming

- Steel frame epoxy coated
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Anti-shock protection
- Cushion in foam - Thickness 8 cm
- PVC coated Jersey



Back wheels Ø 300 mm with central braking system



Front swivel wheels Ø 125 mm
Sliding foot rest



Patient file holder



Polyurethan armrests, can be folded up

Differents colors



ORANGE



GREEN



LAVANDER

Composition : PVC 83 % / Coton 17 %
Thickness : 1,2 mm +/- 0,1 mm
Light resistance : >= 6




FIRE RATING
M1 / AM 18 / EN 1021.1-2-3.

I-MOVE BARIATRIC



ACCESSORIES

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
TSR2-F	Fixed IV pole, 2-stage, 2 hooks
CLORA	Anti-Static system
TROPO	Accessory basket
TROSO2	Holder for oxygen cylinder
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO -EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
PSU	Hooks to hold bag




IV pole, 2-stage,
2 hooks
TSR2



Accessory basket
TROPO



Holder for oxygen
cylinder
TROSO2



Adjustable leg rest
(1 piece)
TRORJ



Coin acceptor
TROMO



Hooks to hold bag
PSU

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Motorized Very easy to use, progressive speed



TRANSFER CHAIR I-MOVE EZ-GO

I-MOVE EZ-GO is a transport chair for your patients' reception and sitting transport. It is equipped with an autonomous motorized wheel, two batteries and an external charger. Thanks to I-MOVE EZ-GO, patients transfer is child's play!




I - M O V E E Z - G O

I-MOVE EZ-GO DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	70 cm	102 cm	H. 55 cm D. 46 cm W. 46 cm	52 Kg

PRODUCT WEIGHT : 52 Kg
SAFE WORKING LOAD : 200 Kg



Motorized wheel including battery

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

Motorized transfer chair I-MOVE EZ-GO for patient welcoming

- Motorized drive system with autonomous wheel
- Switchable battery - 5 km autonomy per battery - Charge time : 5 hours
- Lithium battery without memory effect
- Forward/Backward twist grip with progressive speed acceleration
- Steel frame epoxy coated
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Accessory basket
- Anti-shock protection
- Supplied with 2 batteries and external charger
- Cushion in foam - Thickness 5 cm
- Seamless finish, vinyl coating



Back wheels Ø 300 mm with central braking system



Polyurethan armrests, can be folded up



Sliding footrest
Front swivel wheels Ø 125 mm



Accessory basket

Different colors



FUSCHIA



GREEN




BLUE

Composition : PVC 83 % / Coton 17 %
Thickness: 1,2 mm +/- 0,1 mm
Light resistance: >= 6



FIRE RATING
M1 / AM 18 / EN 1021.1-2-3.




Removable IV pole,
2-stage, 2 hook
TSR2




U-shape anti-theft, fitted with 2 hooks
infusion holder
TROAA



Holder for oxygene
cylinder




Hooks to hold bags
PSU



Adjustable leg rest (1
piece)
TRORJ



Wall fixation
TROMOP



Additional battery
EZOBA



Additional external charger
EZOCH

ACCESSORIES

REFERENCE	DESIGNATION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (incompatble with TSR2)
TROSO2	Holder for oxygen cylinder
PSU	Hooks to hold bags
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO-EP	Coin acceptor for 1€ coin or equivalent
TROMO-EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
EZSBA	Without battery (depreciation)
EZOBA	Additional battery
EZOCH	Additional external charger

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Motorized Very easy to use, progressive speed



TRANSFER CHAIR I-MOVE EZ-GO BARIATRIC

I-MOVE EZ-GO is a transport chair for your patients' reception and sitting transport. It is equipped with an autonomous motorized wheel, two batteries and an external charger.

Thanks to I-MOVE EZ-GO, patients transfer is child's play! It is particularly suitable for the transfer of overweight patients.




I-MOVE EZ-GO BARIATRIC

I-MOVE EZ-GO BARIATRIC DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT	WEIGHT
I-MOVE	105 cm	90 cm	102 cm	H. 55 cm D. 46 cm W. 65 cm	65 Kg

PRODUCT WEIGHT : 65 Kg
SAFE WORKING LOAD: 250 Kg



Motorized wheel including battery

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

Motorized transfer chair I-MOVE EZ-GO for patient welcoming

- Motorized drive system with autonomous wheel
- Switchable battery - 5 km autonomy per battery - Charge time : 5 hours
- Lithium battery without memory effect.
- Forward/Backward twist grip with progressive speed acceleration
- Steel frame epoxy coated RAL 9010
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Accessory basket
- Supplied with 2 batteris and external charger
- Anti-shock protection
- Cushion in foam - Thickness 8 cm



Accessory basket



Polyurethan armrests, can be folded up



Back wheels Ø 300 mm with central braking system



Front swivel wheels Ø 125 mm
Sliding footrest


Différents coloris



ORANGE



GREEN




LAVANDER

Composition : PVC 83 % / Coton 17 %
Thickness: 1,2 mm +/- 0,1 mm
Light resistance: >= 6



FIRE RATING
M1 / AM 18 / EN 1021.1-2-3.

I-MOVE EZ-GO BARIATRIC



Removable IV pole,
2-stage, 2 hook
TSR2



U-shape anti-theft, fitted with 2 hooks
infusion holder
TROAA



Holder for oxygen cylinder
TROSO2




Hooks to hold bags
PSU




Adjustable leg rest
(1 piece)
TRORJ



Wall fixation
TROMOP



Additional battery
EZOBA



Additional external charger
EZOCH

ACCESSORIES

REFERENCE	DESIGNATION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (incompatible with TSR2)
TROSO2	Holder for oxygen cylinder
PSU	Hooks to hold bags
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO-EP	Coin acceptor for 1€ coin or equivalent
TROMO-EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
EZSBA	Without battery (depreciation)
EZOBA	Additional battery
EZOCH	Additional external charger

A medical professional must always rely on his or her own professional clinical judgment when deciding whether to use a product(s) when treating a particular patient.

A healthcare professional anyway, refers to the package insert, product label and/or instructions for use before using any BMB MEDICAL product.

The products depicted are CE marked in accordance with EU Regulations and Directives.

The products are Class 1 Medical Devices.



Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans - CS 11001
01604 TREVOUX Cedex
FRANCE

www.bmb-medical.com