

# EXTRA BED RANGE





## CONTENTS

Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units, as well as any other department where companions may be present.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.



S Y R E E T A



A R E T H A



E L L I O T



C L A S S I Q



C L A S S I Q 2 P L



M A X



M I N



C L E M



C O L O R C H A R T



Comfortable

The most compact chair on the market

Adaptable Various options available

Versatile Chair-bed with three positions: seated/lying down/chaise longue



## SYREETA CONVERTIBLE CHAIR

The SYREETA companion bed is every bit as comfortable as a lounge chair.

This multi-purpose, convertible chair is ideal for visitors, offering a comfy seat and the comfort of a real bed.




# SYREETA

## SYREETA DIMENSIONS

	LENGTH	WIDTH	HEIGHT	MATTRESS	SEAT
SYREETA 70	90 cm (folded) 230 cm (open)	87 cm	90 cm	188 x 70 x 10 cm	H. 50 cm D. 60 cm

Operational reliability load: 160 kg

### SYREETA 70




Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units, as well as any other department where companions may be present. Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

Convertible chair-bed with steel frame. All panels and cushions are retractable for cleaning and repairs.

- Slim armrests
- Sleeping berth: 195 cm long
- Seat made from 35 kg/m<sup>3</sup> high resiliency complex foams (Dacron and HR)
- 4 wheels (ø 80 mm) with independent brakes
- M2 flame retardant cover, 10 color options
- HR foam mattress with a density of 35 kg/m<sup>3</sup>, 10 cm thick
- Bed base: 13 resin slats for better hygiene and easy cleaning

HR35 mattress + Removable flame retardant PU cover



Wheels (ø 80 mm) with independent brakes



Slim armrests



Bed base 13 resin slats

Faux leather covers: VALENCIA M2

### SOOTHING



GRAPHITE  
107-40033



TAUPE  
107-0034



SISAL  
107-1048



ROSÉ  
107-2114

### REFRESHING



TITAN  
107-4042



NIGHT  
107-2108



TEAL  
107-2119

### ENERGIZING



OLIVE  
107-50013



NECTARINE  
107-6003



SAFRAN  
107-6012

Composition : Vinyl / Urethane Topcoat  
Thickness : 1,1 mm  
>300.000 Martindale cycles  
Anti-stain finish  
Resistant to sulfide stains  
Effective against bacteria  
Phthalate free  
Anti-static finish




FIRE RATING  
M2



# SYREETA

## ACCESSORIES

REFERENCE	NAME
FC100	Reading light + USB port
FC101	Light beams
FC102	Fixed headrest
FC107	2 soft cushions
FC110	Manual recliner with adjustable backrest
FC111	Storage compartment




Reading light + USB port  
FC100



Fixed headrest  
FC102



2 soft cushions  
FC107



Light beams  
FC101



Manual recliner with adjustable backrest  
FC110



Storage compartment  
FC111



Comfortable

Adaptable Various options available

The most compact chair on the market

Versatile Chair-bed with three positions: seated and relax



## ARETHA CONVERTIBLE CHAIR

The ARETHA companion bed is every bit as comfortable as a lounge chair.

Ideal for the visitor, this multi-functional convertible armchair offers a pleasant seat and the comfort of a real bed.






# ARETHA

## ARETHA DIMENSIONS

	LENGTH	WIDTH	HEIGHT	MATTRESS	SEAT
ARETHA	100 cm (folded) 240 cm (open)	87 cm	90 cm	188 x 70 x 10 cm	H. 17 cm D. 66 cm

Operational reliability load: 160 kg



Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

Armchair convertible into bed with wooden armrests and steel frame covered with epoxy. All panels and cushions are retractable for cleaning or repairs.

- Independent seat
- Full armrests
- Sleeping berth: 195 cm long
- Seat made from 35 kg/m<sup>3</sup> high resiliency complex foams (Dacron and HR)
- 4 wheels (ø 80 mm) with independent brakes
- M2 flame retardant cover, 10 color options
- HR foam mattress with a density of 35 kg/m<sup>3</sup>, 10 cm thick
- Bed base: 13 resin slats for better hygiene and easy cleaning

HR35 mattress + Removable flame retardant PU cover



Bed base 13 resin slats



Wheels (ø 80 mm) with independent brakes




Full armrests

Faux leather covers: VALENCIA M2


### SOOTHING



### REFRESHING



### ENERGIZING



Composition : Vinyl / Urethane Topcoat  
Thickness : 1,1 mm  
>300.000 Martindale cycles  
Anti-stain finish  
Resistant to sulfide stains  
Effective against bacteria  
Phthalate free  
Anti-static finish



FIRE RATING  
M2

# ARETHA

## ACCESSORIES

REFERENCE	NAME
FC100	Reading light + USB port
FC102	Fixed headrest
FC110	Reclining Backrest from 0° to 15°



Reading light + USB port  
FC100



Fixed headrest  
FC102



Reclining backrest  
FC110



**Adaptable & Functional** Various options

**Very compact armchair**

**Versatile** Chair-bed with three positions: seated and relax



## ELLIOT ARMCHAIR - BED


The ELLIOT companion bed is ideal for the visitor, this multi-functional convertible armchair offers a pleasant seat and the comfort of a real bed.



## ELLIOT DIMENSIONS

	LENGTH	WIDTH	HEIGHT	MATTRESS	SEAT
ELLIOT	80 cm (folded) 210 cm (open)	84 cm	76 cm	190 x 75 x 10 cm	H. 10 cm D. 60 cm

Operational reliability load: 120 kg




Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

Armchair convertible into bed with wooden armrests and steel frame covered with epoxy. All panels and cushions are retractable for cleaning or repairs.

- Independent seat
- Full armrests
- Sleeping mechanism in 75x90
- Seat made from 35 kg/m<sup>3</sup> high resiliency complex foams (Dacron and HR)
- 4 wheels (ø 50 mm)
- M2 flame retardant cover, 10 color options
- Base : tubular frame with welded mesh

Seat made from 35 kg/m<sup>3</sup> high resiliency complex foams (Dacron and HR)



4 wheels (ø 50 mm)




Full armrests

Faux leather covers: VALENCIA M2


### SOOTHING



### REFRESHING



### ENERGIZING



Composition : Vinyl / Urethane Topcoat  
Thickness : 1,1 mm  
>300.000 Martindale cycles  
Anti-stain finish  
Resistant to sulfide stains  
Effective against bacteria  
Phthalate free  
Anti-static finish







## CLASSIQ CHAIR-BED

Easy to use

Comfortable full armrests

Lockable storage compartment

The CLASSIQ companion bed is the narrowest chair-bed (83 cm) in our product range.


Ideal for small rooms in Maternity, Neonatal or Pediatric units.



## CLASSIQ DIMENSIONS

	LENGTH	WIDTH	HEIGHT	BERTH	SEAT
OVERALL	79 cm (folded) 186 cm (open)	83 cm	98 cm	186 x 62 x 10 cm	H. 48 cm D. 49 cm

Operational reliability load: 150 kg




Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

### Convertible chair-bed

- Chair-bed with steel frame
- M1 flame retardant Gingko upholstery
- 4 omnidirectional wheels with brakes ( 75 mm)
- Storage compartment under the seat (key lock closure)
- Comfortable fixed armrests
- HRfoam with 35 kg/m<sup>3</sup> density.
- Plywood internal frame

HR foam with 35 kg/m<sup>3</sup> density.



Wheels (ø 75 mm) with independent brakes

Standard covers: M1 fire retardant Gingko faux leather

For effective care, clean regularly with pH-neutral soap and a soft brush. Always rinse with clean water. Avoid all commercially-available detergents, solvents and stain removers not included in NF French standards, as well as all products not recommended by us. Immediately clean marks such as pen ink, coffee, wine, eosin or colorants from jeans to prevent permanent staining of the coated fabric.

Composition: 85% PVC, 15% Cotton



FIRE RATING  
M1



Color shown:  
Lilas 013 32 008

Comfortable armrests



Storage compartment under the seat (key lock closure)



Citrouille  
013 32 015



Begonia  
013 32 046



Aigue Marine  
013 32 012



Endive  
013 32 032



Menthe  
013 32 042



Taupe  
013 32 043



Tomette  
013 32 014



# CLASSIQ

## ACCESSORIES

REFERENCE	NAME
TENCOUSSIN	Comfortable triangular cushion



Comfortable triangular cushion  
TENCOUSSIN



Easy to use

Comfortable full armrests

Lockable storage compartment

1 or 2 people, as required



## CLASSIQ DOUBLE SOFA-BED

The CLASSIQ double companion bed is the ideal sofa-bed for visitors to the Maternity unit.






# CLASSIQ2PL

## CLASSIQ2PL DIMENSIONS

	LENGTH	WIDTH	HEIGHT	BERTH	SEAT
OVERALL	79 cm (folded) 186 cm (open)	147 cm	98 cm	190 x 123 x 10 cm	H. 48 cm D. 49 cm

Operational reliability load: 150 kg




Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

### Sofa transforms into a bed

- Sofa-bed with steel frame
- M1 flame retardant Ginkgo upholstery
- 4 omnidirectional wheels with brakes ( 75 mm)
- Storage compartment under the seat (key lock closure)
- Comfortable fixed armrests
- HRfoam with 35 kg/m<sup>3</sup> density.
- Plywood internal frame

Wheels (ø 75 mm) with independent brakes



Can be made into a single or double bed

HR foam with 35 kg/m<sup>3</sup> density.



Storage compartment under the seat (key lock closure)


### Standard covers: M1 fire retardant Ginkgo faux leather

For effective care, clean regularly with pH-neutral soap and a soft brush. Always rinse with clean water. Avoid all commercially-available detergents, solvents and stain removers not included in NF French standards, as well as all products not recommended by us. Immediately clean marks such as pen ink, coffee, wine, eosin or colorants from jeans to prevent permanent staining of the coated fabric.

Composition: 85% PVC, 15% Cotton



FIRE RATING  
M1



Color shown:  
Taupe 013 32 043



Citrouille  
013 32 015



Begonia  
013 32 046



Aigue Marine  
013 32 012



Endive  
013 32 032



Menthe  
013 32 042



Lilas  
013 32 008



Tomette  
013 32 014



# CLASSIQ2PL

## ACCESSORIES

REFERENCE	NAME
TENCOUSSIN	Comfortable triangular cushion



Comfortable triangular cushion  
TENCOUSSIN





## MAX MOBILE FOLDING BED

Permissible weight of **200 kg**

**Excellent mobility**

Designed for **intensive use**

Ideal for **facilitating daily decontamination**

The MAX companion bed is ideal for overnight visitors to your facility.

Its PVC mattress support replaces wooden slats or electro-welded mesh, which considerably reduces breakages.


Easy to unfold and designed for intensive use, the Max has been designed and developed to meet the needs of companions as well as nursing staff.



## MAX80 DIMENSIONS

	LENGTH	WIDTH	HEIGHT	MATTRESS
OVERALL	38 cm (folded) 204 cm (open)	81 cm	112 cm (folded) 42 cm (open)	198 x 75 x 10 cm

Operational reliability load: 200 kg



Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

Folding mobile bed with aluminum frame. The entire frame can be disassembled for repairs.

- Folding bed intended for intensive hospital use
- 4 wheels ø 80 mm with independent brakes
- HR foam mattress with a density of 35 kg, 10 cm thick
- M1 fire retardant polyurethane mattress cover
- Aluminum frame
- Blue grey epoxy steel base
- Mattress support made from white synthetic material, ideal for daily cleaning


HR35 mattress + Removable flame retardant PU cover



Easy to open



Wheels (ø 80 mm) with independent brakes



Mattress support made from white synthetic material



Excellent mobility  
Resin slats  
Ease of deployment



## MIN MOBILE FOLDING BED

The MIN companion bed is ideal for the visitor during a stay in your establishment.

Easy to unfold, the MIN has been designed and developed to meet the needs of both carers and nursing staff.





## MIN DIMENSIONS

	LENGTH	WIDTH	HEIGHT	MATTRESS
OVERALL	39 cm (folded) 195 cm (open)	86 cm	111 cm (folded) 42 cm (open)	190 x 80 x 10 cm

Operational reliability load: 200 kg



Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

Folding mobile bed with aluminum frame.

- 4 wheels ø 65 mm
- HR foam mattress with a density of 35 kg, 10 cm thick
- Delivered with its cover of protection mattress in polyurethane not fire
- Steel structure
- 12 Resin slats



HR35 mattress + fireproof PU cover with removable cover



Wheels (ø 65 mm)



Mattress support in resin slats



Delivered with its cover of protection mattress in polyurethane not fire





## CLEM FOLDING LOUNGE CHAIR

The CLEM is a lounge chair designed for companions or custodial staff.

With its small dimensions, the CLEM offers ease of use and saves on space.

Three positions: Seated / Ultra comfort / folded

Space saving

Easy to use

Comfort of an armchair





## CLEM DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	85 cm (chair position) 162 cm (lying position)	70 cm	115 cm (chair position) 80 cm (lying position) 98 cm (folded)
MATRESSES	155 cm	70 cm	

Operational reliability load: 140 kg




Our EXTRA-BED range is suited to Maternity, Pediatric and Neonatal units.

Combining comfort, ergonomics and style, our products are designed to accommodate and meet the needs of companions, parents, and caregivers.

Foldable relaxation chair with steel structure. All the cushions are retractable for cleaning or repairs.

- Manual folding relaxation chair
- Reclining, with ultra-comfort position: places the legs slightly above the heart to provide optimal well-being.
- Steel structure ø 20 mm.
- Foam HR 25 Kg - 7 cm thick
- Fixed black polypropylene armrests
- Retractable headrest
- High stability
- Foldable
- Fireproof M2 coating, 10 colors
- Chair weight : 9 kg

25 kg/m³ HR foam - 7 cm thick  
M2 fire retardant cover



Fixed black polypropylene armrests



Retractable headrest



Faux leather covers : VALENCIA M2

### SOOTHING



GRAPHITE  
107-40033



TAUPE  
107-0034



SISAL  
107-1048



ROSÉ  
107-2114

### REFRESHING



TITAN  
107-4042



NIGHT  
107-2108



TEAL  
107-2119

### ENERGIZING



OLIVE  
107-50013



NECTARINE  
107-6003



SAFRAN  
107-6012

Composition : Vinyl / Urethane Topcoat  
Thickness : 1,1 mm  
>300.000 Martindale cycles  
Anti-stain finish  
Resistant to sulfide stains  
Effective against bacteria  
Phthalate free  
Anti-static finish




FIRE RATING  
M2



# VALENCIA M2


## SOOTHING




GRAPHITE  
107-4003



TAUPE  
107-0034



SISAL  
107-1048




ROSÉ  
107-2114

## REFRESHING



TITAN  
107-4042



NIGHT  
107-2108



TEAL  
107-2119

## ENERGIZING



OLIVE  
107-5001



NECTARINE  
107-6003



SAFRAN  
107-6012

## COLOR CHART



FIRE RATING  
M2

COMPOSITION : Vinyl / Urethane Topcoat  
THICKNESS : 1,1mm

### FINISHINGS AND TREATMENTS



Effective  
against bacteria  
(MRSA)



Antistatic  
finish



Phthalates free



>300000  
Martindale cycles



Stain and sulfide  
resistant finish

A medical professional must always rely on his or her own professional clinical judgment when deciding whether to use a product(s) when treating a particular patient.  
A healthcare professional anyway, refers to the package insert, product label and/or instructions for use before using any BMB MEDICAL product.  
The products depicted are CE marked in accordance with EU Regulations and Directives.  
The products are Class 1 Medical Devices.



Tel: +33 (0)4 74 08 71 71

Fax: +33 (0)4 74 08 71 72

[info@bmb-medical.com](mailto:info@bmb-medical.com)

100 Av. du Formans - CS 11001

01604 Trevoux Cedex

France

[www.bmb-medical.com](http://www.bmb-medical.com)