

HOTEL-STYLE RECLINERS FOR REST & CARE



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units. Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.



CONTENTS

▶ CAREXIA FP

▶ CAREXIA FPE

▶ CAREXIA FPU

▶ CAREXIA FPU E

▶ TIM

▶ VELVET

▶ SWANN

▶ COLOR CHART

Comfortable

Easy to use

Adaptable Various options available

Versatile Adapted to different units



CAREXIA FP FIXED HEIGHT TREATMENT CHAIR

Carexia FP, is a fixed height versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Its pure and harmonious design is the vector of a positive image for your unit.



CAREXIA FP

CAREXIA FP DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	76.4 cm (folded) 170 cm (open)	69.5 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	57 cm	55 cm on seat
INCLINE FOLDER	85 to 60°		

PRODUCT WEIGHT : 65 to 85 Kg depending on the model
SAFE WORKING LOAD: 200 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Versatile treatment chair.

- Adjustable back rest from 80° to 60° with 1 gaz spring
- Fixed height : 55 cm
- Heavy duty design. Steel frame with powder coating. ABS cover
- Integrated headrest
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up
- High resilience foam 35 kg/m3
- Ergonomic shaped backrest and seat - PVC coated jersey



4 self lockable casters with double running wheels

Integrated headrest.



ABS basket


Polyurethan armrests, can be folded up



Faux leather covers : URBAN TEXTILE



ABSINTHE TOMATO



PANSY FUCHSIA PLUM



DEEP OCEAN PETROL GREY

URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 85 % / Coton 15 %
Thickness: 1,2 mm +/- 0,1 mm
Effective against bacteria (MRSA)
Phtalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



ACCESSORIES



4 castors Ø 125 mm with central lock
FPORS1




4 double running wheel castors Ø 150 mm with central lock
FPORS2




Push bar & Paper roll holder
FPOBP & FPOSR



Short legrest synchronized with the backrest. Trend position.
FPORJ



2 holders for IV pole (necessary for TSR2)
FPOSA



Removable IV pole, 2-stage, 2 hooks
TSR2



Comfort pillow (grey)
ACC-CTG



Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOAC3

REFERENCE	DESCRIPTION
FPORS1	4 castors Ø 125 mm with central lock
FPORS2	4 double running wheel castors Ø 150 mm with central lock
FPOBP	Push bar
FPORJ	Short legrest synchronized with the backrest. Trend position.
FPOSA	2 holders for IV pole (necessary for TSR2)
ACC-CTG	Comfort pillow (grey)
FPOAC4	Armrest comfort (pair)
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 2 cm thickness memory foam
SUPTAB	Fixed I-Pad support with antitheft device
TSR2	Removable IV pole, 2-stage, 2 hooks (FPOSA mandatory)
TSR4	Removable IV pole, 2-stage, 4 hooks (FPOSA mandatory)



CAREXIA FPE FIXED HEIGHT TREATMENT CHAIR

Comfortable

Easy to use Electric movements

Adaptable Various options available

Versatile Adapted to different units

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.




CAREXIA FPE

CAREXIA FPE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	76.4 cm (folded) 170 cm (open)	69.5 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	57 cm	55 cm on seat
INCLINE FOLDER	85 to 6°		

PRODUCT WEIGHT : 65 to 85 Kg depending on the model
SAFE WORKING LOAD: 200 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Multi-purpose treatment chair with electric movements.

- Heavy duty design - Steel frame with powder coating - ABS cover
- Integrated headrest
- Electric backrest and legrest
- Trend position
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up
- High resilience foam 35 kg/m³
- Ergonomic shaped backrest and seat - PVC coated jersey




4 self lockable casters with double running wheels



ABS basket



Integrated headrest




Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



ABSINTHE TOMATO



PANSY FUCHSIA PLUM



DEEP OCEAN PETROL GREY

URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 85 % / Coton 15 %
Thickness: 1,2 mm +/- 0,1 mm
Effective against bacteria (MRSA)
Phthalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



CAREXIA FPE



ACCESSORIES

REFERENCE	DESCRIPTION
FPORS1	4 castors Ø 125 mm with central lock
FPORS2	4 double runing wheel castors Ø 150 mm with central lock
FPOBP	Push bar
FPORJL	Long legrest
FPORP	Electric adjustable footrest (only for long legrest) (incompatible with FPOAC4)
FPOCP1	Foot pedal control for relax and trend position
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable IV pole, 2-stage, 2 hooks (FPOSA mandatory)
TSR4	Removable IV pole, 2-stage, 4 hooks (FPOSA mandatory)
ACC-CTG	Comfort pillow (grey)
FPOAC4	Armrest comfort (pair)
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 2 cm thickness memory foam
FPOLAMP	Reading lamp DEL
FPOBA	Battery 2 Ah
SUP-USB	USB Port
PROTEC RJL	Removable Crystal protection for footrest
PROTEC RP	Removable Crystal protection for footrest
SUPTAB	Fixed I-Pad support with antitheft device

ACCESSORIES



4 castors Ø 125 mm with central lock
FPORS1



4 double running wheel castors Ø 150 mm with central lock
FPORS2




Push bar & Paper roll holder
FPOBP & FPOSR




Long legrest & Electric adjustable footrest (only for long legrest)
FPORJL & FPORP



Foot pedal control for relax and trend position
FPOCP1



2 holders for IV pole (necessary for TSR2)
FPOSA



Removable IV pole, 2-stage, 2 hooks (FPOSA mandatory)
TSR2




Comfort pillow (grey)
ACC-CTG




Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) FPOAC3



Reading lamp DEL
FPOLAMP



USB Port
SUP-USB





CAREXIA FPU VARIABLE-HEIGHT CARE CHAIR

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.



Comfortable

Adaptable Various options available

Easy to use variable electric height


Versatile Adapted to different units

CAREXIA FPU

DIMENSIONS CAREXIA FPV

	LENGTH	WIDTH	HEIGHT
OVERALL	76.4 cm (folded) 170 cm (open)	73.2 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	57 cm	55 cm on seat
INCLINE FOLDER	85 to 60°		
VARIABLE HEIGHT			54 to 78 cm

PRODUCT WEIGHT : 65 to 85 Kg depending on the model
SAFE WORKING LOAD: 200 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Height variable electrical care chair

- Electric variable height from 54 to 78 cm
- Adjustable back rest from 80° to 60° with 1 gaz spring
- Heavy duty design - Steel frame with powder coating - ABS cover
- Integrated headrest
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up
- High resilience foam 35 kg/m3
- Ergonomic shaped backrest and seat - PVC coated jersey



4 self lockable casters with double running wheels



ABS basket



Integrated headrest



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 85 % / Coton 15 %
Thickness: 1,2 mm +/- 0,1 mm
Effective against bacteria (MRSA)
Phtalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof




FIRE PROOF
M1

ACCESSORIES



4 castors Ø 125 mm with central lock
FPORS1




4 double runing wheel castors Ø 150 mm with central lock
FPORS2



Push bar & Paper roll holder
FPOBP & FPOSR



Short legrest synchronized with the backrest. Trend position.
FPORJ



2 holders for IV pole (necessary for TSR2)
FPOSA



Removable IV pole, 2-stage, 2 hooks (FPOSA mandatory)
TSR2



Comfort pillow (grey)
ACC-CTG



Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) FPOAC3



USB Port
SUP-USB

REFERENCE	DESCRIPTION
FPORS1	4 castors Ø 125 mm with central lock
FPORS2	4 double runing wheel castors Ø 150 mm with central lock
FPOBP	Push bar
FPORJ	Short legrest synchronized with the backrest. Trend position.
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable IV pole, 2-stage, 2 hooks (FPOSA mandatory)
TSR4	Removable IV pole, 2-stage, 4 hooks (FPOSA mandatory)
ACC-CTG	Comfort pillow (grey)
FPOAC4	Armrest comfort (pair)
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 2 cm thickness memory foam
FPOBA	Battery 2 Ah
SUP-USB	USB Port



ALL ELECTRIC CARE CHAIR CAREXIA FPUE

Comfortable

All-electric care chair

Adaptable Various options available

Easy to use Variable height & electric movements

Versatile Adapted to different units

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.




CAREXIA FPVE

CAREXIA FPVE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	76.4 cm (folded) 177 cm (open) + long leg rest	73.2 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	57 cm	55 cm on seat
INCLINE FOLDER	80 to 6°		
VARIABLE HEIGHT			54 to 78 cm

PRODUCT WEIGHT : 65 to 85 Kg depending on the model
SAFE WORKING LOAD: 200 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Height variable electrical care chair

- Electric variable height from 54 to 78 cm
- Heavy duty design - Steel frame with powder coating - ABS cover
- Integrated headrest
- Electric backrest and legrest
- Trend position
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up
- High resilience foam 35 kg/m³
- Ergonomic shaped backrest and seat - PVC coated jersey



4 self lockable casters with double running wheels




Integrated headrest



ABS basket



Trend position



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 85 % / Coton 15 %
Thickness: 1,2 mm +/- 0,1 mm
Effective against bacteria (MRSA)
Phthalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



FIRE PROOF
M1

CAREXIA FPUE



ACCESSORIES

RÉFÉRENCE	DÉSIGNATION
FPORS1	4 castors Ø 125 mm with central lock
FPORS2	4 double runing wheel castors Ø 150 mm with central lock
FPOBP	Push bar
FPORJL	Long legrest
FPORP	Electric adjustable footrest (only for long legrest) (incompatible with FPOAC4)
FPOCP1	Foot pedal control for relax and trend position
FPOCP2	Foot pedal control for variable height on both sides
FPOSA	2 holders for IV pole (necessary for TSR2)
ACC-CTG	Comfort pillow (grey)
FPOAC4	Armrest comfort (pair)
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 2 cm thickness memory foam
FPOLAMP	Reading lamp DEL
FPOBA	Battery 2 Ah
SUP-USB	USB Port
PROTEC RJL	Removable Crystal protection for footrest
PROTEC RP	Removable Crystal protection for footrest
SUPTAB	Fixed I-Pad support with antitheft device
TSR2	Removable IV pole, 2-stage, 2 hooks (FPOSA mandatory)
TSR4	Removable IV pole, 2-stage, 4 hooks (FPOSA mandatory)

ACCESSORIES



4 castors Ø 125 mm with central lock
FPORS1



4 double runing wheel castors Ø 150 mm with central lock
FPORS2




Push bar & Paper roll holder
FPOBP & FPOSR




Long legrest & Electric adjustable footrest (only for long legrest)
FPORJL & FPORP



FPOCP1 & FPOCP2



2 holders for IV pole (necessary for TSR2)
FPOSA



Removable IV pole, 2-stage, 2 hooks (FPOSA mandatory)
TSR2



Comfort pillow (grey)
ACC-CTG




Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOAC3



Reading lamp DEL
FPOLAMP



USB Port
SUP-USB





TIM MANUAL RECLINER

Versatile: Adapts to different hospital services

Easy to use Manually programmed positions

Adaptable Various options available

Customizable Choose your wheelbases and armrests

The TIM is a manual recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The TIM is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.



TIM DIMENSIONS

	LENGTH	WIDTH	HEIGHT	DEPTH
OVERALL	180 cm	72 cm 80 cm (with aluminum armrests)	118 cm	88 cm
SEAT		55 cm	45 cm	50 cm
LEG REST	40 cm	55 cm		
BACKREST	80 cm	58 cm		
HEADREST		55 cm		

PRODUCT WEIGHT : 35 Kg
SAFE WORKING LOAD: 130 Kg



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Manual recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Assisted synchronized backrest and legrest until the comfort position
- The backrest goes independently of the legrest until the recovery position
- Reinforced steel frame
- Rack adjustment headrest
- Seat on coil springs
- Removable comfort armrests
- Front exit
- Retractable footrest under the seat
- Tubular frame :
 - Front: 2 adjustable support pads
 - Back : 2 self lockable casters (Ø50 mm)
- Seat and back cushions in HR 35 Kg density foam, fireproof, and covered




Tubular frame: Front: 2 adjustable pads, Rear: 2 wheels with brakes (Ø50 mm)

"Leg" section under the seat opens fully



Detachable comfortable armrests



Synchronized backrest and leg rest

T I M

Adjustable headrest using a rack and pinion system



Seat cushions and backrest are made from 35 kg density HR foam




M2 flame retardant coating




Faux leather covers: VALENCIA M2


SOOTHING




REFRESHING



ENERGIZING



Composition : Vinyl / Urethane Topcoat
 Thickness : 1,1 mm
 >300.000 Martindale cycles
 Anti-stain finish
 Resistant to sulfide stains
 Effective against bacteria
 Phthalate free
 Anti-static finish



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

TIM



ACCESSORIES

REFERENCE	NAME
FA007	Reading light + USB port
FA021	Fixed aluminum armrests
FA025	Tubular wheelbase with 4 wheels with independent brakes Ø50 mm
FA030	Star wheelbase: 3 wheels with brakes + 2 fixed pads at the front
FA019	Stainless steel base Ø 690 mm
FA031	Cervical roll cylindrical neck pillow
FA034	Legrest crystal protection



ACCESSORIES



Reading light + USB port
FA007



Fixed aluminum armrests
FA021



Cervical roll cylindrical neck
pillow
FA031



Tubular wheelbase with 4
wheels with independent
brakes Ø50mm
FA025



Stainless steel base
Ø 690 mm
FA019



Tubular wheelbase with 2
rear wheels + 2 front pads
FA030



Legrest crystal protection
FA034





VELVET ELECTRIC RECLINER

Versatile: Adapts to different hospital services

Easy to use Electrically programmed positions

Adaptable Various options available

The VELVET is an electric recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The VELVET is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.



VELVET

VELVET DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	100 cm (folded) 185 cm (open)	82 cm (with ALU armrests) 74 cm (with comfort armrests)	125 cm in seated position
SEAT	50 cm	58 cm	50 cm
LEG REST	40 cm	55 cm	
BACKREST	80 cm	58 cm	
HEADREST		55 cm	

PRODUCT WEIGHT : 40 Kg
SAFE WORKING LOAD: 160 Kg



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Electric recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Electric independent backrest and legrest
- Reinforced steel frame, epoxy coated
- Rack adjustment headrest
- Seat on coil springs
- Foamed fireproof M1 vinyl, with silicone dacron
- Armrests with 7 cm HR foam padding and memory foam
- Front exit
- Retractable footrest under the seat
- Fixed height
- Star base with self lockable casters \varnothing 75 mm
- Remote control holder
- Seat and back cushions in HR 35 Kg density foam, fireproof




Star base on wheels \varnothing 75 mm with independent brakes

"Leg" section under the seat opens fully



Aluminum armrests with strap for remote control




Independent electric backrest and leg rest.

VELVET

Adjustable headrest using a rack and pinion system



Seat cushions and backrest are made from 35 kg density HR foam



M2 flame retardant coating



Faux leather covers: VALENCIA M2

SOOTHING




GRAPHITE
107-4003



TAUPE
107-0034



SISAL
107-1048



ROSE
107-2114

REFRESHING



TITAN
107-4042



NIGHT
107-2108



TEAL
107-2119

ENERGIZING



OLIVE
107-5001



NECTARINE
107-6003



SAFRAN
107-6012

Composition : Vinyl / Urethane Topcoat
Thickness : 1,1 mm
>300.000 Martindale cycles
Anti-stain finish
Resistant to sulfide stains
Effective against bacteria
Phthalate free
Anti-static finish



FIRE RATING
M2

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

VELVET



ACCESSORIES

REFERENCE	DESCRIPTION
FA006	Emergency electric recline function on remote control
FA007	Reading light + USB port
FA008VE	Bilateral holder for IV pole
TSR2	Removable IV pole, 2-stage, 2 hooks (FA008VE mandatory)
TSR4	Removable IV pole, 2-stage, 4 hooks (FA008VE mandatory)
FA011	Examination sheet roll holder
FA015B	Retractable aluminum armrests with 2 positions
FA022	DETACHABLE comfortable "full" armrests (pair) incompatible with FA018V
FA009P	Small Shelf : 42 x 32 cm (FA018V-D or FA018V-G necessary)
FA009G	Large Shelf : 81 x 41 cm (FA018V-D or FA018V-G necessary)
FA018V-D	Support for meal tray on the right side of the patient (compatible with the FA018V-G)
FA018V-G	Support for meal tray on the right side of the patient (compatible with the FA018V-D)
FA025	Tubular wheelbase with 4 wheels with independent brakes Ø50 mm
FA033	Tubular wheelbase with 2 rear wheels + 2 front pads
FA019	Stainless steel base Ø 690 mm
FA031	Cervical roll cylindrical neck pillow
FA034	Legrest crystal protection

ACCESSORIES




Reading light + USB port
FA007



Cervical roll cylindrical neck
pillow
FA031



Emergency electric
recline
on remote control
FA006



Examination sheet roll holder
FA011



Retractable aluminum
armrests with 2 positions
FA015B



DETACHABLE comfortable
"full" armrests (pair) incom-
patible with FA018V
FA022



Tubular wheelbase with 4
wheels with independent
brakes Ø50mm
FA025



Tubular wheelbase with 2
rear wheels + 2 front pads
FA033



Stainless steel base
Ø 690 mm
FA019



Small Shelf : 42 x 32 cm
(FA018V necessary)
FA009P



Large Shelf : 81 x 41 cm
(FA018V necessary)
FA009G



Velvet Shelf support
FA018V



Bilateral holder for IV pole
FA008VE



Legrest crystal protection
FA034





SWANN ELECTRIC RECLINER

Versatile Adapts to different hospital services

Easy to use Electrically programmed positions

Adaptable Various options available



The SWANN is an electric healthcare recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The SWANN is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.




SWANN

SWANN DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	100 cm (folded) 185 cm (open)	80 cm	126 cm in seated position
SEAT	50 cm	58 cm	51 to 78 cm
LEG REST	40 cm	55 cm	
BACKREST	80 cm	58 cm	
HEADREST		55 cm	

PRODUCT WEIGHT : 80 Kg
SAFE WORKING LOAD: 160 Kg



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Electric recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Fixed height tubular frame
- 2 pivoting self lockable back casters ø 125 mm
- 2 fixed front casters ø 125 mm
- Emergency electric recovery position on the remote control
- Aluminium armrest
- Rack adjustment headrest
- 2 holders for IV pole
- Electric backrest and legrest
- Push bar
- PVC coated jersey



2 swiveling rear wheels ø 125 mm with independent brakes

2 fixed front wheels ø 125 mm



Aluminum armrests with strap for remote control



Remote-controlled electric emergency recline

SWANN



Adjustable headrest using a rack and pinion system

Electrically operated backrest and leg rest.



Push handle



M2 flame retardant coating



2 holders for telescopic IV poles



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.


Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Faux leather covers: VALENCIA M2


SOOTHING




REFRESHING



ENERGIZING



Composition : Vinyl / Urethane Topcoat
 Thickness : 1,1 mm
 >300.000 Martindale cycles
 Anti-stain finish
 Resistant to sulfide stains
 Effective against bacteria
 Phthalate free
 Anti-static finish





ACCESSORIES

REFERENCE	DESCRIPTION
FA007	Reading light + USB port
FA008	IV pole
FA010	Battery
FA011	Examination sheet roll holder
FA012	Anatomic shaped armrest - Dimensions 400/200/40 mm (FA012S + FA015 or FA015B required)
FA012S	Support for Anatomic shaped armrest - The pair (FA015 or FA015B required)
FA014	Electrically-powered adjustable height
FA014R	Variable height with lift
FA017	Bilateral foot pedal for reclining (Electric height adjustment necessary FA014)
FA015	Retractable armrests (pair) with 3 positions (FA014 obligatory)
FA015B	Retractable aluminum armrests with 2 positions
FA009P	Small Shelf : 42 x 32 cm (FA018S-D or FA018S-G necessary)
FA009G	Large Shelf : 81 x 41 cm (FA018S-D or FA018S-G necessary)
FA018S-D	Support for meal tray on the right side of the patient (compatible with the FA018S-G)
FA018S-G	Support for meal tray on the right side of the patient (compatible with the FA018S-D)
FA026	Centralized rear wheel locking via a pedal Ø 150 mm with double loading rollers
FA029	2 swiveling double loading rollers to the front with brakes Ø 100 mm
FA027	Floor-facing light beams
FA031	Cervical roll cylindrical neck pillow
FA034	Removable legrest crystal protection

ACCESSORIES



Reading light + USB port
FA007



Cervical roll cylindrical neck
pillow
FA031


Battery
FA010



Small Shelf : 42 x 32 cm
(FA018S necessary)
FA009P



Large Shelf : 81 x 41 cm
(FA018S necessary)
FA009G



SwannShelf support
FA018S



Examination sheet roll holder
FA011



Phlebotomy armrest
+ Stand
FA012 - FA012S



Electrically-powered adjust-
able height
FA014



Variable height with lift
FA014R



Bilateral foot pedal
for reclining
FA017



Retractable armrests (pair)
with 3 positions
(FA014 obligatory)
FA015



Retractable aluminum
armrests with 2 positions
FA015B



2 swiveling double load-
ing rollers to the front
with brakes Ø 100 mm
FA029



Removable legrest crystal protection
FA034



Centralized rear wheel lock-
ing via a pedal Ø 150 mm
with double loading rollers
FA026



URBAN TEXTILE M1


CAREXIA RANGE




ABSINTHE




TOMATO




PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

COLOR CHART




CLASSIFICATION


M1

Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 85 % / Coton 15 %
Thickness: 1,2 mm +/- 0,1 mm



Effective
against bacteria
(MRSA)



Solvent-free



Phtalate-free




Martindale
abrasion 100000
cycles



Waterproof

VALENCIA M2


SOOTHING




GRAPHITE
107-4003



TAUPE
107-0034



SISAL
107-1048




ROSE
107-2114

REFRESHING



TITAN
107-4042



NIGHT
107-2108



TEAL
107-2119

ENERGIZING



OLIVE
107-5001



NECTARINE
107-6003



SAFRAN
107-6012

COLOR CHART



FIRE RATING
M2

COMPOSITION : Vinyl / Urethane Topcoat
THICKNESS : 1,1mm

FINISHINGS AND TREATMENTS



Effective
against bacteria
(MRSA)



Antistatic
finish



Phtalates free



>300000
Martindale cycles



Stain and sulfide
resistant finish

A medical professional must always rely on his or her own professional clinical judgment when deciding whether to use a product(s) when treating a particular patient.
A healthcare professional anyway, refers to the package insert, product label and/or instructions for use before using any BMB MEDICAL product.
The products depicted are CE marked in accordance with EU Regulations and Directives.
The products are Class 1 Medical Devices.



Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans - CS 11001

01604 TREVOUX Cedex

FRANCE

www.bmb-medical.com