



HOTEL-STYLE RECLINERS FOR REST & CARE



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units. Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.



CONTENTS

- ▶ CAREXIA FP
- ▶ CAREXIA FPE
- ▶ CAREXIA FPU
- ▶ CAREXIA FPUe
- ▶ TIM
- ▶ VELVET
- ▶ SWANN
- ▶ COLOR CHART

Comfortable

Easy to use

Adaptable Various options available

Versatile Adapted to different units



CAREXIA FP FIXED HEIGHT TREATMENT CHAIR

Carexia FP, is a fixed height versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Its pure and harmonious design is the vector of a positive image for your unit.



CAREXIA FP

CAREXIA FP DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	77 cm (folded) 170 cm (open)	70 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	55 cm	54 cm on seat
INCLINE FOLDER	80 to 60°		

PRODUCT WEIGHT : 65 to 85 Kg depending on the model
SAFE WORKING LOAD: 200 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Versatile treatment chair.

- Adjustable back rest from 80° to 60° with 1 gaz spring
- Fixed height : 54 cm
- Heavy duty design. Steel frame with powder coating. ABS cover
- Integrated headrest
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up
- High resilience foam 35 kg/m3
- Ergonomic shaped backrest and seat - PVC coated jersey
- Range of colors URBAN



4 self lockable casters with double running wheels



ABS basket



Integrated headrest.



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



URBAN

Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : Mass-processed PVC 86%, Cotton 14%.
Thickness : 1,04 mm +/- 10%
Effective against bacteria (MRSA)
Phthalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



FIRE PROOF
M1



4 castors Ø 125 mm
with central lock
FPORS1



4 double running wheel castors
Ø 150 mm with central lock
FPORS2




Push bar & Paper roll
holder
FPOBP & FPOSR



Short legrest synchron-
ized with the backrest.
Trend position.
FPORJ



2 holders for IV pole
(necessary for TSR2)
FPOSA



Removable IV pole,
2-stage, 2 hooks
TSR2



Comfort pillow (grey)
ACC-CTG



Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in height
(2 positions) and in all directions.
Retractable behind the backrest for lateral
transfer of the patient. (pair)
FPOAC3

REFERENCE	DESCRIPTION
FPORS1	4 castors Ø 125 mm with central lock.
FPORS2	4 double running wheel castors Ø 150 mm with central lock.
FPOBP	Push bar
FPORJ	Short legrest synchronized with the backrest. Trend position (Gas cylinders)
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable 2-stage IV pole, 2 hooks (need the option FPOSA)
TSR4	Removable IV pole, 2-stage, 4 hooks (need the option FPOSA)
ACC-CTG	comfort pillow (grey)
FPOAC4	Armrest comfort
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient.
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 2 cm thickness memory foam.
SUPTAB	Fixed Ipad support with antitheft device



Comfortable

Easy to use Electric movements

Adaptable Various options available

Versatile Adapted to different units



CAREXIA FPE FIXED HEIGHT TREATMENT CHAIR

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.


CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.



CAREXIA FPE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	77 cm (folded) 170 cm (open)	70 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	55 cm	54 cm on seat
INCLINE FOLDER	26 to 83°		

PRODUCT WEIGHT : 65 to 85 Kg depending on the model
SAFE WORKING LOAD: 200 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Multi-purpose treatment chair with electric movements.

- Fixed height : 54 cm
- Heavy duty design - Steel frame with powder coating - ABS cover
- Integrated headrest
- Electric backrest and legrest (Electric cylinders)
- Trend position
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up
- High resilience foam 35 kg/m3
- Ergonomic shaped backrest and seat - PVC coated jersey
- Range of colors URBAN



4 self lockable casters with double running wheels



Polyurethan armrests, can be folded up



ABS basket



Integrated headrest


Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUSCHIA




PLUM



DEEP OCEAN



PETROL



GREY

URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : Mass-processed PVC 86%, Cotton 14%.
Thickness : 1,04 mm +/- 10%
Effective against bacteria (MRSA)
Phthalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



FIRE PROOF
M1

CAREXIA FPE



ACCESSORIES

REFERENCE	DESCRIPTION
FPORS1	4 castors Ø 125 mm with central lock.
FPORS2	4 double runing wheel castors Ø 150 mm with central lock.
FPOBP	Push bar
FPORJL	Long legrest (Electric cylinders)
FPORP	Electric adjustable footrest (only for long legrest) (incompatible with FPOAC4)
FPOCP1	Foot pedal control for relax and trend position
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable 2-stage IV pole, 2 hooks (need the option FPOSA)
TSR4	Removable IV pole, 2-stage, 4 hooks (need the option FPOSA)
ACC-CTG	comfort pillow (grey)
FPOAC4	Armrest comfort
FPOAC3	Anatomic shaped armrest adjustable in height and in all directions. Retractable behind the backrest for lateral transfer of the patient.
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 2 cm thickness memory foam.
FPOLAMP	Reading lamp DEL.
FPOBA	Battery 2 Ah
SUP-USB	USB port
PROTEC RJL	Removable Crystal protection for long legrest.
PROTEC RP	Removable Crystal protection for footrest
SUPTAB	Fixed lpad support (on right or left) with antitheft device

ACCESSORIES



4 castors Ø 125 mm
with central lock
FPORS1



4 double running wheel castors
Ø 150 mm with central lock
FPORS2




Push bar & Paper roll
holder
FPOBP & FPOSR




Long legrest & Electric
adjustable footrest
(only for long legrest)
FPORJL & FPORP



Foot pedal control for relax
and trend position
FPOCP1



2 holders for IV pole
(necessary for TSR2)
FPOSA



Removable IV pole,
2-stage, 2 hooks
(FPOSA mandatory)



Comfort pillow (grey)
ACC-CTG




Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in
height (2 positions) and in all
directions. Retractable behind the backrest
for lateral transfer of the patient. (pair)
FPOAC3



Reading lamp DEL
FPOLAMP



USB Port
SUP-USB





CAREXIA FPU VARIABLE-HEIGHT CARE CHAIR

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.



Comfortable

Adaptable Various options available

Easy to use variable electric height


Versatile Adapted to different units

CAREXIA FPU

DIMENSIONS CAREXIA FPV

	LENGTH	WIDTH	HEIGHT
OVERALL	77 cm (folded) 170 cm (open)	70 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	55 cm	54 cm on seat
INCLINE FOLDER	70 to 83°		
VARIABLE HEIGHT			54 to 78 cm

PRODUCT WEIGHT : 65 to 85 Kg depending on the model
SAFE WORKING LOAD: 200 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Height variable electrical care chair

- Electric variable height from 54 to 78 cm (Electric cylinders)
- Adjustable back rest from 70° to 83° with 1 gaz spring
- Heavy duty design - Steel frame with powder coating - ABS cover
- Integrated headrest
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up
- High resilience foam 35 kg/m3
- Ergonomic shaped backrest and seat - PVC coated jersey
- Range of colors URBAN



4 self lockable casters with double running wheels



ABS basket



Integrated headrest



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE




URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : Mass-processed PVC 86%, Cotton 14%.
Thickness : 1,04 mm +/- 10%
Effective against bacteria (MRSA)
Phthalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



FIRE PROOF
M1



4 castors Ø 125 mm
with central lock
FPORS1



4 double running wheel castors
Ø 150 mm with central lock
FPORS2




Push bar & Paper roll
holder
FPOBP & FPOSR



Short legrest synchro-
nized with the backrest.
Trend position.
FPORJ



2 holders for IV pole
(necessary for TSR2)
FPOSA



Removable IV pole,
2-stage, 2 hooks
(FPOSA mandatory)
TSR2



Comfort pillow (grey)
ACC-CTG



Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in
height (2 positions) and in all
directions. Retractable behind the
backrest for lateral transfer of the pa-
tient. (pair)
FPOAC3



USB Port
SUP-USB

RÉFÉRENCE	DÉSIGNATION
FPORS1	4 castors Ø 125 mm with central lock.
FPORS2	4 double running wheel castors Ø 150 mm with central lock.
FPOBP	Push bar
FPORJ	Short legrest synchronized with the backrest - Trend position (Gas cylinders)
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable 2-stage IV pole, 2 hooks (need the option FPOSA)
TSR4	Removable IV pole, 2-stage, 4 hooks (need the option FPOSA)
ACC-CTG	comfort pillow (grey)
FPOAC4	Armrest comfort
FPOAC3	Anatomic shaped armrest adjustable in height and in all directions. Retractable behind the backrest for lateral transfer of the patient.
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 2 cm thickness memory foam
FPOBA	Battery 2 Ah
SUP-USB	USB port
SUPTAB	Fixed Ipad support with antitheft device



Comfortable

All-electric care chair

Adaptable Various options available

Easy to use Variable height & electric movements

Versatile Adapted to different units



ALL ELECTRIC CARE CHAIR CAREXIA FPUE

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.




CAREXIA FPVE

CAREXIA FPVE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	77 cm (folded) 177 cm (open) + long leg rest	70 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	55 cm	54 cm on seat
INCLINE FOLDER	26 to 83°		
VARIABLE HEIGHT			54 to 78 cm

PRODUCT WEIGHT : 65 to 85 Kg depending on the model
SAFE WORKING LOAD: 200 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Height variable electrical care chair


- Electric variable height from 54 to 78 cm (Electric cylinders)
- Heavy duty design - Steel frame with powder coating - ABS cover
- Integrated headrest
- Electric backrest and legrest (Electric cylinders)
- Trend position
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up
- High resilience foam 35 kg/m3
- Ergonomic shaped backrest and seat - PVC coated jersey
- Range of colors URBAN



4 self lockable casters with double running wheels



ABS basket



Integrated headrest



Trend position



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



URBAN
Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : Mass-processed PVC 86%, Cotton 14%.
Thickness : 1,04 mm +/- 10%
Effective against bacteria (MRSA)
Phtalate-free
Solvent-free
Martindale abrasion 100 000 cycles
Waterproof



FIRE PROOF
M1

CAREXIA FPUE

ACCESSORIES

RÉFÉRENCE	DÉSIGNATION
FPORS1	4 castors Ø 125 mm with central lock.
FPORS2	4 double runing wheel castors Ø 150 mm with central lock.
FPOBP	Push bar
FPORJL	Long legrest (Electric cylinders)
FPORP	Electric adjustable footrest (only for long legrest) (incompatible with FPOAC4)
FPOCP1	Foot pedal control for relax and trend position on both sides (Electric cylinders)
FPOCP2	Foot pedal control for variable height on both sides. (Electric cylinders)
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable 2-stage IV pole, 2 hooks (need the option FPOSA)
TSR4	Removable IV pole, 2-stage, 4 hooks (need the option FPOSA)
ACC-CTG	comfort pillow (grey)
FPOAC4	Armrest comfort
FPOAC3	Anatomic shaped armrest adjustable in height and in all directions. Retractable behind the backrest for lateral transfer of the patient.
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 2 cm thickness memory foam.
FPOLAMP	Reading lamp DEL.
FPOBA	Battery 2 Ah
SUP-USB	USB port
PROTEC RJL	Removable Crystal protection for long legrest.
PROTEC RP	Removable Crystal protection for footrest
SUPTAB	Fixed Ipad support with antitheft device

ACCESSORIES



4 castors Ø 125 mm
with central lock
FPORS1



4 double runing wheel castors
Ø 150 mm with central lock
FPORS2



Push bar & Paper roll
holder
FPOBP & FPOSR



Anatomic shaped armrest adjustable in
height (2 positions) and in all
directions. Retractable behind the backrest
for lateral transfer of the patient. (pair)
FPOAC3



Reading lamp DEL
FPOLAMP




USB Port
SUP-USB




Long legrest & Electric
adjustable footrest
(only for long legrest)
FPORJL & FPORP



FPOCP1 & FPOCP2



2 holders for IV pole
(necessary for TSR2)
FPOSA



Removable IV pole,
2-stage, 2 hooks
(FPOSA mandatory)
TSR2



Comfort pillow (grey)
ACC-CTG



Armrest comfort (pair)
FPOAC4





TIM MANUAL RECLINER

Versatile: Adapts to different hospital services

Easy to use Manually programmed positions

Adaptable Various options available

Customizable Choose your wheelbases and armrests

The TIM is a manual recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The TIM is a versatile recliner with all of the functions required for procedures with medical equipment.


The recliner has been designed for routine use, as well as to simplify regular maintenance.



TIM DIMENSIONS

	LENGTH	WIDTH	HEIGHT	DEPTH
OVERALL	85 cm (chair) 153 cm(extended)	72 cm 80 cm (with aluminum armrests)	115 cm	
SEAT		57 cm	50 cm	55 cm
LEG REST	44 cm	57cm		
BACKREST	76 cm	55 cm		
HEADREST	29 cm	55 cm		16 cm

PRODUCT WEIGHT : 35 Kg
SAFE WORKING LOAD: 130 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Manual recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Assisted synchronized backrest and legrest until the confort position (Gas cylinders)
- The backrest goes independently of the legrest until the recovery position
- Reinforced steel frame
- Rack adjustment headrest
- Seat on coil springs
- Removable comfort armrests
- Front exit
- Retractable footrest under the seat
- Tubular frame :
 - Front: 2 adjustable support pads
 - Back : 2 casters (Ø50 mm)
- Seat and back cushions in HR 35 Kg density foam, fireproof, and covered
- M2 fire-resistant coating Valencia, 10 colours



Tubular frame:
Front: 2 adjustable pads
Rear: 2 wheels with brakes (Ø50 mm)

"Leg" section under the seat
opens fully



Synchronized backrest and leg rest

Detachable comfortable armrests



Adjustable headrest using a rack and pinion system



Seat cushions and backrest are made from 35 kg density HR foam
M2 flame retardant coating

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.



Faux leather covers: VALENCIA M2

SOOTHING



GRAPHITE
107-4003



TAUPE
107-0034



SISAL
107-1048



ROSE
107-2114

REFRESHING



TITAN
107-4042



NIGHT
107-2108



TEAL
107-2119

ENERGIZING



OLIVE
107-5001



NECTARINE
107-6003



SAFRAN
107-6012

Composition : Polyester knitted fabric with fireproof PVC coating on one side
Thickness : 1,1 mm
>300.000 Martindale cycles
Anti-stain finish
Resistant to sulfide stains
Effective against bacteria
Phthalate free
Anti-static finish



TIM

ACCESSORIES

REFERENCE	NAME
FA007	Reading lamp + USB Port
FA021	Fixed aluminium armrests
FA025	4 self lockable casters (Ø50 mm)
FA030	Inox base Ø 690mm
FA019	Star base : 3 back casters + 2 front fixed pads
FA031	Neck cylindrical cushion
FA034	Legrest crystal protection



ACCESSORIES



Reading light + USB port
FA007



Fixed aluminum armrests
FA021



Cervical roll cylindrical
neck pillow
FA031




Tubular wheelbase with 4
wheels with independent
brakes Ø50mm



Stainless steel base
Ø 690 mm
FA019



Tubular wheel base
with 2 rear wheels + 2
front pads
FA030



Legrest crystal protection
FA034





VELVET ELECTRIC RECLINER

The VELVET is an electric recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The VELVET is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.



Versatile: Adapts to different hospital services

Easy to use Electrically programmed positions

Adaptable Various options available



VELVET

VELVET DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	100 cm (folded) 177 cm (open) 185 cm (Trendelenburg)	82 cm (with ALU armrests) 74 cm (with comfort armrests)	125 cm in seated position
SEAT	53 cm	57 cm	51 cm
LEG REST	40 cm	57 cm	
BACKREST	82 cm	57 cm	
HEADREST	32 cm	55 cm	

PRODUCT WEIGHT : 40 Kg
SAFE WORKING LOAD: 160 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Electric recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Electric independant backrest and legrest (Electric cylinders)
- Reinforced steel frame, epoxy coated
- Rack adjustment headrest
- Seat on coil springs
- Foamed fireproof M1 vinyl, with silicone dacron
- Armrests with 7 cm HR foam padding and memory foam
- Front exit
- Retractable footrest under the seat
- Fixed height
- Star base with self lockable casters ø 75 mm
- Remote control holder
- Seat and back cushions in HR 35 Kg density foam, fireproof
- M2 fire-resistant coating Valencia, 10 colours




Star base on wheels ø 75 mm with independent brakes

"Leg" section under the seat opens fully



Independent electric backrest and leg rest.

Aluminum armrests with strap for remote control



VELVET

Adjustable headrest using a rack and pinion system



Seat cushions and backrest are made from 35 kg density HR foam




M2 flame retardant coating



Faux leather covers: VALENCIA M2

SOOTHING



GRAPHITE
107-4003



TAUPE
107-0034



SISAL
107-1048



ROSE
107-2114

REFRESHING



TITAN
107-4042



NIGHT
107-2108



TEAL
107-2119

ENERGIZING



OLIVE
107-5001



NECTARINE
107-6003



SAFRAN
107-6012

Composition : Polyester knitted fabric with fireproof PVC coating on one side
Thickness : 1,1 mm
>300.000 Martindale cycles
Anti-stain finish
Resistant to sulfide stains
Effective against bacteria
Phthalate free
Anti-static finish



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

RÉFÉRENCE	DÉSIGNATION
FA006	Emergency electric recovery position on remote control
FA007	Reading lamp + USB Port
FA0088V	Bilateral holder for IV Pole (incompatible FA022)
FA011	Paper roll dispenser
FA015B	Removable aluminium armrests 2 positions
FA022	Removable comfort armrests - incompatible with FA018V-G/D
FA031	Neck cylindrical cushion
FA009P	Small Shelf : 42 x 32 cm (FA018V-D or FA018V-G necessary)
FA009G	Large Shelf : 81 x 41 cm (FA018V-D or FA018V-G necessary)
FA018V-D	Support for meal tray on the right side of the patient (compatible with the FA018V-G)
FA018V-G	Support for meal tray on the right side of the patient (compatible with the FA018V-D)
FA025	4 self lockable casters (Ø50 mm)
FA033	2 back casters + 2 front fixed pads
FA019	Inox base Ø 690mm
FA034	Legrest crystal protection
TSR2	Removable IV pole, 2-stage, 2 hooks (FA0088VE mandatory)
TSR4	Removable IV pole, 2-stage, 4 hooks (FA0088VE mandatory)

ACCESSORIES



Reading light + USB port
FA007



Cervical roll cylindrical
neck pillow
FA031



Emergency electric
recline
on remote control
FA006



Examination sheet roll
holder
FA011



Retractable aluminum
armrests with 2 posi-
tions FA015B



DETACHABLE comfortable
"full" armrests (pair) incompat-
ible with FA018V
FA022



Tubular wheelbase with 4
wheels with independent
brakes Ø50mm
FA025



Tubular wheelbase with 2
rear wheels + 2 front pads
FA033



Stainless steel base
Ø 690 mm
FA019



Small Shelf : 42 x 32 cm
(FA018V necessary)
FA009P



Large Shelf : 81 x 41 cm
(FA018V necessary)
FA009G



Velvet Shelf support
FA018V



Bilateral holder for IV pole
FA008VE



Legrest crystal protection
FA034





SWANN ELECTRIC RECLINER

The SWANN is an electric healthcare recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The SWANN is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.



Versatile Adapts to different hospital services

Easy to use Electrically programmed positions

Adaptable Various options available



SWANN DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	104 cm (folded) 183 cm (open)	80 cm	131 cm in seated position
SEAT	53 cm	57 cm	54 to 84 cm
LEG REST	40 cm	57 cm	
BACKREST	86 cm	55 cm	
HEADREST	32 cm	55 cm	

PRODUCT WEIGHT : 80 Kg
SAFE WORKING LOAD: 160 Kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Electric recliner with steel frame. All cushions are retractable for cleaning or repairs.


- Fixed height tubular frame
- 2 pivoting self lockable back casters ø 125 mm
- 2 fixed front casters ø 125 mm
- Emergency electric recovery position on the remote control
- Aluminium armrest
- Rack adjustment headrest
- 2 holders for IV pole
- Electric backrest and legrest (Electric cylinders)
- Push bar
- PVC coated jersey
- M2 fire-resistant coating Valencia, 10 colours



2 swiveling rear wheels Ø 125 mm with independent brakes



2 fixed front wheels Ø 125 mm



Remote-controlled electric emergency recline



Aluminum armrests with strap for remote control



Adjustable headrest using a rack and pinion system


Electrically operated backrest and leg rest.



Push handle



M2 flame retardant coating



2 holders for telescopic IV poles




Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.



Faux leather covers : VALENCIA M2

SOOTHING




GRAPHITE
107-4003



TAUPE
107-0034



SISAL
107-1048



ROSE
107-2114

REFRESHING



TITAN
107-4042




NIGHT
107-2108



TEAL
107-2119

ENERGIZING



OLIVE
107-5001



NECTARINE
107-6003



SAFRAN
107-6012

Composition : Polyester knitted fabric with fireproof PVC coating on one side
Thickness : 1,1 mm
>300.000 Martindale cycles
Anti-stain finish
Resistant to sulfide stains
Effective against bacteria
Phthalate free
Anti-static finish

FIRE RATING
M2





ACCESSORIES

RÉFÉRENCE	DÉSIGNATION
FA007	Reading lamp + USB Port
FA008	IV Pole
FA010	Battery
FA011	Paper roll dispenser
FA012	Phlebotomy armrest - Dimensions 400/200/40 mm (FA012S + FA015 or FA015B required)
FA012S	Support for Phlebotomy armrest- The pair (FA015 or FA015B required)
FA014	Electric height adjustment (Electric cylinders)
FA014R	Variable height with lift (Electric cylinders)
FA017	Two-sided foot pedals for Trendelenburg adjustment (Electric height adjustment necessary FA014)
FA015	Retractable armrests 3 positions (FA014 necessary)
FA015B	Retractable aluminium armrests 2 positions
FA009P	Small Shelf : 42 x 32 cm (FA018S-D or FA018S-G necessary)
FA009G	Large Shelf : 81 x 41 cm (FA018S-D or FA018S-G necessary)
FA018S-D	Support for meal tray on the right side of the patient (compatible with the FA018S-G)
FA018S-G	Support for meal tray on the right side of the patient (compatible with the FA018S-D)
FA026	Back double running castors Ø 150 mm with central lock.
FA029	2 pivoting self lockable front casters ø 100 mm
FA027	Beam of light
FA031	Neck cylindrical cushion
FA034	Removable legrest crystal protection

ACCESSORIES



Reading light + USB port
FA007




Cervical roll cylindrical
neck pillow
FA031



Examination sheet roll
holder
FA011



Phlebotomy armrest
+ Stand
FA012 - FA012S



Electrically-powered
adjustable height
FA014




Retractable aluminum
armrests with 2 posi-
tions FA015B




Bilateral foot pedal
for reclining
FA017




Retractable armrests
(pair) with 3 positions
(FA014 obligatory)
FA015



Centralized rear wheel lock-
ing via a pedal Ø 150 mm
with double loading rollers
FA026



2 swiveling double load-
ing rollers to the front
with brakes Ø 100 mm
FA029



Removable legrest crystal protection
FA034



Small Shelf : 42 x 32 cm
(FA018S necessary)
FA009P



Large Shelf : 81 x 41 cm
(FA018S necessary)
FA009G



SwannShelf support
FA018S



Variable height with lift
FA014R



URBAN TEXTILE M1


CAREXIA RANGE




ABSINTHE




TOMATO




PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

COLOR CHART



CLASSIFICATION


M1

Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : Mass-processed PVC 86%, Cotton 14%.
Thickness: 1,04 mm +/- 10%



Effective
against bacteria
(MRSA)



Solvent-free



Phthalate-free




Martindale
abrasion 100000
cycles



Waterproof

VALENCIA M2


SOOTHING




GRAPHITE
107-4003



TAUPE
107-0034



SISAL
107-1048




ROSE
107-2114


REFRESHING



TITAN
107-4042



NIGHT
107-2108




TEAL
107-2119

ENERGIZING



OLIVE
107-5001



NECTARINE
107-6003



SAFRAN
107-6012

COLOR CHART



FIRE RATING
M2

COMPOSITION : Polyester knitted fabric with fireproof PVC coating on one side
THICKNESS: 1,1mm

FINISHINGS AND TREATMENTS



Effective
against bacteria
(MRSA)



Antistatic
finish



Phtalates free



>300000
Martindale cycles



Stain and sulfide
resistant finish

A medical professional must always rely on his or her own professional clinical judgment when deciding whether to use a product(s) when treating a particular patient.

A healthcare professional anyway, refers to the package insert, product label and/or instructions for use before using any BMB MEDICAL product.

The products depicted are CE marked in accordance with EU Regulations and Directives.

The products are Class 1 Medical Devices.



Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans - CS 11001
01604 TREVOUX Cedex
FRANCE

www.bmb-medical.com