

TRANSFER CHAIR H-MOVE





CONTENTS

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position.

Manageable, robust and attractive it will bring for a long time a touch of modernity and color to your department.

A wide range of accessories allows it to adapt to each use.



I - M O U E



I - M O U E O U T D O O R



I - M O U E B A R I A T R I C



I - M O U E E Z - G O



I - M O U E E Z - G O B A R I A T R I C





TRANSFER CHAIR I-MOVE

I-MOVE is a transfer chair, it allows the reception of the patient and his transport in sitting position.

Sturdy

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit


Easy to handle & Stackable for optimal storage



I-MOVE DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT
I-MOVE	105 cm	75 cm	100 cm	H. 55 cm D. 46 cm W. 46 cm

PRODUCT WEIGHT : 35 Kg
SAFE WORKING LOAD : 200 Kg



Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

I-MOVE transfer chair for patient welcoming

- Overall dimensions : 105 (length) x 75 (width) x 100 (height) cm
- Dimensions seat and back: 46 x 46 cm
- Seat height : 55 cm
- Steel frame epoxy coated, white RAL 9003
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Anti-shock protection
- Cushion in foam - Thickness 5 cm
- Seamless finish, vinyl coating



Back wheels Ø 300 mm with central braking system



Polyurethan armrests, can be folded up



Patient file holder



Sliding foot rest
Front swivel wheels Ø 125 mm

Different colors



FUSCHIA

GREEN


BLUE

Composition : Opaque flexible PVC film, fireproofed in the mass
Thickness: 0,5 mm +/-0,1mm
Light resistance: >= 6



FIRE RATING
M1

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR2-F	Fixed IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
CLORA	Anti-static system
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with IV pole TSR2)
TROPO	Accessory basket
TROSO2	Holder for oxygen cylinder
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO - EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
PSU	Hooks to hold bags




Removable IV pole,
2-stage, 2 hooks
TSR2



Accessory basket
TROPO



Holder for oxygen
cylinder
TROSO2



Adjustable leg rest
(1 piece)
TRORJ



Coin acceptor
TROMO



Hooks to hold bag
PSU





TRANSFER CHAIR I-MOVE OUTDOOR

Sturdy It is designed to withstand weather conditions

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Easy to handle & Stackable For optimal storage

I-MOVE is a transfer chair, it allows the reception of the patient and his transport in sitting position, outdoor.

Its coating is weather resistant, so it can be stored outside.




I - M O V E O U T D O O R

I-MOVE OUTDOOR DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT
I-MOVE	105 cm	75 cm	100 cm	H. 55 cm D. 46 cm W. 46 cm

PRODUCT WEIGHT : 35 Kg
SAFE WORKING LOAD : 200 Kg



Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position indoor and outdoor.

I-MOVE OUTDOOR transfer chair for patient welcoming

- Steel frame epoxy coated, white RAL 9003
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Anti-shock protection
- Polypropylene plastic fitting
- Color : Blue



Back wheels Ø 300 mm with central braking system



Polyurethan armrests, can be folded up



Patient file holder




Front swivel wheels Ø 125 mm
Sliding and self-blocking foot rest

Composition : Polypropylene



I-MOVE OUTDOOR

ACCESSORIES



Removable IV pole,
2-stage, 2 hooks
TSR2



U-shape anti-theft, fitted with 2 hooks
infusion holder
TROAA



Accessory basket
TROPO



Holder for oxygen cylinder
TROS02



Wall fixation
TROMOP



Hooks to hold bags
PSU

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR2-F	Fixed IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
CLORA	Anti-static system
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with TSR2)
TROPO	Accessory basket
TROS02	Holder for oxygen cylinder
TROMO - EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
PSU	Hooks to hold bags

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit



TRANSFER CHAIR I-MOVE BARIATRIC

I-MOVE is a bariatric transport chair for your patients' reception and sitting transport.




I-MOVE BARIATRIC

I-MOVE BARIATRIC DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT
I-MOVE	105 cm	90 cm	102 cm	H. 55 cm D. 46 cm W. 65 cm

PRODUCT WEIGHT : 50Kg
SAFE WORKING LOAD : 250 Kg



Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position.

I-MOVE transfer chair for patient welcoming

- Steel frame epoxy coated, white RAL 9003
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Anti-shock protection
- Cushion in foam - Thickness 8 cm
- PVC coated Jersey



Back wheels Ø 300 mm with central braking system



Front swivel wheels Ø 125 mm
Sliding foot rest




Patient file holder



Polyurethan armrests, can be folded up

Differents colors




Composition : Fireproof fabric jersey cotton (15%) coated with PVC (80%)
Thickness: 1,4 mm +/-0,1mm
Light resistance: >= 6

I-MOVE BARIATRIC

ACCESSORIES

REFERENCE	DESCRIPTION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR2-F	Fixed IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
CLORA	Anti-static system
TROPO	Accessory basket
TROSO2	Holder for oxygen cylinder
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO - BP	Coin acceptor for 1€ coin or equivalent
TROMO - BJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
PSU	Hooks to hold bags



IV pole, 2-stage, 2 hooks
TSR2



Accessory basket
TROPO



Holder for oxygen cylinder
TROSO2



Adjustable leg rest (1 piece)
TRORJ



Coin acceptor
TROMO



Hooks to hold bag
PSU



TRANSFER CHAIR I-MOVE EZ-GO

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Motorized Very easy to use, progressive speed



I-MOVE EZ-GO is a transport chair for your patients' reception and sitting transport. It is equipped with an autonomous motorized wheel, two batteries and an external charger.

Thanks to I-MOVE EZ-GO, patients transfer is child's play!




I - M O V E E Z - G O

I-MOVE EZ-GO DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT
I-MOVE	105 cm	70 cm	102 cm	H. 55 cm D. 46 cm W. 46 cm

PRODUCT WEIGHT : 52 Kg
SAFE WORKING LOAD : 200 Kg



Motorized wheel including battery

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

Motorized transfer chair I-MOVE EZ-GO for patient welcoming

- Motorized drive system with autonomous wheel
- Switchable battery - 5 km autonomy per battery - Charge time : 5 hours
- Lithium battery without memory effect
- Forward/Backward twist grip with progressive speed acceleration
- Steel frame epoxy coated
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Accessory basket
- Anti-shock protection
- Supplied with 2 batteries and external charger
- Cushion in foam - Thickness 5 cm
- Seamless finish, vinyl coating



Accessory basket



Back wheels Ø 300 mm with central braking system



Polyurethan armrests, can be folded up




Sliding foot rest
Front swivel wheels Ø 125 mm

Different colors



Composition : Opaque flexible PVC film, fireproofed in the mass
Thickness: 0,5 mm +/-0,1mm
Light resistance: >= 6





Removable IV pole,
2-stage, 2 hook
TSR2



U-shape anti-theft, fitted with 2 hooks
infusion holder
TROAA



Holder for oxygene
cylinder
TROS02



Hooks to hold bags
PSU



Adjustable leg rest (1
piece)
TRORJ



Wall fixation
TROMOP



Batterie supplémentaire
EZ0BA



External wall charger
EZ0CH



Charging station 2/3/4 chargers
(one multi-socket included)
EZ0S2 - EZ0S3 - EZ0S4



Charging station 1 charger
EZ0S1



External charger mobile version
EZ0PO

REFERENCE	DESIGNATION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with TSR2)
TROS02	Holder for oxygen cylinder
PSU	Hooks to hold bags
CLORA	Anti-static system
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO - EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
EZSBA	Without battery (depreciation)
EZOBA	Additional battery
EZOCH	External wall charger
EZOPO	External charger mobile version
EZOS1	Charging station 1 charger
EZOS2	Charging station 2 chargers (one multi-socket included)
EZOS3	Charging station with 3 chargers (one multi-socket included)
EZOS4	Charging station with 4 chargers (one multi-socket included)



TRANSFER CHAIR I-MOVE EZ-GO BARIATRIC

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Motorized Very easy to use, progressive speed



I-MOVE EZ-GO is a transport chair for your patients' reception and sitting transport. It is equipped with an autonomous motorized wheel, two batteries and an external charger.


Thanks to I-MOVE EZ-GO, patients transfer is child's play ! It is particularly suitable for the transfer of overweight patients.

I-MOVE EZ-GO BARIATRIC

I-MOVE EZ-GO BARIATRIC DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT
I-MOVE	105 cm	90 cm	102 cm	H. 55 cm D. 46 cm W. 65 cm

PRODUCT WEIGHT : 65 Kg
SAFE WORKING LOAD: 250 Kg



Motorized wheel including battery

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

Motorized transfer chair I-MOVE EZ-GO for patient welcoming

- Motorized drive system with autonomous wheel
- Switchable battery - 5 km autonomy per battery - Charge time : 5 hours
- Lithium battery without memory effect.
- Forward/Backward twist grip with progressive speed acceleration
- Steel frame epoxy coated RAL 9010
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Accessory basket
- Supplied with 2 batteries and external charger
- Anti-shock protection
- Cushion in foam HR 35kg/m³ - Thickness 8 cm
- Seamless finish, vinyl coating
- Seamless finish, vinyl coating



Accessory basket



Polyurethan armrests, can be folded up



Back wheels Ø 300 mm with central braking system



Front swivel wheels Ø 125 mm
Sliding footrest

Different colors



ORANGE



GREEN



LAVANDER

Composition : Fireproof fabric jersey cotton (15%) coated with PVC (80%)
Thickness: 1,4 mm +/-0,1mm
Light resistance: >= 6





Removable IV pole,
2-stage, 2 hook
TSR2




U-shape anti-theft, fitted with 2 hooks
infusion holder
TROAA



Holder for oxygen cylinder
TROSO2



Hooks to hold bags
PSU



Adjustable leg rest
(1 piece)
TRORJ



Wall fixation
TROMOP



Batterie supplémentaire
EZOBA



External wall charger
EZOCH



Charging station 2/3/4 chargers
(one multi-socket included)
EZOS2 - EZOS3 - EZOS4



Charging station 1 charger
EZOS1



External charger mobile version
EZOPO

REFERENCE	DESIGNATION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with TSR2)
TROSO2	Holder for oxygen cylinder
PSU	Hooks to hold bags
CLORA	Anti-static system
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO - EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
EZSBA	Without battery (depreciation)
EZOBA	Additional battery
EZOCH	External wall charger
EZOPO	External charger mobile version
EZOS1	Charging station 1 charger
EZOS2	Charging station 2 chargers (one multi-socket included)
EZOS3	Charging station with 3 chargers (one multi-socket included)
EZOS4	Charging station with 4 chargers (one multi-socket included)

A medical professional must always rely on his or her own professional clinical judgment when deciding whether to use a product(s) when treating a particular patient.

A healthcare professional anyway, refers to the package insert, product label and/or instructions for use before using any BMB MEDICAL product.

The products depicted are CE marked in accordance with EU Regulations and Directives.

The products are Class 1 Medical Devices.



Tél : +33 (0)4 74 08 71 71

Fax : +33 (0)4 74 08 71 72

info@bmb-medical.com

100 Av. du Formans - CS 11001

01604 TREVOUX Cedex

FRANCE

www.bmb-medical.com